

US ENVIRONMENTAL PROTECTION AGENCY  
OFFICE OF PESTICIDES PROGRAMS  
REGISTRATION DIVISION (FS-767)  
WASHINGTON, DC 20460

EPA REGISTRATION NO.

46851-9

DATE OF ISSUANCE

AUG 2 1995

TERM OF ISSUANCE

NAME OF PESTICIDE PRODUCT

Omnicide Plus

NOTICE OF PESTICIDE:  REGISTRATION  
 REREGISTRATION

(Under the Federal Insecticide, Fungicide,  
and Rodenticide Act, as amended)

NAME AND ADDRESS OF REGISTRANT (Include ZIP code)

ProChem Company  
A Wholly Owned Subsidiary of Cottrell, Ltd.  
7399 South Tucson Way  
Englewood, CO 80112

NOTE: Changes in labeling formula differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above U.S. EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby Registered/Reregistered under the Federal Insecticide, Fungicide, and Rodenticide Act.

A copy of the labeling accepted in connection with this Registration/Reregistration is returned herewith.

Registration is in no way to be construed as an indorsement or approval of this product by this Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA sec. 3(c)(7)(A) provided that you:

1. Submit/cite all data required for registration/reregistration of your product under FIFRA section 3(c)(5) when the Agency requires all registrants of similar products to submit such data; and submit acceptable responses required for reregistration of your product under FIFRA section 4.

2. Make the labeling changes listed below before you release the product for shipment.

3. On the Container Labeling:


a. Add the phrase "EPA Registration No. 46851-9".

b. The "Virucidal" claim must be keyed by a symbol to the paragraph listing the specific tested viruses.

c. Specify the major areas in which the product is recommended for use (e.g. homes, school, hospitals).

ATTACHMENT IS APPLICABLE

SIGNATURE OF APPROVING OFFICIAL



DATE

8/2/95

1081

26, 22  
0

- d. Delete the term "ProCide Plus" wherever it may appear on the label. This is not your product name.
- 4. On the package insert labeling, change the product name "ProCide Plus" to "Omnicide Plus".
- 5. You are reminded that the container labeling must meet EPA's current labeling requirements.

A stamped copy of the labeling is enclosed for your records. Submit one (1) copy of the final printed label prior to release of the product for shipment.

If you have any questions concerning this letter, please contact Martha Terry at 703-305-6982.

Sincerely,



Marion J. Johnson, Jr.  
Product Manager (31)  
Antimicrobial Program Branch  
Registration Division (7505C)

Enclosure

CONCURRENCES

SYMBOL							
SURNAME							
DATE							



4 57 22

**DIRECTIONS FOR USE**

IT IS VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING

**Activation:** Activate the ProCide plus solution by adding the entire contents of the Activator Bottle (which is attached to the ProCide plus solution container) to the container. Place cap on container and shake well. On activation the solution immediately changes color to blue, thereby indicating solution is ready to use. Record the date of activation on the indicated space below, in a log book or a label affixed to any secondary container used for activated solution. See package insert for additional instructions and information regarding activated solution.

**Cleaning/Decontamination:** Blood and other body fluids must be thoroughly cleaned from the surfaces, lumens, and objects before application of the disinfectant or sterilant. Blood and other body fluids should be autoclaved and disposed of according to all applicable Federal, State and local regulations for infectious waste disposal.

For complete disinfection or sterilization of medical instruments and equipment, thoroughly clean, rinse and rough dry objects before immersing in ProCide plus solution. Cleanse and rinse the lumens of hollow instruments before filling with ProCide plus solution. See package insert for additional cleaning/decontamination instructions.

**Sterilization:** Immerse medical equipment/devices completely in ProCide plus Solution for a minimum of ten hours at 68°F (20°C) to eliminate all microorganisms,

including *Clostridium sporogenes* and *Bacillus subtilis* spores. Remove equipment from the solution using sterile technique and rinse thoroughly with sterile water. See package insert for complete instructions/information on sterilization.

**High Level Disinfection:** Immerse medical equipment/devices completely in ProCide plus solution for a minimum of 45 minutes at 68°F (20°C) to destroy all pathogenic microorganisms, except for large numbers of bacterial endospores, but including *Mycobacterium bovis* (Quantitative TB Method).

Remove devices and equipment from the solution and rinse thoroughly with sterile water or potable water. The quality of rinse water used is dependent upon the intended use of the instrument/equipment. See package insert for complete instructions/information on high level disinfection.

**Intermediate Level Disinfection:** Immerse medical instruments/equipment completely in ProCide plus for a minimum of 10 minutes at 68°F (20°C) to destroy all vegetative bacteria, specified viruses and fungi, and 99.99% of *Mycobacterium* strains (Quantitative TB) as represented by *bovis* and *terae*. See package insert for complete instructions/information on intermediate level disinfection.

Refer to package insert for more detailed usage/product data.

Note Dates Here  
Activation Date

Expiration Date For  
Activated Solution

Check solution upon activation and rinse to each use with Test kit 4. Chloramphenicol Concentration Indicator. This product must be discarded after 28 days, even if the Test kit 4, Chloramphenicol Concentration Indicator indicates "Pass".

Control Number and Expiration Date for Unactivated Solution on Stocker

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:

AUG 2 1985

Under the Federal Insecticide,  
Fungicide, and Rodenticide Act  
as amended, for the pesticide  
registered under EPA Reg. No.

46831-9

BEST COPY AVAILABLE

intermediate level disinfectant when used according to the Instructions for Use.

1. **Germicide Level of Activity:** Omnicide™Plus can be used at the following germicide levels of activity:

**Sterilant:** Omnicide™Plus is a sterilant when used or reused, according to the Directions For Use, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of **ten hours** for a use period not to exceed 28 days.

**High Level Disinfectant:** Omnicide™Plus is a high level disinfectant when used or reused, according to the Directions for Use, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of 45 minutes for a use period not to exceed 28 days.

**Intermediate Level Disinfectant:** Omnicide™Plus is an intermediate level disinfectant when used or reused, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of ten minutes, according to Directions for Use in Section E, for a use period not to exceed 28 days.

A ten minute immersion at 68°F (20°C) will destroy all vegetative bacteria including *S. aureus*, *S. choleraesuis*, *P. aeruginosa*, *E. coli*, all pathogenic fungi, representative viruses and 99.99% of *Mycobacterium* strains (Quantitative TB) as represented by *bovis* and *terrae*.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:  
AUG 2 1995  
Under the Federal Insecticide,  
Fungicide, and Rodenticide Act  
as amended, for the pesticide  
registered under EPA Reg. No.  
176851-9

2. **Reuse Period:** Omnicide™ Plus has demonstrated efficacy in the presence of 5 percent organic soil contamination and a simulated amount of microbiological burden during reuse. **OMNICIDE™ PLUS SOLUTION MAY BE REUSED ONLY WHILE THE MINIMUM EFFECTIVE CONCENTRATION (MEC) AS DETERMINED BY THE PROCHEK G GLUTARALDEHYDE CONCENTRATION INDICATOR TEST STRIP, PH AND TEMPERATURE MEET THE REQUIREMENTS BASED UPON MONITORING AS DESCRIBED IN INDICATOR DIRECTIONS FOR USE.** Efficacy of Omnicide™ Plus solution during its use-life must be verified by the ProCheK G Glutaraldehyde Concentration Indicator Test Strip to determine that at least the MEC as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip is present. The product must be discarded after 28 days.

3. **General Information:** Choose a germicide with the level of microbial activity that is appropriate for the reusable medical device or equipment surface. Follow the reusable device labeling and standard institutional practices. In the absence of complete instructions, use the following guidance:

First, for patient contacting devices, determine whether the reusable device to be processed is a critical, semi-critical, or non-critical device.

- A critical device routinely penetrates the skin or mucous membranes during use or are otherwise used in normally sterile tissues of the body.
- A semi-critical device makes contact with mucous membranes but does not ordinarily penetrate sterile areas of the body.
- A non-critical device contacts only intact skin during routine use.

Second, determine if sterilization, high level or intermediate level disinfection

Critical device: Sterilization is required.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:  
AUG 2 1995  
Under the F... an insecticide,  
Fungicide, and Rodenticide Act  
as amended, for the pesticide  
registered under EPA Reg. No.  
46857-9



Lipid Medium Virus

10 minutes

- Herpes simplex
- HIV-1 (AIDS Virus)

60 seconds at full strength

\* Testing was performed using Omnicide™Plus solution which had been diluted to 1.5 percent including 5 percent bovine calf serum.

5. **Material Compatibility:** Omnicide™Plus solution is recommended for usage with medical devices made from the materials shown below.

- |                    |                          |
|--------------------|--------------------------|
| polypropylene      | vinyl and Tygon tubing** |
| ABS                | nickel plating*          |
| polyethylene*      | acrylic bar*             |
| polycarbonate      | polyethylene tubing*     |
| black oxide steel* | PVC*                     |
| Mylar*             |                          |

\*represents four weeks of exposure

\*\*represents 13 days of exposure

Following sterilization or disinfection, the sterilized or disinfected medical device should be rinsed according to the Instruction for Use, Rinsing (Section E.4), and dried according to manufacturers instructions.

6. **Pre-cleaning Agent Compatibility:** Omnicide™Plus is compatible with enzymatic detergents which are neutral in pH, low foaming and easily rinsed from equipment. Detergents that are either highly acidic or alkaline are contraindicated as precleaning agents since improper rinsing could effect the efficacy of the Omnicide™Plus solution by altering its pH.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:  
AUG 2 1995  
Under the Federal Insecticide,  
Fungicide, and Rodenticide Act  
as amended, for the pesticide  
registered under EPA Reg. No.  
46851-9

B. CONTRAINDICATIONS

1. **Sterilant Usage:** Routine biological monitoring is not possible with Omnicide™Plus solution and therefore Omnicide™Plus solution should NOT be used to sterilize reusable medical devices that are compatible with other sterilization processes that can be biologically monitored.

Omnicide™Plus solution should not be used for sterilization of critical devices intended for single use (e.g. catheters).

2. **High Level Disinfectant Usage:** Omnicide™Plus solution should not be used to high level disinfect a semi-critical device when sterilization is practical.

3. **Endoscope Usage:** Omnicide™Plus solution is not the method of choice for sterilization of rigid endoscopes which the device manufacturer indicates are compatible with steam sterilization.

C. WARNINGS

OMNICIDE™PLUS ACTIVATED DIALDEHYDE SOLUTION IS HAZARDOUS TO HUMANS AND DOMESTIC ANIMALS.

DANGER: Keep Out of Reach of Children  
Contains Glutaraldehyde

- 1. Direct contact is corrosive to exposed tissue, causing eye damage and skin irritation/damage. Do not get into eyes, on skin or on clothing.
- 2. Avoid contamination of food.
- 3. Use in well ventilated area in closed containers.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:  
AUG 2 1995

Under the...  
Fungicide...  
as an...  
registered...  
46857-9

10 9 22

In case of contact, immediately flush eyes or skin with copious amounts of water for at least 15 minutes. For eyes, get medical attention.

Harmful if swallowed. Drink large quantities of water and call a physician immediately.

Probable mucosal damage from oral exposure may contraindicate the use of gastric lavage.

Emergency, safety, or technical information about Omnicide™Plus solution can be obtained from Cottrell, Ltd. Hotline at 1-800-843-3343, or by contacting Chemtrac at 800-535-5053 or by contacting your local Cottrell, Ltd. representative.

**D. PRECAUTIONS**

- 1. Appropriate hand eye, and face protection as well as liquid proof gowns should be worn when cleaning and sterilizing/disinfecting soiled devices and equipment.

More detailed information regarding the handling of the products along with compatible materials is included in the MSDS sheet attached to the product container.

- 2. Contaminated, reusable devices **MUST BE THOROUGHLY CLEANED** prior to disinfection or sterilization, since residual contamination will decrease effectiveness of the germicide.
- 3. The user **MUST** adhere to the **Directions for Use** since any modification will affect the safety and effectiveness of the germicide.
- 4. The reusable device manufacturer should provide the user with a validated reprocessing procedure for that device using Omnicide™Plus solution.

ACCEPTED  
 with COMMENTS  
 in EPA letter Dated:  
 AUG 2 1995  
 Under the Federal Insecticide,  
 Fungicide, and Rodenticide Act  
 as amended, for the pesticide  
 registered under EPA Reg. No.  
 46857-9



1.2 9-22  
6

### 3. Usage

It is a violation of the Federal Law to use this product in a manner inconsistent with its labeling.

#### a. Sterilization (Bucket/Tray Manual System)

Prior to immersing medical equipment/devices, test the solution with a ProChek G Glutaraldehyde Indicator Test Strip to assure that the glutaraldehyde concentration is above its MEC.

Immerse medical equipment/devices completely in Omnicide™Plus Solution for a minimum of ten hours at 68°F (20°C) to eliminate all microorganisms including Clostridium sporogenes and Bacillus subtilis spores. Remove equipment from the solution using sterile technique and rinse thoroughly with sterile water following the rinsing instructions below.

#### b. High Level Disinfection (Bucket/Tray Manual System)

Prior to immersing medical equipment/devices, test the solution with a ProChek G Glutaraldehyde Indicator Test Strip to assure that the glutaraldehyde concentration is above its MEC.

Immerse medical equipment/devices completely in Omnicide™Plus solution for a minimum of 45 minutes at 68°F (20°C) to destroy all pathogenic microorganisms, except for large numbers of bacterial endospores, but including Mycobacterium tuberculosis (Quantitative TB Method). Remove devices and equipment from the solution and rinse thoroughly following the rinsing instructions below.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:  
AUG 2 1995  
Under the Federal Insecticide,  
Fungicide, and Rodenticide Act  
as amended, for the pesticide  
registered under EPA Reg. No.  
46851-9

12

13 9 22

c. Intermediate Level Disinfection

Prior to immersing medical equipment/devices, test the solution with a ProChek G Glutaraldehyde Indicator Test Strip to assure that the glutaraldehyde concentration is above its MEC.

Omnicide™Plus is an intermediate level disinfectant when used or reused, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of ten minutes for a use period not to exceed 28 days.

4. **Rinsing Instructions**

Following immersion in Omnicide™Plus solution thoroughly rinse the equipment or medical device by immersing in two gallons of water. Repeat this procedure a second time with a fresh two gallon volume of water.

For endoscopic instruments with lumens, a minimum of 500 ml of water should be flushed through lumens during each separate rinse unless otherwise noted by the device or equipment manufacturer. Use fresh volumes of water for each rinse. Discard the water following each rinse. Do not reuse the water for rinsing or any other purpose as it will become contaminated with glutaraldehyde.

Refer to the reusable device/equipment manufacturers labeling for rinsing instructions.

**Sterile Water Rinse:**

Critical devices which are sterilized with Omnicide™Plus must be rinsed with sterile water.

**Potable Water Rinse:**

A sterile water rinse is recommended when practical, for all devices. Alternatively a high quality potable water (one that meets Federal Clean Water Standards at point of use) may be used.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:

AUG 2 1995

Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No.

146857-9

The use of potable water for rinsing, increases the risk of contaminating the device or medical equipment with Pseudomonad and atypical (fast growing) Mycobacteria that are often present in potable water supplies. The devices (e.g. colonoscope) need to be completely dried, because any moisture remaining provides an ideal situation for rapid colonization of bacteria. Additionally, mycobacteria are highly resistant to drying, therefore, rapid drying will avoid possible colonization but may not result in a device free from atypical mycobacteria. A final rinse using a 70 percent isopropyl alcohol solution should be used to speed the drying process and reduce the numbers of any organism present as a result of rinsing with potable water.

F. REUSE

Omnicide™Plus solution has demonstrated efficacy in the presence of 5 percent organic soil contamination and a simulated amount of microbiological burden during reuse. The glutaraldehyde concentration of this product during its use-life must be verified by the ProChek G Glutaraldehyde Concentration Indicator Test Strip to determine the solution is above the Minimum Effective Concentration requirement (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip is present. This solution may be used and reused within the limitations indicated above for up to 28 days after activation. Omnicide™ must be discarded after 28 days, even if the ProChek G Glutaraldehyde Indicator Test Strip indicates pass.

G. MONITORING OF GERMICIDE TO ENSURE SPECIFICATIONS ARE MET

It is recommended that the Omnicide™Plus solution be tested with the ProChek G glutaraldehyde test strip prior to each usage. This is to insure that the appropriate concentration of glutaraldehyde is present and to guard against a dilution which may lower the concentration of the glutaraldehyde below its MEC. During the use of Omnicide™Plus as a high level disinfectant and/or sterilant, it is also highly recommended that a thermometer and timer be used to ensure that optimum conditions are met. The pH of the activated solution may be periodically checked to verify that the pH of the solution is between 8.0 and 9.0.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:

AUG 2 1995

United States Environmental Protection Agency,  
Environmental and Systems Research Institute,  
as provided for the products  
registered under EPA Reg. No.  
46851-7

15 9 22

H. POST-PROCESSING HANDLING AND STORAGE OF REUSABLE DEVICES

Sterilized or disinfected reusable devices are either to be immediately used or stored in a manner to minimize contamination. Refer to reusable device equipment manufacturers labeling for additional storage and/or handling instructions.

I. STORAGE CONDITIONS AND EXPIRATION DATE

1. Prior to activation, Omnicide™Plus solution should be stored in its original sealed container at controlled room temperature 15°-30° C (59°-86° F).
2. The expiration date of the unactivated Omnicide™Plus solution and activator will be found on the side of the immediate container.
3. The use period for activated Omnicide™Plus is for no longer than as indicated by ProCheK G Glutaraldehyde Concentration Indicator Test Strip or 28 days following activation. Once activated, the solution requires no further dilution prior to its usage.

J. SAFETY INFORMATION

Emergency, safety, or technical information about Omnicide™Plus solution can be obtained from Cottrell, Ltd. at 1-800-843-3343, Infotrac at 800-535-5053, or by contacting your Cottrell, Ltd. representative.

K. USER TRAINING

The user should be adequately trained in the decontamination and disinfection or sterilization of medical devices and the handling of liquid chemical germicides. Additional information about Omnicide™Plus solution can be obtained from Cottrell, Ltd. at 1-800-843-3343, or by contacting your local Cottrell, Ltd. representative.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:

U.S. Environmental Protection Agency  
Washington, D.C. 20460

46851-9

L. DISPOSAL INFORMATION

0.946 L (1 quart), 3785 L (1 gallon), and 9462 L (2.5 gallon) size container must be triple rinsed and disposed of in accordance with local or state regulations.

M. REORDER INFORMATION

Reorder	Description	Case Contains
PC1032	0.946 L (1 quart)	4 x 0.946 L (4 qts/case)
PC1128	3.785 L (1 gallon)	4 x 3785 L (4 gals/case)
PCG660 Strips	ProCheK G Concentration Indicator	60 strips/bottle (6 btl/case)
PCG615	ProCheK G Concentration	15 strips/pkg

ACCEPTED  
with COMMENTS  
in EPA Letter Dated  
AUG 2 1995

United States Insecticide,  
Fungicide and Rodenticide Act  
as amended for the pesticide  
registered under EPA Reg. No.  
46851-9

17 9 22

OMNICIDE™PLUS SOLUTION  
PACKAGE INSERT  
ENDOSCOPE REPROCESSING  
DRAFT LABELING

A. INTENDED USE/INSTRUCTIONS FOR USE

Omnicide™Plus solution is a liquid chemical sterilant and high level disinfectant for flexible endoscopes when used according to the Instructions For Use. (See Section E).

- 1. Germicide Level of Activity: Omnicide™Plus can be used at the following germicide levels of activity:

Flexible Endoscopes, when expected to penetrate the skin or mucous membranes or are used in otherwise normally sterile tissues of the body during use are critical devices and therefore, are required to be sterile.

**Sterilant:** Omnicide™Plus is a sterilant for flexible endoscopes when used or reused, according to the Directions For Use, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F(20°C) with an immersion time of ten hours for a use period not to exceed 28 days.

Flexible Endoscopes when expected to come in contact without penetration of mucous membranes are semi-critical devices and therefore may be high level disinfected.

**High Level Disinfectant:** Omnicide™Plus is a high level disinfectant for flexible endoscopes when used or reused, according to the Directions for Use, at or above its Minimum Effective Concentration (MEC) as determined by Prochek G Glutaraldehyde Concentration Test Strip at 68°F(20°C) with an immersion time of 45 minutes according to the instructions for use in Section E for a use period not to exceed 28 days.

1500 10  
MAY 1995  
in EPA Reg. No. 40857-9

AUG 2 1995

U.S. Environmental Protection Agency,  
EPA and the Federal Food, Drug, and Cosmetic Act  
as amended, for the pesticide  
registered under EPA Reg. No.  
40857-9

# General Procedure for High Level Disinfection of Flexible Endoscopes

(This procedure is recommended in the absence of specific directions from the device manufacturer)

## Trained Personnel

- Personnel involved in the reprocessing of endoscopes should have the ability to read, understand, and implement instructions from manufactures and regulatory agencies as they relate to endoscopic disinfection.
- The person(s) to whom the job of reprocessing endoscopes is given should have the opportunity to become completely familiar with the mechanical aspects of the endoscopic equipment.
- Training should include familiarization with OSHA regulations and in-house policies on how to appropriately and safely handle liquid chemical germicides.
- Training should also include information on the safe handling of instruments that are contaminated with body fluids after use. This should include familiarization with universal precautions.

## Cleaning of flexible endoscopes

- **Cleaning at the Examination Room**

Reflux of body fluids from the patient may occur in any of the standard channels. Cleaning of endoscopes and accessories should be performed promptly after removing the endoscope from the patient to prevent drying of secretions.

1. Personnel should don all personnel protective equipment.
2. Prepare an enzyme detergent (e.g., Pro EZ or Pro EZ Plus) or one recommended by the scope manufacturer.

ACCEPTED  
with COMMENTS  
in EPA Letter dated

AUG 2 - 1995

Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No.

46091-9

3. Gently wipe all debris from the insertion tube with a moistened gauze or the like.
4. Place the distal end of the flexible endoscope into the water and enzyme detergent solution and aspirate through the biopsy/suction channel for 5-10 seconds or until the solution is visibly clean. Alternate aspiration of the detergent solution and air several times. Finish by suctioning air.
5. Flush or blow out air and water channels in accordance with the endoscope manufacturers instructions.
6. Transport the endoscope to the reprocessing area.

Cleaning at the Reprocessing Area

1. Attach any necessary water-tight caps to the electrical portions of the umbilicus.
2. Before proceeding with any further cleaning steps, the flexible endoscope should be leak tested. (Refer to manufacturers leakage test instructions). Follow the manufacturers instructions if the instrument appears damaged.
3. Fill a sink or basin with a freshly made enzyme (e.g., Pro EZ or Pro EZ Plus) detergent solution.
4. Immerse the endoscope. All channels should be irrigated with copious amounts of detergent and tap water to soften, moisten, and dilute the organic debris. All detachable parts (e.g., hoods and suction valves) should be removed and soaked in the detergent solution. The insertion tube should be washed with detergent solution and rinsed.
5. Use a small soft brush to scrub all detachable parts.
6. Use a brush to clean under the suction valve, air/water valve and biopsy port openings.

**ACCEPTED**  
**with COMMENTS**  
**in EPA Letter Dated:**  
**AUG 2 . 1995**

Under the Federal Insecticide,  
 Fungicide, and Rodenticide Act  
 as amended, for the pesticide  
 registered under EPA Reg. No.  
46857-9

- 7. Brush the entire suction/biopsy system including the body, the insertion tube, and the umbilicus of the endoscope in accordance with the manufacturer's instructions.
- 8. Accessible channel(s) should be brushed to remove particulate matter, and the detergent solution must be suctioned or pumped through all channels to remove dislodged material. Channel irrigators may be useful for this step.  
Fill all channels with detergent solution and soak 2-5 minutes in accordance to the Pro EZ instructions.

Rinse after cleaning:

- 9. Rinse the endoscope and all detachable parts in clean water.
- 10. Rinse all channels well with water to remove debris and detergent.
- 11. Purge water from all channels and wipe dry the exterior of the endoscope with a soft clean cloth to prevent dilution of the Omnicide™Plus disinfectant used in subsequent steps.

Manual Disinfection

- 12. Activate the Omnicide™Plus Long Life activated dialdehyde by adding the entire content of the Activator bottle which is attached to the Omnicide™Plus solution container. Recap the container and shake well. On activation, the solution immediately changes color to blue, thereby indicating the solution is ready. Use ProChek G Glutaraldehyde Concentration Indicator Test Strips after activation to determine that the solution is above the Minimum Effective Concentration (MEC) before use.

Test the activated Omnicide™Plus solution using the ProChek G Glutaraldehyde Concentration Indicator at 68°F (20°C), before each use.

- 13. Attach channel irrigators/adapters and cover the biopsy port in accordance with the manufacturer's instructions.
- 14. Pour the activated Omnicide™Plus into an appropriate

ACCEPTED  
 46857-9  
 AUG 2 - 1995  
 Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No. 46857-9

21 0, 22  
/

15. Completely immerse the endoscope in the basin of Omnicide™Plus. Note: in order to prevent damage to the endoscope, DO NOT soak any other accessory equipment with the endoscope.
16. Inject the Omnicide™Plus solution into all channels of the endoscope until it can be seen exiting the opposite end of each channel. Assure that all channels are filled with disinfectant and that no air pockets remain within the channels.
17. Cover the disinfectant soaking basin with a tight fitting lid to minimize chemical vapor exposure.
18. Soak the endoscope for 45 minutes. Use a timer to ensure adequate soaking time.
19. Before completely removing the endoscope from the disinfectant, flush all channels with air to remove disinfectant.

Rinse After Manual Disinfection

20. Rinse 1: Fill a basin with two gallons of water (preferably sterile water). Place the endoscope into the basin and thoroughly rinse the exterior of the scope. Attach channel irrigators/adapters to the endoscope and flush with 500 ml of water through the channel irrigator. Empty basin.
21. Rinse 2: Fill a basin with two gallons of water (preferably sterile water). Place the endoscope into the basin and thoroughly rinse the exterior of the scope and flush with 500 ml of water through the channel irrigator.
22. Purge all channels with air.
23. Flush all channels with 70% alcohol until the alcohol can be seen exiting the opposite end of each channel.
24. Purge all channels with air.
25. Remove all adapters and devices.

ACCEPTED  
WITH COMMENTS  
in EPA 5-Step Data  
AUG 2 1005  
Under the provisions of the  
Freedom of Information Act  
as amended, for the purpose  
of release under E.O. 13526  
46857-9

13

22 22

Storage

- 26. Dry the exterior of the endoscope with a soft clean cloth. Do not attach detachable parts to the endoscope prior to storage. Storage of endoscopes with the removable parts detached lowers the risk of trapping liquid inside the instrument and facilitates continued drying of the channels and channel openings. To prevent the growth of waterborne organisms, the endoscope and all detached parts should be thoroughly dried prior to storage.
- 27. Hang the endoscope vertically with the distal tip hanging freely in a well ventilated, dust-free cabinet.

References:

- 1. ASTM:F 1518-94, Standard for Cleaning and Disinfection of Flexible Fiberoptic and Video Endoscopes Used in the Examination of the Hollow Visera, Current edition approved May 15, 1994, Published July 1994.
- 2. Martin, M.A., MD, Reichelderfer, M., APIC guideline for infection prevention and control in flexible endoscopy, Association for Professionals in Infection Control and Epidemiology, Inc., AJIC Am J Infect Control 1994: 22: 19-38.
- 3. Vesley, D. et. al., Significant factors in the disinfection and sterilization of flexible endoscopes, AJIC, December 1992, pg. 292.
- 4. Axon, A.T.R., Bond, B., Bottill, F.M., Cowen, A.E., Fleisher, D.E. and Tandon, R.K., Endoscopic disinfection, Working Party Reports, Blackwell Scientific Publications, 1990, 46-50.

ACCEPTED  
 WITH COMMENTS  
 In EPA Letter Dated:  
 AUG 2 1995  
 Under the Clean Air Act  
 Federal Government  
 for the purpose of the  
 regulation under EPA  
 46851-2

85