

---

**LIST OF TABLES**

|                   |  |
|-------------------|--|
| Table A-1         | Data Quality Objectives  |
| Table A-2         | Decision Criteria used for Initial Disposal Classification of Sediments under RCRA and TSCA Disposal Rules |
| Table B-1         | Example Sample ID and Horizontal Coordinates   |
| Table B-2         | River Section 1 Program Summary  |
| Table B-3         | River Section 2 Program Summary  |
| Table B-4         | River Section 3 Program Summary  |
| Table B-5         | Sample Container and Preservation Requirements   |
| Table B-6a – B-6j | Reference Limit and Evaluation Tables for Analytical Methods   |
| Table B-7a – B-7n | Measurement Performance Criteria Tables for Analytical Methods   |
| Table B-8         | Data Collected During Sediment Core Collection   |
| Table B-9         | Data Collected During Sample Processing in the Field Lab   |
| Table B-10        | Valid Values for PCBs  |
| Table C-1         | Summary of Analyses for Initial PE Acceptance Criteria Development   |
| Table C-2         | Summary of Analyses for Inter-Laboratory Comparison Study  |
| Table D-1         | Format of an Environmental Standards Quality Assurance Review  |

**APPENDICES**

|             |   |
|-------------|---|
| Appendix 1  | SOP for Sediment Core Collection  |
| Appendix 2  | SOP for Bathymetric Survey  |
| Appendix 3  | SOP for Sub-Bottom Acoustic and Electromagnetic Surveying Equipment   |
| Appendix 4  | SOP for Sediment Probing  |
| Appendix 5  | SOP for the Analysis of PCBs by SW-846 Method 8082 (GEHR8082)   |
| Appendix 6  | SOP for the Extraction and Cleanup of Sediment/Solid Samples for PCB Analysis Using the Pressurized Fluid Extraction by SW-846 Method 3545 (GEHR3545) |
| Appendix 7  | SOP for the Extraction and Cleanup of Sediment/Soil Samples for PCB Analysis Using the Soxhlet Extraction by SW-846 Method 3540C (GEHR3540C)          |
| Appendix 8  | SOP for Analysis of PCB Homologs by EPA Method 680 (GEHR680)  |
| Appendix 9  | SOP for Data Package Deliverable (DPSOP)  |
| Appendix 10 | SOP for Grain Size  |
| Appendix 11 | SOP for Atterberg Limit   |
| Appendix 12 | SOP for Specific Gravity  |
| Appendix 13 | SOP for Bulk Density  |

---

**APPENDICES (Cont.)**

- Appendix 14 *(This appendix is no longer necessary. It has been left for convenience for potential future use).*
- Appendix 15 SOP for Total Organic Carbon
- Appendix 16 SOP for USCS Classification
- Appendix 17 SOP for Side Scan Survey Procedures
- Appendix 18 SOP for Core Processing
- Appendix 19 *(This appendix is no longer necessary. It has been left for convenience for potential future use).*
- Appendix 20 SOP for <sup>137</sup>Cs
- Appendix 21 SOP for TCLP Preparation by 1311
- Appendix 22 SOP for VOAs With TCLP Preparation
- Appendix 23 SOP for SVOAs With TCLP Preparation
- Appendix 24 SOP for Pesticides With TCLP Preparation
- Appendix 25 SOP for Herbicides With TCLP Preparation
- Appendix 26 SOP for Preparation of Metals and Mercury and Analysis of Mercury in Leachate  
(Refer to Appendix 29 for the Analysis of Metals by ICP)
- Appendix 27 SOP for Ignitability
- Appendix 28 SOP for PCDD/PCDF
- Appendix 29 SOP for Preparation and Analyses of Metals and Mercury in Sediment
- Appendix 30 Performance and Reporting of Field Audits
- Appendix 31 Performance and Reporting of Analytical Laboratory Audits
- Appendix 32 SOP for Data Validation of VOA Data (DV8260B)
- Appendix 33 SOP for Data Validation of SVOA Data (DV8270C)
- Appendix 34 SOP for Data Validation of Pesticide Data (DV8081A)
- Appendix 35 SOP for Data Validation of Herbicide Data (DV8151A)
- Appendix 36 SOP for Data Validation of PCBs (DV8082)
- Appendix 37 SOP for Data Validation of PCB (Homolog) Data (by GEHR680)
- Appendix 38 SOP for Data Validation of PCDD and PCDF Data (DV1613B)
- Appendix 39 SOP for Data Validation of ICP Metals Data (DV6010B)
- Appendix 40 SOP for Data Validation of Mercury Data (DV74707471)
- Appendix 41 SOP for Data Validation of TOC Data (DVTOC)
- Appendix 42 EDD Specifications

## **APPENDICES**

## **APPENDIX 21**

Controlled Copy  
Copy No. \_\_\_\_\_  
Implementation Date: \_\_\_\_\_

SOP No. CORP-IP-0004  
Revision No. 1  
Revision Date: 02/01/00  
Page: 1 of 55


**STL<sup>®</sup> STANDARD OPERATING PROCEDURE**


**TITLE: TOXICITY CHARACTERISTIC LEACHING PROCEDURE  
AND SYNTHETIC PRECIPITATION LEACHING PROCEDURE**


**(SUPERSEDES: REVISION 0)**

Prepared by:   
Kurt Ill

Reviewed by:   
Technology Standardization Committee

Approved by:   
Director, Quality Assurance

Approved by:   
Director, Environmental Health and Safety

Approved by:   
Management

Proprietary Information Statement:

This document has been prepared by and remains the sole property of STL<sup>®</sup> Incorporated. It is submitted to a client or government agency solely for its use in evaluating STL's<sup>®</sup> qualifications in connection with the particular project, certification, or approval for which it was prepared and is to be held proprietary to STL<sup>®</sup>.

The user agrees by its acceptance or use of this document to return it upon STL's<sup>®</sup> request and not to reproduce, copy, lend, or otherwise disclose or dispose of the contents, directly or indirectly, and not to use it for any purpose other than that for which it was specifically furnished. The user also agrees that where consultants or others outside of the user's organization are involved in the evaluation process, access to these documents shall not be given to those parties unless those parties also specifically agree to these conditions.

**Table of Contents**

|  |    |
|--|----|
| 1. SCOPE AND APPLICATION .....   | 3  |
| 2. SUMMARY OF METHOD .....   | 4  |
| 3. DEFINITIONS .....   | 5  |
| 4. INTERFERENCES .....   | 5  |
| 5. SAFETY .....  | 6  |
| 6. EQUIPMENT AND SUPPLIES .....  | 7  |
| 7. REAGENTS AND STANDARDS .....  | 8  |
| 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE .....   | 10 |
| TABLE 1 - HOLDING TIMES (DAYS) .....   | 11 |
| 9. QUALITY CONTROL .....   | 12 |
| 10. CALIBRATION AND STANDARDIZATION .....  | 13 |
| 11. PROCEDURE .....  | 13 |
| 11.1. GENERAL COMMENTS .....   | 13 |
| 11.2. PRELIMINARY EVALUATIONS .....  | 13 |
| TABLE 2 - MINIMUM REQUIRED LEACHATE VOLUMES .....  | 14 |
| 11.3. PROCEDURE: NON-VOLATILE CONSTITUENTS: SEMI-VOLATILES, PESTICIDES, HERBICIDES,<br>METALS (BOTTLE LEACH) ..... | 21 |
| 11.4. PROCEDURE: VOLATILE CONSTITUENTS (ZHE) .....   | 24 |
| 12. DATA ANALYSIS AND CALCULATIONS .....   | 28 |
| 13. METHOD PERFORMANCE .....   | 31 |
| 14. POLLUTION PREVENTION .....   | 31 |
| 15. WASTE MANAGEMENT .....   | 31 |
| 16. REFERENCES .....   | 31 |
| 17. MISCELLANEOUS (METHOD INTERPRETATIONS, RECORD MANAGEMENT, ETC..) .....   | 32 |

**LIST OF APPENDICES**

|                                |    |
|--------------------------------|----|
| APPENDIX A - TABLES .....      | 35 |
| APPENDIX B - FIGURES .....     | 37 |
| APPENDIX C - WORKSHEETS .....  | 41 |
| APPENDIX D - FLOW CHARTS ..... | 48 |

## 1. SCOPE AND APPLICATION

- 1.1. This SOP describes the application of the Toxicity Characteristic Leaching Procedure (TCLP), SW846 Method 1311. The Toxicity Characteristic (TC) of a waste material is established by determining the levels of 8 metals and 31 organic chemicals in the aqueous leachate of a waste. The TC is one of four criteria in 40 CFR Part 261 to determine whether a solid waste is classified as a hazardous waste. The other three are corrosivity, reactivity and ignitability. The TC Rule utilizes the TCLP method to generate the leachate under controlled conditions which were designed to simulate leaching through a landfill. EPA's "worst case" waste disposal model assumes mismanaged wastes will be exposed to leaching by the acidic fluids generated in municipal landfills. The EPA's model also assumes the acid/base characteristics of the waste will be dominated by the landfill fluids. The TCLP procedure directs the testing laboratory to use a more acidic leaching fluid if the sample is an alkaline waste, again in keeping with the model's assumption that the acid fluids will dominate leaching chemistry over time.
- 1.2. The specific list of TC analytes and regulatory limits may be found in Appendix A.  
**Note:** The list in Appendix A does not include the December 1994 EPA rule for Universal Treatment Standards for Land Disposal Restrictions. Those requirements include 216 specific metallic and organic compounds and, in some cases, lower detection limit requirements (see 40 CFR 268.40). TCLP leachates are part of the new Universal Treatment Standards, but the conventional analytical methods will not necessarily meet the new regulatory limits. Consult with the client and with STL<sup>®</sup> Technical Specialists before establishing the instrumental methods for these regulations.
- 1.3. This SOP also describes the application of the Synthetic Precipitation Leaching Procedure (SPLP) which was designed to simulate the leaching that would occur if a waste was disposed in a landfill and exposed only to percolating rain water. The procedure is based on SW846 Method 1312. The list of analytes for SPLP may extend beyond the toxicity characteristic compounds shown in Appendix A. With the exception of the use of a modified extraction fluid, the SPLP and TCLP protocols are essentially equivalent. Where slight differences may exist between the SPLP and TCLP they are distinguished within this SOP.
- 1.4. The procedure is applicable to liquid, solid, and multiphase wastes.
- 1.5. The results obtained are highly dependent on the pH of the extracting solution, the length of time that the sample is exposed to the extracting solution, the temperature

during extraction, and the particle size/surface area of the sample. These parameters must be carefully controlled.

- 1.6. The reporting limits are based on the individual samples as well as the individual analysis techniques. However, the sample is determined to be hazardous if it contains any analyte at levels greater than or equal to the regulatory limits.
- 1.7. If a total analysis of the waste demonstrates that individual analytes are not present in the waste, or that they are present but at such low concentrations that the appropriate regulatory levels could not possibly be exceeded, the procedure need not be run. If the total analysis results indicate that TCLP is not required, the decision to cease TCLP analysis should be remanded to the client.
- 1.8. If an analysis of any one of the liquid fractions of the procedure leachate indicates that a regulated compound is present at such a high concentration that, even after accounting for dilution from the other fractions of the leachate, the concentration would be equal to or above the regulatory level for that compound, then the waste is hazardous and it may not be necessary to analyze the remaining fractions of the leachate. However, the remaining analyses should not be terminated without the approval of the client.
- 1.9. Volatile organic analysis of the leachate obtained using a bottle extraction, normally used for extractable organics and metals, can be used to demonstrate that a waste is hazardous, but only the ZHE option can be used to demonstrate that the concentration of volatile organic compounds is below regulatory limits.

## **2. SUMMARY OF METHOD**

- 2.1. For liquid wastes that contain less than 0.5% dry solid material, the waste, after filtration through 0.6 to 0.8  $\mu\text{m}$  glass fiber filter, is defined as the TCLP leachate.
- 2.2. For wastes containing greater than or equal to 0.5% solids, the liquid, if any, is separated from the solids and stored for later analysis. The particle size of the remaining solid phase is reduced, if necessary. The solid phase is extracted with an amount of extraction fluid equal to 20 times the weight of the solid phase. For TCLP, the extraction fluid employed for extraction of non-volatile analytes is a function of the alkalinity of the solid phase of the waste. For SPLP, the extraction fluid employed is a function of the region of the country where the sample site is located if the sample is a soil. If the sample is a waste or wastewater the extraction fluid employed is a pH 4.2 solution. Two leachates may be generated: a) one for analysis of non-volatile constituents (semi-volatile organics, pesticides, herbicides and metals and/or b) one from a Zero Headspace Extractor (ZHE) for analysis of volatile organic constituents.

Following extraction, the liquid leachate is separated from the solid phase by filtration through a 0.6 to 0.8  $\mu\text{m}$  fiber filter.

- 2.3. If compatible (i.e., multiple phases will not form on combination), the initial liquid phase of the waste is added to the liquid leachate and these are prepared and analyzed together. If incompatible, the liquids are analyzed separately and the results are mathematically combined to yield a volume-weighted average concentration.

### **3. DEFINITIONS**

- 3.1. "Leachate" is used to refer to the TCLP solution generated from this procedure.
- 3.2. "Percent Wet Solids" is that fraction of a waste sample (as a percentage of the total sample) from which no liquid may be forced out by an applied pressure.

### **4. INTERFERENCES**

- 4.1. Oily wastes may present unusual filtration and drying problems. As recommended by EPA (see Figure 3), oily wastes will be assumed to be 100% liquid and analysis for total concentrations of contaminants will be performed. This applies specifically to samples containing viscous non-aqueous liquids that would be difficult to filter.
- 4.2. Wastes containing free organic liquids (i.e., those with separable non-aqueous liquid phases) will be assumed to be 100% liquid and totals analysis will be performed to determine if the oil exceeds TCLP limits.
- 4.3. Solvents, reagents, glassware and other sample processing hardware may yield artifacts and/or interferences to sample analysis. All these materials must be demonstrated to be free from interferences under the conditions of the analysis by analyzing method blanks as described in the Section 9.0 and the individual determinative SOPs.
- 4.4. Glassware and equipment contamination may result in analyte degradation. Soap residue on glassware and equipment may contribute to this. All glassware and equipment should be rinsed very carefully to avoid this problem.
- 4.5. Phthalates may be eliminated by proper glassware cleanup and by avoiding plastics. Only glass, Teflon or Type 316 stainless steel tumblers may be used for leachates to be analyzed for organics. Plastic tumblers may be used for leachates to be analyzed for the metals.

- 4.6. Overexposure of the sample to the environment will result in the loss of volatile components.
- 4.7. Potential interferences that may be encountered during analysis are discussed in the individual analytical methods.

## 5. SAFETY

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL<sup>®</sup> associates.
- 5.2. Eye protection that satisfies ANSI Z87.1 (as per the Chemical Hygiene Plan), laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:
  - 5.3.1. Chemicals that have been classified as **carcinogens**, or **potential carcinogens**, under OSHA include:

Methylene chloride
  - 5.3.2. Chemicals known to be **flammable** are:

Methanol
  - 5.3.3. The following materials are known to be **corrosive**:

Hydrochloric acid, nitric acid, sulfuric acid, acetic acid, sodium hydroxide

The following materials are known to be **oxidizing agents**:

Nitric Acid.
- 5.4. Gas pressurized equipment is employed in this procedure. Be sure all valves and gauges are operating properly and that none of the equipment, especially tubing, is over-pressurized. CAUTION: Do not open equipment that has been pressurized until it has returned to ambient pressure.

- 5.5. A rotary agitation apparatus is used in this procedure. Certain samples may break the glass jars used in the procedure. For these samples, extra caution, including plastic or polyethylene overwraps of the glass jar, may be necessary.
- 5.6. Secure tumbler and extraction apparatus before starting rotary agitation apparatus.
- 5.7. During sample rotation, pressure may build up inside the bottle. Periodic venting of the bottle will relieve pressure.
- 5.8. Exposure to chemicals must be maintained as low as reasonably achievable, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.9. The preparation of standards and reagents and glassware cleaning procedures that involve solvents such as methylene chloride will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.10. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL<sup>®</sup> associate. The situation must be reported **immediately** to a laboratory supervisor.
- 5.11. Due to the potential for ignition and/or flammability, do not attempt to dry non-aqueous liquid samples in an oven.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Extraction vessels
  - 6.1.1. For volatile analytes - zero-headspace extraction (ZHE) vessel, gas-pressure actuated, Millipore YT3009OHW or equivalent (see Figure 2).
  - 6.1.2. For metals - either borosilicate glass jars (1/2 - 1 gallons, with Teflon lid inserts) or 1 L HDPE (Nalgene or equivalent) bottles may be used.
  - 6.1.3. For non-volatile organics - only borosilicate glass may be used.
- 6.2. Vacuum filtration apparatus, capable of 0 - 50 psi and stainless steel pressure filtration apparatus (142 mm diameter), capable of 0 - 50 psi.
- 6.3. Borosilicate glass fiber filters, 0.6 - 0.8  $\mu\text{m}$  (Whatman GF/F 14.2 cm, 0.7  $\mu\text{m}$  or equivalent). When analyzing for metals, wash the filters with 1 N nitric acid and de-ionized water prior to use. As an alternative, certified pre-washed filters may be used. Glass fiber filters are fragile and should be handled with care.

- 6.4. Rotary agitation apparatus, multiple-vessel, Associated Design and Manufacturing Company 3740-6 or equivalent (see Figure 1). The apparatus must be capable of rotating the extraction vessel in an end-over-end fashion at  $30 \pm 2$  rpm.
- 6.5. ZHE Extract Collection Devices are used to collect the initial liquid phase and the final extract of the waste from the ZHE device, either of the following may be used:
  - 6.5.1. Gas-tight syringes, 100 mL capacity, Hamilton 0158330 or equivalent, or
  - 6.5.2. Tedlar bags
- 6.6. Top loading balance, capable of  $0 - 4000 \pm 0.01$ g (all measurements are to be within  $\pm 0.1$  grams).
- 6.7. pH meter and probe capable of reading to the nearest 0.01 unit, and with automatic temperature compensation.
- 6.8. pH probes.
- 6.9. Magnetic stirrer/hotplate and stirring bars.
- 6.10. VOA vials, 20 mL, with caps and septa.
- 6.11. Glass jars, 1/2 - 1 gallon, with Teflon lid-inserts.
- 6.12. Nalgene plastic bottles, 1 liter.
- 6.13. Miscellaneous laboratory glassware and equipment.

## 7. REAGENTS AND STANDARDS

- 7.1. Reagent water for non-volatile constituents must be produced by a Millipore DI system or equivalent. For volatile constituents, water must be passed through an activated carbon filter bed (Milli-Q or tap water passed through activated carbon). Reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks.

- 7.2. Hydrochloric acid, 1 N: Carefully add 83 mL concentrated reagent grade HCl to 800 mL reagent water, cool and dilute to 1 liter with reagent water. Cap and shake to mix well.
- 7.3. Nitric acid, 1 N: Carefully add 64 mL concentrated reagent grade HNO<sub>3</sub> to 800 mL of reagent water, cool and dilute to 1 liter with reagent water. Cap and shake to mix well.
- 7.4. Sodium hydroxide, 1 N: Carefully add 40 g reagent grade NaOH pellets to 800 mL reagent water, stir until the pellets are completely dissolved, cool and dilute to 1 liter with reagent water.

**CAUTION:** Heat is generated during this process.

- 7.5. Acetic acid, glacial: concentrated, reagent grade liquid (HOAc).
- 7.6. pH calibration solutions: buffered to a pH of 4, 7, and 10. Commercially available. Fresh buffer solution must be used each day of analysis.
- 7.7. TCLP Leaching Fluids
  - 7.7.1. General Comments
    - 7.7.1.1. The pH of both solutions listed below should be monitored daily and the pH probes are to be calibrated prior to use.
    - 7.7.1.2. The leaching fluids **MUST** be prepared correctly. If the desired pH range is not achieved and maintained, the TCLP may yield erroneous results due to improper leaching. If the pH is not within the specifications, the fluid must be discarded and fresh extraction fluid prepared.
    - 7.7.1.3. Additional volumes of extraction fluids listed above may be prepared by multiplying the amounts of acetic acid and NaOH by the number of liters of extraction fluid required.
  - 7.7.2. TCLP Fluid #1: Carefully add 5.7 mL glacial acetic acid and 64.3 mL of 1 N NaOH to 500 mL reagent water in a 1 liter volumetric flask. Dilute to a final volume of 1 L with reagent water, cap and shake to mix well. When correctly prepared, the pH of this solution is  $4.93 \pm 0.05$ .
  - 7.7.3. TCLP Fluid #2: Carefully add 5.7 mL glacial acetic acid to 500 mL reagent water in a 1 liter volumetric flask. Dilute to a final volume of 1 L with

reagent water, cap and shake to mix well. When correctly prepared, the pH of this solution is  $2.88 \pm 0.05$ .

- 7.8. Nitric acid, 50% solution: Slowly and carefully add 500 mL concentrated  $\text{HNO}_3$  to 500 mL reagent water. Cap and shake to mix well.
- 7.9. Sulfuric acid / nitric acid (60/40 weight percent mixture)  $\text{H}_2\text{SO}_4/\text{HNO}_3$ . Cautiously mix 60 g of concentrated sulfuric acid with 40 g of concentrated nitric acid.
- 7.10. SPLP Leaching fluids
  - 7.10.1. SPLP solutions are unbuffered and exact pH may not be attained. The pH of TCLP and SPLP fluids should be checked prior to use. If not within specifications, the fluid should be discarded and fresh fluid prepared.
  - 7.10.2. SPLP fluid #1: Add 60/40 weight percent mixture of sulfuric and nitric acids to reagent water until the pH is  $4.20 \pm 0.05$ . This fluid is used for soils from a site that is east of the Mississippi River and for wastes and wastewaters.
  - 7.10.3. SPLP fluid #2: Add 60/40 weight percent mixture of sulfuric and nitric acids to reagent water until the pH is  $5.00 \pm 0.05$ . This fluid is used for soils from a site that is west of the Mississippi River.
  - 7.10.4. SPLP fluid #3: This fluid is reagent water and is used for leaching of volatiles. Additionally, any cyanide-containing waste or soil is leached with fluid #3 because leaching of cyanide containing samples under acidic conditions may result in the formation of hydrogen cyanide gas.
- 7.11. Methanol and methylene chloride - used to aid in cleaning oil contaminated equipment.

## **8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**

- 8.1. Samples being analyzed for non-volatile organic compounds should be collected and stored in glass containers with Teflon lid liners. Chemical preservatives shall NOT be added UNTIL AFTER leachate generation.
- 8.2. Samples being analyzed for metals only can be collected in either glass or polyethylene containers.
- 8.3. When the waste is to be evaluated for volatile analytes, care should be taken to minimize the loss of volatiles. Samples shall be collected and stored in a manner intended to prevent the loss of volatile analytes (e.g., samples should be collected in

Teflon lined septum capped vials with minimal headspace and stored at  $4 \pm 2$  °C). Samples should be opened only immediately prior to extraction.

- 8.4. Samples should be refrigerated to  $4 \pm 2$  °C unless refrigeration results in irreversible physical changes to the waste. If precipitation occurs, the entire sample (including precipitate) should be extracted.
- 8.5. The minimum TCLP sample collection size is determined by the physical state or states of the waste and the analytes of concern. The amount of waste required varies with the percent solids. The lower the percent solids, the more waste will be required for preliminary and final testing. For aqueous samples containing between 0.5 and 10% solids, several kilograms of sample are required to complete the analyses. The general minimal requirements when the samples are 100% solids include: 1 - 32 oz jar for semi-volatile organic analysis and metals, and 1 - 4 oz jar for volatile organic analysis. Low density sample materials, such as rags or vegetation, will require larger volumes of sample. For liquid samples (less than 50% solids), minimum requirements are 2 - 32 oz jars for semi-volatile organic analysis and metals, and 2 - 8 oz jars for volatile organic analysis. If volatile organic analysis is the only requested parameter, 2 separate jars are required. If matrix spike or duplicate control samples are requested, additional sample volume is required. If sufficient sample volumes were not received, analyses cannot be started and the client should be notified as soon as possible.
- 8.6. TCLP leachates should be prepared for analysis and analyzed as soon as possible following extraction. Leachates or portions of leachates for metallic analyte determinations must be acidified with nitric acid to a pH less than 2, unless precipitation occurs. If precipitation occurs upon addition of nitric acid to a small aliquot of the leachate, then the remaining portion of the leachate shall not be acidified and the leachate shall be analyzed as soon as possible. All other leachates should be stored under refrigeration ( $4 \pm 2$  °C) until analyzed. ZHE leachates must be stored in VOA vials filled to eliminate all headspace.
- 8.7. Samples are subject to appropriate treatment within the following time periods:

**Table 1 - Holding Times (days)**

| Parameter       | Collection to Start of TCLP Leach | End of TCLP Tumble to Preparation | Start of TCLP Leach or Semi-volatile Prep Extraction to Analysis | Total Elapsed Time |
|-----------------|-----------------------------------|-----------------------------------|--|--------------------|
| Volatiles:      | 14                                | N/A                               | 14   | 28                 |
| Semi-volatiles: | 14                                | 7                                 | 40   | 61                 |
| Mercury:        | 28                                | N/A                               | 28   | 56                 |
| Other Metals:   | 180                               | N/A                               | 180  | 360                |

**NOTE:** The initial holding time is measured from date of collection to date TCLP extraction started. (This should be the TCLP extraction date in QuanTims.) Semi-volatile method prep holding time is measured from the day tumbling is complete to the start of method extraction. Subsequent analysis holding times are measured from the date extraction (TCLP or method prep) starts. If sample holding times are exceeded, the values obtained will be considered minimal concentrations. Exceeding holding times is not acceptable in establishing that a waste does not exceed the regulatory level. Exceeding the holding time will not invalidate characterization if the waste exceeds the regulatory limit. The Total Elapsed Time is to be used as guidance. If preps are initiated at the last possible moment of a holding time, the elapsed times may be exceeded.

## 9. QUALITY CONTROL

- 9.1. Quality Control Batch (QC Batch) - QA-003 defines a QC Batch as a set of up to 20 field samples of similar matrix that behave similarly and are processed using the same procedures, reagents and standards within the same time period. The same lot of reagents must be used within a batch. A minimum of one TCLP extraction blank (Method Blank), one Laboratory Control Sample (LCS), one Matrix Spike (MS), and one Matrix Spike Duplicate (MSD) will be prepared with each TCLP leachate batch.
- 9.2. Batching Samples - Groups of samples with visibly different bulk matrices (e.g., petroleum sludge and soil samples) must be batched separately for QC testing purposes.
- 9.3. TCLP Extraction Blanks - A minimum of one blank (using the same extraction fluid as used for the samples) must be prepared and analyzed for every batch of samples extracted in a particular vessel type. The blanks are generated in the same way as the samples (i.e., blanks will be tumbled and filtered with the samples). Extraction vessels will be uniquely numbered. Each time a new batch is set up the blank should be rotated sequentially to the next vessel to ensure all vessels are periodically checked. Consult the STL<sup>®</sup> QC Program and the individual analysis SOPs for blank acceptance criteria.
- 9.4. Laboratory Control Sample (LCS) - A LCS is required with each batch of 20 or fewer samples. The LCS shall be generated after a batch of TCLP leachates have been generated (i.e., at the time of the preparative digestion or extraction) by spiking an aliquot of the appropriate extraction fluid used for that batch. Consult the individual analysis SOPs for additional LCS guidance (i.e., spike amounts, spike levels, recovery criteria, etc.).
- 9.5. Matrix Spike (MS/MSD) - Matrix spikes are used to monitor the performance of the analytical methods on the matrix and to assess the presence of interferences. A MS/MSD pair are required with each batch of 20 or fewer samples.
  - 9.5.1. Matrix spikes are to be added after filtration of the TCLP leachate. Spikes are not to be added prior to the TCLP leaching. For metals, matrix spikes are to be added before preservation with nitric acid.

- 9.5.2. The use of internal calibration or alternate methods may be needed when the recovery of the matrix spike is below the expected performance (see Section 9.6.2).
- 9.5.3. Consult the individual analysis SOPs for additional guidance on spike compounds and levels.

#### 9.6. Corrective Actions

- 9.6.1. Consult the STL<sup>®</sup> QC Program and individual analysis SOPs for corrective action for blanks and LCS
- 9.6.2. Method of Standard Additions (MSA) shall be used for metals if all of the following conditions are met:
- Recovery of the analyte in matrix spike is not at least 50%,
  - the concentration of the analyte does not exceed the regulatory level, and
  - the concentration of the analyte measured in the sample is within 20% of the appropriate regulatory level.

If the matrix spike recovery is 5% or less due to dilution or matrix interference, contact the project manager and client for guidance. The client should also be contacted prior to initiation of any MSA steps. Refer to the individual analysis SOPs for details on how to perform MSA analysis.

### 10. CALIBRATION AND STANDARDIZATION

- 10.1. Refer to appropriate analysis SOPs.

### 11. PROCEDURE

#### 11.1. GENERAL COMMENTS

- 11.1.1. One-time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented on a Nonconformance Memo kept in the project file and described in the final report. The variation must be approved by a project manager, Technical Specialist and QA Manager. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

- 11.2. PRELIMINARY SAMPLE EVALUATIONS (Refer to Flow Chart #1, Appendix D)

- 11.2.1. Preliminary TCLP evaluations (percent solids, particle size, selection of extraction fluid, and fluid/leachate compatibility) are required to be done using a minimum of a 100 gram aliquot of waste. This aliquot may also undergo the actual TCLP or SPLP extraction for Non-volatiles ONLY IF it has NOT been oven dried. If the solid portion is oven dried, a separate aliquot must be used for the actual leaching procedure.
- 11.2.2. Consult the holding times for the appropriate tests (Section 8.7) and prioritize extractions such that holding times are not exceeded.
- 11.2.3. Determine the total volume of TCLP leachate (solid phase leachate + liquid filtrate) that needs to be generated for analysis according to the following:

**Table 2. Minimum Required Leachate Volume**

| Analysis       | Required Volume (mL) |
|----------------|----------------------|
| Volatiles      | 3 x 20               |
| Semi-volatiles | 400                  |
| Pesticides     | 200                  |
| Herbicides     | 200                  |
| Metals         | 300                  |

- 11.2.3.1. SPLP - similar volumes are required for volatiles and metals. If semivolatiles, pesticides or herbicides are required, a full 1 L volume must be prepared for each test requested.
- 11.2.3.2. For TCLP and SPLP samples used for matrix spike and matrix spike duplicate analysis, three times the listed volumes are required.
- 11.2.4. **Sample Description** (enter data on Worksheet 1)
  - 11.2.4.1. Record the number of phases observed in the sample. It is common that when more than one container of multi-phasic materials is received from the field, each container will show different amounts of each phase.
    - 11.2.4.1.1. If the sample has multiple phases and is received in more than 1 bottle then the contents of each bottle should be combined in a single larger container prior to processing the sample further. If this is not

possible, then the alternate procedure described in the following section should be used.

11.2.4.1.2. Properly record the relative amounts of each phase by measuring the depth of the layers in each container after the contents have been allowed to settle. Determine the combined volume of each phase for all containers. Then mark the phase composition on a single container, mix thoroughly to obtain a representative subsample, and accurately measure the phase composition according to the following procedure. The two sets of values (combined volumes per phase and phase composition for one container) are used to determine the correct volume/mass adjustments on the final result. This procedure is not appropriate if testing will be done for volatile organic compounds.

11.2.4.2. Solid - record the visible presence of a solid material heavier than water. If the sample contains more than one solid phase (e.g., wood and sediment mixed with water), attach additional documentation to Worksheet 1, line A1.

11.2.4.3. Liquid - record the number of liquid phases observed in the sample according to apparent density. It may be impossible to distinguish apparent density if only one liquid phase is observed and there is no indication on the COC form. If this is the case, record it as aqueous material and let the subsequent analytical record show if the liquid is organic.

11.2.5. **Percent Solid Phase** (enter data on Worksheet No. 1)

11.2.5.1. Percent Solids and ZHE Extractions - The ZHE filtration apparatus cannot accurately determine percent solids less than 5%. If an extraction is to be performed solely for volatile organic compounds and the percent solids concentration is apparently greater than 5%, proceed to Section 11.4 (Procedure: ZHE Extraction Procedure, Volatile Constituents). Otherwise, continue with the steps in this section. The aliquot of sample used here cannot be used again for the ZHE extraction.

11.2.5.2. Determine Type of Filtration Apparatus Needed

- 11.2.5.2.1. If the waste will obviously yield no free liquid when subjected to pressure filtration (i.e., it is 100% solid), then proceed to Section 11.2.6 (Particle-size Reduction).
- 11.2.5.2.2. If the sample is mostly a non-viscous liquid (water or non-viscous organic liquid) of low solids content (<10%) or a highly granular, liquid containing waste vacuum filtration may be used.
- 11.2.5.2.3. If the sample is viscous (sludge or has high solids content), use pressure filtration.
- 11.2.5.3. Weight of filter - Measure and record this value before loading the filter into the filter holder (Worksheet 1, line B1).
- 11.2.5.4. Weight of subsample and filtrate for percent solids measurement
  - 11.2.5.4.1. Assemble the filtration apparatus (use blunt forceps to handle the 0.6 to 0.8  $\mu\text{m}$  filter membrane).
  - 11.2.5.4.2. Homogenize the waste, transfer a minimum of a 100 g subsample to the weighing vessel. Measure and record the gross weight (Worksheet 1, line B2a).
  - 11.2.5.4.3. Measure and record the tare weight of the filtration vessel (Worksheet 1, line B3b).
  - 11.2.5.4.4. Transfer the sample to the filtration device attempting to spread the waste sample evenly over the surface of the filter. Measure and record the tare weight of the empty weighing vessel and any residual sample (Worksheet 1, line B2b).
  - 11.2.5.4.5. Calculate and record the net weight of sample used for testing (Worksheet 1, line B2c).
- 11.2.5.5. Filtration for percent solids
  - 11.2.5.5.1. Slowly apply gentle pressure or vacuum of 10 psi to the filtration apparatus. Allow the sample to filter until no SIGNIFICANT additional liquid has passed through the filter during a 2 minute period.

- 11.2.5.5.2. Repeat previous step by increasing the pressure in 10 psi increments until a maximum of 50 psi is reached. Stop the filtration when no additional filtrate is generated within a 2 minute period.

**Note:** Some samples will contain liquid material that does not filter (e.g., oil). Do not attempt to filter the sample again by exchanging filters. Viscous oils or any wastes which does not pass through the filter is classified as a solid.

- 11.2.5.5.3. Remove the filtrate collection vessel, weigh and record the gross weight (Worksheet 1, line B3a).

- 11.2.5.5.4. Calculate and record the net weight of filtrate (Worksheet 1, line B3c). This result will be used in the percent solids calculation.

- 11.2.5.5.5. Pour the filtrate into a graduated cylinder. Measure and record the volume of the aqueous phase (Worksheet 1, line B7). Measure and record the volume of any organic phase (Worksheet 1, line B8). If more than one organic phase is observed, enter “see below” and provide a description at the bottom of Worksheet 1. These results will be used in the final sample calculations (Worksheets 5 & 6).

- 11.2.5.5.6. Retain the filtrate for use in Section 11.2.8 (Determination of Filtrate/Extraction Fluid Compatibility), and for possible recombination with the filtrate obtained in Section 11.3.

#### 11.2.5.6. Percent of Wet Solids

- 11.2.5.6.1. Calculate the total weight of wet solids using Equation 0 on Worksheet 1 and record the result on line B4.

- 11.2.5.6.2. Calculate the weight percent of wet solids using Equation 1 on Worksheet 1 and record the result on line B5.

- 11.2.5.6.3. If the percent wet solids result is  $\geq 0.5\%$  and  $< 5.0\%$ , and it is noticed that a small amount of the aqueous filtrate is entrained in the wetting of the filter, proceed

to Section 11.2.5.7 to complete the percent solids measurement on a dry-weight basis.

**Note:** If obviously oily (non-aqueous) material is entrained on the filter, do not dry the filter; proceed to Section 11.2.6 (Particle-Size Reduction).

11.2.5.6.4. If the percent wet solids result is greater than 5.0%, proceed to Section 11.2.6 (Particle-Size Reduction) and mark “ $\geq 0.5\%$ ” on Worksheet 1, line B6c.

11.2.5.6.5. If the percent wet solids result is less than 0.5%, discard the solid phase. No leaching will be necessary; the filtrate is equivalent to the final leachate.

11.2.5.7. Weight percent of dry solids (skip this step for oily samples).

**Note:** These steps are required only if it is noticed that a small amount of the filtrate is entrained in wetting of the filter and the percent wet solids content is  $\geq 0.5\%$  and  $< 5.0\%$ .

11.2.5.7.1. Remove the filter with the wet solids from the filtration apparatus.

11.2.5.7.2. Dry the filter and solid phase at  $100 \pm 20^\circ \text{C}$ .

11.2.5.7.3. Remove the filter from the oven and allow to cool in a desiccator.

11.2.5.7.4. Weigh and record the gross dry weight (Worksheet 1, line B6a).

11.2.5.7.5. Repeat the drying step. Weigh and record the second gross dry weight (Worksheet 1, line B6b). If the two weightings do not agree within 1%, perform additional drying and weighing until successive weightings agree within 1%.

11.2.5.7.6. Calculate the weight percent of dry solids using Equation 2 on Worksheet 1 and record the result on line B6c.

11.2.5.7.7. If the dry solids result is  $\geq 0.5\%$  and the sample will be extracted for non-volatile constituents, proceed to

Section 11.2.6 (Particle Size Reduction) using a fresh wet portion of waste.

11.2.5.7.8. If the percent solids result is less than 0.5%, discard the solid phase. No leaching will be necessary; the filtrate is the TCLP leachate. Proceed to Section 11.2.8 (Determination of Filtrate/Leachate Compatibility) to determine whether or not the material is a non-aqueous, immiscible liquid.

**11.2.6. Particle-size Reduction for Fluid Selection** (enter data on Worksheet 2)

11.2.6.1. The subsample used for fluid selection must consist of particles less than 1 mm in diameter (versus the less than 1 cm requirement for the material used for the actual extraction). The method requires a smaller particle size to partially compensate for the shorter duration of contact time with the leachate solution as compared to the full extraction. Inappropriate use of coarser materials could result in the selection of the wrong fluid type.

11.2.6.2. Surface area exclusion - size reduction is not required if the sample surface area is greater than or equal to 3.1 cm<sup>2</sup> per gram. Enter "No" on Worksheet 2, line C1.

11.2.6.3. If the sample contains particles greater than 1 mm in diameter, crush, cut, or grind the solids to the required size. Enter "Yes" on Worksheet 2, line C1.

11.2.6.4. Consult your supervisor or manager when dealing with unusual sample matrices (e.g., wood, cloth, metal, brick).

**11.2.7. Determination of Appropriate Extraction Fluid** (Worksheet 2)

11.2.7.1. If the solid content is greater than or equal to 0.5%, and if the sample is being analyzed for metals or nonvolatile organic compounds, the type of leaching solution must be determined.

11.2.7.2. Follow times, temperature, and particle size specified in this section as closely as possible. If reaction time between the acid solution and solid waste is too short or too long, the procedure may produce false pH readings.

11.2.7.3. For SPLP, refer to Section 7.10 for fluid selection. Matrix type must be specified by the client. Check special instructions or see

the project manager, then put a check mark by the fluid type selected (Worksheet 2, D).

- 11.2.7.4. The TCLP leaching fluid for all volatiles is Fluid #1.
- 11.2.7.5. For TCLP leach fluid determination for non-volatile analytes, continue with the following steps.
- 11.2.7.6. Calibrate the pH meter with fresh buffer solution in accordance with the pH SOP.
- 11.2.7.7. Weigh out a  $5.0 \pm 0.1$  g subsample (less than 1 mm particle size) of the solid phase into a 250-mL beaker, and enter a “✓” on Worksheet 2, line C2. If 5.0 grams not used, enter the actual weight on line C2.
- 11.2.7.8. Add  $96.5 \pm 1.0$  mL of reagent water, cover with a watchglass, and stir for 5 minutes on a stirrer, and enter a “✓” on Worksheet 2, line C3. If a different volume used, enter the actual volume on line C3.
- 11.2.7.9. Measure and record the sample pH (Worksheet 2, line C4).  
**Note:** To avoid damaging the pH probe when organic liquid is present, use narrow range pH indicator paper.
- 11.2.7.10. If the pH is less than or equal to 5.0, use Fluid #1 and proceed to Section 11.2.8 (Fluid Compatibility).
- 11.2.7.11. If the fluid pH is greater than 5.0, add 3.5 mL 1 N HCl, cover with a watchglass. Slurry the sample briefly then heat at  $50^{\circ}\text{C}$  for 10 minutes. Enter a “✓” on Worksheet 2, line C5 and C6).  
**Note:** The heating cycle is a critical step. If the solid waste does not remain in contact with the acidic solution under specified time and temperature conditions, an erroneous pH may be measured.
- 11.2.7.12. Cool to room temperature.
- 11.2.7.13. Measure and record the pH immediately after the sample has reached room temperature (Worksheet 2, line C7).
  - 11.2.7.13.1. If the pH is less than or equal to 5.0, use Fluid #1. Enter a “✓” on Worksheet 2, line D1.

11.2.7.13.2. If the pH is greater than 5.0, use Fluid #2. Enter a “✓” on Worksheet 2, line D2.

11.2.8. **Determination of Filtrate/Extraction Fluid Compatibility** (skip this step for SPLP extractions)

11.2.8.1. Place 5 mL of the appropriate leaching fluid (determined in the previous step) into a 20-25 mL vial.

**Note:** Use fluid type # 1 if simply testing the filtrate for a sample with less than 0.5% solids.

11.2.8.2. Add 5 mL of the initial filtrate, cap and shake.

11.2.8.3. If the phases are miscible, the initial filtrate and solid phase leachate will be physically recombined upon completion of the leachate generation. Enter a “✓” on Worksheet 2, line D3 for TCLP, or on line D1 for SPLP.

11.2.8.4. If the phases are NOT miscible, enter “No” on Worksheet 2, line D3. The initial filtrate and the solid phase leachate will be prepared and analyzed separately and the results mathematically combined (see Section 12.1.4).

11.2.9. For samples requiring analysis for semi-volatile organics, pesticides, herbicides or metals proceed to Section 11.3.

11.2.10. For samples requiring analysis for volatile organics (ZHE), proceed to Section 11.4.

11.3. **BOTTLE EXTRACTION PROCEDURE: NON-VOLATILE CONSTITUENTS: SEMI-VOLATILES, PESTICIDES, HERBICIDES, METALS** (Refer to Flow Chart #2, Appendix D)

11.3.1. All masses should be recorded to the nearest 0.1 g.

11.3.2. The aliquot used in the Preliminary Evaluation MAY be used for this procedure ONLY if it was not oven dried. If the sample is 100% solid or if the preliminary aliquot was not oven dried proceed directly to Section 11.3.7 (Particle Size Reduction). If the Preliminary Evaluation aliquot was oven dried then, using a fresh aliquot of sample, continue as described in Sections 11.3.3 through 11.3.6.

- 11.3.3. Examine the sample and determine the type of filtration to employ per Section 11.2.5.2.
- 11.3.4. Repeat the steps outlined in Sections 11.2.5.3 through 11.2.5.5.3.
- 11.3.5. Determine and record the volume (mass) of the initial filtrate. Cover with aluminum foil and retain for use as defined in Section 11.3.18.
- 11.3.6. Determine and record the “solid” phase mass by subtracting the mass of the liquid filtrate from the mass of the subsample.
- 11.3.7. Evaluate the solid portion of the waste for particle size. If it contains particles greater than 1 cm in size, prepare the solid portion of the waste for leaching by crushing, cutting, or grinding such that all particles are less than 1 cm in size (i.e., capable of passing through a 9.5 mm, 0.375 inch, standard sieve). Size reduction is not required if the sample surface area is greater than or equal to 3.1 cm<sup>2</sup> per gram. If particle size reduction was required, record this on Worksheet 3, line E1.
  - 11.3.7.1. Consult your supervisor or manager when dealing with unusual sample matrices (e.g., wood, cloth, metal, brick). Scissors or shears may be used to cut cloth, plastic or sheet metal. Saws may be used for wood or solid metal. Bricks, rocks, or other solids amenable to grinding should be subcontracted out for particle size reduction. (Contact PA or PM.) Note that size reduction to fine powder is not appropriate, and could invalidate results. If necessary, consult client for guidance.
- 11.3.8. Determine the minimum total volume of solid phase leachate that needs to be generated. Refer to Section 11.2.3.
- 11.3.9. Divide the total volume of solid phase leachate required by 20 to determine the mass of solid phase required for leaching. Round this mass UP to the nearest 5g.
- 11.3.10. Weigh the required mass of solid phase into an appropriate bottle (plastic for metals only, glass for all others) and **slowly** add 20 times its mass of appropriate leaching fluid as determined under Section 11.2.7 (e.g., 20 g of sample would require 400 g of leaching fluid). Record the weight of the sample aliquoted for the extraction on Worksheet 3, line E2 and the amount of extraction fluid added on line F1.
- 11.3.11. Ensure any effervescence has stopped before capping the bottle tightly. Secure in a rotary agitator and rotate end-over-end at 28-32 rpm for 16-20

hours. The temperature of the room should be  $23 \pm 2^{\circ}\text{C}$ . The room temperature and time should be checked at both the start and end of the extraction and recorded on Worksheet 3, lines G1 through G3.

**NOTE:** As agitation continues, pressure may build up within the bottle for some types of wastes. To relieve excessive pressure, the bottle may be removed and opened periodically in a properly vented hood to relieve any built-up pressure.

- 11.3.12. Remove the bottle and filter the sample using vacuum or pressure filtration by filtering through a new glass fiber filter as discussed in Sections 11.2.5.5.1 - 11.2.5.5.2. For final filtration of the TCLP leachate, the glass fiber filter may be changed, if necessary, to facilitate filtration. Filters must be acid washed if metals are to be determined (see Section 6.3). The entire sample need not be filtered; however, sufficient volume should be generated to support the required analyses. Record the date and time the filtration is completed on Worksheet 3, Line G4.
- 11.3.13. If the waste contained no initial filtrate, this solution from 11.3.12 is defined as the TCLP leachate.
- 11.3.14. If the waste did yield an initial filtrate, consult the worksheet for initial filtrate/leachate compatibility. If they are compatible, they are to be combined in the correct proportions (see Section 12.1.4) and mixed well. This combined solution is defined as the TCLP leachate.
- 11.3.15. If the individual phases are NOT compatible, they are to be prepared and analyzed separately and the results combined mathematically. See Section 12.1.5.
- 11.3.16. Measure and record the pH of the TCLP leachate on Worksheet 3, Line 5. (Do not attempt to measure the pH of oily samples as the probe may be rendered inoperable.)
- 11.3.17. Prepare subsamples for metals for MS/MSD quality control testing using the appropriate TCLP spiking solution (do not spike for organics). Refer to the appropriate determinative SOPs for further guidance on the spike components, levels and action criteria.
- 11.3.18. Immediately preserve the leachate as follows:

|            |  |
|------------|--|
| Metals     | pH < 2 w/50% HNO <sub>3</sub> for non-oils (do not acidify oils) |
| All others | Refrigerate to $4 \pm 2^{\circ}\text{C}$                         |

**Note:** Refer to Section 8.6 if precipitation occurs upon preservation.

11.3.19. Label each sample with the appropriate information and submit to the appropriate analytical groups for prep and analysis with copies of the TCLP preparation worksheets.

11.4. **ZHE EXTRACTION PROCEDURE: VOLATILE CONSTITUENTS** (Refer to Flow Chart #3, Appendix D)

11.4.1. Use the ZHE device to obtain a TCLP leachate for analysis of volatile compounds only. Leachate resulting from the use of the ZHE shall NOT be used to evaluate the mobility of non-volatile analytes (e.g., metals, pesticides, etc...).

11.4.2. Due to some shortcomings of the method, losses of volatile compounds may occur. Extra care should be observed during the ZHE procedure to ensure that such losses are minimized. Charge the ZHE with sample only once and do not open the device until the final extract has been collected. Do not allow the waste, the initial liquid phase or the extract to be exposed to the atmosphere any longer than necessary.

11.4.3. If the TCLP extraction is for volatile components only, refer to Section 11.2.5.1 before proceeding.

11.4.4. All masses should be recorded to the nearest 0.1 g.

11.4.5. Assemble the ZHE apparatus. Test for leakage by closing all valves except the gas inlet/outlet valve and pressurizing to 50 psi. Allow to stand for 15 minutes and check the pressure on the built-in gauge to make sure it is not leaking. If the pressure is NOT 50 psi, consult your supervisor.

11.4.6. Adjust the ZHE piston in the ZHE body to the appropriate height (slightly moisten the O-rings with leaching fluid if necessary).

11.4.7. Consult the worksheet and examine the sample. If the sample appears to be different from the preliminary information found on the worksheet, consult your supervisor.

11.4.8. If the preliminary evaluations indicated the need for particle size reduction, homogenize the waste, weigh out a sufficient size subsample and prepare for leaching by crushing, cutting, or grinding such that all particles are less than 1 cm in size as measured with a ruler (Do NOT sieve the sample). Size reduction is not required if the sample surface area is greater than or equal to

**TABLE II. ICP and FLAA Soil Matrix Spike and Aqueous LCS Levels**

| ELEMENT     | Working LCS/MS Standard (mg/L) | Aqueous LCS/MS Level* (ug/L) | Soil MS Level ** (mg/Kg) |
|-------------|--------------------------------|------------------------------|--------------------------|
| Aluminum    | 200                            | 2000                         | 200                      |
| Antimony    | 50                             | 500                          | 50                       |
| Arsenic     | 200                            | 2000                         | 200                      |
| Barium      | 200                            | 2000                         | 200                      |
| Beryllium   | 5                              | 50                           | 5                        |
| Cadmium     | 5                              | 50                           | 5                        |
| Calcium     | 5000                           | 50000                        | 5000                     |
| Chromium    | 20                             | 200                          | 20                       |
| Cobalt      | 50                             | 500                          | 50                       |
| Copper      | 25                             | 250                          | 25                       |
| Iron        | 100                            | 1000                         | 100                      |
| Lead        | 50                             | 500                          | 50                       |
| Lithium     | 100                            | 1000                         | 100                      |
| Magnesium   | 5000                           | 50000                        | 5000                     |
| Manganese   | 50                             | 500                          | 50                       |
| Molybdenum  | 100                            | 1000                         | 100                      |
| Nickel      | 50                             | 500                          | 50                       |
| Phosphorous | 1000                           | 10000                        | 1000                     |
| Potassium   | 5000                           | 50000                        | 5000                     |
| Selenium    | 200                            | 2000                         | 200                      |
| Silver      | 5                              | 50                           | 5                        |
| Sodium      | 5000                           | 50000                        | 5000                     |
| Strontium   | 100                            | 1000                         | 100                      |
| Thallium    | 200                            | 2000                         | 200                      |
| Vanadium    | 50                             | 500                          | 50                       |
| Zinc        | 50                             | 500                          | 50                       |
| Boron       | 100                            | 1000                         | 100                      |
| Silica      | 1000                           | 10000                        | 1000                     |
| Tin         | 200                            | 2000                         | 200                      |
| Titanium    | 100                            | 1000                         | 100                      |

\* Levels shown indicate the spike concentration in the final digestate of the aqueous LCS or matrix spike based on the addition of 1.0 mL working spike (7.3) to 100 mL of sample.

\*\* Final soil spike concentration based on the addition of 1.0 mL working spike (7.3) to 1.0 g of sample/100 mL final volume (assumes 100% solids).

**TABLE III. GFAA Soil Matrix Spike and Aqueous LCS Spike Levels**

3.1 cm<sup>2</sup> per gram If particle size reduction was required record this on Worksheet 4, Line H1.

**Note:** To minimize loss of volatiles, samples for volatiles that require particle size reduction should be kept in sample storage (at 4 °C) until immediately before size reduction. Aggressive reduction which would generate heat should be avoided and exposure of the waste to the atmosphere should be avoided to the extent possible. Size reduction to a fine powder is not appropriate. Also see Section 11.3.11.

11.4.8.1. Consult your supervisor or manager when dealing with unusual sample matrices (e.g., wood, cloth, metal, brick). Scissors or shears may be used to cut cloth, plastic or sheet metal. Saws may be used for wood or solid metal. Bricks, rocks, or other solids amenable to grinding should be subcontracted out for particle size reduction (Contact PA or PM).

11.4.9. Place the ZHE apparatus on the balance and tare the balance.

11.4.10. Determine the appropriate size subsample to weigh using the percent solids information from Section 11.2.5 and record the weight used on Worksheet 4, Line H2.

11.4.10.1. For wastes that are 100% solids, a 25 g sample is used.

11.4.10.2. For wastes containing < 0.5% solids the liquid portion of the waste, after filtration, is defined as the TCLP leachate. Filter enough of the sample to support all of the volatile analyses required.

11.4.10.3. For wastes containing  $\geq 0.5\%$  and  $< 5.0\%$  solids, a 500 g subsample of waste is recommended.

11.4.10.4. If the sample has  $\geq 5.0\%$  solids, the appropriate sample size should be determined using the equation in Section 12.1.2.

**Note:** For wastes containing greater than 0.5% wet or dry solids (Section 11.2.5), the “solids” value from the ZHE filtration process may be used to determine the volume of fluid to load into the ZHE. This approach is recommended since the solids value from Section 11.2.5 may differ from the filtration solids due to sample variability or differences in the filtration apparatus.

- 11.4.11. Homogenize and transfer an appropriate size subsample of the waste into the ZHE and record the mass on Worksheet 4, Line I1.
- 11.4.12. Carefully place the glass fiber filter between the support screens and secure to the ZHE. Tighten all the fittings.
- 11.4.13. Place the ZHE in a vertical position; open both the gas AND liquid inlet/outlet valves. Attach a gas line to the gas inlet/outlet valve.
- 11.4.14. If the waste is 100% solid, slowly increase the pressure to a maximum of 50 psi to force out as much headspace as possible and proceed to Section 11.4.18.
- 11.4.15. If the waste is < 100% solids, carefully apply gentle pressure of 10 psi (or more, if necessary) to force all headspace slowly out of the ZHE. At the FIRST appearance of liquid from the liquid inlet/outlet valve, quickly close the valve and discontinue gas pressure.
- 11.4.16. Assemble a syringe and place the plunger in all the way. Adjust the tension on the plunger to provide slight drag. Attach the pre-weighed syringe or Tedlar bag to the liquid inlet/outlet valve and open the valve. Record the tare weight of the collection device on Worksheet 4, Line 3b.
- 11.4.17. Carefully apply gas pressure of no more than 10 psi to force out the liquid phase. Allow the sample to filter until no SIGNIFICANT additional filtrate has passed in a 2 minute period.

**Note:** If the capacity of the syringe is reached, close the liquid inlet/outlet valve, discontinue gas pressure, remove the syringe and return to Section 11.4.15.

- 11.4.18. Repeat previous step increasing the pressure in 10 PSI increments until 50 psi is reached and no significant liquid has passed in a 2 minute period. Remove the collection device and record the total weight of the collection device with filtrate on Worksheet 4, Line 3a. Close the valve and discontinue gas pressure. Transfer the filtrate to VOA vials and label appropriately. Calculate the weight of filtrate collected and record on Worksheet 4, Line I3c.

**Notes:** If the original waste contained less than 0.5% solids (Section 11.2.5), this filtrate is defined as the TCLP leachate and you may proceed to Section 11.4.28. Otherwise, save the vials by storing at 4 C under minimal headspace conditions, for recombination as in Section 11.4.27.

The material remaining in the ZHE is defined to be the “solid” phase. Calculate the weight of the solid phase using Equation 4.

- 11.4.19. Based on the information from Sections 11.2.5 and 11.4.11 and using the formula in 12.1.3, determine the weight of fluid to load into the ZHE on the “solid” phase. The ZHE device has approximately a 500-mL capacity. Based on the need to add an amount of extraction fluid equal to 20 times the mass of the “solid” phase, the ZHE can therefore accommodate a maximum of 25 grams of “solid”.

**Note:** The TCLP ZHE prep uses only TCLP fluid #1; the SPLP ZHE prep uses only SPLP fluid #1.

- 11.4.20. Load the fluid transfer reservoir with an excess of Fluid #1 and preflush the transfer line to eliminate air pockets. Be sure the required volume remains.
- 11.4.21. Attach the transfer line to the liquid inlet/outlet valve and open the valve. Carefully pump the required volume into the ZHE and close the valve. Disconnect the transfer line.
- 11.4.22. Check the ZHE to make sure all the valves are closed and manually rotate the ZHE (end-over-end) 2 or 3 times. Reposition the ZHE in the vertical position.
- 11.4.23. Pressurize the ZHE to 5-10 psi. Allow to stand for 10 minutes, and then recheck the pressure. If the ZHE appears to be leaking, follow the corrective action protocols recommended by the manufacturer and repeat the analysis.
- 11.4.24. Slowly open the liquid inlet/outlet valve to bleed out any headspace that may have been introduced during the introduction of the Fluid. Upon the first sign of liquid from the valve, close the valve.
- 11.4.25. Repressurize the ZHE to 5-10 psi and place in the rotary agitator. Rotate at 28-32 rpm for 16-20 hours. Room temperature should be  $23 \pm 2$  °C. The room temperature and time should be checked at both the start and end of the extraction and recorded on Worksheet 4, lines J1, J3 and J4.
- 11.4.26. Confirm that the pressure of 5-10 psi was maintained throughout the leaching. If it was NOT maintained, return to Section 11.4.1 and repeat the leachate with a new aliquot of sample.
- 11.4.27. Attach a syringe or Tedlar bag and open the liquid inlet/outlet valve to collect the aqueous leachate and proceed as outlined in 11.4.19 - 11.4.20.

Record the volume/mass of the leachate and any oil phase on Worksheet 4, Lines J7 and J7a. Record the date and time the filtration is completed on Worksheet 4, Line J6.

**Notes:** If the waste contained an initial liquid phase, the liquid may be filtered directly into the same collection device holding the initial liquid phase of the waste.

A separate filtrate collection container must be used if combination would create multiple phases or there is not enough volume left within the filtrate collection container.

11.4.28. If the waste contained an initial filtrate (Section 11.4.18) that is miscible with the solid phase leachate (as determined in Section 11.2.8), the solid phase leachate and the initial filtrate are directly recombined in the correct proportions (see Section 12.1.4). If the individual phases are NOT compatible, they are to be collected, prepped and analyzed separately.

**Note:** Chill the filtrate and receiving vessels before recombining.

11.4.29. Following collection, store the TCLP leachate in 3 20-mL VOA vials with minimal headspace at  $4 \pm 2$  °C and prepare for analysis as soon as possible using the appropriate organic extraction procedure (see Section 16.3).

11.4.30. If the individual phases are analyzed separately, combine the results mathematically by using the recombination calculation in Section 12.1.5.

## 12. DATA ANALYSIS AND CALCULATIONS

### 12.1. Calculations

12.1.1. Calculation of Percent Wet Solids:

$$\text{Percent Wet Solids} = 100 \left( \frac{\text{Mass, "solid" phase}}{\text{Mass, initial subsample}} \right)$$

12.1.2. Calculation of weight of waste to charge to ZHE:

$$\text{Weight of waste to charge to ZHE} = 100 \left( \frac{25}{\% \text{wet solids}} \right)$$

12.1.3. Calculation of weight of extraction fluid to use:

$$\text{Weight of Extraction fluid} = \frac{20x \% \text{ wet solids} \times \text{weight of waste to be extracted}}{100}$$

- 12.1.4. Calculation of volume of initial filtrate phase to recombine with solid phase leachate:

$$\text{Volume of filtrate for recombination} = \left( \frac{\text{Weight of solids leached}}{\text{Total weight of solids}} \right) \left( \frac{\text{Leachate recovered}}{\text{Fluid added}} \right) (\text{Volume of initial aqueous filtrate})$$

- 12.1.5. Mathematical recombination of analytical results:

$$\text{Final Analyte Concentration} = \frac{(V_1 \times C_1) + (V_2 \times C_2)}{V_1 + V_2}$$

$V_1$  = total volume of the initial filtrate phase (L).

$C_1$  = analyte concentration in initial filtrate phase (mg/L).

$V_2$  = volume of the theoretical solid phase leachate (L).

$C_2$  = analyte concentration in solid phase leachate (mg/L).

## 12.2. REPORTING REQUIREMENTS

- 12.2.1. Follow these reporting conventions for multi-phase samples:

12.2.1.1. If both phases have positive results, use the values from each phase to calculate the recombined result. Use the reporting limit for each phase to calculate the recombined reporting limit.

12.2.1.2. If both phases are “ND,” not detected, the recombined result is “ND,” and the reporting limit is calculated from the reporting limit for each phase.

12.2.1.3. If one phase is “ND” and the other phase has a positive result, use the reporting limit for the “ND” phase and the positive value for the other phase to calculate the combined result. The combined reporting limit is based on the reporting limit for both phases. If the combined result is less than the combined reporting limit, then supply a footnote to indicate that “a positive result was detected below the calculated detection limit.”

- 12.2.2. Units - regardless of the nature of the sample, all TCLP and SPLP results are reported in units of mg/L.
- 12.2.3. For limits and significant figures, consult the appropriate analytical methods (Section 16.3).
- 12.2.4. Anomalies - all anomalies observed during the leach procedure must be noted on the worksheet or an anomaly form. Some examples of such anomalies are:
  - 12.2.4.1. Sample was monolithic - subsample was obtained by crushing, cutting, grinding, sawing, etc.
  - 12.2.4.2. Insufficient sample - less than the required 100 g minimum was available.
  - 12.2.4.3. Multiple phases - "X" phases were present.
  - 12.2.4.4. Sample was oil - single phase.
  - 12.2.4.5. Sample contained liquid which did not filter under test conditions.

### 12.3. REVIEW REQUIREMENTS

- 12.3.1. Review all applicable holding times. If a holding time was exceeded, confirm that a holding time violation form was properly documented and routed.
- 12.3.2. If Total analysis results are available, those results may be compared with the TCLP analysis results according to the following:

$$Total \geq 20 \times TCLP$$

**NOTE:** Assumes the sample is 100% Solids.

- 12.3.3. Total constituent analysis results can be used to demonstrate the TCLP protocol is unnecessary. In performing a TCLP analysis, there is a 20:1 dilution of the original sample with the leaching solution. Thus, if the "total constituent" result is less than 20 times the TC level, it is impossible for the leachate to "fail" and the TCLP does not need to be performed. For example, the TC level for lead is 5.0 mg/L (ppm). Therefore, if a sample of lead-contaminated soil contains less than 100 ppm total lead, a TCLP test need not be run for lead.

### **13. METHOD PERFORMANCE**

13.1. Refer to individual analysis SOPs.

13.2. Training Qualification:

The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

### **14. POLLUTION PREVENTION**

14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

### **15. WASTE MANAGEMENT**

15.1. Waste generated in this procedure must be segregated and disposed according to the facility's hazardous wastes procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

### **16. REFERENCES**

16.1. Method 1311, Toxicity Characteristic Leaching Procedure, Revision 0, July 1992, SW-846 Final Update I.

16.2. Method 1312, Synthetic Precipitation Leaching Procedure, Revision 0, November 1992, SW-846 Proposed Update II.

16.3. Related Documents

16.3.1. Toxicity Characteristic: Corrections to Final Rule. Method 1311, Federal Register, Vol. 55, No. 126, Friday, June 29, 1990.

16.3.2. Toxicity Characteristic: Final Rule. Method 1311, Federal Register, Vol. 55, No. 61, Thursday, March 29, 1990.

16.3.3. Technical Background Document and Response To Comments, Method 1311, Toxicity Characteristic Leaching Procedure, USEPA/OSW, April, 1989.

16.3.4. QA-003, STL<sup>®</sup> QC Program

16.3.5. CORP-IP-0003: Acid Digestion of Aqueous Samples by SW846 and MCAWW 200 Series Methods.

- 16.3.6. CORP-MT-0001: Inductively Coupled Plasma-Atomic Emission Spectroscopy, Spectrometric Method for Trace Element Analysis, Method 6010A and Method 200.7.
- 16.3.7. CORP-MT-0003: Graphite Furnace Atomic Absorption Spectroscopy, SW846 Method 7000A and MCAWW 200 Series Methods.
- 16.3.8. CORP-MT-0005: Preparation and Analysis of Mercury in Aqueous Samples by Cold Vapor Atomic Absorption, SW-846 7470A and MCAWW 245.1.
- 16.3.9. CORP-IP-0003: Acid Digestion of Aqueous Samples by SW846 and MCAWW 200 Series Methods.
- 16.3.10. CORP-MS-0002: Determination of Volatile Organics by GC/MS based on Methods 8240B and 8260A.
- 16.3.11. CORP-MS-0001 : GC/MS Analysis Based on Method 8270B, SW846.
- 16.3.12. CORP-GC-0001: Gas Chromatographic Analysis Based on Methods 8000A, 8010B, 8020A, 8080A and 8150B, SW846.
- 16.3.13. CORP-OP-0001: Extraction and Cleanup of Organic Compounds from Waters and Soils, Based on SW846 3500 Series, 3600 Series, 8150 and 600 Series Methods.

## 17. MISCELLANEOUS

### 17.1. Modifications/Interpretations from Reference Methods

- 17.1.1. Section 8: Preliminary Evaluations. Section 7.1 of the source method states that the sample aliquot used for the preliminary evaluation "...may not actually undergo TCLP extraction." Section 7.1.5 of the source method indicates that the portion used for the preliminary evaluation may be used for either the ZHE or non-volatile extraction if the sample was 100% solid. Section 7.1.5 further indicates that if the sample was subjected to filtration (i.e., < 100% solid) that this aliquot may be used for the non-volatile extraction procedure only as long as sufficient sample is available (minimum 100 g). Samples which have been subjected to the oven drying step may not be used for TCLP extraction because solid phase degradation may result upon heating.
- 17.1.2. Section 11.2.5.6.3: Percent Solids Determination. Section 7.1.2 of the source method indicates that "if the percent wet solids is  $\geq 0.5\%$  and it is noticed that a small amount of the filtrate is entrained in wetting of the filter"

that the filter should be oven dried to determine percent dry solids “. Drying of oil or organic matrices can both be hazardous and inappropriate. Additionally, it may be impossible to achieve a constant weight when performing this step. Due to safety concerns, if obviously oily or heavy organic matrices are entrained on the filter, the filter is not oven dried.

- 17.1.3. Section 11.2.8: Preliminary Determination of Filtrate/Extraction Fluid Compatibility. Section 7.2.13 of the source method provides no guidance as to how to make this determination. As a result, the procedure herein was developed and incorporated into the Preliminary Determinations section.
- 17.1.4. Section 9.2: TCLP Extraction Blanks. Section 8.1 of the source method states that a minimum of one blank for every 20 extractions “...that have been conducted in an extraction vessel.” STL<sup>®</sup> has interpreted this to mean one blank per twenty samples leached per TYPE of leaching vessel (i.e., Bottle or ZHE) per leach fluid used.
- 17.1.5. Section 11.2.7.9: Determination of Appropriate Extraction Fluid. Method 1311 does not address the appropriate approach to take if the pH equals 5.0. This SOP requires that Fluid #1 must be used if the pH is less than or equal to 5.0.
- 17.1.6. Section 9.4: QA/QC - Matrix Spikes. Section 8.2 of the source method states “A matrix spike shall be performed for each waste type...” and “A minimum of one matrix spike must be analyzed for each analytical batch.” Further, Section 8.2.3 of the source method also states “The purpose of the matrix spike is to monitor the performance of the analytical methods used, and to determine whether matrix interferences exist.” The standard STL<sup>®</sup> QAPP is designed to address the performance monitoring of analytical methodology through the LCS program. A minimum of one MS and MSD will be prepared for each TCLP leachate batch. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, the MS/MSD results have immediate bearing only on the specific sample spiked and not all samples in the batch.
- 17.1.7. Section 8.2.2 of the source method states that “In most cases, matrix spikes should be added at a concentration equivalent to the corresponding regulatory level.” The method also states “If the analyte concentration is less than one half the regulatory level, the spike concentration may be as low as one half of the analyte concentration but may not be less than five times the method detection limit”. For several analytes, spiking at the regulatory level is inappropriate to the range of analysis afforded by the determinative

methods. Due to the wide range in these levels, STL<sup>®</sup> spikes at the levels specified in the determinative SOPs.

17.2. Modifications from Previous SOP

None

17.3. Facility Specific SOPs

Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none. Refer to the Appendices for any facility specific information required to support this SOP.

17.4. Documentation and Record Management

The following documentation comprises a complete TCLP preparation raw data package:

- Completed worksheets (Appendix C).
- Non-conformance summary (if applicable).
- Anomaly documentation (if applicable).

**APPENDIX A**  
**TABLES**

**Table 3 - Toxicity Characteristic Analytes and Regulatory Levels (Final Rule)**

| Contaminant                                  | mg/L  |
|--|-------|
| Arsenic                                      | 5.0   |
| Barium                                       | 100.0 |
| Benzene                                      | 0.5   |
| Cadmium                                      | 1.0   |
| Carbon tetrachloride                         | 0.5   |
| Chlordane                                    | 0.03  |
| Chlorobenzene                                | 100.0 |
| Chloroform                                   | 6.0   |
| Chromium                                     | 5.0   |
| o-Cresols                                    | 200.0 |
| m-Cresols                                    | 200.0 |
| p-Cresols                                    | 200.0 |
| Total Cresols (used if isomers not resolved) | 200.0 |
| 2,4-D  | 10.0  |
| 1,4-Dichlorobenzene                          | 7.5   |
| 1,2-Dichloroethane                           | 0.5   |
| 2,4-Dinitrotoluene                           | 0.13  |
| 1,1-Dichloroethylene                         | 0.7   |
| Endrin                                       | 0.02  |
| Heptachlor (& epoxide)                       | 0.008 |
| Hexachlorobenzene                            | 0.13  |
| Hexachlorobutadiene                          | 0.5   |
| Hexachloroethane                             | 3.0   |
| Lead   | 5.0   |
| Lindane                                      | 0.4   |
| Mercury                                      | 0.2   |
| Methoxychlor                                 | 10.0  |
| Methyl ethyl ketone                          | 200.0 |
| Nitrobenzene                                 | 2.0   |
| Pentachlorophenol                            | 100.0 |
| Pyridine                                     | 5.0   |
| Selenium                                     | 1.0   |
| Silver                                       | 5.0   |
| Tetrachloroethylene                          | 0.7   |
| Toxaphene                                    | 0.5   |
| Trichloroethylene                            | 0.5   |
| 2,4,5-Trichlorophenol                        | 400.0 |
| 2,4,6-Trichlorophenol                        | 2.0   |
| 2,4,5-TP (Silvex)                            | 1.0   |
| Vinyl chloride                               | 0.2   |

**APPENDIX B**  
**FIGURES**

**Figure 1 & 2 - Rotary Agitation Apparatus and Zero Headspace Extraction Vessel (ZHE)**





Figure 1. Rotary Agitation Apparatus



**Figure 3 - US Environmental Protection Agency Memorandum #35, Page 1**



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY  
WASHINGTON, D.C. 20460

OFFICE OF  
SOLID WASTE AND EMERGENCY RESPONSE

**MEMORANDUM # 35**

**DATE:** June 12, 1992  
**SUBJECT:** Notes on RCRA Methods and QA Activities  
**From:** Gail Hansen, Chief *Gail Hansen*  
Methods Section (OS-331)

This memo addresses the following topics:

- o 1992 Symposium on Waste Testing and Quality Assurance
- o SW-846 Update
  - Final Rule for January 23, 1989 Proposed Rule
  - Notice, Proposed Rulemaking for the Second Update to the Third Edition
- o Chlorofluorocarbon 113 (CFC-113) Solvent Replacement Update
- o Environmental Monitoring Methods Index (EMMI)
- o Sampling Work Group Formation
- o MICE Update
- o Oily Waste Analysis
- o Electronic SW-846 Availability.

Figure 3 - US Environmental Protection Agency Memorandum #35, Page 10

Oily Waste Analysis

One of the most frequently asked questions on the MICE Service concerns the application of the TCLP, Method 1311, to oily wastes. Many callers request technical guidance on the extraction of oily wastes due to the difficulty in the filtration on these types of waste. In many cases, an oily waste does not filter completely due to premature clogging of the glass fiber filter. This can result in the retention of standing liquid on the glass fiber filter. Material that do not pass through the glass fiber filter at the conclusion of the filtration step is defined by the method as the solid phase of the waste. The solid phase is then subjected to the leaching procedure of the TCLP. For oily wastes, clogging of the glass fiber filter can result in an overestimation of the amount of solid material available for leaching.

To solve this problem, the Agency recommends a conservative approach, one that probably will overestimate the amount of leaching. Rather than performing the TCLP extraction on the unfiltered portion of the oily waste, assume the waste is 100% liquid (e.g., will pass through the glass fiber filter) and perform a totals analysis on the oily waste to determine if the oil exceeds the appropriate regulatory level.

Filterable waste oil generated during the TCLP must be analyzed for a variety of organic and inorganic analytes. The OSW recognizes the difficulty in achieving acceptable performance for the analysis of waste oil using methods currently provided in SW-846. As a result, the Agency will provide several new methods for the preparation and analysis of oil samples to the Organic Methods Workgroup in July. In addition, a microwave assisted digestion procedure should improve the analysis of metals and will be proposed as part of the Second Update of the Third Edition of SW-846. Brief descriptions of these techniques are provided below, for additional information on the organic procedures contact Barry Lesnik at (202) 260-7459. For additional information on microwave digestion contact Ollie Fordham (202) 260-4778.

The use of purge-and-trap (Method 5030) for volatiles in oil generally results in severe contamination of analytical instrumentation. Traps, transfer lines and chromatography columns may become contaminated with oil. This leads to elevated baselines, hydrocarbon background in subsequent analyses, and cross-contamination. Headspace (Method 3810) is currently allowed only as a screening procedure in SW-846. The Agency is evaluating the use of headspace in conjunction with isotope dilution mass spectrometry for the quantitative analysis of volatiles in oil. Headspace reduces interference problems encountered with purge-and-trap. However, headspace quantitation can be questionable because the distribution of analytes is not

**APPENDIX C**  
**WORKSHEETS**

## APPENDIX C - WORKSHEETS

**TCLP/SPLP Worksheet 1 - Sample Description**

|   |  |  |  |  |  |
|---|--|--|--|--|--|
| Lot Number                                    |  |  |  |  |  |
| Laboratory Sample No.                         |  |  |  |  |  |
| <b>A. Sample Description</b>                  |  |  |  |  |  |
| Number of phases                              |  |  |  |  |  |
| 1. Solid                                      |  |  |  |  |  |
| 2. Liquid                                     |  |  |  |  |  |
| a. lighter than water                         |  |  |  |  |  |
| b. water                                      |  |  |  |  |  |
| c. heavier than water                         |  |  |  |  |  |
| <b>B. Percent Solid Phase</b>                 |  |  |  |  |  |
| 1. weight of filter (g)                       |  |  |  |  |  |
| 2. weight of subsample                        |  |  |  |  |  |
| a. gross weight (g)                           |  |  |  |  |  |
| b. tare weight (g)                            |  |  |  |  |  |
| c. net weight (g)                             |  |  |  |  |  |
| 3. weight of filtrate                         |  |  |  |  |  |
| a. gross weight (g)                           |  |  |  |  |  |
| b. tare weight (g)                            |  |  |  |  |  |
| c. net weight (g)                             |  |  |  |  |  |
| 4. total weight wet solids <sup>0</sup>       |  |  |  |  |  |
| 5. weight percent wet solids <sup>1</sup> (%) |  |  |  |  |  |
| 6. weight percent dry solids (%)              |  |  |  |  |  |
| a. gross dry weight (g)                       |  |  |  |  |  |
| b. final gross dry weight (g)                 |  |  |  |  |  |
| c. percent dry solids <sup>2</sup> (%)        |  |  |  |  |  |
| 7. volume of initial aqueous filtrate (mL)    |  |  |  |  |  |
| 8. volume of initial organic filtrate (mL)    |  |  |  |  |  |

$$^0 \text{ Total Weight of wet solids} = (B2c - B3c)$$

$$^1 \text{ Weight percent wet solids} = 100 \left( \frac{\text{Total weight wet solids, B4}}{\text{Weight of subsample, B2c}} \right)$$

$$^2 \text{ Weight percent dry solids} = 100 \left( \frac{\left( \text{Gross dry weight, B6b} \right) - \left( \text{weight of filter, B1} \right)}{\text{Weight of subsample, B2c}} \right)$$

Comments: \_\_\_\_\_

Analyst: \_\_\_\_\_ Date: \_\_\_\_\_

## APPENDIX C - WORKSHEETS

**TCLP/SPLP Worksheet 2 - Selection of Extraction Fluid****TCLP**

|  |  |  |  |  |  |
|--|--|--|--|--|--|
| Lot Number.  |  |  |  |  |  |
| Laboratory Sample No.  |  |  |  |  |  |
| <b>C. Extraction Fluid Determination—does not apply to determination of volatile organic components.</b> |  |  |  |  |  |
| 1. particle size reduction?<br>yes/no (<1 mm)  |  |  |  |  |  |
| 2. sample weight, ✓ if $5.0 \pm 0.1$ g   |  |  |  |  |  |
| 3. volume of water, ✓ if $96.5 \pm 1.0$ mL   |  |  |  |  |  |
| 4. initial pH (after 5 min. mixing time)   |  |  |  |  |  |
| 5. if pH >5.0, ✓ if 3.5 mL 1N HCL  |  |  |  |  |  |
| 6. ✓ if heated and held at 50°C for ten minutes  |  |  |  |  |  |
| 7. secondary pH (at room temperature)  |  |  |  |  |  |
| <b>D. Selection of Extraction Fluid</b>  |  |  |  |  |  |
| 1. ✓ if pH from § C (4) or § C (7) is <5.0, use extraction fluid No. 1.                                  |  |  |  |  |  |
| 2. ✓ if pH from § C (7) IS >5.0, use extraction fluid No. 2  |  |  |  |  |  |
| 3. ✓ if filtrate and fluid are miscible  |  |  |  |  |  |

**SPLP**

|   |  |  |  |  |  |
|---|--|--|--|--|--|
| Laboratory Sample No.   |  |  |  |  |  |
| Field Sample No.  |  |  |  |  |  |
| <b>D. Selection of Extraction Fluid (✓ one)</b>                       |  |  |  |  |  |
| Fluid 1: Soils—East of the Mississippi River; Wastes; or Wastewaters. |  |  |  |  |  |
| Fluid 2: Soils—West of Mississippi River                              |  |  |  |  |  |
| Fluid 3: If VOCs or Cyanide containing wastes.                        |  |  |  |  |  |
| 1. ✓ if filtrate and fluid are miscible                               |  |  |  |  |  |

Comments: \_\_\_\_\_  
 \_\_\_\_\_

Analyst: \_\_\_\_\_ Date: \_\_\_\_\_



APPENDIX C - WORKSHEETS

**TCLP/SPLP Worksheet 4 - Zero Headspace Extraction (ZHE)**

|   |  |  |  |  |  |
|---|--|--|--|--|--|
| Lot Number  |  |  |  |  |  |
| Laboratory Sample No.   |  |  |  |  |  |
| <b>H. Determination of Sample Size for Leach Testing—maximum 25 grams</b>         |  |  |  |  |  |
| 1. particle size reduction?<br>yes/no (<9.5 mm)                                   |  |  |  |  |  |
| 2. weight of wet solid (g)  |  |  |  |  |  |
| <b>I. Determination of Amount of Extraction Fluid No. 1</b>                       |  |  |  |  |  |
| 1. weight of material added to ZHE (g)  |  |  |  |  |  |
| 2. volume of filtrate in syringe (mL)   |  |  |  |  |  |
| 3. weight of filtrate in syringe  |  |  |  |  |  |
| a. gross weight (g)   |  |  |  |  |  |
| b. tare weight (g)  |  |  |  |  |  |
| c. net weight (g)   |  |  |  |  |  |
| 4. wet solids in ZHE <sup>4</sup> (g)   |  |  |  |  |  |
| 5. weight of fluid to add <sup>5</sup> (g)  |  |  |  |  |  |
| <b>J. Record of ZHE Extraction Text—the extraction period is as 18 ± 2 hours.</b> |  |  |  |  |  |
| 1. room temperature initial (°C)<br>final (°C)                                    |  |  |  |  |  |
| 2. ZHE vessel number  |  |  |  |  |  |
| 3. extraction start date & time   |  |  |  |  |  |
| 4. extraction stop date & time  |  |  |  |  |  |
| 5. ✓ if still under positive pressure   |  |  |  |  |  |
| 6. filtration completed date & time   |  |  |  |  |  |
| 7. volume of filtrate recovered after leaching (mL)                               |  |  |  |  |  |
| a. volume of oil recovered after leaching (mL)                                    |  |  |  |  |  |
| 8. volume of initial aqueous filtrate for recombination <sup>6</sup> (mL)         |  |  |  |  |  |
| 9. combined initial aqueous + leachate (mL)                                       |  |  |  |  |  |

$$^4 \text{Solids remaining in ZHE} = (\text{Material added, I1}) - (\text{Weight of filtrate in syringe, I3}) \qquad ^5 \text{Weight of fluid to add} = 20(\text{Wet solids in ZHE, I4})$$

$$^6 \text{Volume of filtrate used for recombination} = \left( \frac{\text{Weight of solids leached, I4}}{\text{Total weight of solids, B4}} \right) \left( \frac{\text{Leachate recovered, J7}}{\text{Fluid added, I5}} \right) (\text{Volume initial aqueous filtrate, B7})$$

Comments: \_\_\_\_\_  
\_\_\_\_\_

Analyst: \_\_\_\_\_ Date: \_\_\_\_\_

**TCLP/SPLP Worksheet 5 - Organic Results**






**APPENDIX D**  
**FLOW CHARTS**


APPENDIX D - FLOW CHARTS

**Flow Chart 1. Preliminary Sample Evaluation  
(Section 11.2)**




## APPENDIX D - FLOW CHARTS

**Flow Chart 2. Bottle Extraction, Non-Volatile Constituents  
(Section 11.3)**




## APPENDIX D - FLOW CHARTS

**Flow Chart 3. ZHE Extraction, Volatile Constituents  
(Section 11.4)**



**Flow Chart 3. ZHE Extraction  
(Continued)**



APPENDIX D – Bench Sheets



**TCLP (Method 1311)**

| Analyst   | Room Temp.          | Date                | Checked By | Date                      | pH Instrument             | Probe     | Balance         |  |            |                      |                          |                         |
|---|---------------------|---------------------|------------|---------------------------|---------------------------|-----------|-----------------|--|------------|----------------------|--------------------------|-------------------------|
| Lot #   | Client Codes        |                     |            | Solution # 1 - Log book # | Solution # 2 - Log book # |           |                 |  |            |                      |                          |                         |
| Sample ID   | Fluid Determination |                     |            | TCLP Extraction           |                           |           |                 | Conc. HNO <sub>3</sub> to Acidity (for Metals) |            |                      |                          |                         |
|   | Wgt/Vol             | Init. pH            | Final pH   | Tumbler RPM               | Extract fluid #1 or #2    | Sample ID | Wgt (g) Wet/Dry |  | Fluid Vol. | Vessel glass/plastic | Start time of filtration | Final pH after Tumbling |
| 1.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 2.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 3.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 4.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 5.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 6.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 7.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 8.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 9.  |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 10.   |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 11.   |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 12.   |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 13.   |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 14.   |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| 15.   |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| Extract(s)  |                     | Extract(s) Received |            | Extract(s) Retinquished   |                           |           |                 |  |            |                      |                          |                         |
| (Record line number from above)   |                     | Date                | Time       | Date                      | Time                      | Date      | Time            | Date   | Time       | Analyst              | Analyst                  | Location                |
|   |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| <p>Comments</p> <p>* - Sample determined to have free liquid, % solid determination was performed.</p> <p>&lt;5 = Extraction fluid # 1 5.7mL Glacial acetic acid dil 500mL + 64.3mL HaOH dil to 1L (pH 4.93 ± 0.05)</p> <p>&gt;5 = Extraction fluid #2 5.7mL Glacial acetic acid dil to 1L (pH 2.88 ± 0.05)</p> |                     |                     |            |                           |                           |           |                 |  |            |                      |                          |                         |
| Date/Time On  |                     |                     |            | Date/Time Off             |                           |           |                 | Agitation Apparatus RPM's 30±2 checked by      |            |                      |                          |                         |
|   |                     |                     |            |                           |                           |           |                 | (7.9) (10.0)                                   |            |                      |                          |                         |
| Date  |                     |                     |            | Date                      |                           |           |                 | pH calibration (4.0) (7.9)                     |            |                      |                          |                         |
| Date  |                     |                     |            | Date                      |                           |           |                 | Filter Lot #                                   |            |                      |                          |                         |



APPENDIX D – Bench Sheets

**STL - Pittsburgh**      **(SPLP) LEACHATE LOGBOOK**      **Logbook ID: OP91**

| Analyst                    | Date     | Received by                  | Date                       | pH instrument                | Probe        | Balance                |                         |  |          |
|----------------------------|----------|------------------------------|----------------------------|------------------------------|--------------|------------------------|-------------------------|--|----------|
| Lot Number(s)              |          | Solution #1 - Logbook Number |                            | Solution #2 - Logbook Number |              |                        |                         |  |          |
| Work Order #               | Comments | Int pH                       | Extract Fluid #1, #2 or #3 | Weight (gm) Wt. / Dry        | Fluid Volume | Vessel Glass / Plastic | Final pH after sampling | Conc HNO <sub>3</sub> to Acidity (or Metals) |          |
| 1.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 2.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 3.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 4.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 5.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 6.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 7.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 8.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 9.                         |          |                              |                            |                              |              |                        |                         |  |          |
| 10.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 11.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 12.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 13.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 14.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 15.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 16.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 17.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 18.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 19.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 20.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 21.                        |          |                              |                            |                              |              |                        |                         |  |          |
| 22.                        |          |                              |                            |                              |              |                        |                         |  |          |
| Extract(s)                 | Data     | Time                         | Received                   | Analyst                      | Location     | Date                   | Time                    | Analyst                                      | Location |
| (record time # from above) |          |                              |                            |                              |              |                        |                         |  |          |

\* = Sample determined to have free liquid, % solid determination was performed  
 Extraction Fluid 1 = Samples that are EAST of the Mississippi River. Prepared by adding an aliquot of a 60/40 % by weight mixture of H<sub>2</sub>SO<sub>4</sub> and HNO<sub>3</sub> to reagent water (pH 4.20 ± 0.05)  
 Extraction Fluid 2 = Samples that are WEST of the Mississippi River. Prepared by adding an aliquot of a 60/40 % by weight mixture of H<sub>2</sub>SO<sub>4</sub> and HNO<sub>3</sub> to reagent water (pH 5.00 ± 0.05)  
 Extraction Fluid 3 = This is a reagent water and is used to determine cyanide leachability

|              |               |                     |
|--------------|---------------|---------------------|
| Date/Time on | Date/Time off | Agitation Apparatus |
|              | (4.0)         | RPVA 30 + 2         |
| Date         | (7.0)         | (10.0)              |
| Date         |               |                     |