

Q4.5-5.0PB-4.5

DISINFECTANT - DEODORIZER - FUNGICIDE - VIRUCIDE**
WHIRLPOOL CLEANER DESIGNED FOR HOSPITAL,
NURSING HOME, AND INSTITUTIONAL USE

Active Ingredients:

n-Alkyl (60% C_{14} , 30% C_{16} , 5% C_{12} , 5% C_{18}) dimethyl benzyl ammonium chlorides.....	2.25%
n-Alkyl (68% C_{12} , 32% C_{14}) dimethyl ethylbenzyl ammonium chlorides.....	2.25%
Inert Ingredients:.....	95.50%

KEEP OUT OF REACH OF CHILDREN

DANGER

SEE LEFT PANEL FOR ADDITIONAL
PRECAUTIONARY STATEMENTS



CENTRAL SOLUTIONS, INC.

3130 BRINKERHOFF ROAD
KANSAS CITY, KANSAS 66115

EPA REG. NO. 211-40
EPA EST. NO. 211-KS-1

NET CONTENTS

GAL.

ACCEPTED
NOV 22 1994

Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 211-40

(Left Panel)

Precautionary Statements
Hazards to Humans and Domestic Animals

DANGER

Keep Out of Reach of Children. Corrosive. Causes severe eye and skin damage. Do not get in eyes, on skin, or on clothing. Harmful if swallowed. Wear goggles or face shield and rubber gloves when handling. Avoid contamination of food. Remove and wash contaminated clothing.

STATEMENT OF PRACTICAL TREATMENT

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink large quantities of water. Call a physician.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression, and convulsion may be needed.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage or disposal.

STORAGE - Store in securely closed original container. Store in area that will prevent cross contamination of other pesticides, fertilizer, food, and feed. Avoid storage at temperature extremes or in sunlight.

PESTICIDE DISPOSAL - Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

CONTAINER DISPOSAL - Triple rinse (or equivalent), then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration, or if allowed by state and local authorities, by burning. If burned, stay out of smoke.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

ACCEPTED
NOV 22 1994
Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under 211-40
EPA Reg. No.

BEST AVAILABLE COPY

PRODUCT COMPOSITION

Q4.5-5.OPB-4.5 is a disinfectant, deodorizer, fungicide, and virucide* cleaner designed for hospital, nursing home, and institutional use. This formulation is designed to clean and disinfect inanimate hard surfaces such as Central Solutions' fiberglass, stainless steel, chrome plated and/or acrylic bathing and physical therapy equipment such as whirlpools and patient transfer equipment.

Q4.5-5.OPB-4.5 contains two highly active twin chain quaternary ammonium compounds (listed as active ingredients) for bactericidal efficacy. Inert ingredients include a balance of alkaline builders, non-ionic surfactants, and odor counteragents to enhance the cleaning and odor control characteristics. This formulation is designed to clean and disinfect bathing and therapy equipment. When used as directed, **Q4.5-5.OPB-4.5** has been found effective for use on all fiberglass, stainless steel, and chrome plated surfaces of Central Solution's bath, therapy, and patient transfer equipment.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

DISINFECTION: To disinfect hard, non-porous, inanimate surfaces (such as fiberglass and stainless steel tubs and chair surfaces, chrome plated intakes and lifts, etc.), apply properly diluted **Q4.5-5.OPB-4.5** with mop, cloth sponge, or hand pump trigger sprayer so as to wet all surfaces thoroughly. For routine disinfection, proper dilution is 2 ounces of **Q4.5-5.OPB-4.5** per gallon of water (1:64). Allow to remain wet for 10 minutes and then let air dry. For heavily soiled areas, a pre-cleaning step is required. Prepare a fresh solution for each use.

FOR CLEANING CENTRAL SOLUTIONS' BATH AND THERAPY EQUIPMENT: To remove body oils, dead tissue, hard water, and all other buildups of organic matter on inanimate surfaces after using the whirlpool unit, drain the water and refill with fresh water to just cover the intake valve. Add 2 ounces of **Q4.5-5.OPB-4.5** for each gallon of water in the unit at this point. Briefly start the pump to circulate the solutions. Turn off pump. Wash down the unit sides, seat of the chairlift, and any/all related equipment with a clean swab or sponge. After the unit has been thoroughly cleaned, drain solutions from the unit and rinse any/all cleaned surfaces with fresh water. The unit is ready for re-use. Product to surface contact time must be at least 10 minutes for proper disinfection.

When used as directed, **Q4.5-5.OPB-4.5** has been found effective for use on all fiberglass, stainless steel, and chrome plated surfaces of bath, therapy, and patient transfer equipment.

ACCEPTED
NOV 22 1994
U.S. Environmental Protection Agency
EPA Office of Pesticide Regulation
Office of Registration Services
EPA Reg. No. 211-40

BEST AVAILABLE COPY

PERFORMANCE DATA

Efficacy tests have demonstrated that Q4.5-5.0PB-4.5 is an effective bactericide, fungicide, and virucide* on hard, non-porous, inanimate surfaces in the presence of moderate amounts of organic soil (5% blood serum). All tests have been performed using A.O.A.C. protocol testing as required by the EPA at a use-dilution of 1:64 dilution (two ounces per gallon of water) in the presence of 5% blood serum (used to simulate organic soil load), at 20 degrees Celsius with an exposure time of 10 minutes.

GERMICIDAL

The broad spectrum activity of Q4.5-5.0PB-4.5 has been confirmed against the following organisms:

- | | |
|--|-------------------------|
| Pseudomonas aeruginosa | Staphylococcus aureus |
| Escherichia coli | Klebsiella pneumoniae |
| Serratia marcescens | Salmonella choleraesuis |
| Proteus vulgaris | Candida albicans |
| Methicillin resistant Staphylococcus aureus (MRSA) | |

FUNGICIDAL

Q4.5-5.0PB-4.5's efficacy as a fungicide has been confirmed through testing using the AOAC Fungicidal Test against the pathogenic fungi, Trichophyton mentagrophytes.

VIRUCIDAL

*Testing of Q4.5-5.0PB-4.5 has confirmed its virucidal activity against many lipophilic viruses including:

- Influenza A2-Asian.....(representative of the common flu virus)
- Herpes simplex Type 1...(causative agent of fever blisters)
- Adeno virus type 2.....(causative agent of upper respiratory infections)
- Vaccinia virus.....(representative of the pox virus group)

KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of inanimate surfaces/objects likely to be soiled with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type-1 (HIV-1) (associated with AIDS). Contact time ten (10) minutes.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST AIDS (HIV-1 VIRUS) ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS. with ten (10) minutes Contact time.

PERSONAL PROTECTION: Disposable gloves should be worn when cleaning blood, body fluids, or contaminated items. Hands should be washed immediately and thoroughly after gloves are removed, or if hands are contaminated with blood or body fluids. Masks and protective eyewear or faceshields should be worn during cleaning procedures involving blood or body fluids. Gowns or aprons should be worn during cleaning procedures for protection against splashes of blood or other body fluids.

CLEANING PROCEDURE: Blood and other fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant. After cleaning is complete, follow disinfection described previously.

DISPOSAL OF INFECTIOUS MATERIALS: Contaminated waste and blood and other body fluids should be autoclaved and disposed according to local regulations for infectious waste disposal.

ACCEPTED

NOV 22 1994

U.S. Department of Health and Human Services
 Center for Disease Control and Prevention
 Division of Field Epidemiology
 as required for the pesticide
 registered under
 EPA Reg. No. 211-40

BEST AVAILABLE COPY