

PM 30 6:15 15

12/4

**FRONT LABEL**

**PART OF YOUR INFECTION CONTROL PROGRAM**

**LEHN & FINK  
INSTRUMENT GERMICIDE**

**Instrument Disinfectant Presoak  
(Contains Rust Inhibitor)**

**Phenolic-Based • Ready-to-Use  
Instrument Disinfectant Presoak**

|                    |   |
|--------------------|---|
| <b>GERMICIDAL*</b> | Staphylocidal • Pseudomonacidal • Tuberculocidal  |
| <b>FUNGICIDAL*</b> | Against pathogenic fungi  |
| <b>VIRUCIDAL*</b>  | Against Influenza A <sub>2</sub> (Japan), Herpes Simplex Type 1,<br>Adenovirus Type 2, Vaccinia, and HIV-1 (AIDS Virus) |

\*On environmental surfaces.

|                            |                              |         |
|----------------------------|------------------------------|---------|
| <b>ACTIVE INGREDIENTS:</b> | o-Phenylphenol.....          | 0.518%  |
|                            | o-Benzyl-p-chlorophenol..... | 0.250%  |
| <b>INERT INGREDIENTS:</b>  | .....                        | 99.232% |

Contains sodium nitrate.

**KEEP OUT OF REACH OF CHILDREN**

**CAUTION:**

Harmful if swallowed. Avoid contact with skin and eyes. Avoid contamination of food. See side/back panel for statement of practical treatment and additional precautionary statements.

EPA Reg. No. 675-15  
EPA Est. No. 34797-IA-03  
EPA Est. No. 777-OH-1

**NATIONAL LABORATORIES (logo)  
L&F Products, Inc • Montvale, NJ 07645**

**Net Contents:  
1 U.S. Gallon (3.79 Liters)**

|  |
|--|
| <b>ACCEPTED</b>  |
| <b>DEC 21 1994</b>   |
| <small>Under the Federal Insecticide, Fungicide, and<br/>Rodenticide Act, as amended, for the<br/>pesticide registered under<br/>EPA Reg. No. 675-15</small> |

**HOW TO USE LEHN & FINK INSTRUMENT GERMICIDE**

**DIRECTIONS FOR USE:** It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Lehn & Fink Instrument Germicide is intended to help prevent protein matter from hardening on instruments and to help reduce exposure to pathogens during the decontamination process.

This product is not to be used as a terminal sterilant/high level disinfectant on a surface or instrument that (1) is introduced directly into the human body, either or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical devices prior to sterilization or high level disinfection.

Lehn & Fink Instrument Germicide is recommended for use in hospitals, nursing homes, dental and doctors' offices, and other institutional facilities.

**CLEANING SURGICAL AND DENTAL INSTRUMENTS:** Instruments must be cleaned prior to disinfection to remove all adhering blood and serous exudates, body fluids, and human wastes.

**DISINFECTING SURGICAL AND DENTAL INSTRUMENTS:** Immerse instruments in Lehn & Fink Instrument Germicide for 15 minutes to wet all surfaces thoroughly. Flush and fill the lumen of hollow instruments/equipment, avoiding entrapment of bubbles. If instruments are to remain in Lehn & Fink Instrument Germicide for substantially longer periods of time, add 1 gram of sodium bicarbonate to each liter of solution to maintain clarity; stir to dissolve completely. Fresh solution must be prepared at least daily or more often if the solution becomes diluted or soiled.

**FOR TUBERCULOCIDAL EFFECTIVENESS:** Thoroughly clean instruments to remove excess organic debris. Rinse and rough dry before immersing in Lehn & Fink Instrument Germicide for 10 minutes at 20°C (68°F).

**HYPODERMIC NEEDLES AND SYRINGES:** Flush with Lehn & Fink Instrument Germicide, immerse for 15 minutes, then flush with sterile distilled water.

**RUBBER GLOVES, TUBES AND OTHER RUBBER ARTICLES:** Immerse for 15 minutes. Remove, then immerse in 70% alcohol for one minute.

If rapid drying is required, immerse for a few seconds in alcohol after removing from the Lehn & Fink Instrument Germicide.

**NOTE:** Before using Lehn & Fink Instrument Germicide, consult the instrument manufacturer's instructions for care and disinfection of their instrument.

Lehn & Fink Instrument Germicide may cause damage to certain types of plastic trays. Refer to manufacturer's instructions for proper use.

**KILLS HUMAN IMMUNODEFICIENCY VIRUS TYPE 1 (HIV-1) (AIDS VIRUS) ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of (HIV-1) (associated with AIDS).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.**

**PERSONAL PROTECTION:** When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks or eye coverings.

**CLEANING PROCEDURES:** Blood or other body fluids must be thoroughly cleaned from surfaces and objects before application of Lehn & Fink Instrument Germicide.

**CONTACT TIME:** Allow surface to remain wet for 10 minutes.

**DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

LEHN & FINK

**PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

**KEEP OUT OF REACH OF CHILDREN**

**CAUTION:** Harmful if swallowed. May cause eye and skin irritation. Avoid contact with eyes or skin. Avoid contamination of food.

**STATEMENT OF PRACTICAL TREATMENT:**

**In case of eye contact, IMMEDIATELY** flush eyes thoroughly with plenty of water. Remove any contact lenses and continue to flush eyes with plenty of water for at least 15 minutes. Get medical attention. **If on skin**, wash with plenty of soap and water. If irritation persists, get medical attention. **If swallowed**, do not induce vomiting. Drink a large amount of water. Call a physician.

**STORAGE AND DISPOSAL**

**STORAGE:** Store in original container in areas inaccessible to small children. Keep securely closed.

**DISPOSAL:** Do not reuse empty container. Wrap and discard in trash.

Please encourage your state and local authorities to establish a recycling program so that this bottle can be recycled.

Lehn & Fink Instrument Germicide is part of a system of infection control products provided by National Laboratories. People respect the LYSOL Brand name for providing high quality and trusted surface disinfection, and gain the added benefit of a clean, sanitary environment. We invite you to call us at 1-800-677-9218, Monday through Friday, 8:45 a.m. to 4:30 p.m. (Eastern time), if you have any questions or comments on our products.

