

To increase detergency: Increase solution strength to 1%-2%.

For heat sterilization: Boil for 5 minutes in 1/2%-1% solution.  
Boiling in 1% STAPHENE destroys B. cereus spores in less than 5 minutes.

For cold instrument disinfection: Soak for 15 minutes in 1/2% Staphene.  
Increase to 1%-2% if instruments are soiled.

When used as directed, STAPHENE is harmless to fabrics, resilient flooring, tile, terrazzo, painted and varnished surfaces.



CAUTION: KEEP OUT OF REACH OF CHILDREN. Harmful if swallowed. Avoid prolonged skin contact and avoid splashing in eyes. Wash eyes with large amount of flowing water if exposed.

Manufactured By

VESTAL LABORATORIES

Division of W. R. Grace & Co.

St. Louis, Mo.

**VESTAL**

**VES-PHENE-O<sup>®</sup>**  
germicidal detergent

**ACTIVE INGREDIENTS:** o-benzyl-p-chlorophenol 4.96%; o-phenylphenol 4.62%; p-tertiary amyphenol 4.46%; 2,2 methyl-enebis (3,4,6-trichlorophenol) 0.40%.

**INERT INGREDIENTS:** 85.56%.

<b>CONTENTS:</b> <b>ONE U. S. GALLON</b>	<b>CLEANS</b> <b>DISINFECTS</b> <b>DEODORIZES</b>
---	---

**DANGER: KEEP OUT OF REACH OF CHILDREN**  
See back panel for additional precautions

6743

**VES-P**

E.P.A. REGISTR

**GERM**  
**TUBERC**  
**FUNC**  
**VIRU**  
**PSEUDOM**

VES-PHENE-O is a concentrated phenolic germicidal detergent system which contains a special grade of biodegradable vegetable oil soap as its cleansing agent. It shows a broad range of germicidal effectiveness against pathogenic bacteria, fungi and many lipophilic viruses including S. aureus, S. typhosa, Ps. aeruginosa, M. tuberculosis, T. interdigitale, E. floccosum, and Asian Influenza Virus. In addition to its germicidal properties VES-PHENE-O is an excellent multi-purpose cleaner, particularly adapted to maintaining the patina of terrazzo surfaces. VES-PHENE-O is a free rinsing product that does not leave an insulating film and is recommended for use on conductive flooring in hospital operating rooms.


Read carefully the following the contents of this container.

**DANGER:** Corrosive to eye and skin damage. If not get in eyes, on skin or rubber gloves when handling. Do not reuse empty container.

**FIRST AID:** If in eyes, flush with water for 15 minutes. If on skin, wash with plenty of water. If inhaled, get fresh air and drink water. If swallowed, do not induce vomiting.

M. M.

VESTAL LABORATORIES • Division of C



**ACCEPTED**

Nov 9 1971

UNDER THE FEDERAL INSECTICIDE  
FUNGICIDE AND RODENTICIDE ACT  
FOR ECONOMIC PESTICIDE REGISTER-  
ED UNDER NO. 1049-10. SUBJECT  
TO APPLICABLE REGULATIONS.

**VESTAL**

# VES-PHENE-O<sup>®</sup>

## germicidal detergent

ACTIVE INGREDIENTS: o-benzyl-p-chlorophenol 4.96%; o-phenylphenol 4.62%; p-tertiary amylphenol 4.46%; 2,2 methyl-enebis (3,4,6-trichlorophenol) 0.40%.

INERT INGREDIENTS: 85.56%.

<b>CONTENTS:</b> <b>ONE U. S. GALLON</b>	<b>CLEANS</b> <b>DISINFECTS</b> <b>DEODORIZES</b>
---	---

**DANGER: KEEP OUT OF REACH OF CHILDREN**  
See back panel for additional precautions

6393

VES-PHENE-O is a concentrated phenolic germicidal detergent system which contains a special grade of biodegradable vegetable oil soap as its cleansing agent. It shows a broad range of germicidal effectiveness against pathogenic bacteria, fungi and many lipophilic viruses including S. aureus, S. typhosa, Ps. aeruginosa, M. tuberculosis, T. interdigitale, E. floccosum, and Asian Influenza Virus. In addition to its germicidal properties VES-PHENE-O is an excellent multi-purpose cleaner, particularly adapted to maintaining the patina of terrazzo surfaces. VES-PHENE-O is a free rinsing product that does not leave an insulating film and is recommended for use on conductive flooring in hospital operating rooms.

**VES-PH**  
E.P.A. REGISTRAT  
**GERMI**  
**TUBERCU**  
**FUNGI**  
**VIRUS**  
**PSEUDOM**

Read carefully the following the contents of this container

**DANGER:** Corrosive to tissue, eye and skin damage. Harmful if not get in eyes, on skin, or on rubber gloves when handling. Do not reuse empty container.

**FIRST AID:** In case of contact with eyes or skin with plenty of water. If in eyes, call a physician. Remove clothing before reuse. If swallowed, drink large quantities of water immediately.

Manufa

VESTAL LABORATORIES • Division of Ch