

**Precautionary Statements  
HAZARD TO HUMANS AND DOMESTIC ANIMALS**

**CAUTION**

Harmful if swallowed. Avoid prolonged skin contact and avoid splashing in eyes. Wash eyes with large amounts of flowing water if exposed.

**DIRECTIONS FOR USE**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

All instruments and devices must be thoroughly cleaned to remove excess organic debris, rinsed and then rough dried before immersion. Clean and rinse the lumens of hollow instruments before filling with solution. Use a 1:32 dilution with water for cold disinfection of precleaned equipment. The removal of heavy soils and the cleaning of all surfaces prior to application of this product is required.

**DILUTION TABLE**


QUANTITY OF AMERSE INSTRUMENT GERMICIDE TO MAKE A 1:32 AQUEOUS USE-DILUTION				
Water	1 Pint	1 Quart	1/2 Gallon	1 Gallon
AMERSE	(1 stroke) 1/2 fl. oz.	(2 strokes) 1 fl. oz.	(4 strokes) 2 fl. oz.	(8 strokes) 4 fl. oz.

To remove residue of previous disinfectant, thoroughly clean and rinse instruments, trays or containers prior to first use of this product. Prepare a new solution daily or when solution turns cloudy or dirty for disinfecting precleaned items.

**FOR COLD INSTRUMENT DISINFECTION:** After cleaning, immerse articles and wet all surfaces thoroughly in a 1:32 solution for 10 minutes. Rinse and dry.

**FOR METAL INSTRUMENT STORAGE—**A 1:32 solution is non-corrosive and does not rust, dull or otherwise stain or attack quality instruments such as carbide dental burs. The solution may be used for prolonged storage of all but lens instruments, plastic or rubber items. This is a complete product. Do not add other chemicals. Use only as directed. If frozen, thaw and remix before use.

**ACCEPTED**  
 JUN 6 1969  
 Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 1043-36



**BEST AVAILABLE COPY**