

CENTER PANEL

amerse

Hospital Instrument Germicide

EPA Reg. No. 1043-36
EPA Est. No. 1043-MO-1

Germicide • Tuberculocide • Fungicide
FOR DISINFECTION AND STORAGE OF SURGICAL
AND DENTAL INSTRUMENTS

Active ingredients:

o-Benzyl-o-chlorophenol	5.25%
o-Phenylphenol	100%

Inert ingredients 93.75%

KEEP OUT OF THE REACH OF CHILDREN


CAUTION

STATEMENT OF PRACTICAL TREATMENT

In case of eye contact, flush with water for 15 minutes. Get prompt medical attention.
For skin contact, flush with water. See back panel for additional precautionary statements and directions for use.

Net Contents: 64 fl. oz. (1/2 U.S. Gallon)

ACCEPTED
NOV 06 1989
Under the Federal Insecticide
Fungicide and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 1043-36



[Faint, illegible text at the bottom left of the page]

FILE COPY

Amerse Instrument Germicide is a concentrated disinfectant to be used as a 1:32 dilution with water for cold disinfection of precleaned equipment, thermometers, surgical and dental instruments. Recommended for use in hospitals and dental offices or other facilities requiring disinfection of instruments and equipment. The 1:32 dilution kills pathogenic bacteria and fungi such as *S. aureus*, *S. choleraesuis*, *Ps. aeruginosa*, *M. tuberculosis* (effective in 10 minutes at 20°C), and *C. albicans*. Amerse is virucidal against Influenza A2 (Hong Kong), Herpes simplex I and II, Type 2 Adenovirus and Vaccinia on hard, inanimate, surfaces in 10 minutes at 20°C.

Minimum Use-Dilution Confirmation (A.O.A.C.-10 Min.-20°C) *S. choleraesuis*-1:32, *S. aureus*-1:32, *Ps. aeruginosa*-1:32

A 1:32 dilution is effective against *M. tuberculosis* on hard, inanimate surfaces in 10 minutes; at 20°C. Remove heavy soil or gross filth and thoroughly clean surfaces.

STORAGE AND DISPOSAL

Prohibitions: Do not contaminate water, food, or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container.

Pesticide Disposal: This germicide, its solutions or rinsings from empty containers, should be disposed of in a toilet or service sink served by a sanitary sewer or in a landfill approved for pesticides.

Container Disposal: Triple rinse and deposit in a waste container for incineration or burial in a land fill approved for pesticide containers.

General: Consult Federal, State, or Local disposal authorities for approved alternative procedures such as limited open burning.

6401-06A



DIVISION OF CALGON CORPORATION
SUBSIDIARY OF MERCK & CO. INC.
St. Louis, Missouri 63123

ACCEPTED
NOV 06 1969
Under the Federal Insecticide,
Fungicide, and Rodenticide Act
as amended, for the pesticide
registered under
EPA Reg. No. 1047-36

9530

ALL COPY

**Precautionary Statements
HAZARD TO HUMANS AND DOMESTIC ANIMALS**

CAUTION

Harmful if swallowed. Avoid prolonged skin contact and avoid splashing in eyes. Wash eyes with large amounts of flowing water if exposed.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

All instruments and devices must be thoroughly cleaned to remove excess organic debris, rinsed and then rough dried before immersion. Clean and rinse the lumens of hollow instruments before filling with solution.

Use a 1:32 dilution with water for cold disinfection of precleaned equipment. The removal of heavy soils and the cleaning of all surfaces prior to application of this product is required.

DILUTION TABLE


QUANTITY OF AMERSE INSTRUMENT GENMICIDE TO MAKE A 1:32 AQUEOUS USE-DILUTION				
Water	1 Pint	1 Quart	1/2 Gallon	1 Gallon
AMERSE	(1 stroke) 1/2 fl. oz.	(2 strokes) 1 fl. oz.	(4 strokes) 2 fl. oz.	(8 strokes) 4 fl. oz.

To remove residue of previous disinfectant, thoroughly clean and rinse instruments, trays or containers prior to first use of this product. Prepare a new solution daily or when solution turns cloudy or dirty for disinfecting precleaned items.

FOR COLD INSTRUMENT DISINFECTION: After cleaning, immerse articles and wet all surfaces thoroughly in a 1:32 solution for 10 minutes. Rinse and dry.

FOR METAL INSTRUMENT STORAGE—A 1:32 solution is non-corrosive and does not rust, dull or otherwise stain or attack quality instruments such as carbide dental burs. The solution may be used for prolonged storage of all but lens instruments, plastic or rubber items. This is a complete product. Do not add other chemicals. Use only as directed. If frozen, thaw and remix before use.

ACCEPTED
 JUL 20 1969
 Under the Federal Insecticide,
 Fungicide, and Rodenticide Act
 as amended, for the pesticide
 registered under
 EPA Reg. No. 10. 3-76



CENTER PANEL

amerse

Hospital Instrument Germicide

EPA Reg. No. 1043-36
EPA Est. No. 1043-MO-1

Germicide • Tuberculocide • Fungicide
FOR DISINFECTION AND STORAGE OF SURGICAL
AND DENTAL INSTRUMENTS

Active ingredients:

o-Benzyl-p-chlorophenol
o-Phenylphenol

5.25%
1.00%

Inert ingredients

93.75%

KEEP OUT OF THE REACH OF CHILDREN

CAUTION

STATEMENT OF PRACTICAL TREATMENT

In case of eye contact, flush with water for 15 minutes. Get prompt medical attention.
For skin contact, flush with water. See back panel for additional precautionary
statements and directions for use

Net Contents: 64 fl. oz. (1/2 U.S. Gallon)

ACCEPTED
NOV 06 1989
Under the Federal Insecticide,
Fungicide and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 1043-36



Amerse Instrument Germicide is a concentrated disinfectant to be used as a 1:32 dilution with water for cold disinfection of precleaned equipment, thermometers, surgical and dental instruments. Recommended for use in hospitals and dental offices or other facilities requiring disinfection of instruments and equipment. The 1:32 dilution kills pathogenic bacteria and fungi such as S. aureus, S. choleraesuis, Ps. aeruginosa, M. tuberculosis (effective in 10 minutes at 20°C), and C. albicans. Amerse is virucidal against Influenza A2 (Hong Kong), Herpes simplex I and II, Type 2 Adenovirus and Vaccinia on hard, inanimate, surfaces in 10 minutes at 20°C.

Minimum Use-Dilution Confirmation (A.O.A.C.-10 Min.-20°C) S. choleraesuis-1:32, S. aureus-1:32, Ps. aeruginosa-1:32

A 1:32 dilution is effective against M. tuberculosis on hard, inanimate surfaces in 10 minutes at 20°C. Remove heavy soil or gross filth and thoroughly clean surfaces.

STORAGE AND DISPOSAL

Prohibitions: Do not contaminate water, food, or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container.

Pesticide Disposal: This germicide, its solutions or rinsings from empty containers, should be disposed of in a toilet or service sink served by a sanitary sewer or in a landfill approved for pesticides.

Container Disposal: Triple rinse and deposit in a waste container for incineration or burial in a land fill approved for pesticide containers.


General: Consult Federal, State, or Local disposal authorities for approved alternative procedures such as limited open burning.

6401-06A



DIVISION OF CALGON CORPORATION
SUBSIDIARY OF MERCK & CO. INC.
St. Louis, Missouri 63122

ACCEPTED
NOV 06 1989
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 1043-36



CENTER PANEL

amerse®

Hospital Instrument Germicide

EPA Reg. No. 1043-36

EPA Est. No. 1043-MO-1

Germicide • Tuberculocide • Fungicide
FOR DISINFECTION AND STORAGE OF SURGICAL
AND DENTAL INSTRUMENTS

Active ingredients:

o-Benzyl-p-chlorophenol
o-Phenyphenol

5.25%
100%

Inert ingredients

93.75%


KEEP OUT OF THE REACH OF CHILDREN

CAUTION

STATEMENT OF PRACTICAL TREATMENT

In case of eye contact flush with water for 15 minutes. Get prompt medical attention.
For skin contact flush with water. See back panel for additional precautionary
statements and directions for use.

Net Contents: 64 fl. oz. (1/2 U.S. Gallon)



Amerse — Right Panel

Amerse Instrument Germicide is a concentrated disinfectant to be used as a 1:32 dilution with water for cold disinfection of pre-cleaned equipment, thermometers, surgical and dental instruments. Recommended for use in hospitals and dental offices or other facilities requiring disinfection of instruments and equipment. The 1:32 dilution kills pathogenic bacteria and fungi such as *S. aureus*, *S. choleraesuis*, *Ps. aeruginosa*, *M. tuberculosis* (effective in 10 minutes at 20°C), and *C. albicans*. Amerse is virucidal against Influenza A₂ (Hong Kong), Herpes simplex I and II, Type 2 Adenovirus and Vaccinia on hard, inanimate, surfaces in 10 minutes at 20°C.

**Minimum Use-Dilution Confirmation (A.O.A.C.-10 Min.-20°C) S
choleraesuis-1:32, S. aureus-1:32, Ps. aeruginosa-1:32**

A 1:32 dilution is effective against *M. tuberculosis* on hard, inanimate surfaces in 10 minutes at 20°C. Remove heavy soil or gross filth and thoroughly clean surfaces.

STORAGE AND DISPOSAL

Prohibitions: Do not contaminate water, food, or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container

Pesticide Disposal: This germicide, its solutions or rinsings from empty containers, should be disposed of in a toilet or service sink served by a sanitary sewer or in a landfill approved for pesticides.


Container Disposal: Triple rinse and deposit in a waste container for incineration or burial in a land fill approved for pesticide containers.

General: Consult Federal, State, or Local disposal authorities for approved alternative procedures such as limited open burning.

6401-06A



DIVISION OF CALGON CORPORATION
SUBSIDIARY OF MERCK & CO. INC.
St. Louis, Missouri 63133



**Precautionary Statements
HAZARD TO HUMANS AND DOMESTIC ANIMALS**

Harmful if swallowed. Avoid prolonged skin contact and avoid splashing in eyes. Wash eyes with large amounts of flowing water if exposed.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

All instruments and devices must be thoroughly cleaned to remove excess organic debris, rinsed and then rough dried before immersion. Clean and rinse the lumens of hollow instruments before filling with solution. Use a 1:32 dilution with water for cold disinfection of precleaned equipment. The removal of heavy soils and the cleaning of all surfaces prior to application of this product is required.

DILUTION TABLE

QUANTITY OF AMERSE INSTRUMENT GERMICIDE TO MAKE A 1:32 AQUEOUS USE-DILUTION				
Water	1 Pint	1 Quart	½ Gallon	1 Gallon
AMERSE	(1 stroke) ½ fl. oz.	(2 strokes) 1 fl. oz.	(4 strokes) 2 fl. oz.	(8 strokes) 4 fl. oz.

To remove residue of previous disinfectant, thoroughly clean and rinse instruments, trays or containers prior to first use of this product. Prepare a new solution daily or when solution turns cloudy or dirty for disinfecting precleaned items.

FOR COLIFORM BACTERIA DISINFECTION: After cleaning, immerse articles and wet all surfaces thoroughly in a 1:32 solution for 10 minutes. Rinse and dry.

FOR METAL INSTRUMENT STORAGE—A 1:32 solution is non-corrosive and does not rust, dull or otherwise stain or attack quality instruments such as carbide dental burs. The solution may be used for prolonged storage of all but lens instruments, plastic or rubber items. This is a complete product. Do not add other chemicals. Use only as directed. If frozen, thaw and remix before use.

