

HAZARD TO HUMANS AND DOMESTIC ANIMALS Concentrate corrosive. Causes eye and skin damage. Do not get in eyes, on skin or on clothing. Wear goggles or face shield and rubber gloves when handling. Harmful if swallowed.

STATEMENT OF PRACTICAL TREATMENT: In case of eye contact, immediately flush eyes with water for 15 minutes and get prompt medical attention. For skin contact, flush thoroughly with water. Remove and wash contaminated clothing before reuse. If swallowed, drink large amounts of water. Do not induce vomiting. Call a physician or poison control center immediately. **NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

ENVIRONMENTAL HAZARD STATEMENT: (For use on containers 5 gallons and larger) Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

CLEANING AND DISINFECTING WASHABLE, HARD, NON-PORUS SURFACES:

To make solution, add one-half fl. oz (15 cc) to each pre-measured gallon of water and gently mix until solution is uniform. **NOTE:** Always add this product (concentrate) to pre-measured water. Apply solution with a cloth, sponge, mop or brush using normal cleaning methods. Allow treated surfaces to soak for 10 minutes. Then remove excess solution with a wrung-out applicator. As solution becomes dirty, discard and replace with a fresh one. A properly prepared solution of Vesphene 15c intended for use as a hard surface disinfectant has a shelf life of 14 days (use life) when stored in a closed container such as a spray bottle. Before cleaning begins, always sweep or dust mop area to be cleaned thoroughly, removing any gross filth. This is a complete product. Do not add other chemicals. If frozen, thaw and re-mix before use. Use only as directed.

Directions for use with Twin-Neck Dispensing Bottle: To make one gallon (1.256) solution, loosen cap from dispensing chamber. Gently squeeze bottle to fill measuring chamber to the 1/2 fluid ounce line. Remove cap from dispensing chamber. Pour measured product (concentrate) from the chamber into one gallon (128 ounces) of pre-measured cold water and mix well. Apply as directed above.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (HUMAN IMMUNODEFICIENCY VIRUS OR AIDS VIRUS) OF SURFACES / OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Wear appropriate barrier protection such as latex gloves, gowns, masks or eye coverings.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of a 1.256 use-solution (one-half ounce per 1 gallon). Prepare and apply solution as directed in paragraph above.

Contact Time: While the HIV-1 virus is inactivated in 1 minute, use a 10-minute contact time for disinfection of all organisms on this label.

Infectious Materials Disposal: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious disposal.

CLEANING AND DISINFECTING FOOD PREPARATION, PROCESSING FACILITIES AND EQUIPMENT:

Cover or remove all food and packaging materials. Remove all gross soils. Saturate all surfaces with the use-solution (one-half ounce per gallon prepared as directed above). Scrub to loosen all soils. Allow to soak for 10 minutes, then pick up excess solution with a wrung-out applicator. Thoroughly rinse all wetted and cleaned food contact surfaces with potable water.

STORAGE AND DISPOSAL

PROHIBITIONS: Do not contaminate water, food, or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container.

PESTICIDE DISPOSAL: This germicide, its solutions or rinsings from empty containers, should be disposed of in a toilet or service sink served by a sanitary sewer, or in a landfill approved for pesticides.

CONTAINER DISPOSAL: Triple rinse and deposit in a waste container for incineration or in a landfill approved for pesticide containers.

GENERAL: Consult Federal, State, or Local disposal authorities for approved alternative procedures such as limited open burning.

LpH[®] se

ONE-STEP GERMICIDAL DETERGENT

Germicide • Fungicide • Virucide • Tuberculocide
• Deodorizer • Cleaner • Detergent

When diluted with soft or hard water (400 ppm CaCO₃) even in the presence of organic soil (5% serum)

EPA REG. NO. 1043-92

FOR INDUSTRIAL AND INSTITUTIONAL USE ONLY

Active Ingredients:

o-phenylphenol.....7.7%

p-tertiary amylphenol.....7.6%

Inert Ingredients:.....84.7%

KEEP OUT OF REACH OF CHILDREN

DANGER

STATEMENT OF PRACTICAL TREATMENT

In case of eye contact, flush with water for 15 minutes.

Get prompt medical attention. For skin contact, flush with water immediately. See side panel for additional

precautionary statements.

ACCEPTED

FEB - 7 2001

NET CONTENTS:

(Logo) STERIS Corporation

5035 Manchester Avenue
St. Louis, Mo 63110 U.S.A.

1043-92

This non-alkaline phenolic germicidal detergent is a state-of-the-art, concentrated scientifically developed product effective in hard or soft water even in the presence of organic soil (5% blood serum). The use-solution of this product when diluted one-half ounce per gallon of water is specially formulated to clean, disinfect and deodorize in a simple one-step process. With its non-alkaline base, it will remove soap scum and hard water salts which can harbor bacteria. Use this product on washable, hard, non-porous environmental surfaces such as floors, walls, counters, tools, carts, and other equipment in food processing preparation and service facilities, veterinary clinics, animal research facilities and in other areas where antibacterial control measures are necessary.

This product is not for use in the neonatal (newborn) units in hospitals specifically isolettes and other infant contact surfaces.

Use for one-step cleaning, disinfecting and deodorizing of stainless or galvanized steel and chrome, porcelain, plastics, glass, tile, washable painted or varnished surfaces, concrete, as well as resilient vinyl, asphalt, linoleum, rubber, terrazzo and other combination type floors.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

GERMICIDAL: When tested against *S. aureus*, *S. choleraesuis*, *Ps. aeruginosa*, it passes the A.O.A.C. Germicidal Use-Dilution Test when diluted with 400 ppm A.O.A.C. hard water to make a 1.256 (one half ounce per 1 gallon) solution, in the presence of 5% organic soil (serum), 10 minutes at 20°C.

BROAD SPECTRUM DATA: The following additional organisms pass the A.O.A.C. Use-Dilution Test in 400 ppm A.O.A.C. hard water at a dilution of 1:256 in the presence of 5% organic soil (serum), 10 minutes at 20°C.

<i>Acinetobacter calcoaceticus</i> , ATCC 19606	<i>Salmonella typhi</i> , ATCC 6539
<i>Candida albicans</i> , Clinical Isolate	<i>Salmonella typhimurium</i> , ATCC 14028
<i>Candida parapsilosis</i> , Clinical Isolate	<i>Serratia marcescens</i> , ATCC 8100
<i>Citrobacter freundii</i> , ATCC 8090	<i>Shigella flexneri</i> , ATCC 12022
<i>Enterobacter aerogenes</i> , ATCC 13048	<i>Shigella sonnei</i> , ATCC 25931
<i>Enterobacter cloacae</i> , ATCC 23355	<i>Staphylococcus aureus</i> , ATCC 25923
<i>Escherichia coli</i> , ATCC 25922	<i>Staphylococcus aureus</i> , (MRSA), Multiply (Methicillin)-Resistant Clinical Isolate
<i>Klebsiella pneumoniae</i> , ATCC 13883	<i>Staphylococcus epidermidis</i> , ATCC 12228
<i>Proteus mirabilis</i> , Clinical Isolate	<i>Streptococcus faecalis</i> , ATCC 19433
<i>Proteus vulgaris</i> , ATCC 13315	<i>Streptococcus pyogenes</i> , ATCC 19615
<i>Pseudomonas aeruginosa</i> , ATCC 27853	

FUNGICIDAL: When tested against *T. mentagrophytes*, it passes the A.O.A.C. Fungicidal Test when diluted with 400 ppm A.O.A.C. hard water to make a 1.256 (one-half ounce per 1 gallon) solution, in the presence of 5% organic soil (serum), 10 minutes at 20°C.

TUBERCULOCIDAL: When tested against *Mycobacterium tuberculosis* var. bovis, BCG, it passes the A.O.A.C. Tuberculocidal Test when diluted with 400 ppm A.O.A.C. hard water to make a 1.256 (one-half ounce per 1 gallon) solution, in the presence of 5% organic soil (serum), 10 minutes at 20°C.

***VIRUCIDAL:** When tested against Influenza A₁ (Japan), Herpes simplex Type 2, Vaccinia and Adenovirus Type 2, it passes the Virucidal Assay (EPA Proposed Method) when diluted with 400 ppm A.O.A.C. hard water to make a 1.256 (one-half ounce per 1 gallon) solution in the presence of 5% organic soil (serum), 10 minutes at 20°C.

When tested by an EPA-approved Dilution Method, the HIV-1 (AIDS) virus, with added 5% organic soil (serum), was completely inactivated by a 1.256 (one-half ounce per 1 gallon) solution in 400 ppm A.O.A.C. hard water in 60 seconds at 20°C. Although efficacy at 1 minute contact time has been shown to be adequate against HIV-1, this time would not be sufficient for other organisms. Use a 10-minute contact time for disinfection against all of the organisms claimed.

ODOR CONTROL: It eliminates many odors by killing odor-causing bacteria while simultaneously chemically neutralizing their odors and leaving a pleasant fragrance.

1043-92

2-7-2001

41