



U.S. ENVIRONMENTAL PROTECTION AGENCY
 Office of Pesticide Programs
 Antimicrobials Division (7510P)
 1200 Pennsylvania Ave., N.W.
 Washington, D.C. 20460

EPA Reg. Number:

1043-138

Date of Issuance:

7/19/22

NOTICE OF PESTICIDE:

Registration
 Reregistration
 (under FIFRA, as amended)

Term of Issuance:

Conditional

Name of Pesticide Product:

Vesta-Syde SQ 64 st Ready-To-Use
 Disinfectant Wipes

Name and Address of Registrant (include ZIP Code):

Kent Ediger
 STERIS Corporation
 7501 Page Avenue
 St. Louis, MO 63133 USA

Electronic Transmittal: Kent_Ediger@steris.com

Note: Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Antimicrobials Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby registered under the Federal Insecticide, Fungicide and Rodenticide Act (FIFRA).

Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA section 3(c)(7)(A). You must comply with the following conditions:

1. Submit and/or cite all data required for registration/reregistration/registration review of your product under FIFRA when the Agency requires all registrants of similar products to submit such data.

Signature of Approving Official:

Heather A. Garvie, Senior Regulatory Advisor
 RMB II, Antimicrobials Division (7510P)

Date:

7/19/22

2. You are required to comply with the data requirements described in the DCI identified below:

a. DDAC GDCI-069149-30869

You must comply with all of the data requirements within the established deadlines. If you have questions about the Generic DCI listed above, you may contact the Reevaluation Team Leader (Team 36): <http://www2.epa.gov/pesticide-contacts/contacts-office-pesticide-programs-antimicrobial-division>

3. The data requirements for storage stability and corrosion characteristics (Guidelines 830.6317 and 830.6320) are not satisfied. A one year study is required to satisfy these data requirements. You have 18 months from the date of registration to provide these data.

4. Make the following label changes before you release the product for shipment:

- Revise the EPA Registration Number to read, "EPA Reg. No. 1043-138."

5. Submit one copy of the final printed label for the record before you release the product for shipment.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under FIFRA and is subject to review by the Agency. See FIFRA section 2(p)(2). If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) lists examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process, FIFRA section 12(a)(1)(B). Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Assurance.

If you fail to satisfy these data requirements, EPA will consider appropriate regulatory action including, among other things, cancellation under FIFRA section 6(e). Your release for shipment of the product constitutes acceptance of these conditions. A stamped copy of the label is enclosed for your records. Please also note that the record for this product currently contains the following CSFs:

- Basic CSF dated 5/25/2021

If you have any questions, please contact Aidan Fife by phone at 202-566-0822, or via email at fife.aidan@epa.gov.


Enclosure

Vesta-Syde® SQ 64 st Ready-To-Use Disinfectant Wipes

Bactericidal, Fungicidal, Virucidal*

EPA REG. NO. 1043-138

EPA Est. No. 1043-M0-2 - EPA Est. No. 92998-NC-1



FOR INDUSTRIAL AND INSTITUTIONAL USE

Active Ingredients:

Didecyl dimethyl ammonium chloride	0.159%
Isopropyl alcohol	10.000%

Other Ingredients:	<u>89.841%</u>
---------------------------	----------------

Total:	100.000%
---------------	----------

KEEP OUT OF REACH OF CHILDREN (KOROC)
CAUTION

NOT TO BE USED AS A BABY WIPE

FIRST AID

In case of emergency, call a poison control center or doctor for treatment advice. Have product container or label with you when calling a Poison Control Center or doctor for treatment advice.

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a Poison Control Center or doctor for treatment advice.

IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a Poison Control Center or doctor for treatment advice.

For chemical emergency, spill, leak, fire, exposure and accident, call Chemtrec, day or night (800) 424-9300, (703) 527-3887. For general information on product, use, etc., call the National Pesticides Information Center at 1-800-858-7378, 8:030 AM to 12:00 PM Pacific time (PT), Monday – Friday; email: npic@ace.orst.edu; or website: <http://npic.orst.edu/>. You may also contact the poison control center at 1-800-222-1222 for emergency medical treatment information.

.REORDER NUMBER: () NET CONTENTS: ()

Manufactured by:
STERIS Corporation
7501 Page Avenue
St. Louis, MO 63133
PRODUCT MADE IN THE U.S.A.

See (side panel) (back panel) (insert) for additional (precautionary statements) (directions for use).

PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION. Causes moderate eye irritation. Harmful if absorbed through the skin. Avoid contact with skin, eyes, or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the restroom. Remove and wash contaminated clothing before reuse.

PHYSICAL OR CHEMICAL HAZARDS

Combustible. Do not use or store near heat or open flame.



This ready-to-use product is intended for use on washable hard, non-porous environmental surfaces such as sealed floors, walls, counters, tools, carts, and other non-critical and semi non-critical medical equipment surfaces used in pharmaceutical and cosmetic manufacturing facilities, and in food processing, preparation and service facilities, in veterinary clinics and animal research facilities.

Use this product for cleaning and disinfection of hard, non-porous surfaces such as metal surfaces made of stainless or galvanized steel and chrome, glazed porcelain, plastics, glass, glazed tile, washable painted or varnished surfaces, sealed concrete, as well as resilient vinyl, linoleum, rubber or terrazzo.

DISINFECTANT: This ready-to-use product is effective in 10 minutes at 20°C on washable hard, non-porous surfaces, in the presence of 5% organic soil against the following: **Bacteria:** *Pseudomonas aeruginosa* (ATCC 15442), *Staphylococcus aureus* (ATCC 6538), *Staphylococcus epidermidis* (ATCC 12228), *Burkholderia cepacia* (ATCC 25416), *Enterobacter cloacae* (ATCC 13047), *Streptococcus pyogenes* (ATCC 19615), *Listeria monocytogenes* (ATCC 19111), *Micrococcus yunnanensis* (ATCC 7468), *Serratia marcescens* (ATCC 14756). **Fungi:** *Trichophyton interdigitale* (ATCC 9533), *Candida albicans* (ATCC 10231), *Aspergillus niger* (ATCC 6275). ***Viruses:** Influenza A2 Virus (A/2/Japan/305/57) (H2N2)*, Human Immunodeficiency Virus Type 1 (HIV-1)*, Avian Influenza (H9N2) Virus*, and Human Rotavirus*, Strain: WA (ATCC VR-2018).

ODOR CONTROL: This product eliminates many odors by killing odor-causing bacteria. Follow "**CLEANING AND DISINFECTING HARD, NON-POROUS SURFACES**" directions.

***KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type HIV-1 (associated with AIDS).

DIRECTIONS FOR USE


It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

CLEANING AND DISINFECTING HARD, NON-POROUS SURFACES:


Remove visible filth mechanically by sweeping or other appropriate cleaning techniques. Roll outer pouch to release liquid onto the wipe material. Peel back the label. Remove wetted wipes from pouch. Use as many wipes as needed to thoroughly wet surfaces for a contact time of 10 minutes. The surfaces must remain visibly wet for the duration of the contact time to ensure clarity for end users. After use, discard pouch and any unused wipes.

This product requires saturation of the dry wipes contained in the package at point of use.


PREPARATION OF WIPES:



1) Ensure package opening is facing up. Use two hands to roll one end of the pouch and apply pressure until you feel the inner packet containing the **Vesta-Syde SQ 64 st Ready-To-Use Disinfectant Wipes** release the disinfectant solution.



2) Invert pouch until wipes are fully saturated.



3) Using the resealable cover, remove wipes. Reseal once desired number of wipes have been removed.

CLEANING AND DISINFECTING WASHABLE HARD NON-POROUS FOOD PREPARATION SURFACES AND PROCESSING FACILITIES AND EQUIPMENT: Cover or remove all food and packaging materials. Remove all visible soil. Use product as directed in paragraphs above. Surfaces must remain visibly wet for the duration of the 10 minute contact time to ensure clarity for end users. Thoroughly rinse all wetted and cleaned surfaces with potable water.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (HUMAN IMMUNODEFICIENCY VIRUS OR AIDS VIRUS) ON HARD, NON-POROUS SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks, and eye coverings.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application. Use product as directed in paragraphs above.

Contact Time: Allow surface to remain visibly wet for 10 minutes.

Infectious Materials Disposal: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious disposal.

STORAGE AND DISPOSAL

Do not contaminate food, feed or water by storage or disposal.

Pesticide Storage: Store in original container in a cool, well-ventilated area. Do not store near heat or open flame. If frozen, thaw and remix before use.

Pesticide Disposal: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your state Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

Container Handling:

Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available.

Wipes Disposal:

Nonreusable one time use container. Discard used and unused wipes in trash. Do not flush wipes down toilet.

Optional Label Text is enclosed in parentheses.

(STERILE UNLESS PACKAGING IS COMPROMISED)