

ACCEPTED

SEP 1 1971

UNDER THE FEDERAL INSECTICIDE
FUNGICIDE AND RODENTICIDE ACT
FOR ECONOMIC POISON REGISTERED
UNDER NO. 10607-1

CHLOROPHENYL

DIRECTIONS FOR USE: CHLOROPHENYL is a powerful disinfecting solution. It is non-injurious to metallic instruments, and non-selective when used as directed. Use full strength. Although the solution is inherently stable it should be periodically changed in practice to avoid pollution or dilution. Prolonged exposure to light will darken solution but will not affect its properties or germicidal potency.

Instruments must be thoroughly cleaned and free from blood before immersion in the solution. Laboratory tests show that common vegetative bacteria (except tubercle bacilli) are destroyed in 5 minutes, but instruments should be immersed for not less than fifteen minutes for adequate safety. When blood or other body fluids have been allowed to dry on instruments before cleaning, long immersion is highly desirable.

The solution is not intended for the destruction of pathogenic spores or the etiologic agent of viral hepatitis. It is not recommended for the disinfection of hypodermic syringes, needles, or blood lancets. It should not be relied upon to disinfect corroded instruments or those with defective plating, deep, narrow crevices, hinged or otherwise inaccessible surfaces. Covered containers such as the Bard-Parker No. 300 are recommended.

WARNING: AVOID CONTACT WITH SKIN AND EYES. IN CASE OF CONTACT FLUSH WITH PLENTY OF WATER. FOR EYES, GET MEDICAL ATTENTION. AVOID CONTAMINATION OF FOOD.

NOTICE: If spilled solution is allowed to dry out it may leave a yellow stain. Wipe off or wash out spilled solution before it dries. Due to possible deterioration of anticorrosive and rust-inhibiting compounds, we do not recommend using solution after date stamped on bottom of container.

B. F. B. P. Bard Parker & CHLOROPHENYL ARE TRADEMARKS. U.S.D.A. REG. NO. 10607-1. MADE AND PRINTED IN U.S.A.

BACK PANEL

