

**TRIADINE(R) 10 BACTERICIDE-FUNGICIDE  
INDUSTRIAL MICROBIOSTAT**

**EPA REG. # 1258-990**

**ACTIVE INGREDIENTS:**

**SODIUM 2-PYRIDINETHIOL-1-OXIDE.....6.4%**

**HEXAHYDRO-1,3,5-TRIS(2-HYDROXYETHYL)-S-TRIAZINE.....63.6%**

**INERT INGREDIENTS:.....30.0%**

**DIRECTIONS FOR USE:**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**FOR THE PRESERVATION OF AQUEOUS ANALYTICAL AND DIAGNOSTIC REAGENTS  
USED IN CHEMICAL AND CLINICAL ANALYSIS.**

Addition of up to 0.1% (by volume) of this product can inhibit microbial growth (bacteria and fungi) in aqueous reagent solutions used in chemical and clinical analysis.

For example, isotonic diluents used in blood diagnostic analyzers can be preserved by addition of 0.03% TRIADINE 10 (approximately 30 mls of the product into 100 L of reagent solution, or 4 oz. per 100 gallons of solution.)

Approval from the FDA must be secured before offering the end product for FDA regulated use.


**DANGER**

**KEEP OUT OF REACH OF CHILDREN**

**SEE PRINCIPLE LABEL FOR COMPLETE PRECAUTIONARY INFORMATION AND  
STORAGE AND HANDLING INSTRUCTIONS**

All applicable directions, restrictions, and precautions on the principle registered label are to be followed.

(t10diag)



**OLIN CHEMICALS  
OLIN CORPORATION  
120 LONG RIDGE ROAD  
STAMFORD, CONNECTICUT 06904**

Reg # 1258-990  
PM-31

Technical Data Sheet

191

TRIADINE(R) 10 BACTERICIDE-FUNGICIDE  
INDUSTRIAL MICROBIOSTAT

EPA REG. # 1258-990

ACTIVE INGREDIENTS:

SODIUM 2-PYRIDINETHIOL-1-OXIDE.....6.4%

HEXAHYDRO-1,3,5-TRIS(2-HYDROXYETHYL)-S-TRIAZINE.....63.6%

INERT INGREDIENTS:.....30.0%

DIRECTIONS FOR USE:

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

FOR THE PRESERVATION OF AQUEOUS ANALYTICAL AND DIAGNOSTIC REAGENTS  
USED IN CHEMICAL AND CLINICAL ANALYSIS.

Addition of up to 0.1% (by volume) of this product can inhibit microbial growth (bacteria and fungi) in aqueous reagent solutions used in chemical and clinical analysis.

For example, isotonic diluents used in blood diagnostic analyzers can be preserved by addition of 0.03% TRIADINE 10 (approximately 30 mls of the product into 100 L of reagent solution, or 4 oz. per 100 gallons of solution.)


Approval from the FDA must be secured before offering the end product for FDA regulated use.

DANGER

KEEP OUT OF REACH OF CHILDREN

SEE PRINCIPLE LABEL FOR COMPLETE PRECAUTIONARY INFORMATION AND  
STORAGE AND HANDLING INSTRUCTIONS

All applicable directions, restrictions, and precautions on the principle registered label are to be followed.



OLIN CHEMICALS  
OLIN CORPORATION  
120 LONG RIDGE ROAD  
STAMFORD, CONNECTICUT 06904