

PM 32 1677-39

RECEIVED  
FEB 13 1989  
1677-39

PLATE 3

DISINFECTANT  
PSEUDOMONACIDAL • TUBERCULOCIDAL  
MICROBICIDAL • HOSPITALS • NURSING HOMES • SCHOOLS

ECOLAB

Net Contents  
1 Gal. (3.78 Liters)

OF HEAT SENSITIVE INSTRUMENTS: Soak for 15 minutes in a solution of 2 ounces MIKRO-BAC per gallon of water. Rinse with water.

GLASS AND PLASTIC GOODS: Soak catheters, tubes, gloves, masks, and no more in a 2 ounce MIKRO-BAC per gallon of water solution. When rinsing, catheters should be soaked for one minute in 70% alcohol.

CONVEYORS: Rinse well. Install and adjust automatic feeding of MIKRO-BAC per gallon of diluted lubricant solution. Apply to conveyor central distribution and existing application points. Allow a minimum 10 minute contact time. Following the 10 minute contact time, the MIKRO-BAC feed can be discontinued. Rinse conveyors following the MIKRO-BAC application.

For specialized cleaning and disinfecting operations, contact your distributor.

For more information, call 1-800-EL-CLEAN (352-5326).

Mikro-Bac®

ACTIVE INGREDIENTS:

Potassium phenylphenate	4.5%
Potassium o-benzyl p-chlorophenolate	5.7%
Potassium p-1-aminophenolate	0.8%
INERT INGREDIENTS:	88.9%

\*Active as cleaning agents.

KEEP OUT OF REACH OF CHILDREN

DANGER

See Statement of Practical Treatment and additional precautionary statements on back panel.

15750

23 1/2 PICAS

MAIN PANEL 27 PICAS

18 1/2"

BACK UP PLATE

RECOMMENDED MINIMUM TYPE  
IS: 8 PT. REGULAR FOR POS TYPE  
AND 8 PT. BOLD FOR REVERSE PRINTING.

USE 18 AND 24 PT LOGOS

COLOR: PMS

ura.

TYPE SMALLER THAN 8 PT. MAY RESULT  
IN BROKEN, OR FILLED IN CHARACTERS.

10/11

CENTER | PANEL 4

Contents  
3.78 Fl. Oz.

**DIRECTIONS FOR USE**

It is the policy of Federal law to sell this product in a wrapper conformant with its labeling.

**COMBINATION DISINFECTION AND CLEANING:** MICRO BAC is effective against *S. aureus*, *S. pneumoniae*, *S. typhimurium*, *S. choleraesuis*, and *Pseudomonas aeruginosa* at 1:200 in hard water (400 ppm as CaCO<sub>3</sub>) and 5% blood serum. Remove heavy soil or gross filth prior to disinfection. Damp mop or dust walls, floors, equipment with 1 ounce MICRO BAC per gallon of water (1:200). Allow surfaces to remain wet for 10 minutes, then wipe dry. Use MICRO BAC in surgical and obstetrical suites (suitable for conductive floors), housekeeping services, physical therapy departments, nursing services, autopsy facilities.

**FUNGICIDAL:** Kills *Trichophyton mentagrophytes* at 1 ounce per gallon of water (1:200).

**TUBERCULOCIDAL:** MICRO BAC passes the rigid requirements as a tuberculocide at 1 ounce per gallon of water (1:200) in the presence of 5% blood serum on porous and nonporous surfaces at 20° C. Remove heavy soil or gross filth prior to disinfection. Allow surfaces to remain wet for 10 minutes, then wipe dry.

**VIROCIDAL:** At 1 ounce per gallon of water (1:200), MICRO BAC passes virucidal modified methods against Herpes simplex Type 1, Adenovirus type 3, and Molluscum virus in the presence of 5% serum on hard, inanimate surfaces.

EPA Reg. No. 1677-39  
EPA Est. No. 1677 IL 2 (J), 1677 NJ 1 (W), 1677 TX 1 (D), 1677 GA 1 (M),  
1677 MN 1 (E), 1677 CA 1 (S), 1677 PR 1 (B), 6574 CA 1 (L)

The superscript refers to the first letter of the date code

71677/04M 1988 Printed in U.S.A. © 1981 Ecolab Inc.

5750

ART/COPY MAX - 3 1/2" (20-9 PICAS)

3 5/8 CYLINDER BACKUP

2.3 PICAS

COLOR: PMS \_\_\_\_\_

1/4" ALLOWANCE ALL SIDES FOR CYLINDER BACKUP (EXCEPT THIS RIGHT SIDE - COLOR WAVE & BACKUP CYLINDER MATCH AS SHOWN)

Handwritten mark

CENTER PANEL

PELIGRO: Si no se lee estrictamente las instrucciones de uso apropiadas antes de trabajar con este producto.

**PRECAUTIONARY STATEMENTS**

**HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

**DANGER:** Contains eye irritant and skin irritant. Do not get in eyes or on skin or clothes. Protect skin and eyes when handling. Harmful if swallowed. Avoid contamination of food.

**PHYSICAL AND CHEMICAL HAZARDS:** Do not use pump sprayer or spray near heat or open flame. Do not mix with anything but water.

**STATEMENT OF PRACTICAL TREATMENT**

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. If swallowed, call a physician. Remove and wash contaminated clothing before reuse.

If swallowed, drink promptly large quantities of water. Avoid alcohol.

**CALL A POISON CONTROL CENTER OR PHYSICIAN IMMEDIATELY FOR EMERGENCY MEDICAL INFORMATION. CALL TOLL FREE 1 800 328 0026**

**NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage. Measures against respiratory shock, respiratory depression and convulsion may be needed.

**STORAGE AND DISPOSAL**

**DO NOT CONTAMINATE WATER, FOOD OR FEED BY STORAGE OR DISPOSAL**

**PESTICIDE DISPOSAL:** Pesticide spray mixture or rinse water that cannot be used according to label instructions must be disposed of according to applicable federal, state or local procedures.

**CONTAINER DISPOSAL:** Do not reuse empty container. Wrap container and put in trash.

**ECOLAB** Institutional Products Division  
Ecolab Center, St. Paul, MN 55102

23 1/2 PICS

1 gallon polyethylene specification with back up (cylinder) plate.  
SIZE: 18 1/2" x 3 5/8" (back up cylinder) 18 7/16" x 3 1/2" copy/art print or  
QUALITY PLATE SPECIFICATION, For plants: Joliet, Midway, Muni  
Dean Ackerman(Quality Platemaking) & Weg(Ecolab) April 1988 (Revised)



MAX HEIGHT 33" MAX DIAMETER 17.8" DIA

PELIGRO: NO SE DEBE USAR SI SE USAN APROPRIADAS ANTES DE

PRECAUTIONARY STATEMENTS  
HAZARDS TO HUMANS AND DO  
DANGER: Contains Gases or  
Protect skin and eyes when handling  
PHYSICAL AND CHEMICAL HA  
Mix with anything but water


STATEMENT OF PRACTICAL TR  
In case of contact, immediately  
Call a physician, remove and wash

If swallowed, drink promptly and  
CALL A POISON CONTROL CENTER  
FOR EMERGENCY MEDICAL IN

NOTE TO PHYSICIAN: Probable mu  
against respiratory tract - respirator

STORAGE AND DISPOS  
DO NOT CONTAMINATE WATER  
PESTICIDE DISPOSAL: Pesticide  
instructions must be disposed of  
CONTAINER DISPOSAL: Do not

**ECOLAB**



1 gallon polyethy  
SIZE: 18 1/2" x 3 5/8"  
QUALITY PLATE  
Dean Ackerman (Quality f