

ionox

DISINFECTANT CLEANER

ACTIVE INGREDIENT:

IODINE * 1.75

INERT INGREDIENTS: 98.25

TOTAL INGREDIENTS: 100.00%

*From alpha-(p-Nonylphenyl)-omega-Hydroxypoly (oxyethylene)-iodine complex

EPA REG. NO. 2296-93

EPA EST. NO. 2296-PA-1

**KEEP OUT OF REACH OF CHILDREN
DANGER**

**SEE BACK PANEL FOR ADDITIONAL
PRECAUTIONARY STATEMENTS**

Mfg. by

National Chemical Laboratories

401 North 10th Street
PHILADELPHIA, PA 19123

of PA, Inc.

Net Contents: 5 gallons / 18.93 Liters

Recommended for cleaning, disinfecting and sanitizing in hospitals, institutions, food plants and commercial buildings. A fresh solution should be prepared daily or more often if the solution becomes diluted or soiled. This product can be used in hard water with up to 400 ppm CaCo3 and cold water as low as 41° F (5°C). *This product is tuberculocidal when used on pre-cleaned, hard, non-porous surfaces at a dilution of 6 ounces to 5 gallons of water (providing 150 ppm titratable iodine) for a contact time of 10 minutes at 20° C."

DIRECTIONS FOR USE

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING

For General Cleaning and Disinfecting in One operation: Use 3 ounces of disinfectant cleaner to 5 gallons of water (provides 75 ppm titratable iodine).

For Cleaning and Disinfecting in Hospitals: Use 6 ounces of disinfectant cleaner to 5 gallons of water (provides 150 ppm titratable iodine).

For Sanitizing Previously Cleaned Food Contact Surfaces: Use 1 ounce disinfectant cleaner to 5 gallons of water (provides 25 ppm titratable iodine).

Hospitals: Dilute 6 ounces of disinfectant cleaner in 5 gallons of water for use in operating areas, patient rooms, floors, walls, counter tops, rest rooms, shower areas, refuse containers, bed frames, morgues. Laboratories, cafeterias, commercial freezers, bathroom fixtures, (stainless steel, chrome, porcelain and other non porous surfaces). Remove gross filth or heavy soil prior to use of the disinfecting solution. Thoroughly wet surfaces and then scrub with a mop, sponge, brush, etc.. Allow 10 minutes contact time and drain dry without rinsing.

Medical / Surgical Instruments: Before applying disinfecting solution thoroughly clean interior and exterior surfaces of thermometers, scalpels, tubing, catheters and other, instruments. Completely immerse small objects and instruments in the disinfecting solution. Flush and fill tubing, catheters and the lumens of hollow instruments with the disinfecting solution. Allow at least 10 minutes of contact time then drain dry without rinsing

Commercial Buildings, Institutions and General Housekeeping: Dilute 3 ounces of disinfectant in 5 gallons of water for use in cleaning and disinfecting walls, floors, lavatories, shower rooms, refuse containers, kitchen areas and lunch rooms. Allow 10 minutes contact time. Rinse all food handling equipment with potable water prior to use.

STATEMENT OF PRACTICAL TREATMENT

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink promptly a large quantity of water, do not induce vomiting, avoid alcohol. Get Medical Attention.

Food Plant Sanitization: Sanitize previously cleaned and rinsed hard non porous surfaces such as glass, metal, plastic, and porcelain. Use 1 ounce of disinfectant to 5 gallons of water. Spray or immerse equipment with sanitizing solution and allow 1 minute contact time. Allow to drain and air dry but do not rinse. When used at this dilution (1:640) food processing equipment does not require a final rinse with potable water in accordance with federal food additive regulation 178.1010.

PRECAUTIONARY STATEMENTS

Hazard to humans and domestic animals.

DANGER

Keep Out Of Reach Of Children . Corrosive. Causes severe eye and skin damage. Do not get in eyes, skin, or on clothing. Harmful if swallowed. Wear goggles or face shield and rubber gloves when handling. Avoid contamination of food. Remove and wash contaminated clothing.

NOTE TO PHYSICIAN

Probable mucosal danger may contraindicate the use of gastric lavage.

STORAGE AND DISPOSAL

Do not Contaminate Water, Food Or Feed By Storage Or Disposal.

Storage: Store in a dry place between 50°F and 120°F.

Pesticide Disposal: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of according to label instructions contact your State Pesticide or Environmental Control Agency or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

Container Disposal: Do not reuse empty container. Triple rinse empty container with water. Rinse metal drum, then offer for reconditioning or puncture and dispose of in a sanitary landfill, or by other procedures approved by State or Local authorities. Plastic containers may be disposed of in a sanitary landfill, incinerated or, if allowed by local authorities, by burning. If burned, stay out of smoke.

Environmental Hazards: Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

INX-0225

MAY 22 1996

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 2296-93