

Cetylcide®

An instrument and hard surface disinfectant for use in hospitals and dental/medical operatories.

Active Ingredients:

Cetyl Dimethyl Ethyl Ammonium Bromide	6.5%
Alkyl (50% C ₁₂ , 30% C ₁₄ , 17% C ₁₆ , 3% C ₁₈) Dimethyl Benzyl Ammonium Chloride	6.5%
Isopropyl Alcohol	13.0%
Inert Ingredients	74.0%
Total Ingredients	100.0%

Formulation Contains Sodium Nitrite

Keep out of reach of children

DANGER

See back panel for additional Precautionary Statements.



**CETYLITE
INDUSTRIES, INC.**
9051 River Road
Pennsauken, N.J. 08110

Net Contents: 1 Quart (946 ml.)

32 Oz.

ACCEPTED

JUL 07 1987

Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 3150-1

Directions for Use

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Per Gallon: 1 oz. (30 ml) Cetylcide to one gallon (3.785 liters) tap water. In known hard water areas, use distilled water. Clean instruments thoroughly ... rinse free of all soap before immersion, and allow them to remain in Cetylcide solution 15 minutes. Germicides containing a Quaternary Ammonium derivative should not be relied on to destroy spore forming organisms, *Mycobacterium tuberculosis* or the etiologic agent of viral hepatitis. Thus, instruments suspected of such contamination should be sterilized by heat. Needles, hypodermic syringes, corroded instruments or those with deep narrow crevices, as well as hinged or defective plated instruments should be sterilized by heat. To maintain disinfection, they may then be immersed in Cetylcide solution. Discard the used solution daily.

Dilution directions for use against *Pseudomonas Aeruginosa*.

Per Gallon: 2 oz. (60 ml) Cetylcide to one gallon (3.785 liters) tap water.

Per Quart: ½ oz. (15 ml) Cetylcide to one quart (946 ml) tap water.

Storage and Disposal

Do not reuse empty container.

Rinse empty container thoroughly and discard in trash.

Cetylcide ... when diluted as directed, is a potent germicidal solution combining bactericidal and bacteriostatic action against many kinds of pathogenic bacteria. Cetylcide is stable in both concentrated and diluted form. When diluted for use as directed, Cetylcide makes a RUST-INHIBITING ODORLESS COLORLESS, efficient germicidal solution for the chemical disinfection of all METAL, RUBBER, and PLASTIC SURGICAL and DENTAL INSTRUMENTS and APPLIANCES.

Precautionary Statements

Hazards to Humans and Domestic Animals

Corrosive. Causes eye and skin damage. Do not get in eyes, on skin, or on clothing. Wear goggles or face shield and rubber gloves when handling. Harmful or fatal if swallowed. Avoid contamination of food.

Statement of Practical Treatment

Skin and Eyes: In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse.

Ingestion: If swallowed, drink PROMPTLY a large quantity of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water.

Avoid alcohol. Call a physician immediately.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage.

EPA Reg. No. 3150-1, EPA Est. No. 3150-NJ-1

Made in U.S.A. Rev. 4/87