PM 31

3150-6

1/2

PROPOSED FRONT LABEL

Net contents 12.8 oz. (380 mL)

Cetylcide-G®

DISINFECTING AND STERILIZING SOLUTION

For use in Hospitals, Medical and Dental Offices and Other Health Care Facilities

Bactericidal, Virucidal*, Tuberculocidal, Fungicidal, Sporicidal.

Cetylcide-G^o must be used only with use-dilution (1:64) Cetylcide³

Active Ingredients:

Keep out of reach of children

Danger

Precautionary Statements
Wear goggles or safety glasses when handling.

Hazards to Humans and Domestic Animals

Causes severe eye and skin irritation; harmful if swallowed. Avoid contact with eyes, skin and clothing. Avoid contamination of food.

Statment of Practical Treatment

In cases of eye contact, immediately flush with water for 15 minutes, then get medical attention. For skin, flush with water. Wash clothing before re-use. <u>IF SWALLOWED</u>: drink promptly a large quantity of water. Do not induce vomiting. Avoid alcohol. Get Medical attention. <u>NOTE TO PHYSICIAN</u>: Probable mucosal damage may contraindicate the use of gastric lavage.

EPA No. 3150-4 EPA Establishment No. 4150-RU4

Made in U.S.A.

ACCEPTED

MAY 3 1 1990

Directions for use:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Do not use this product on food or feed.

Mix use-dilution (1:64) Cetylcide according to directions on bottle: ie., 60 ml. Cetylcide to one gallon tap water. In known hard water areas, use softened, deionized or distilled water.

Per quart of use-dilution Cetylcide, add two pre-measured portions 63 ml. Cetylcide-G. To measure portion, remove cap and seal (if unopened bottle) from Cetycide-G bottle, gently squeeze lower portion of indented bottle allowing solution to fill the top section of bottle to black fill line. Pour contents of top section into quart of use-dilution Cetylcide. The resulting mixture will be green. To one-gallon of use-dilution Cetylcide, add eight pre-measured portions, 252 ml. Cetylcide-G.

Label the Cetylcide-G use-dilution with the label(s) provided.

Clean instruments thoroughly. Hinse free of all soap before immersion. Clean and rinse the lumens of hollow instruments before filling with solution.

For use as a disinfectant, immerse instruments in Cetylcide-G use-dilution for 20 minutes at room temperature. This will destroy bacteria in 10 minutes. Fungi in 10 minutes, Viruses" including: Herpes simplex Type I, Herpes simplex Type II, Poliovirus Type II, Poliovirus Type II, Poliovirus Type III, Coxsackie Virus BI, Vaccinia Virus, Adenovirus, in 10 minutes on hard, inanimate surfaces. Mycobacterium tuberculosis in 20 minutes at 0°C. Can be re-used as a disinfectant use-dilution for up to 28 days.

For use as a sterilant, immerse instruments in Cetylcide-G use-dilution for 3 hours at room temperature. Remove instruments from the solution using a sterile technique and rinse thoroughly with sterile water. If sterile water is not available, use 70% isopropyl alcohol for rinse. This sterilization procedure destroys pathogenic spores. Use-dilution can be re-used as a sterilant for up to 28 days.

Storage and Disposal

Do not contaminate water, food or feed by storage disposal. Pesticide wastes are toxic. Improper disposal of excess pesticide, spray use mixture or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label directions, contact your State Pesticide or your Environmental Control Agency, or the Hazardous Waste representative at the nearest EFA Regional Office for guidance. Waste resulting from the use of this product may be disposed of on site.

Storage

This product should be stored at room temperature, preferably in a cool place.

Container Disposal

Do not reuse empty container. Triple rinse (or equivalent) top section of the indented bottle, then punctual top section and dispose of entire container in a sanitary landfill or by incirrigation if allowed by state and local authorities.

MAY 3 1 1990

MAY 3 1 1990

MAY 3 1 1990

MANUAL PROPERTY OF THE SERVICE OF THE S

BEST AVAILABLE COPY