

May 17-1969

Will not rust or corrode surgical steel, or other metals. 513541

May be used for disinfection of metal, plastic, rubber, or glass instruments and apparatus including catheters, gloves, hypodermic equipment, and tubing.

Not to be relied upon to destroy bacterial spores, Mycobacterium tuberculosis or the etiological agent of viral hepatitis.

NET CONTENTS—1 GALLON

ERLEN[®]

UNICIDE

GERMICIDAL SOLUTION

BRAND OF BENZALKONIUM CHLORIDE
AQUEOUS SOLUTION
FOR DISINFECTION OF DENTAL AND
SURGICAL INSTRUMENTS

ACTIVE INGREDIENTS:

N-ALKYL (C₁₂—40%, C₁₄—50%, C₁₆—10%) DIMETHYL
BENZYL AMMONIUM CHLORIDES 0.27%

SODIUM CARBONATE (ANHYDROUS) 0.42%

INERT INGREDIENTS 99.31%

CAUTION: KEEP OUT OF REACH OF CHILDREN
AVOID EYE CONTACT. IN CASE OF EYE CONTACT, FLUSH
WITH WATER. USE PRECAUTIONS.



WATER SERVICES BURBANK, CALIFORNIA 91502

UNIVERSAL OIL PRODUCTS COMPANY

DIRECTIONS

The solution should be used full strength. Clean instruments well and rinse. Submerge article in solution for at least 10 minutes.

Instruments may be used without further rinsing or drying.

Steam autoclaving is mandatory for destruction of spore forming bacteria, after which instruments or needles may be stored in a sterile covered pan in this solution until used.