



TINTED

Brand of 12.9% benzalkonium chloride solution for the chemical disinfection of surgical or dental instruments.

ACTIVE INGREDIENTS: Alkyl (C12-67%, C14-25%, C15-7%, C16, C17 and C18-1%) dimethyl benzyl ammonium chloride 12.9%.

INERT INGREDIENTS: 87.1% (including sodium nitrite 15.9%)

DIRECTIONS: Dissolve contents in one gallon of either tap water or d.stilled water. Completely immerse thoroughly cleaned rubber or metal instruments for 15 minutes or longer. Solution may be used for storage of heat sterilized instruments. See package insert for complete directions. Incorporated in each individual package is sufficient anti-rust material for one gallon of water. Do not use anti-rust tablets with BENZ-ALL TINTED as they are not needed.

DANGER: Keep Out of Reach of Children. Corrosive. Causes eye damage and severe skin irritation. Do not get in eyes, on skin or on clothing. Protect eyes and skin when handling. Harmful if swallowed.

FIRST AID: In case of contact immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol. Call a physician immediately.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

CONTENTS: 15—1.35 Fl. Oz. (40cc) Packets
Each packet makes a 1 gallon solution of 1:750 solution of Benzalkonium Chloride



XTTRIUM LABORATORIES
Chicago, Illinois 60609

EPA,Reg. No. 5464-1 EPA Est. 5464-IL-1

Stock No. 300-40 T



TINTED



15—1.35 Fl. Oz. (40cc) Packets

Stock No. 300-40 T



TINTED

Brand of 12.9% benzalkonium chloride solution for the chemical disinfection of surgical or dental instruments.

ACTIVE INGREDIENTS: Alkyl (C12-67%, C14-25%, C15-7%, C16, C17 and C18-1%) dimethyl benzyl ammonium chloride 12.9%

INERT INGREDIENTS: 87.1% (including sodium nitrite 15.9%)

DIRECTIONS: Dissolve contents in one quart of either tap water or d.stilled water. Completely immerse thoroughly cleaned rubber or metal instruments for 15 minutes or longer. Solution may be used for storage of heat sterilized instruments. See package insert for complete directions. Incorporated in each individual package is sufficient anti-rust material for one quart of water. Do not use anti-rust tablets with BENZ-ALL TINTED as they are not needed.

DANGER: Keep Out of Reach of Children. Corrosive. Causes eye damage and severe skin irritation. Do not get in eyes, on skin or on clothing. Protect eyes and skin when handling. Harmful if swallowed.

FIRST AID: In case of contact immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol. Call a physician immediately.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

CONTENTS: 15—0.34 Fl. Oz. (10cc) Packets
Each packet makes a 1 quart solution of 1:750 solution of Benzalkonium Chloride



XTTRIUM LABORATORIES
Chicago, Illinois 60609

EPA Reg. No. 5464-1 EPA Est. 5464-IL-1

Stock No. 300-10T




TINTED



15—0.34 Fl. Oz. (10cc) Packets

Stock No. 300-10 T





PROD. NO. 01265

# D-TROL

## DISINFECTANT, SANITIZER AND ALGICIDE

**ACTIVE INGREDIENTS:**

n-alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> )	5.25%
dimethyl benzyl ammonium chlorides	5.25%
n-alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethyl benzyl ammonium chlorides	10.50%

**INERT INGREDIENTS:**

Water, Food Color and Ethyl Alcohol	89.50%
EPA Registration No. 5739- AA	100.00%

**USING DIRECTIONS:**

**For Final Sanitizing Rinsing in Food Plants:** Clean equipment with a suitable detergent and rinse. For final sanitizing rinsing use ¼ ounce of D-TROL (200 ppm Quat) per one gallon of water.

D-TROL complies with Paragraph 121.2547, Sanitizing Solution, Part 121, Food Additive Law for final sanitizing rinsing of food processing equipment and on food contact surfaces in bars and restaurants. No rinse required in dairies unless specified by area health authorities. This product fulfills the criteria of Appendix F of the Grade A Pasteurized Milk Ordinance 1965 Recommendations of the U.S. Public Health Service in waters up to 750 ppm of hardness calculated as CaCO<sub>3</sub> when tested by the AOAC Germicidal and Detergent Sanitizers Official Method.

**For Sanitizing Equipment in Meat and Poultry Plants:** Clean equipment with a suitable detergent and rinse. Apply ¼ ounce of D-TROL in 1 gallon of water (200 ppm Quat) as a final sanitizing rinse, provided that surfaces are adequately drained before food contact. Listed as a Code G, sanitizer and Code N, compound for treating closed systems in the "List of Chemical Compounds authorized for use under USDA Poultry, Meat, Rabbit and Egg Product Inspection Program."

**For Disinfecting in Hospitals and Rest Homes:** Clean equipment with a suitable detergent and rinse. For disinfecting hard surfaces such as floors, walls, furnishings and carts use ¾ ounces of D-TROL per 1 gallon

of water (600 ppm Quat). For disinfecting porous surfaces use 1 ½ ounces of D-TROL per 1 gallon of water (1200 ppm Quat).

See separate D-TROL Technical Data Sheets for additional information.

**APPLICATIONS:** After cleaning surfaces with suitable detergent and rinsing, the sanitizing solution of D-TROL may be applied by spraying with suitable spray equipment, a clean mop or brush.

**DANGER!**  
KEEP OUT OF REACH OF CHILDREN

**CORROSIVE:** Causes eye damage and skin irritation. Do not use on skin, or on clothing. Harmful if swallowed. Avoid contact with food. Protect eyes and skin when handling.

**First Aid:** In case of contact flush eyes or skin with plenty of water for at least 15 minutes. For eyes call a physician. Remove contaminated clothing before reuse. If swallowed, drink plenty of water. Do not induce vomiting. If necessary, give 1-2 glasses of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol. Call a physician immediately. **NOTE TO PHYSICIAN:** Probable mucosal irritation. Contraindicate the use of gastric lavage. Measures against shock, respiratory depression and convulsion may be needed. Do not reuse container. Destroy when empty.

Do not contaminate water by the cleaning of equipment, floors, or wastes.

This product does not contain phosphorus.

**KEEP FROM FREEZING**  
NET CONTENTS SHOWN ELSEWHERE ON CONTAINER