

FEB 17 1989

National Detergents, Inc.
2819 West Lake Street
Chicago, IL 60612

Attention: Brian Egan

Gentlemen:

Subject: S.A.I.D. II
EPA Registration No. 6621-53
Your Submission Dated January 22, 1987

The submission referred to above, submitted in connection with the Label Improvement Program Notice for Antimicrobial Products Used as Food Contact Surface Sanitizers, is acceptable provided that you:

1. Make the labeling changes listed below before you release the product for shipment bearing the amended labeling.
 - a. On the left panel, above the word "DANGER," add the following two headings:

Precautionary Statements
Hazards to Humans and Domestic Animals

- b. Revise the statement "Protect eyes and skin when handling concentrate" to read "Wear goggles or face shield and rubber gloves when handling."
- c. Refer to PR Notice 83-3 for appropriate pesticide disposal instructions and for appropriate storage and disposal instructions for container sizes greater than 1 gallon.
- d. Revise the "First Aid" heading to read "Statement of Practical Treatment."
- e. Delete the word "milk."
- f. The statement "Keep Out of Reach of Children" must be placed directly above the signal word "DANGER" (center panel only).

50260:I:DeLaney:L-15:KENCO:02/01/89:03/16/89:CL:EK:AS

CONCURRENCES

SYMBOL								
SURNAME								
DATE								

- g. Revise the referral statement to read:

See left side panel for additional
precautionary statements.

- h. Revise the "Use Directions" heading to "Directions For Use"
and add the following misuse statement directly below to
read:

It is a violation of Federal law to use this
product in a manner inconsistent with its
labeling.


- i. In the "Directions For Use" instructions, change "a pre-
cleaning step is recommended" to read "a pre-cleaning step
is required" (two places).
- j. Include a 10-minute contact time for the disinfection section.
- k. Include a 1-minute contact time for the sanitization section.
2. Submit five (5) copies of final printed labeling (deleting the
words "Catalog Label" from across the front panel) before you
release the product for shipment.

Your release for shipment of the product bearing the amended labeling
constitutes acceptance of these conditions.

A stamped copy of the labeling is enclosed for your records.

If you have any questions concerning this letter, please contact
Martha DeLaney at (703) 557-6982.

Sincerely yours,

John H. Lee 
Product Manager (31)
Antimicrobial Program Branch
Registration Division (TS-767C)

Enclosure

S.A.I.D. II DEODORIZER

DISINFECTANT SANITIZER

Cold Sterilization

For sterilization of previously cleaned surgical instruments, barber tools and dental equipment, simply submerge instruments into solution containing 1% (one) part S.A.I.D. II per gallon of water for 10 minutes.

FOR ROOM FOGGING

Directions:

To minimize the danger of cross contamination from environments, surfaces, and solutions of S.A.I.D. II as a fog is described below, before applying standard procedure for room cleaning and disinfecting practices: 1. Remove all human, animal and plant life from room. 2. Open closet doors and drawers. 3. Set up fog receiving platform in center of room. 4. Mount fogging MASTER MIXTURE (6208) or other fogging device delivering equivalent spray at 3 RPM. 5. Fill spray reservoir with 1.5 oz. of S.A.I.D. II in 1 gallon of water to produce 1150 ppm active ingredients. 6. Set sprayer mechanism to deliver 1 gallon of solution. 7. Fog for 15 minutes for an average 2150 cu. ft. room. For different room sizes, vary spray time proportionately to volume. 8. After setting of exposed surfaces. 9. Wait 2 hours before entering room after fogging. 10. Remove fogging equipment thoroughly with clear water following use.

DANGER

KEEP OUT OF REACH OF CHILDREN

Caution: Causes severe eye and skin damage. Do not get in eyes, on face, or on clothing. Protect eyes and skin when handling concentrate. Avoid contact with skin. Avoid contact with food. **NEVER EMPTY CONTAINER THOROUGHLY WITH WATER. USE AND WASH CONTAINER THOROUGHLY TO REMOVE RESIDUALS.**

FIRST AID

In case of eye contact, flush eyes with water for at least 15 minutes. If on face, flush with water. If on skin, wash with soap and water. If inhaled, drink milk, egg whites, or other liquid if available. If severe symptoms of water, call a physician.

NOTE TO PHYSICIAN

Proteinuria and renal damage may contraindicate the use of quaternary ammonium compounds. Consult your physician. Hypertension, depression and confusion may be needed. Rinse empty container thoroughly with water and discard.

ACTIVE INGREDIENTS:

n-Alkyl (60% C12, 30% C16, 5% C12, 0% C18) dimethyl benzyl ammonium chlorides 5%
n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides 5%

INERT INGREDIENTS:

90%
100%

S.A.I.D. II is a concentrate and should be diluted before using

ACCEPTED
with COMMENTS
in EPA Letter Dated:

FEB 17 1989

Under the Federal Insecticide, Fungicide, and Rodenticide Act is amended for the pesticide registered under EPA Reg. No. 6621-53

DANGER

KEEP OUT OF REACH OF CHILDREN

(See other precautionary and first aid statements on left side panel)

EPA Est No. 6621-II-1
EPA Registration No. 6621-53

WARRANTY

Our recommendations for use of this product are based upon tests believed to be reliable. The use of this product being beyond the control of the manufacturer, no guarantee, expressed or implied, is made as to the effects of such or the results to be obtained if not used in accordance with directions or established safe practice. The buyer must assume all responsibility, including injury or damage, resulting from its misuse as such or in combination with other materials. Statements concerning the use of these products are not to be construed as recommending the infringement of any patent, and no liability for infringement arising out of any such use is assumed.

USE OF
New use
For use
Shower
Sephro
for use
Happo
For use
Sme. Top
per 5 gal
effective
temperature
Food Pro
For use
dual 10
names. A
Dance
To use
use the
S.A.I.D. II
current f
Reserve
Ester
S.A.I.D. II
The pro
5 gallons
and V