

**STAF
HOSPITAL SPRAY DISINFECTANT
PRECAUTIONARY STATEMENTS
HAZARDS TO HUMANS AND
DOMESTIC ANIMALS**

WARNING

Avoid getting spray concentrate in eyes. In case of contact, flush with water. If irritation persists, get medical attention.

PHYSICAL OR CHEMICAL HAZARDS

Contents under pressure. Do not use or store near heat or open flame. Exposure to temperatures above 130° F may cause bursting. Do not puncture or incinerate container.

DISPOSAL

Do not re-use empty container. Wrap container and put in trash collection.

DIRECTIONS FOR USE

General Classification

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

For use in Hospitals, Nursing Homes, Institutions, Offices, Schools, Motels and Hotels.

Kills Staphylococcus Aureus, Salmonella Choleraesuis, Pseudomonas Aeruginosa, and Mycobacterium Tuberculosis.

For use on telephones, door knobs, light switches, bed springs, diaper pails, toilet seats, sink tops, bed pans, garbage pails, hampers, linen carts, and desk tops.

TO DISINFECT. 1. First remove gross filth and heavy soil deposits from surface to be treated.

2. Hold can six inches from surface and spray until completely wet, allowing the surface to remain wet for 10 minutes.

To deodorize, air refreshen and sanitize by temporarily reducing the number of air borne bacteria, point spray upwards to center of room and spray for three (3) seconds.

EPA Reg No 7056 67 EPA Est 7056-TX 1

**CHEM SPRAY
STAF
HOSPITAL
SPRAY
DISINFECTANT**

**Eliminates Odors
Kills Most Germs
Pleasantly Scented**

ACTIVE INGREDIENTS:

n Alkyl (60% C ₁₂ , 30% C ₁₄ , 5% C ₁₆)	
5% C ₁₂ Dimethyl Benzyl Ammonium Chloride	00 10%
n Alkyl (60% C ₁₂ , 32% C ₁₄) Dimethyl Ethylbenzyl Ammonium Chloride	00 10%
Isopropanol	52 00%
Essential Oils	2 00%
Triethylene glycol	8 00%
INERT INGREDIENTS:	39 80%

Contains Sodium Nitrite as a corrosion inhibitor.

**Keep Out Of Reach Of Children
WARNING**

See Additional Cautions on Back Panel

Mfg by
CSA Ltd., Inc.
Houston TX 77090
Net. Wt. 11 oz.

7056-67