

Reorder No.
2253

REUSABLE

Sterilizing and Disinfecting Solution

Rapid-acting / Non-staining / Non-rusting

Sporicidal

Tuberculocidal

Virucidal†

Bactericidal

Fungicidal

Pseudomonacidal

Cidex*

activated dialdehyde solution

EPA REG. No. 7078-1
EPA EST. No. 30126-PR-1

Active Ingredient
Glutaraldehyde 2.0%
Inert Ingredients 98.0%
Total 100.0%

NOTE — Attached vial contents *must* be added to solution before this product is effective.
See Directions for Activating.

Keep out of reach of children.
DANGER: See right side panel for precautions

† See Directions for Use.
* TRADEMARK

ACCEPTED
AUG 22 1989
Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 7078-1

SURGIKOS
a Johnson & Johnson company

0013

0443

PRECAUTIONARY STATEMENTS

DANGER:

HAZARD TO HUMANS

FOLLOW ALL LABEL DIRECTIONS

Eyes — Contact with eyes causes damage. In case of contact, flush with water immediately and get medical attention.

Skin — Avoid skin contact, as possibility of sensitization exists. May cause skin irritation. For skin contact, flush thoroughly with water.

Ingestion — Avoid contamination of food. Harmful if swallowed.

Inhalation — Use in ventilated area.

STORAGE AND DISPOSAL:

Keep in a cool place. Triple rinse empty container with water and dispose in an incinerator or landfill approved for pesticide containers.

Especially recommended for sterilization or disinfection of . . .

Lensed instruments

Plastic objects

Anesthesia equipment

Sharp instruments

Respiratory therapy equipment

Thermometers

Rubber objects

Flexible fiberoptic endoscopes

CONTENTS: 2½ GALS. (320 FLUID OZ.)

MADE IN U.S.A.

CONTROL NUMBER ON BOTTOM


2253 C1-3

SURGIKOS

a Johnson & Johnson company

00110

W443



For complete sterilization or disinfection of medical instruments and equipment, *thoroughly clean, rinse and rough dry* objects before immersing in CIDEX® Activated Diacetyde Solution. Cleanse and rinse the lumens of hollow instruments before filling with CIDEX Solution.

DIRECTIONS FOR ACTIVATING:

Add contents of attached vial to this container of CIDEX Solution. Shake, solution is ready to use. *When activated, solution changes color to green.* In addition to the standard bioburden challenge, this product has also demonstrated efficacy in the presence of 2% organic soil contamination during reuse. This solution may be used and reused for up to 14 days after activation. Do not use activated solution beyond 14 days.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

For repeated use in manual (bucket and tray) systems.


Disinfection: Immerse completely for a minimum of 10 minutes to destroy vegetative organisms including *Pseudomonas aeruginosa*, pathogenic fungi, and viruses (Poliovirus Type 1; HIV-1; Adenovirus Type 2; Herpes simplex Type 1, 2; Influenza Type A [WSN/33]; Vaccinia; Coronavirus, Cytomegalovirus; Rhinovirus Type 14) on inanimate surfaces. For *Mycobacterium tuberculosis* CIDEX Solution destroys 99.8% in 10 minutes and 100% in 45 minutes at 25°C on inanimate surfaces. Remove equipment from CIDEX Solution and rinse thoroughly.

Sterilization: Immerse equipment completely for a minimum of 10 hours at 25°C, to destroy pathogenic spores as represented by *Clostridium sporogenes* and *Bacillus subtilis*. Remove equipment from CIDEX Solution using sterile technique and rinse thoroughly with sterile water.

NOTE DATES HERE

MIXING DATE _____

EXPIRATION DATE _____



MARKETED BY:



SURGIKOS

a Johnson & Johnson company

SURGIKOS, INC., ARLINGTON, TEXAS 76010

* TRADEMARK

© SURGIKOS 1986



Personal Protection

Disposable latex gloves, eye protection, and face masks should be worn when cleaning and disinfecting soiled equipment.

Cleaning Procedure

Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant or sterilant.

Contact Time

Immerse completely for a minimum of 5 minutes to destroy HIV-1 (AIDS virus). However, a 10 minute immersion is still required for all other vegetative organisms.

Disposal of
Infections
Materials

Blood and other body fluids should be disposed of according to approved hospital procedures.

ACCEPTED
AUG 22 1989
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-1

25
25

Cidex*

activated dialdehyde solution

Activator Vial

For Use With 2 1/2 Gallon
CIDEX Solution Only

© TRADEMARK

DANGER: Contact with eyes causes damage. In case of contact with eyes, flush with water immediately and get medical attention. For skin contact, flush thoroughly with water.

DISPOSAL: Rinse empty container thoroughly with water prior to discarding.

NET CONT. WT. 51 GRAMS
MADE IN U.S.A.
CONTROL NUMBER ON BOTTOM

DISTRIBUTED BY:

SURGIKOS

a Johnson & Johnson company
SURGIKOS, INC., ARLINGTON, TEXAS 76010
© SURGIKOS, INC.

ACCEPTED

AUG 22 1989

Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under 7078-1
EPA Reg. No.

