

✓ Sterilizing and Disinfecting Solution

Rapid-acting / Non-staining / Non-rusting

Sporicidal

Tuberculocidal

Virucidal †

Bactericidal

Fungicidal

Pseudomonacidal

Cidex3X

CONCENTRATE *

long-life activated dialdehyde

Product
CC-850

Active Ingredient
Glutaraldehyde 6.0
Inert Ingredients 94.0
Total 100.0%

Contains Sodium Formate
as a corrosion inhibitor

Note Attached vial
contents *must* be added
to solution before this
product is effective.
See directions.

Rinse empty container
thoroughly with water
prior to discarding.

ACCEPTED
MAY 28 1985
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-6

DANGER: See Back Label for Precautions

Cidex3X
long-life activated dialdehyde

Activator Vial

for use with Cidex3X Solution

Arbrook inc.

Directions: AFTER addition of CIDE X 3X
Concentrate, slowly add contents of this
vial to 3 gallons of CIDE X Fluid. Stir well.
Activated solution will be green.

DANGER: Contact with eyes causes
damage. In case of contact with eyes, flush
with water immediately and get medical aid
at once. For skin contact, flush thoroughly
with water.

TRADE MARK ARBROOK 1975 PATENT PENDING

Arbrook

DIRECTIONS FOR DILUTION and Activation:

Add contents of this container of CIDE X 3X Concentrate to 2 gallons of
distilled or deionized water. Then slowly add contents of attached vial.
Stir well. Solution is ready to use when salts are completely dissolved.

Directions for Use:

For complete sterilization or disinfection of medical instruments and
appliances, *thoroughly clean, rinse and rough dry* objects before immersing
in CIDE X Solution. Cleanse and rinse the lumens of hollow instruments
before filling with CIDE X Solution.

Disinfection and Terminal Decontamination: Immerse completely for a
minimum of 10 minutes to destroy vegetative pathogens including
Mycobacterium tuberculosis, Pseudomonas aeruginosa and viruses
(Poliovirus Type 1, 2, 3, Influenza Type A₁, A₂ (Hong Kong), B, Vaccinia,
Coxsackievirus B-1, Reovirus Type 3, Herpes simplex) on inanimate
surfaces.

Sterilization: Immerse completely for a minimum of 10 hours to destroy
resistant pathogenic spores including **Clostridium sporogenes** and
Clostridium tetani. Remove instruments from CIDE X Solution using
sterile technique and rinse thoroughly with sterile water.

When activated, solution changes color to green

Do not use activated stock solution beyond 28 days after mixing

Discard used solution after each day's use.

NOTE DATES HERE
MIXING DATE

E. P. A. Reg. No. 7078-6

E. P. A. Reg. No. 7078 TX 1

ACCEPTED

EXPIRATION DATE

MAY 28 1985

Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-6

DANGER:

FOLLOW DIRECTIONS FOR USE

HARMFUL IF SWALLOWED

Skin — Avoid skin contact, as possibility of sensitization exists.
May cause skin irritation. For skin contact flush thoroughly with water.

Eyes — Contact with eyes causes damage. In case of contact
flush with water immediately and get medical attention.

Food — Avoid contamination of food.

Inhalation — Use in ventilated area.


Instruments — CIDE X Solution will not prevent or remove oxidation
caused by poor or improper plating, cleaning, or care of instruments.

KEEP IN A COOL PLACE

Especially recommended for terminal decontamination procedures
as well as sterilization or disinfection of:

Lensed instruments
Anesthesia equipment
Respiratory therapy equipment
Rubber objects

Plastic objects
Sharp instruments
Thermometers



CONTENTS ONE GAL (128 FLUID OZ)
MADE IN U.S.A.
PAT. 3,016,328
CONTROL NUMBER ON BOTTOM

Arbrook inc.
Arlington, Texas 76010