


Cidex*
industrial disinfecting
solution concentrate

ACTIVATOR
for use with
CIDEX Industrial Fluid only

*Trademark

DANGER: Contact with eyes causes damage. In case of contact with eyes, flush with water immediately and get medical attention. For skin contact, flush thoroughly with water.

Rinse empty container thoroughly with water prior to discarding.



DISTRIBUTED BY:

SURGIKOS

SURGIKOS, INC., ARLINGTON, TEXAS 76010

*Trademark

SURGIKOS 1982

PRECAUTIONARY STATEMENTS

DANGER:

HAZARD TO HUMANS AND DOMESTIC ANIMALS

FOLLOW ALL LABEL DIRECTIONS

Prevent contact with skin or eyes by using eye protection and rubber gloves when handling or pouring.

STATEMENT OF PRACTICAL TREATMENT

EYES — Corrosive, causes eye damage. Do not get in eyes. In case of contact, flush with water immediately and get medical attention.

SKIN — Avoid skin contact, as possibility of sensitization exists. Can cause skin irritation. For skin contact, flush thoroughly with water. Get medical attention if irritation occurs.

FOOD — Avoid contamination of food.

INHALATION — Avoid spillage. Use in a well-ventilated area.

HARMFUL IF SWALLOWED

CIDEX Industrial Solution has been shown to be an effective disinfectant (virucide, fungicide, bactericide, pseudomonacide) when tested by AOAC and EPA test methods.

STORAGE AND DISPOSAL
KEEP IN A COOL PLACE.


PROHIBITIONS — DO NOT CONTAMINATE WATER, FOOD OR FEED BY STORAGE OR DISPOSAL. OPEN DUMPING IS PROHIBITED.

PESTICIDE DISPOSAL — Pesticide spray mixture, or rinse water that cannot be used according to label instructions must be disposed of according to Federal or approved state procedures under Subtitle C of the Resource Conservation and Recovery Act.

CONTAINER DISPOSAL — Triple rinse (or equivalent). Then offer for recycling or reconditioning, or dispose of in a sanitary landfill, or by incineration if allowed by state and local authorities.

GENERAL — CONSULT FEDERAL, STATE OR LOCAL DISPOSAL AUTHORITIES FOR APPROVED ALTERNATE PROCEDURES.

CONTENT 108 FL. OZ. (3.8 LITERS)
MADE IN U.S.A.
CONTROL NUMBER: 15-B-11110
112020C11




DISTRIBUTED BY
SURGIKOS
SURGIKOS, INC., ARLINGTON, TEXAS 76010

Disinfecting Solution
 Rapid-Acting/Nonrusting/Nonstaining
 Pseudomonacidal Bactericidal
 Virucidal† Fungicidal

Cidex*

industrial disinfecting
 solution concentrate



Active ingredient
 Glutaraldehyde 26%
 Inert ingredients 74%
 Totals 100%

NOTE: Activator must be added to
 fluid before this product is
 effective. See directions.

EPA Reg. No. 7078-11
 EPA Est. No. 36126-PR-1

Reorder Number
112020

Keep out of reach of children.

DANGER: See back panel for additional precautionary statements.
 †See Directions for use.

SURGIKOS

*Trademark

GENERAL CLASSIFICATION

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Intended for use to disinfect mixing and storage tanks, filling equipment and accessory lines in industries such as cosmetic, pharmaceutical and medical device manufacturers.

DIRECTIONS FOR ACTIVATION AND USE:

Metal, glass, rubber or plastic equipment or lines should be thoroughly cleaned and rinsed by using acceptable standard industry techniques prior to using CIDEX Industrial Disinfecting Solution. CIDEX Industrial Disinfecting Solution is prepared by mixing, at a 50:1 volume ratio, 50 parts of CIDEX Industrial Disinfecting Concentrate with 1 part of activator. Use CIDEX Industrial Disinfecting Solution immediately after activation. Dilute CIDEX Industrial Disinfecting Solution Concentrate with 34 equivalent parts of less of distilled, dechlorinated or softened water.

DILUTION CHART

CIDEX Industrial Disinfecting Concentrate (oz)	Activator (cc)	Water (gal)	Product Produced (gal)
32	20	8.5	8.75
64	40	17.0	17.5
128 (1 gal)	80	34.0	35.0
256 (2 gal)	160	68.0	70.0
512 (4 gal)	320	136.0	140.0

DISINFECTION:

Depending on equipment, CIDEX Industrial Disinfecting Solution may be applied by sponge or spray. Surfaces must be completely freed. Insure that all surfaces are thoroughly wet. CIDEX Industrial Disinfecting Solution must remain in contact with equipment for a minimum of 15 minutes. Within 15 minutes contact at 20 C, CIDEX Industrial Disinfecting Solution will destroy vegetative pathogens including *Pseudomonas aeruginosa* and hydrophilic and lipophilic viruses (Poliovirus Type 1; Coxsackievirus B-1; Influenza Type A, (Hong Kong), Herpes Simplex) on inanimate surfaces. Drain CIDEX Industrial Disinfecting Solution and thoroughly rinse with water. To assure complete rinsing, the final rinse water may be tested for residual CIDEX Industrial Disinfecting Solution by using the CIDEX Industrial Disinfecting Residual Test Kit.

This is a single-use product only. Do not reuse.

*Trademark SURGIKOS 1982