

Reorder No.
2785

Steriliz

Rapid-Acting
Sporicidal
Virucidal
Fungicidal

ESPECIALLY

Ci

28 da

EPA REG. No. 7078-14
EPA EST. No. 36126-PR-1

Active Ingredient
Glutaraldehyde
Inert Ingredients

Keep out of
DANGER

† See Directions for

*TRADEMARK

tion or disinfection of medical instruments and equipment, *thoroughly clean, rinse and rough dry* objects
CIDEX PLUS® 28 Day Solution. Cleanse and rinse the lumens of hollow instruments before filling with this

DIRECTIONS FOR ACTIVATING:

attached vial to this container. Shake, solution is ready to use. *When activated, solution changes color to green.*
the standard bioburden challenge, this product has also demonstrated efficacy in the presence of 2% organic soil
on during reuse. This solution may be used and reused for up to 28 days after activation. Do not use activated solution
days. Efficacy of this product during its use-life can be verified by the CIDEX PLUS 28 Day Solution Effectiveness Test Kit.

DIRECTIONS FOR USE:


Classification: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
repeated use in manual (bucket and tray) systems.

Disinfection: Immerse completely for a minimum of 10 minutes at 20°C or higher to destroy vegetative organisms including
Pseudomonas aeruginosa, pathogenic fungi, and viruses (Poliovirus Type 1, Adenovirus Type 2; Herpes simplex Type 1, 2; Influenza
Type A [WSN/33]; Vaccinia; Coronavirus, Cytomegalovirus; Rhinovirus Type 14; Coxsackievirus B-1) on inanimate surfaces. To
destroy Mycobacterium tuberculosis on inanimate surfaces, check and insure that solution temperature is 25°C, before
immersing completely for a minimum of 20 minutes. Remove equipment from the solution and rinse thoroughly.

Sterilization: Immerse equipment completely for a minimum of 10 hours at 20°C or higher to destroy resistant pathogenic spores
including Clostridium species and Bacillus subtilis. Remove equipment from the solution using sterile technique and rinse
thoroughly with sterile water.

NOTE DATES HERE

MIXING DATE _____ EXPIRATION DATE _____



DISTRIBUTED BY:

SURGIKOS

a Johnson & Johnson company

SURGIKOS, INC., ARLINGTON, TEXAS 76010
*TRADEMARK ©SURGIKOS 1986

Cidex
28 day solution

Activator

For Use With On
Container **CIDEX**
28 Day Solution

M/
Under the
Fungicide
as named
registered
EPA Reg.

*TRADEMARK



RETAKE

BECAUSE OF SUSPECTED OPERATOR OR
EQUIPMENT ERROR, THE PRECEDING PAGE
WAS REFILMED AND APPEARS AS THE NEXT
IMAGE.

RETAKE

Add contents of attached vial to this container. Shake, solution is ready to use. When activated, solution changes color to green. In addition to the standard bioburden challenge, this product has also demonstrated efficacy in the presence of 2% organic soil contamination during reuse. This solution may be used and reused for up to 28 days after activation. Do not use activated solution beyond 28 days. Efficacy of this product during its use-life can be verified by the CIOEX PLUS 28 Day Solution Effectiveness Test Kit.

DIRECTIONS FOR USE:

General Classification: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

For repeated use in manual (bucket and tray) systems.

Disinfection: Immerse completely for a minimum of 10 minutes at 20°C or higher to destroy vegetative organisms including *Pseudomonas aeruginosa*, pathogenic fungi, and viruses (Poliovirus Type 1, Adenovirus Type 2; Herpes simplex Type 1, 2; Influenza Type A [WSN/33]; Vaccinia; Coronavirus, Cytomegalovirus; Rhinovirus Type 14; Coxsackievirus B-1) on inanimate surfaces. To destroy *Mycobacterium tuberculosis* on inanimate surfaces, check and insure that solution temperature is 25°C, before immersing completely for a minimum of 20 minutes. Remove equipment from the solution and rinse thoroughly.


Sterilization: Immerse equipment completely for a minimum of 10 hours at 20°C or higher to destroy resistant pathogenic spores including *Clostridium sporogenes* and *Bacillus subtilis*. Remove equipment from the solution using sterile technique and rinse thoroughly with sterile water.

EPA REG. No. 7078-14
EPA EST. No. 36126-P

NOTE DATES HERE

MIXING DATE _____ EXPIRATION DATE _____

(K)



DISTRIBUTED BY:

SURGIKOS

a Johnson & Johnson company

SURGIKOS, INC., ARLINGTON, TEXAS 76010

*TRADEMARK

©SURGIKOS 1926

Reorder No.
2785

Sterilizing and Disinfecting Solution

Rapid-Acting / Non-Staining / Non-Rusting
Sporicidal Tuberculocidal
Virucidal† Bactericidal
Fungicidal Pseudomonacidal

ACCEPTED
MAY 12 1986
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-14

ESPECIALLY FORMULATED FOR REUSE

Cidexplus*

28 day solution

Active Ingredient
Glutaraldehyde 3.2%
Inert Ingredients 96.8%
Total 100.0%

NOTE — Attached vial contents
must be added to solution before
this product is effective.
See Directions for Activating

Keep out of reach of children.

DANGER: See right side panel for precautions

† See Directions for Use

*TRADEMARK

SURGIKOS

**PRECAUTIONARY
DANGER:**
HAZARD TO HUMANS
FOLLOW ALL LABEL DIRECTIONS

Prevent contact with eyes by using eye protection.
Eyes — Corrosive, causes eye damage. Do not touch.
Immediately and get medical attention.
Ingestion — Harmful if swallowed. Avoid contact with food.
Inhalation — Avoid spillage. Use in ventilated area.
Skin — Avoid skin contact, as possibility of irritation. In case of
skin contact, wash thoroughly with water. Get medical attention.

STORAGE AND DISPOSAL:
Keep in a cool place. Triple rinse empty containers.
Landfill approved for pesticide containers.

Especially recommended for sterilization of rubber, plastic, glass or metallic objects or therapy equipment, nonlensed surgical instruments.

CONTENTS: ONE GAL (128 FLUID OZ.)
MADE IN U.S.A.
PAT. NO. 4,436,754
CONTROL NUMBER ON BOTTOM
2785 C1-0

Cidexplus*

28 day solution

DANGER: Contact with eyes causes damage. In case of contact with eyes, flush with water immediately and get medical attention. For skin contact, flush thoroughly with water.

DISPOSAL: Rinse empty container thoroughly with water and discard it.

Activator Vial

For Use With One Gallon
Container CIDEX PLUS
28 Day Solution Only

NET CONTENTS: 158 GRAMS
MADE IN U.S.A.
PAT. NO. 4,436,754
CONTROL NUMBER ON BOTTOM

ACCEPTED
MAY 12 1986
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-14

*TRADEMARK

SURGIKOS

DISTRIBUTED BY:
Johnson & Johnson company
SURGIKOS, INC., ARLINGTON, TEXAS 76010
2785V1-0 ©SURGIKOS 1986