

PM 31 7078-14

1974

Reorder No.
2785

Sterilizing and Disinfecting Solution

Rapid-Acting / Non-Staining / Non-Rusting

Sporicidal	Tuberculocidal
Virucidal†	Bactericidal
Fungicidal	Pseudomonacidal

ESPECIALLY FORMULATED FOR REUSE

Cidexplus*

28 day solution

EPA REG. No. 7078-14
EPA EST. No. 36128-PP-1

Active ingredient	
Glutaraldehyde	3.2%
Inert ingredients	96.8%
Total	100.0%

NOTE — Attached vial contents must be added to solution before this product is effective.
See Directions for Activating.

Keep out of reach of children.

DANGER: See right side panel for precautions

† See Directions for Use

*TRADEMARK

SURGIKOS

ACCEPTED
MAY 24 1989

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 7078-14

BEST AVAILABLE COPY

204

For complete sterilization or disinfection of medical instruments and equipment, thoroughly clean, rinse and rough dry objects before immersing in CIDEX PLUS* 28 Day Solution. Cleanse and rinse the lumens of hollow instruments before filling with this solution.

DIRECTIONS FOR ACTIVATING:

Add contents of attached vial to this container. Shake, solution is ready to use. When activated, solution changes color to green. In addition to the standard bioburden challenge, this product has also demonstrated efficacy in the presence of 2% organic soil contamination during reuse. This solution may be used and reused for up to 28 days after activation. Do not use activated solution beyond 28 days. Efficacy of this product during its use-life can be verified by the CIDEX PLUS 28 Day Solution Effectiveness Test Kit.

DIRECTIONS FOR USE:

General Classification: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. For repeated use in manual (bucket and tray) systems.

Disinfection: Immerse completely for a minimum of 10 minutes at 20°C or higher to destroy vegetative organisms including Pseudomonas aeruginosa, pathogenic fungi, and viruses (Poliovirus Type 1, Adenovirus Type 2; Herpes simplex Type 1, 2; Influenza Type A [WSN33]; Vaccinia; Coronavirus, Cytomegalovirus; Rhinovirus Type 14; Coxsackievirus B-1) on inanimate surfaces. To destroy Mycobacterium tuberculosis on inanimate surfaces, check and insure that solution temperature is 25°C, before immersing completely for a minimum of 20 minutes. Remove equipment from the solution and rinse thoroughly.

Sterilization: Immerse equipment completely for a minimum of 10 hours at 20°C or higher to destroy resistant pathogenic spores including Clostridium sporogenes and Bacillus subtilis. Remove equipment from the solution using sterile technique and rinse thoroughly with sterile water.

NOTE DATES HERE

MIXING DATE _____ EXPIRATION DATE _____



DISTRIBUTED BY:

SURGIKOS

a Johnson & Johnson company

SURGIKOS, INC., ARLINGTON, TEXAS 76010

*TRADEMARK

© SURGIKOS 1988

BEST AVAILABLE COPY

ACCEPTED
MAY 24 1989
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-14

3074

PRECAUTIONARY STATEMENTS

DANGER:

HAZARD TO HUMANS

FOLLOW ALL LABEL DIRECTIONS

Prevent contact with eyes by using eye protection when handling or pouring.

Eyes — Corrosive, causes eye damage. Do not get in eyes. In case of contact, flush with water immediately and get medical attention.

Ingestion — Harmful if swallowed. Avoid contamination of food.

Inhalation — Avoid spillage. Use in ventilated area.

Skin — Avoid skin contact, as possibility of sensitization exists. Can cause skin irritation. For skin contact, wash thoroughly with water. Get medical attention if irritation occurs.

STORAGE AND DISPOSAL:

Keep in a cool place. Triple rinse empty container with water and dispose in an incinerator or landfill approved for pesticide containers.

Especially recommended for sterilization or disinfection of...

Rubber, plastic, glass or metallic objects or instruments such as: anesthesia and respiratory therapy equipment, nonreusable surgical instruments, dental instruments, mirrors, burs and thermometers

CONTENTS: ONE GAL. (128 FLUID OZ.)
MADE IN U.S.A.
PAT. NO. 4,438,775
CONTROL NUMBER ON BOTTOM
2785C1-0

SURGIKOS

ACCEPTED
MAY 24 1989
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-14

BEST AVAILABLE COPY

4074

Cidexplus*

28 day solution

DANGER: Contact with eyes causes damage. In case of contact with eyes, flush with water immediately and get medical attention. For skin contact, flush thoroughly with water.

DISPOSAL: Rinse empty container thoroughly with water and discard it.

Activator Vial

For Use With One Gallon Container CIDEX PLUS
28 Day

NET CONTENTS: 158 GRAMS
MADE IN U.S.A
PAT. NO. 4,436,754
CONTROL NUMBER ON BOTTOM

DISTRIBUTED BY:

SURGIKOS

a Johnson & Johnson company

SURGIKOS, INC., ARLINGTON, TEXAS 76010
2785 V1-0 © SURGIKOS 1986

*TRADEMARK

ACCEPTED
MAY 24 1988
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-14

BEST AVAILABLE COPY

481

7008-14

1000

Reorder No.
2785

Sterilizing and Disinfecting Solution

Rapid-Acting / Non-Staining / Non-Rusting

Sporicidal
Virucidal
Fungicidal

Tuberculocidal
Bactericidal
Pseudomonacidal

ESPECIALLY FORMULATED FOR REUSE

Cidexplus*

28 day solution

EPA REG. No. 7078-14
EPA EST. No. 38128-PR-1

Active Ingredient	
Glutaraldehyde	3.2%
Inert Ingredients	96.8%
Total	100.0%

NOTE — Attached vial contents must be added to solution before this product is effective. See Directions for Advertising.

Keep out of reach of children.

DANGER: See right side panel for precautions

† See Directions for Use

*TRADEMARK

SURGIKOS

ACCEPTED
MAY 24 1980

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. **7078-14**

5000

hbc

10-8-80

For complete sterilization or disinfection of medical instruments and equipment, thoroughly clean, rinse and rough dry objects before immersing in CIDEX PLUS® 28 Day Solution. Cleanse and rinse the lumens of hollow instruments before filling with this solution.

DIRECTIONS FOR ACTIVATING:

Add contents of attached vial to this container. Shake, solution is ready to use. When activated, solution changes color to green. In addition to the standard bioburden challenge, this product has also demonstrated efficacy in the presence of 2% organic soil contamination during reuse. This solution may be used and reused for up to 28 days after activation. Do not use activated solution beyond 28 days. Efficacy of this product during its use-life can be verified by the CIDEX PLUS 28 Day Solution Effectiveness Test Kit.

DIRECTIONS FOR USE:

General Classification: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. For repeated use in manual (bucket and tray) systems.

Disinfection: Immerse completely for a minimum of 10 minutes at 20°C or higher to destroy vegetative organisms including *Pseudomonas aeruginosa*, pathogenic fungi, and viruses (Poliovirus Type 1, Adenovirus Type 2; Herpes simplex Type 1, 2, Influenza Type A [WSN/33]; Vaccinia; Coronavirus, Cytomegalovirus; Rhinovirus Type 14; Coxsackievirus B-1) on inanimate surfaces. To destroy *Mycobacterium tuberculosis* on inanimate surfaces, check and insure that solution temperature is 25°C, before immersing completely for a minimum of 20 minutes. Remove equipment from the solution and rinse thoroughly.

Sterilization: Immerse equipment completely for a minimum of 10 hours at 20°C or higher to destroy resistant pathogenic spores including *Clostridium sporogenes* and *Bacillus subtilis*. Remove equipment from the solution using sterile technique and rinse thoroughly with sterile water.

NOTE DATES HERE

MIXING DATE _____ EXPIRATION DATE _____



DISTRIBUTED BY:

SURGIKOS

a Johnson & Johnson company

SURGIKOS, INC., ARLINGTON, TEXAS 76010

TRADEMARK

© SURGIKOS 1989

ACCEPTED
MAY 24 1989
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-14

1.11

384

7078-14

PRECAUTIONARY STATEMENTS

DANGER:

HAZARD TO HUMANS

FOLLOW ALL LABEL DIRECTIONS

Prevent contact with eyes by using eye protection when handling or pouring.

Eyes — Careful, avoid eye damage. Do not get in eyes. In case of contact, flush with water immediately and get medical attention.

Ingestion — Hazard if swallowed. Avoid contamination of food.

Inhalation — Avoid spillage. Use in ventilated area.

Skin — Avoid skin contact, as possibility of sensitization exists. Can cause skin irritation. For skin contact, wash thoroughly with water. Get medical attention if irritation occurs.

STORAGE AND DISPOSAL:

Keep in a cool place. Triple rinse empty container with water and dispose in an incinerator or landfill approved for pesticide containers.

Especially recommended for sterilization or disinfection of ...

Rubber, plastic, glass or metallic objects or instruments such as: anesthesia and respiratory therapy equipment, nonresected surgical instruments, dental instruments, mirrors, burs and thermometers

CONTENTS: ONE GAL (128 FLUID OZ.)
MADE IN U.S.A.
PAT. NO. 4,438,754
CONTROL NUMBER ON BOTTOM
2788 C1-8

SURGIKOS

ACCEPTED
MAY 24 1989
Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 7078-14

PM 31

484
41-2102

Cidexplus*

28 day solution

DANGER: Contact with eyes causes damage. In case of contact with eyes, flush with water immediately and get medical attention. For skin contact, flush thoroughly with water.

DISPOSAL: Rinse empty container thoroughly with water and discard it.

Activator Vial

For Use With One Gallon
Container CIDEX PLUS
28 Day

NET CONTENTS: 150 GRAMS
MADE IN U.S.A.
PAT. NO. 4,436,754
CONTROL NUMBER ON BOTTOM

DISTRIBUTED BY:

SURGIKOS

a Johnson & Johnson company

SURGIKOS, INC., ARLINGTON, TEXAS 76010

2785 V10

© SURGIKOS 1988

*TRADEMARK

