



Form No  
SIV-1-82

**Matheson**  
A SLOAN COMPANY

East Rutherford, N.J. 07073 / Morrow, Ohio 44660  
Greenville, Massachusetts 01533 / Lowell, Massachusetts 01854  
La Porte, Texas 77571 / Covington, Louisiana 70429  
Newark, California 94568 / Liberty, Maryland 21527  
Bridgport, New Jersey 08001 / Go Park, Louisiana 70737  
Dayton, Ohio 45424 / Wilmington, Ohio 45397

**STERILIZING  
GAS MIXTURE**

FOR INDUSTRIAL AND MEDICAL USE ONLY  
IN COMMERCIAL GAS STERILIZER BY MEDICAL  
PROFESSIONALS OR APPROPRIATELY TRAINED  
TECHNICIANS.

ACTIVE INGREDIENT 20%  
ETHYLENE OXIDE

INERT INGREDIENT 80%  
CARBON DIOXIDE  
EPA REG. NO. 7256-2

Gross Wt. \_\_\_\_\_ lbs.

Tare Wt. \_\_\_\_\_ lbs.

Net Wt. \_\_\_\_\_ lbs.

Filling Date \_\_\_\_\_



**DANGER  
GAS AND LIQUID  
UNDER PRESSURE  
DANGER**

**STORAGE AND DISPOSAL**

Use and store in cool, dry, well ventilated place. Keep cylinder temperature from exceeding 125°F. Close valve when not in use. Never refill cylinder. Handle cylinder with care. Discharge in normal vertical position. After discharge close cylinder valve to prevent suck-back of contaminants. Ethylene Oxide is subject to polymerization, use within 2 months after receipt.

Avoid breathing vapors. Avoid contact with skin & eyes. In case of contact, flush thoroughly with water for at least 15 minutes. For eyes, get medical attention. Avoid contact with oxidizing materials.

Do not reuse empty cylinder. Return cylinder to originating shipping point.

**READ MATHESONS STERILIZING  
GAS MIXTURE BROCHURE  
BEFORE USING**

RECEIVED  
1982  
7256-2

BEST DOCUMENT AVAILABLE

BEST EQUIPMENT AVAILABLE

# STERILIZING GAS MIXTURES

For Industrial And Medical Use  
Only In Commercial Gas Sterilizer

RECEIVED  
6-1-62

Li A. Eng. No. 7256-24

From

**Matheson**  
Securities, Inc.

141 East 60th St., N.Y.C. 22, N.Y. 10022  
Dayton, OH 45424  
Dallas, TX 75221  
Denver, CO 80202  
Detroit, MI 48201  
Houston, TX 77001  
Los Angeles, CA 90001  
Miami, FL 33101  
New York, NY 10022  
Philadelphia, PA 19101  
Pittsburgh, PA 15201  
Portland, ME 04101  
San Francisco, CA 94101  
Seattle, WA 98101  
St. Louis, MO 63101  
Washington, DC 20001

## CONTENTS

	PAGE
Description of Sterilizing Gas Mixtures .....	3
Types of Sterilizing Gases Available .....	3
Advantages of Sterilizing Gas Mixtures .....	3
Gas Sterilization Procedures and Conditions .....	4
Handling and Storage of Cylinders .....	6
Toxicity .....	7
First Aid .....	7
Economics .....	8
References .....	8

BEST COPY AVAILABLE

## STERILIZING GAS MIXTURES

**DESCRIPTION** — Sterilizing gas mixtures as provided by Mallinckrodt are composed of ethylene oxide, which is the active sterilizing agent, plus inert gases such as Carbon Dioxide, or a halogenated hydrocarbon, or combination of halogenated hydrocarbons.

Ethylene Oxide is a true bactericidal agent. It acts as a germicide on both bacterial and fungal molds. The addition of proper diluents acts to increase the volatility of the sterilizing gas and decrease the flammability associated with the handling of Ethylene Oxide.

### Penetrative Properties

Sterilizing gases containing Ethylene Oxide have high penetrative qualities. They are capable of readily diffusing through paper, clothing and fabrics. Penetration of powders, dried drugs, ointments, gums depends on the method of packaging and the density of the material itself. Sterilization of prepackaged material is possible through proper use of packaging materials. The gas is known to penetrate polyethylene and polyfilm.

### Types of Sterilizing Gases Available

The following typical Ethylene Oxide mixtures are available from all branches of Mallinckrodt:

- a. 10% Ethylene Oxide — 90% Carbon Dioxide
- b. 20% Ethylene Oxide — 80% Carbon Dioxide
- 12% Ethylene Oxide — 88% Freon-12<sup>™</sup>

Mixtures containing 10% by weight or less of Ethylene Oxide in Carbon Dioxide, or 12% by weight or less of Ethylene Oxide in Freon-12 are completely non-flammable in proportions with air. The use of mixtures containing 20% Ethylene Oxide in the central area of a container is recommended for the sterilization of special handling powders and equipment drags.

The cylinder pressures of sterilizing mixtures will vary with the concentrations of the various constituents. A higher proportion of diluents leads to lower vapor pressures which will cause the total mixture to have

a higher vapor pressure. The vapor pressures at 70°F. (cylinder pressures) of the above mixtures are as follows:

MIXTURE	VAPOR PRESSURE (p.s.i.g.)
10% Ethylene Oxide — 90% Carbon Dioxide	720
20% Ethylene Oxide — 80% Carbon Dioxide	640
12% Ethylene Oxide — 88% Freon-12 <sup>™</sup>	52

Vapor pressure or cylinder pressure should not be confused with sterilizer pressure. Sterilizer pressure is dependent on the amount of sterilizing mixture that appears to be in the gas phase in the sterilizer; however, cylinder pressure is of value in determining whether or not enough liquid mixture is available to have the sterilizer maintain an atmosphere in a so-called sterilizer equipment.

### ADVANTAGES OF STERILIZING GAS MIXTURES

These mixtures have been used to sterilize heat and moisture sensitive materials. Gas sterilization can substantially reduce the time or temperature normally subjected to steam or liquid temperature sterilization processes.

Among the materials that can be sterilized and transported by Ethylene Oxide mixtures are the following:

**PHARMACEUTICAL**  
 Plastic vials, bottles, and droppers  
 Antibiotics and other suitable drugs  
 Rubber materials

MEDICAL	OTHER
Chloroform	Soil
Starches	Straw
Surgical instruments	Bands
Drugs and ointments	Leather
Food packaging	Paper
Chemicals	Electrical equipment
Culture media	Motor oil
	Vaseline

BEST DOCUMENT AVAILABLE

## STERILIZING GAS MIXTURES

Below are given the characteristics of sterilization. The information presented is for general information only. It is suggested that all factors be weighed carefully prior to their adoption.

**GAS STERILIZATION**—Sterilization with Ethylene Oxide mixtures has been practiced for some time. Because of its effectiveness and its advantages it is now commonly used by the surgical dressing and pharmaceutical industries. Some of the most important advantages are as follows:

1. They are not harmful to most materials.
2. They have high penetrative properties.
3. They are readily removed by aeration.
4. They destroy all forms of micro-organisms.
5. They are readily available.
6. They are convenient to store and handle.
7. They reduce the mutagenicity hazard inherent in the use of pure Ethylene Oxide.

### Conditions of Sterilization

Sterilization by Ethylene Oxide mixtures is essentially a function of time of exposure, concentration of Ethylene Oxide in the sterilizer space, temperature and humidity. Other conditions affecting sterilization are the moisture content of the microbial cell at the time of exposure and the type of organism to be destroyed.

### A. TIME OF EXPOSURE

In general, non-spore forming bacteria are killed with relative ease, while a longer exposure is required for a 100% kill of spore-forming bacteria. Generally, exposure time is a function of the Ethylene Oxide concentration in the sterilizer space. The following table illustrates relative exposure times to kill various types of organisms:

SPECIES OF BACTERIA	RELATIVE TIME (X UNIT TIME)
Bacillus globigii spores	5X
Staphylococcus aureus	3X
Mycobacterium phlei	3X
Gaika tetragina	3X
Serratia marcescens	2X
Fibrobacter	X
Klebsiella pneumoniae	X
Escherichia coli	X
Clostridium	0.5X
Bacillus subtilis	0.5X

### B. CONCENTRATION

For sterilization purposes, Ethylene Oxide concentrations of 450-600 mg/l. of chamber space are required to effect a 100% kill of micro-organisms. The higher the concentration of Ethylene Oxide, the shorter the time of exposure required for sterilization, usually doubling the concentration reduces the exposure time by one-half. The "Etron" mixtures have the added advantage of giving a higher volume percent of Ethylene Oxide in the gas phase which still remaining non-flammable, compared with Carbon Dioxide mixtures of similar weight percent, thus offering more rapid sterilization at low chamber pressures.

Fig. 1 shows the variation of chamber pressure for three typical gas mixtures at a few commonly used temperatures, as a function of Ethylene Oxide concentration in the sterilizer space. These curves have been developed assuming the mixtures behave as ideal gases. This graph assumes complete evacuation of the sterilizer before addition of the sterilizing gas. Initial sterilizer absolute pressure should be added to the graph value to determine final sterilizer pressure. The use of the chart is shown by the following example:

What sterilizer pressure is required to utilize a 20% Ethylene Oxide - 80% Carbon Dioxide mixture at a concentration of 600 mg/l. of Ethylene Oxide and 1200 mg/l. of Carbon Dioxide, assuming initial sterilizer pressure of 2 psia?

Follow vertically from 600 on the horizontal scale to the line representing 130 F. in the group of lines representing 20% Ethylene Oxide - 80% Carbon Dioxide mixture. At this intersection, follow horizontally to the left to the pressure scale reading the pressure of approximately 27 psia. Add to this the initial pressure of 2 psia = total of 29 psia. Subtracting 14.7 psia gives the sterilizer final gauge pressure of 14.3 psig.


The following equation can be used to calculate the lbs. of sterilizing gas mixture needed for a sterilizer.

$$\text{lbs. mixture needed} = \frac{0.00524CV}{L}$$

C = Concentration of Ethylene Oxide in Sterilizer space, mg/l.

L = Percent Ethylene Oxide concentration in sterilizing mixture.

V = Sterilizer volume, cu ft.



BEST COPY AVAILABLE

STERILIZING GAS MIXTURES

100  
150  
200

1  
7

## STERILIZING GAS MIXTURES

For example, how many pounds of 25% Ethylene Oxide in helium should be contained in a cylinder required for a 10 cubic foot sterilization cycle, assuming a final Ethylene Oxide concentration of 800 mg/l in the sterilize space?

$$\text{lbs. mixture needed} = \frac{0.025 \times 100 \times 10}{11} = 4.5 \text{ lbs.}$$

This assumes neither leak nor absorption. From a practical standpoint some excess should be used.

### C. TEMPERATURE

While sterilization can be effected with Ethylene Oxide at 20°C, rather low exposure periods are required. As temperature is raised, the rate of sterilization is increased. Generally, for every 30°C rise in temperature, the time required for sterilization is halved. Phillips has demonstrated the temperature effect in his studies on the sterilization of bacterial spores at 30°C, 25°C, and 17°C, in which the increase in temperature afforded a considerable reduction in exposure time. He showed that the effect of Ethylene Oxide to sterilize is increased by a factor of 2.1 for each 10°C rise in temperature, with the attainment of a complete kill of the bacterial spores. Temperatures in the range of 150-170°F are usually employed.

### D. HUMIDITY

Humidity is definitely related to the sterilization efficiency for Ethylene Oxide. It has been shown that relative humidity in the range of 40-70% is optimal for maximum sterilization efficiency. It is important to avoid the gas exposure of dry surfaces, especially in a dry state. If the relative humidity is too low, the gas will not penetrate to very high concentrations in the sterilization space, and the sterilization will be less assured. In practice, humidity control is not necessary, but the cycle should be a minimum of one hour.

### E. RESIDUAL GAS

All materials should be dried prior to use. The release of residual Ethylene Oxide is dependent on a number of factors, including the type of packaging material. The amount of Ethylene Oxide and relative humidity during the sterilization cycle will reduce the amount of residual Ethylene Oxide more rapidly as compared to the type of material used.

room conditions. A thorough analysis of the distribution of Ethylene Oxide in the sterilization cycle is essential to ensure that order kinetics. In all cases, the user should verify that his procedure reduces residual to the acceptable level.

A number of analytical procedures for the analysis of Ethylene Oxide residues are available and can be directed to the literature for the method most suited to individual needs.

The standardization of maximum acceptable limits for Ethylene Oxide residues is being prepared and will be issued by the Federal regulatory agencies.

See Section on Toxicity for more detailed information regarding the effects of Ethylene Oxide.

### Designing of the Gas for Sterilization

In designing the gas system, either a complete sterilization cycle or a minimum 20°C exposure period should be used. The gas mixture should be contained in the cylinder as a liquid, and inverted over small cylinders or containers. The gas should be filled with the Ethylene Oxide mixture, and the cylinder should be sealed. The cylinder should be sterilized by the warm water method, or by the dry heat method. The cylinder should be cooled to the sterilization temperature before use. The gas should be used in a sterile environment, and the use of an inverted cylinder should be avoided. The gas should be used in a sterile environment, and the use of an inverted cylinder should be avoided.

### References

A number of references are given in this section for the user's information. The references are given in alphabetical order. The references are given in alphabetical order. The references are given in alphabetical order. The references are given in alphabetical order. The references are given in alphabetical order.

LIBRARY

## STERILIZING GAS MIXTURES

### HANDLING AND STORAGE OF CYLINDERS

The following general rules should be followed in the handling and storage of sterilizing gas mixtures:

1. Never drop cylinders or permit them to strike each other violently.

2. Cylinders should be assigned a definite area for storage. The area should be dry, cool, well ventilated and preferably free of dirt, grease, oil, and other contaminants. From excessive temperature rise existing in the cylinders from radiators or other sources of heat, storage conditions should comply with local and State regulations.

3. The valve protection cap should be kept in place until the cylinder has been secured against a wall or bench, or placed in a cylinder stand and is ready to be used.

4. Avoid dragging, pushing or pulling cylinders, even for a short distance. They should be moved by using a suitable hand truck.

5. Never tamper with safety devices on valves or cylinders.

6. As soon as a cylinder has been discharged, close the cylinder valve. Replace any valve cap or protective cap or plug supplied with the cylinder. Mark or label the cylinder EX-119.

7. No part of the cylinder should be subjected to a temperature higher than 125°F. Temperature increases of this figure may cause the cylinder to become so full of gas that it will burst if not immediately vented. Never permit a flame to come in contact with any part of a compressed gas cylinder.

8. Contact the gas supplier if a part of the cylinder or valve develops leakage that cannot be completely repaired, for example, by tightening of a parking nut to prevent leakage through valve packing.

9. Federal law prohibits the owner's disposal of cylinder ending in 119's. The owner may either return empty cylinder to original supplier or

**TOXICITY**—The toxic effects of sterilizing gas mixtures stem mainly from the Ethylene Oxide. The method of use described in the present manual is considered to be as simple and safe as that of spraying a toxin present in large amounts. Thus, the amount of gas that is inhaled to the lungs. Certain U. S. Department of Health and Safety Administration regulations require 10 ppm as the maximum allowable concentration of Ethylene Oxide for daily exposure to man. The American Conference of Governmental Hygienists has recommended

recommended 10 ppm. However, in the U. S. of evidence concerning Ethylene Oxide, a cancer center in Montreal has reported a case of cancer in a rat. The A.C.G.I.H. has also reported a case in a rat. The first warning of Ethylene Oxide is the irritation of the eye, and more precautionary must be taken as warning of a high concentration. To avoid serious irritation, avoid high concentrations of vapor over extended periods of time. A result in dizziness, nausea and vomiting. Exposure to low concentrations, if continuous, will result in a numbness of the sense of smell, and irritation of the skin. Irritation without warning. High concentrations of Ethylene Oxide can produce edema of the lungs and irritation of the eyes. Ethylene Oxide is not a flammable gas, is not cumulative to the body and therefore is not a poison due to Ethylene Oxide. However, Ethylene Oxide can act as a vesicant. It is a blistering agent if its vapors are in contact with the skin for a period of time. The pure liquid, if spilled on the skin, can cause frostbite to cause an ulcer. But if it is in contact with the skin, blisters will appear. Rubber and leather articles, in particular, should be stored in a dry place after being exposed to a spill of high concentration of Ethylene Oxide vapors.

High Ethylene Oxide concentrations may cause severe eyeburn.

Protective clothing should be worn when working with the gas. The gas is not a fire hazard. Each cylinder should be checked for leaks. When exposed to the liquid, the skin will be irritated. Burns occur only as a delayed reaction. When working with the gas, no special precautions are required other than to insure proper ventilation.

**PRECAUTIONS**—Because of the nature and use of Ethylene Oxide sterilization mixtures, the possibility of being exposed to excess concentrations of the gas is of the mixture is not high. However, the following first aid suggestions are given for those who are severely exposed. This is not a standard first aid procedure for Ethylene Oxide.

If concentrations of the gas are believed to be common products of industry, then the amount in any particular case should be determined by a competent analytical laboratory. Some of the physical and chemical properties of Ethylene Oxide are given in the following table. It is a colorless, odorless, non-flammable gas. It is a highly flammable gas. It is a highly flammable gas. It is a highly flammable gas. It is a highly flammable gas.

BEST DOCUMENT AVAILABLE



## STERILIZING GAS MIXTURES

### Inhalation

In case breathing has stopped, artificial respiration should be started at once. After a severe exposure, pure Oxygen should be administered as soon as possible by trained personnel. The patient should breathe pure Oxygen under positive exhalation pressure for one-half hour periods every hour for at least 24 hours. If there are no signs of breathing at the end of this period, artificial breathing is necessary and the color of the skin and mucous membranes is good. Oxygen inhalation may be discontinued if there is no improvement; the patients should be kept comfortably warm, but not hot.

### Skin Contact

All clothing contaminated with liquid Ethylene Oxide or exposed to high concentrations of Ethylene Oxide vapors should be removed at once. Clothing, including shoes, should not be worn again until thoroughly aerated. Shoes can be decontaminated, and it is suggested that they be not worn again. All affected areas of skin should be thoroughly washed with soap and water.

### Eyes

Ethylene Oxide, as liquid or vapor, is capable of penetrating

during eye contact. Should it reach the eyes, the eyes should be irrigated copiously with water for fifteen minutes.

The eyes should be held under running water for fifteen minutes to insure contact of the water with all tissues of the eye, coat the eyes and face with a liberal amount of petrolatum jelly, and the eyes should be irrigated for a second period of fifteen minutes. If possible, an eye specialist should be called in attendance.

### ECONOMICS

Once the concentration of Ethylene Oxide is established to a high degree of accuracy by calibration (1) will then supply the quantity of the quantity of any sterilizing gas mixture which is required for any sterilization process. This quantity and process may be determined by the following methods:

All Matheson branches are equipped to supply all types of sterilizing gas mixtures and will be very happy to meet any special requirements of your business meeting your requirements.

### REFERENCES

1. Gabel, H. J. (1954) J. Biol. Chem. 101: 1-10.
2. H. J. Gabel, C. E. H. (1954) J. Biol. Chem. 101: 11-15.
3. K. J. Gabel, C. E. H. (1954) J. Biol. Chem. 101: 16-20.
4. J. R. P. (1954) J. Biol. Chem. 101: 21-25.
5. J. R. P. (1954) J. Biol. Chem. 101: 26-30.
6. J. R. P. (1954) J. Biol. Chem. 101: 31-35.
7. J. R. P. (1954) J. Biol. Chem. 101: 36-40.
8. J. R. P. (1954) J. Biol. Chem. 101: 41-45.
9. J. R. P. (1954) J. Biol. Chem. 101: 46-50.
10. J. R. P. (1954) J. Biol. Chem. 101: 51-55.
11. J. R. P. (1954) J. Biol. Chem. 101: 56-60.
12. J. R. P. (1954) J. Biol. Chem. 101: 61-65.
13. J. R. P. (1954) J. Biol. Chem. 101: 66-70.
14. J. R. P. (1954) J. Biol. Chem. 101: 71-75.
15. J. R. P. (1954) J. Biol. Chem. 101: 76-80.
16. J. R. P. (1954) J. Biol. Chem. 101: 81-85.
17. J. R. P. (1954) J. Biol. Chem. 101: 86-90.
18. J. R. P. (1954) J. Biol. Chem. 101: 91-95.
19. J. R. P. (1954) J. Biol. Chem. 101: 96-100.
20. J. R. P. (1954) J. Biol. Chem. 101: 101-105.

Matheson provides researchers and industry with their specialty gas requirements.

We take pride in our capability to provide gases of high purity in quantities ranging from a few grams to a tube trailer. We make gas cylinders ranging from two cubic feet to ten cubic feet, and also up to 20 cubic feet in size, in a wide variety of materials.

We carry in stock a full line of high quality equipment for regulating, measuring, monitoring and controlling gases.

Our catalog is yours for the asking at any of our branch plants.

BEST DOCUMENT AVAILABLE

100182

**MATHESON**  
A FARMER COMPANY

**STERILIZING GAS MIXTURE**

FOR INDUSTRIAL AND MEDICAL USE ONLY IN  
COMMERCIAL GAS STERILIZER BY MEDICAL  
PROFESSIONALS OR APPROPRIATELY  
TRAINED TECHNICIANS

ACTIVE INGREDIENT ETHYLENE OXIDE 20%  
INERT INGREDIENT CARBON DIOXIDE 80%

**DANGER**


Flammable compressed gas in form explosive mixture  
with air. Keep away from fire, heat and open flames.  
See other precautions on the back of this tag and small,  
cylinder valve tag.

EPA EST 7256 CA-1-2 GA-1 IL-1 MA-1 NJ-1 TX-1

THIS CYLINDER IS  
THE PROPERTY OF  
MATHESON.

EPA REG NO 72562  
THIS CYLINDER IS  
EQUIPPED WITH  
VALVE CHECK  
CGA 350

DO NOT  
REFILL



DO NOT  
REMOVE THIS TAG

If this tag is not  
present on the  
refill of this cylinder  
to Matheson

BEST AVAILABLE

**STERILIZING GAS MIXTURE**  
FOR INDUSTRIAL AND MEDICAL USE ONLY IN  
COMMERCIAL GAS STERILIZER BY MEDICAL  
PROFESSIONALS OR APPROPRIATELY  
TRAINED TECHNICIANS.

**DANGER: GAS AND LIQUID  
UNDER PRESSURE**

EPA REG NO 70762

**STORAGE AND DISPOSAL**

Use and store in cool, dry, well-ventilated area. Keep cylinder temperature from exceeding 125°F. Close valve when not in use. Never refill cylinder. Handle cylinder with care. Discharge in vent, inverted position. After discharge close cylinder valve to prevent suck back of contaminants. Ethylene oxide is subject to polymerization, use within 2 months after receipt.

**DO NOT REUSE EMPTY CYLINDERS  
RETURN TO ORIGINATING SHIPPING POINT.  
READ MATHESON'S STERILIZING GAS MIXTURES  
BROCHURE BEFORE USING.**

**BEST DOCUMENT AVAILABLE**

SHIP TO: \_\_\_\_\_

SHIP TO:

7256-21

RETURN TO MATHESON