

### CAUTION

Not to be used for sterilization purposes, not for disinfecting of instruments contaminated with sporebearing pathogens and rusty instruments.

Scrub instruments with soap and water until free of organic matter, and rinse thoroughly.

Contaminated instruments and hypodermic needles should always be sterilized by boiling or autoclaving.

For technical use of hospitals, medical, dental and veterinary professions.

Rinse empty container thoroughly with water and discard it.

EPA Est. No. 07601CA01

# GERMICIDAL SOLUTION

ALCOHOLIC  
RUST PREVENTING

NON-CORROSIVE TO INSTRUMENTS OR INANIMATE SURFACES

### ALCOHOLIC SOLUTION 1:750

#### ACTIVE INGREDIENTS:

Isopropanol ..... 44.000% by wt  
Alkyl (C14, 50%; C12, 40%; C16, 10%) Dimethyl Benzyl Ammonium Chloride ..... 0.133% by wt

INERT INGREDIENTS ..... 55.867% by wt

Contains Sodium Nitrite

### CAUTION

### KEEP OUT OF REACH OF CHILDREN

#### Not To Be Taken Internally

Harmful if swallowed. Avoid getting the concentrate in eyes. In case of contact, flush with plenty of water. If irritation persists, get medical attention. Avoid contamination of foods.

DO NOT USE, POUR, SPILL OR STORE NEAR HEAT OR OPEN FLAME

NET CONTENTS 1 GALLON

LIBBY LABORATORIES INC.  
2214 Sixth Street  
Berkeley, Ca. 94710

### DIRECTIONS

Instruments must be thoroughly cleaned of contaminating materials such as all adhering blood and serous exudates prior to immersion in the disinfecting solution. Immersion time must be at least 15 minutes or longer.

Hollow instruments should be slightly tilted in the solution to displace the trapped air pockets. Upon removal from the solution and before storing instruments, wipe dry with sterile towel.

Do not use with soap or other cleaning agents.

EPA Reg. No. 7601-11

**ACCEPTED**  
  
MAY 24 1976  
  
Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the purpose of EPA Reg. No. 7601-11 registered under