

MAR-V-CIDE
DISINFECTANT AND GERMICIDE

ACTIVE INGREDIENTS:
 n-alkyl (C12-40%, C14-50%, C16-10%) dimethyl benzyl ammonium chlorides...10%
 Isopropyl Alcohol.....20%

INERT INGREDIENTS: 70%

Phenol Coefficients: 50 (Salmonella typhosa) 75 (Staphylococcus aureus) (Std. A. O. A. C. Method)

Non-flammable, non-volatile, odorless, non-corrosive, non-irritating in use dilution.

Does not contain carbolic acid, iodine, available chlorine, mercury or other metals.

Copy "A"

FOR PROFESSIONAL USE ONLY
WILLIAM MARVY CO.
 1540 ST. CLAIR AVE.
 ST. PAUL, MINN. 55105
 NET CONTENTS: 1 PINT

size: 1-7/8 x 5 1/4

Copy "A"
CAUTION (10 point - caps) Keep out of reach of children (6 point)
 (See other warnings on back panel)

USDA REGISTRATION NO. 8296-1

ACCEPTED
 10-6-67

UNDER THE FEDERAL INSECTICIDE FUNGICIDE AND RODENTICIDE ACT FOR ECONOMIC PEST CONTROL REGISTERED UNDER NO. 8296-1

DIRECTIONS
DISINFECTION

Mix 1 tablespoon Mar-V-Cide in one quart of water or make stock solution by mixing 2 ounces (1/4 cup) Mar-V-Cide in one gallon of water.

(Makes solution containing 1600 ppm active quaternary)

Completely immerse cleaned brushes, combs, razors, clipper heads, scissors, manicure instruments, etc. for 10 minutes (or as required by local health authorities). Change Mar-V-Cide solution every two weeks. For a 2340 ppm solution, use 1 1/2 tablespoons Mar-V-Cide to one quart of water.

COLD DISINFECTION OF SURGICAL INSTRUMENTS

Clean thoroughly, rinse, and immerse in a solution of one ounce of Mar-V-Cide to one gallon of water for ten minutes.

Mar-V-Cide is a concentrated disinfectant. Do not take internally. Avoid contamination of food. Avoid getting concentrate in eyes; in case of contact, immediately flush eyes thoroughly with plenty of water for at least 15 minutes, and get medical attention.

Do not mix with soaps or synthetic detergents or other disinfectants.