

Continuous Residual Activity  
up to 6 Months\*



**Sporicidin**<sup>®</sup>  
BRAND

**DISINFECTANT  
SOLUTION**

ALCOHOL  
FREE

**CLEANS • DISINFECTS • DEODORIZES**

- BACTERICIDAL
- VIRUCIDAL\*\*
- KILLS MOLD and MILDEW\*\*\*
- FUNGICIDAL
- TUBERCULOCIDAL
- ELIMINATES ODORS

**Active Ingredients**

Phenol .....1.56%  
 Sodium Phenate .....0.06%  
 Inert Ingredients .....98.38%  
 Total.....100%

EPA Reg. No. 8383-3  
 EPA Est. No. 047179-IL-001  
 Net Contents: 1 U.S. Gal. (3.79 ltr.)

Keep Out of Reach of Children

**CAUTION**

See additional precautions on back panel

**Non-corrosive • Non-flammable • Non-staining**

**ACCEPTED**  
 FEB 20 2003  
 Under the Federal Insecticide, Fungicide, and  
 Rodenticide Act as amended, for the  
 pesticide, registered under  
 EPA Reg. No. 8383-3



# PRESOAK DISINFECT

An effective cleaner, disinfectant and active deodorant for hospitals, clinics, medical and veterinary offices, laboratories, industrial clean rooms, homes, nursing homes, ambulances, hotels, restaurants, schools, airplanes, trains, boats, autos, buses, health spas and toilets.

**DISINFECTANT AND DEODORIZER FOR HARD, NON-POROUS SURFACES** in medical treatment and patient rooms, operating rooms, emergency rooms, recovery rooms and bathrooms. Use on medical and dental equipment such as beds, surgical carts, counter tops, manikins, telephones, furniture, wheelchairs, walkers, sinks, floors, walls, light switches, linen hampers, shower stalls, toilet seats, trash containers and damp and musty areas, including animal enclosures. Reduces the hazards of cross-contamination.

**ACTIVE DEODORIZER\***: Provides continuous residual bacteriostatic activity against odor-causing organisms for *up to 6 months*, if not removed from surfaces including plastics, latex, vinyl, glass, wood, metal, porcelain, tile and paint.

- Kills Trichophyton mentagrophytes (Virus that causes Athlete's Foot) on bathroom, shower and pool surfaces
- Penetrates cracks
- Non-staining
- Contains no harmful abrasives
- Self-sanitizes surfaces against odor-causing bacteria.

**SPORICIDIN Brand DISINFECTANT SOLUTION** is germicidal when used as directed for disinfection but differs in composition from other Sporidicin products.

**PASSES AOAC EFFICACY STANDARDS FOR INSTITUTIONAL AND HOSPITAL DISINFECTION**

- Kills HIV-1 (AIDS virus) *in one minute*. (See Sporidicin Bulletin No. 303)
- Kills Staphylococcus aureus, Salmonella choleraesuis, Pseudomonas aeruginosa, Methicillin resistant Staphylococcus aureus (MRSA), Vancomycin resistant Enterococcus faecium (VRE) and Trichophyton mentagrophytes (Athlete's Foot fungus) *in three minutes*.
- Kills vegetative organisms including Streptococcus pyogenes, Streptococcus salivarius, Escherichia coli, \*\*Herpes simplex types 1/F and 2/G (oral, ocular, and genital), Influenza A<sub>2</sub> (Japan 305/57 Asian strain), Pellicularia filamentosa, Penicillium variable, Vaccinia, Canine parvovirus, Cytomegalovirus, Coronavirus and Polio type 1 viruses and Mycobacterium tuberculosis *in ten minutes* at 20°C/68°F or above.

**DIRECTIONS FOR USE:** It is a violation of U.S. Federal Law to use this product in a manner inconsistent with its labeling.

PRECLEAN surfaces using Sporidicin Brand Disinfectant Solution to remove soil and filth. Wipe dry with paper towel, cloth or sponge.

**TO DISINFECT AND DEODORIZE:** Thoroughly wet pre-cleaned surface with disinfectant. Use contact times as stated above at room temperature (20°C/68°F or above) to kill listed organisms. NOTE: See Sporidicin Bulletin No. 303 for HIV-1 (AIDS virus) instructions.

**\*\*\*TO KILL, CLEAN AND PREVENT MOLD AND MILDEW:** Thoroughly wet pre-cleaned contaminated area with sprayer or saturated Sporidicin towelette. Repeat application if new growth appears. Effective against *Stachybotrys chartarum*, *Chaetomium globosum* and *Aspergillus niger*.

**TO CLEAN AIR DUCTS:** Apply to hard, nonporous, surfaces.

**TO CLEAN, DEODORIZE AND DECONTAMINATE CARPETS:** For manual cleaning, spray product onto carpet and allow to dry. For machine cleaning, apply product in accordance with the machine manufacturer's guidelines.

**FOR FOGGING CLEAN ROOMS:** See Sporidicin Bulletin No. 301.

**PRECAUTIONARY STATEMENTS:** Hazard to Humans and Domestic Animals.

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

**FIRST AID:** If in eyes, hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have product container or label with you when calling a poison control center or doctor, or going for treatment.

**STORAGE AND DISPOSAL:** Store at room temperature. Discard empty container in regular trash.

**NOTE:** This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.




Rockville, MD 20852

# Sporicidin®

BULLETIN NOS.  
303 and 403

---

**SPECIAL INSTRUCTIONS  
FOR  
CLEANING AND  
DECONTAMINATION OF  
SURFACES/OBJECTS  
SOILED WITH  
BLOOD/BODY FLUIDS  
AGAINST  
HIV-1 (AIDS VIRUS)**



**SPORICIDIN BULLETIN NO. 303:**

**Sporicidin Brand  
Disinfectant Solution and Towelette  
(EPA Reg. Nos. 8383-3 and 8383-7)**

- a. Personal Protection:**  
Persons handling items soiled with blood or body fluids likely to be contaminated with HIV-1 (AIDS virus) must protect themselves by use of disposable latex gloves, gowns, masks and eye covering.
- b. Decontamination of Contaminated Objects:**  
Immediately after use, items contaminated with blood, body fluids, or other soils that may contain HIV-1 (AIDS virus) must be thoroughly wetted with Sporicidin brand disinfectant solution or towelette for at least 1 minute before further processing.\*
- c. Cleaning Procedure:**  
Following decontamination, objects and surfaces are moved to the cleaning and preparation areas where they are washed; rinsed and rough dried to remove all blood and body fluids before application of the disinfectant.
- d. Contact Time for Disinfection:**  
Thoroughly wet entire object or surface with solution or towelette. HIV-1 (AIDS virus) is killed within 1 minute contact time. \*Note: Use a 10 minute contact time for other organisms. This product should not be used for disinfection of critical or semi-critical items.
- e. Disposal of Infectious Materials:**  
Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

**SPORICIDIN BULLETIN NO. 403:**

**Sporicidin Brand  
Disinfectant Spray  
(EPA Reg. No. 8383-4)**

- a. Personal Protection:**  
Persons handling items soiled with blood or body fluids likely to be contaminated with HIV-1 (AIDS virus) must protect themselves by use of disposable latex gloves, gowns, masks and eye covering.
- b. Decontamination of Contaminated Objects:**  
Immediately after use, items contaminated with blood, body fluids, or other soils that may contain HIV-1 (AIDS virus) must be sprayed to thoroughly wet with disinfectant for 1 minute before further processing.\*
- c. Cleaning Procedure:**  
Following decontamination, objects and surfaces are moved to the cleaning and preparation area where they are washed; rinsed and rough dried to remove all blood and body fluids before application of the disinfectant.
- d. Contact Time for Disinfection:**  
Spray to thoroughly wet entire object or surface. HIV-1 (AIDS virus) is killed within 1 minute contact time. \*Note: Use a 10 minute contact time for other organisms. This product should not be used for disinfection of critical or semi-critical items.
- e. Disposal of Infectious Materials:**  
Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.