


5. FIRST AID

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothes before reuse.

6. If swallowed, drink PROMPTLY a large quantity of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol. CALL a physician immediately.

7. NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock as well as oxygen and measures to support breathing manually or mechanically may be needed. If persistent, convulsions may be controlled by the cautious intravenous injection of a short-acting barbiturate drug.



# HYA-CIDE

## CONCENTRATED GERMICIDAL SOLUTION

— WITH —

## RUST-INHIBITOR

INCLUDED

NOT INJURIOUS TO CARBIDE BURS



CONTENTS: 1 QUART

ACTIVE INGREDIENTS: —

n-alkyl (C<sub>14</sub> 50%, C<sub>12</sub> 40%, C<sub>16</sub> 10%) dimethyl  
benzyl ammonium chloride .....17.4%

ETHANOL ..... 4.3%

INERT INGREDIENTS .....78.3%

Propylene Glycol, Wetting Agent  
Detergent, Sodium Nitrite, Water

## DANGER

KEEP out of reach of children.

Corrosive. Causes eye damage and severe skin irritation. DO NOT get in eyes, on skin or on clothing. Protect eyes and skin when handling concentrate. HARMFUL or FATAL if swallowed. AVOID contamination of FOOD.

EPA Reg. No. 8701-1-AA

RINSE EMPTY CONTAINER WITH  
WATER BEFORE DISCARDING

See Side Panel for First Aid

## FOR THE CHEMICAL DISINFECTION OF DENTAL & SURGICAL INSTRUMENTS

### INSTRUCTIONS:

1. Dilute one oz. to one gallon of water.
2. Scrub all instruments clean of blood and pus before immersing in solution. Instruments should remain in solution a minimum of 15 minutes or longer.
3. DISCARD SOLUTION AT THE END OF EACH DAY. SCRUB CONTAINER CLEAN BEFORE ADDING NEW SOLUTION. BE SURE NO SOAP REMAINS.
4. GERMICIDAL solutions should not be relied upon to disinfect hinged instruments, those with narrow or deep crevices, HYPODERMIC NEEDLES or SYRINGES. No solution should be relied on to kill spores or mycobacterium tuberculosis. THESE instruments should be BOILED or AUTOCLAVED.

Nonflammable, nonvolatile, odorless

Do not mix with soap or anionic detergents.

HYA-CIDE, Los Angeles, Calif.

This product will not destroy the etiologic agent for viral Hepatitis.