

9480-4

04/24/2009

1/17

APR 24 2009

Christina Griffin
Agent for PDI
NCH Corporation
12510 Prosperity Drive, Suite 160
Silver Spring, MD 20904

Subject: Sani-Cloth Germicidal Wipes Germicidal Disposable Cloth
EPA Reg. No. 9480-4
Application Dated: March 25, 2009
Receipt Date: March 26, 2009

Dear Christina Griffin

The following notification submitted in connection with registration under the provisions of PR Notice 98-10, Federal Insecticide, Fungicide, and Rodenticide Act (FIFRA) section 3(c)9. is acceptable.

Proposed Notification

- Revise Ingredient Statement by adding Does Not Include the weight of the Cloth..... or Wipe" for clarity.
- Updating the Net Contents
- Revised establishment numbers
- Revising the Web address section
- Adding the Statement "made in USA"
- Adding Marketing Language "Suitable for use on equipment requiring alcohol-based Products.

Comments


Based on a review of the material submitted, the following comments apply: This application for notification to revise the label per the above mentioned comments is acceptable. A copy has been placed in our records for future reference.

CONCURRENCES							
SYMBOL							
SURNAME							
DATE							

2/17

Should you have any questions or comments concerning this letter, please contact Velma Noble at (703) 308-6233.

Sincerely,



Velma Noble
Product Manager (31)
Regulatory Management Branch1
Antimicrobial Division (7510P)



United States
Environmental Protection Agency
 Washington, DC 20460

Registration
 Amendment
 Other

OPP Identifier Number

NOTIF

Application for Pesticide - Section I

1. Company/Product Number 9480-7	2. EPA Product Manager Velma Noble	3. Proposed Classification <input type="checkbox"/> None <input type="checkbox"/> Restricted
4. Company/Product (Name) Professional Disposables International, Inc. (PDI) / Sani-Wipe	PM# 31	
5. Name and Address of Applicant (Include ZIP Code) Professional Disposables International, Inc. (PDI) c/o Delta Analytical Corp. 12510 Prosperity Drive, Suite 160 Silver Spring, MD 20904 <input type="checkbox"/> Check if this is a new address		6. Expedited Review. In accordance with FIFRA Section 3(c)(3)(b)(i), my product is similar or identical in composition and labeling to: EPA Reg. No. _____ Product Name _____

Section - II

Amendment - Explain below
 Resubmission in response to Agency letter dated _____
 Notification - Explain below.

Final printed labels in response to Agency letter dated _____
 "Me Too" Application.
 Other - explain below.

Explanation: Use additional page(s) if necessary. (For section I and Section II.)

Notification per PR Notice 98-10 to add 1 alternate brand name and minor label changes.

This notification is consistent with the provisions of PR Notice 98-10 and EPA regulations at 40 CFR 152.46, and no other changes have been made to the labeling or the confidential statement of formula of this product. I understand that it is a violation of 18 U.S.C. Sec. 1001 to willfully make any false statement to EPA. I further understand that if this notification is not consistent with the terms of PR Notice 98-10 and 40 CFR 152.46, this product may be in violation of FIFRA and I may be subject to enforcement action and penalties under sections 12 and 14 of FIFRA.

Section - III

1. Material this Product will be Packaged in:

Child-Resistant Packaging <input type="checkbox"/> Yes* <input type="checkbox"/> No	Unit Packaging <input type="checkbox"/> Yes <input type="checkbox"/> No	Water Soluble Packaging <input type="checkbox"/> Yes <input type="checkbox"/> No	2. Type of Container <input type="checkbox"/> Metal <input type="checkbox"/> Plastic <input type="checkbox"/> Glass <input type="checkbox"/> Paper <input type="checkbox"/> Other (Specify) _____
* Certification must be submitted.		If "Yes," Unit Package wgt. No. per container	

3. Location of Net Contents Information
 Label Container

4. Size(s) of Retail Container

5. Location of Label Direction:
 On Label
 On Labeling accompanying product

6. Manner In Which Label Is Affixed to Product
 Lithograph
 Paper glued
 Stenciled
 Other _____

Section - IV

1. Contact Point (Complete items directly below for identification of individual to be contacted, if necessary, to process this application.)

Name Cristina Griffin	Title Agent for Professional Disposables International, Inc. (PDI)	Telephone No. (Include Area Code) 301-680-7971
---------------------------------	--	--

Certification

I certify that the statements I have made on this form and all attachments thereto are true, accurate and complete. I acknowledge that any kind of knowingly false or misleading statement may be punishable by fine or imprisonment or both under applicable law.

2. Signature 	3. Title Agent for Professional Disposables International, Inc. (PDI)	6. Date Application Received (Stamped)
4. Typed Name Cristina Griffin	5. Date March 25, 2009	

Demetrius

9/17



March 25, 2009

Document Processing Desk (NOTIF)
 Office of Pesticide Programs (7504P)
 U.S. Environmental Protection Agency
 One Potomac Yard
 2777 S. Crystal Drive, Room S-4900
 Arlington VA 22202

Attn: Velma Noble, PM 31

RE: Notification per PR Notice 98-10 to make minor changes to first page of label
Product: Sani-Cloth Germicidal Wipes, EPA Reg. No. 9480-4
Company: Professional Disposables International, Inc. (PDI)

Dear Ms. Noble:

On behalf of PDI, I am submitting a notification per PR Notice 98-10 for the product Sani-Cloth Germicidal Wipes, EPA Reg. No. 9480-4.

The minor label changes on page 1:

- Under the ingredient statement we are adding "[Does not include the weight of the [cloth][towelette][or][wipe]]" to clarify the ingredient statement.
- Updating Net Contents to add [and/or Net Wt.] _____oz (g) [or lbs. (kg)]
- Revising establishment numbers
- Deleting the website address, and adding [For information call: [phone number to be inserted]]
- Added optional [Made in USA]

Change on page 3:

- In the first paragraph we are adding the minor marketing language, "Suitable for use on equipment requiring alcohol-based products."

No other changes to the label.

Enclosures

- EPA form 8570-1
- 1 copy of the revised label (additions highlighted)

If you have any questions regarding this submission, please contact me at 301-680-7971 or cgriffin@delta-ac.com.

Sincerely,

A handwritten signature in cursive script that reads "Cristina Griffin".

Cristina Griffin
 Agent for PDI

Enclosures

cc: David Jones, PDI

5/17

SANI-CLOTH® GERMICIDAL WIPES GERMICIDAL DISPOSABLE CLOTH

ACTIVE INGREDIENTS:

n-Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzyl ammonium chlorides.....	0.25%
n-Alkyl (60% C ₁₄ , 30% C ₁₆ , 5% C ₁₂ , 5% C ₁₈) dimethyl benzyl ammonium chlorides.....	0.25%
OTHER INGREDIENTS	99.50%
TOTAL	100.00%

[Does not include the weight of the [cloth] [towelette] [or] [wipe]]

KEEP OUT OF REACH OF CHILDREN
DANGER

See back [side] panel for additional precautionary statements

NET CONTENTS [and/or Net Wt.] _____ oz (g) [or lbs. (kg)]

____ # OF CLOTHS [SIZE of each cloth]

EPA Reg. No.: 9480-4

EPA Est. No. : 9480-NY-1

or

EPA Est. No. A=9480-NY-1, B=72956-WI-1, C=72956-AR-1

Alpha character will precede batch code on product


[Manufactured by [or for]]:

Professional Disposables International, Inc.

Two Nice-Pak Park, Orangeburg, NY 10962-1376

[For information call: [phone number to be inserted]]

[Made in USA]



First Aid:

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment

If in eyes:

-Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing.

If on skin:

-Take off contaminated clothing.

-Rinse skin immediately with plenty of water for 15-20 minutes.

[For additional information in case of emergency call toll free [phone number to be inserted]].

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

6/17

Precautionary Statements

Hazards to Humans and Domestic Animals

Danger. Causes irreversible eye damage. Harmful if absorbed through skin. Do not get in eyes or on clothing. Wear suitable hand protection (e.g. gloves) when dispensing and using this product. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using restroom. Remove and wash contaminated clothing before reuse.

Physical or chemical hazard

Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

Bactericidal

Tuberculocidal

Virucidal*

[Bactericidal, Tuberculocidal and Virucidal* in 2 minutes]

*Adenovirus Type 5 [ATCC VR-5], Herpes Simplex Type 2 [ATCC VR-734], Influenza A2/Hong Kong [ATCC VR-544], Influenza A (H1N1) virus [ATCC VR-98] Strain A/Malaya/302/54, Human Coronavirus [ATCC VR-740] Strain 229E, Respiratory Syncytial Virus [RSV], Rhinovirus [ATCC VR-1110], Rotavirus Strain WA, Vaccinia virus [ATCC VR-1354], HIV-1[†], Hepatitis B Virus[†], Hepatitis C Virus[†]

[*]Virucidal against Adenovirus Type 5, Herpes Simplex Type 2, Influenza A2/Hong Kong, Influenza A (H1N1) virus, Human Coronavirus, Respiratory Syncytial Virus [RSV], Rhinovirus, Rotavirus, Vaccinia virus, HIV-1[†], Hepatitis B Virus[†], Hepatitis C Virus[†]

[Product name] kills bacteria listed on this label

- Acinetobacter baumannii [ATCC 19606]
 - Burkholderia cepacia [ATCC 25416]
 - Campylobacter jejuni [ATCC 29428]
 - Candida albicans [ATCC 14053]
 - Escherichia coli (E. coli) [ATCC 11229]
 - Escherichia coli O157:H7[ATCC 35150]
 - ESBL Producing Escherichia coli (E. coli) [ATCC BAA-196]
 - Klebsiella pneumoniae [ATCC 4352]
 - Methicillin Resistant Staphylococcus aureus (MRSA) [ATCC 33592]
 - Mycobacterium bovis BCG (Tuberculosis) [ATCC 35743]
 - Pseudomonas aeruginosa [ATCC 15442]
 - Salmonella enterica [ATCC 10708]
 - Staphylococcus aureus [ATCC 6538]
 - Vancomycin Resistant Enterococcus faecalis (VRE) [ATCC 51299]
- in the presence of organic soil (5% blood serum)

[*][Product name] is virucidal against Adenovirus Type 5, Herpes Simplex Type 2 (causative agent of genital herpes), Influenza A2/Hong Kong (representative of the common flu virus), Influenza A (H1N1) virus, Human Coronavirus, Respiratory Syncytial Virus [RSV], Rhinovirus, Rotavirus, and Vaccinia virus in the presence of organic soil (5% blood serum).

[†]KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68°- 77°F) or (20°- 25°C) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[Product name] is a premoistened, non-woven durable wipe containing a quaternary/alcohol based solution. Recommended for use in hospitals and other critical care areas where the control of the hazards of cross-contamination between treated surfaces is required. Suitable for use on equipment requiring alcohol-based products. Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica® and glass. Some organisms are removed from the surface by thoroughly wiping the surface with the wipe. Most remaining organisms are killed within two (2) minutes by exposure to the liquid in the wipe.

Areas of Use:

Hospital and Healthcare Settings: Ambulances, Ambulatory Surgical Centers (ASC), Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctors Offices, Donation Centers [blood][plasma][semen][milk][apheresis], Emergency Medical Settings, Emergency Rooms, E.R., Emergency Vehicles, Eye Surgical Centers, Healthcare Settings, Health Care Facilities, Home Health Care, Hospices, Hospitals, ICU, Isolation Areas, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating Rooms (OR), Ophthalmic Offices, Orthopedics, Out-Patient Surgical Centers (OPSC), Patient Care Areas, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory therapy, Surgery, Surgical Centers, Transport vehicles and X-Ray.

and/or

Critical Care Areas: CCU, Emergency Rooms, E.R., ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU) and Surgery.

and

Hospital, Healthcare, and Critical Care Use Sites: May be used on hard non-porous surfaces of: ambulance equipment surfaces; bed railings; cabinets; carts; chairs; counters; dental chairs; dental countertops; dental unit instrument trays; endodontic equipment such as apex locators; exam tables; gurneys; infant incubators [interior and exterior surfaces of]; infant warmers [interior and exterior surfaces of]; isolettes; IV poles; laboratory equipment and surfaces; loupes; operating room tables and lights; operator light switches; oxygen hoods; physical therapy (PT) equipment surfaces; slit lamps; spine back boards; stethoscopes; stretchers; **stools**; tables; telephones; toilet seats; and hard non-porous outside surfaces of: amalgamators and dental curing lights; anesthesia machines and respiratory therapy equipment; apheresis machines, diagnostic equipment, dialysis machines; patient monitoring equipment, patient support and delivery equipment, pulp testers and motors; toilets; and ultrasound transducers and probes.

Grocery Stores or Shopping Centers or Supermarkets or Delicatessens: Disinfects grocery cart handles and grocery cart child seats. Disinfects shopping cart handles and child seats.

Nursery or Daycare Centers: Bassinets [interior and exterior surfaces of], changing tables, cribs, diaper changing stations, diaper pails, infant/child care equipment surfaces, trash cans, toys and other hard, non-porous surfaces listed on this label.

Bathroom[s] or Public Restrooms: Basins, [bath]tubs, bathroom fixtures, cabinets, countertops, doorknobs, faucet[s], handles, hampers, handrails, shower[s] [stalls], showers, sinks, glazed tiled walls, glazed tiles, toilet [handle] [rim] [seats] [tank] [and other exterior surfaces of], urinals [exterior surfaces of] and vanity tops, and other hard non-porous surfaces listed on this label.

Health club facilities and Gyms: Gym equipment, weight-lifting surfaces and other hard non-porous surfaces listed on this label.

Schools or Classrooms or Office Buildings or Public Facilities or Manufacturing Facilities or Libraries or Public Telephones: Cabinets, chairs, computers, computer monitors, computer mouses, counters, countertops, desks, desktops, doorknobs, fax machines, filing cabinets, filing doors, filing handles, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, telephones, toys, trash cans, work stations, and other hard non-porous surfaces listed on this label.

Transport vehicles, school busses, public busses, fire trucks, police cars: Hard non porous surfaces of vehicles including steering wheel, hard plastic or metal flooring, seats, and walls.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Restaurants: Hard non porous, non food contact surfaces only including chairs, display cabinets, garbage cans, booster seats, floor, sinks, bathroom surfaces, glazed tiled walls, telephones, cash registers, and doorknobs.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

[Use one of the following as applicable to type of container:]

[To Dispense TOWELETTES] [or WIPES] [or TOWELS] [or CLOTHS]:

Remove lid [and discard inner seal] [from container or canister.] Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the [container or] canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep [center cap of] lid closed to prevent moisture loss.

TRAVEL PACK:

Peel back label slowly. Pull out wipes as needed and reseal label. Unfold wipe.

BOX:


Open along perforations as indicated. Remove packets as needed.

PACKET/PACKAGE:

Tear open packet [package]. Unfold wipe and use as directed. Discard after single use.

SOLO® PAK:

Open package lid. Completely remove inner seal. Pull out wipe. (The next wipe pops up automatically.) Snap lid closed after use to retain moisture.



DO NOT FLUSH • STORE AT ROOM TEMPERATURE
DISPOSE OF USED WIPES IN TRASH RECEPTACLE • KEEP SOLO®
LID CLOSED WHEN NOT IN USE
SOLO® is a registered trademark of Nice-Pak Products, Inc.

REFILL PACKET

This refill package is made for use only in a [product name] [container or] canister; do not use this product in a [container or] canister for another product. Do not use the refill wipes in a container that is cracked or broken. Discard the cracked or broken container in accordance with the container disposal instructions on the label.

[Option 1]

Remove [entire] lid from [container or] canister. Tear open refill [package] at [notch] [mark] [arrow] and place [wipe] roll in empty [container or] canister. Follow information on [container or] canister to dispense wipes.

[Option 2]

Remove [entire] lid from [container or] canister. Tear open refill [package] at [notch] [mark] [arrow]. Place [entire] refill package of wipes in empty [container or] canister. Follow information on [container or] canister to dispense wipes.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

9/17

GENERIC:

Open lid. Pull out wipe. Close lid tightly.

Bucket [or pail or container]

1. Remove the lid from the [pail or bucket or container.]
2. Press open the dispenser at center of the lid.
3. Open the refill pack.
4. Locate the first wipe at the center of the roll.
5. Feed the first wipe through the dispensing orifice.
6. Reapply the lid to the [pail or bucket or container.]
7. Dispense the wipe.

For the refills:

1. Remove the lid from the [pail or bucket or container.]
2. Remove the empty refill pack material from the [pail or bucket or container.]
3. Press open the dispenser at center of the lid.
4. Open the refill pack.
5. Locate the first wipe at the center of the roll.
6. Feed the first wipe through the dispensing orifice.
7. Reapply the lid to the pail.
8. Dispense the wipe.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

TO DISINFECT [AND DEODORIZE]: To disinfect nonfood contact surfaces only: Use a wipe to remove heavy soil. Unfold a clean wipe and thoroughly wet surface. Treated surface must remain visibly wet for a full two (2) minutes. Use additional wipe(s) if needed to assure continuous two (2) minute wet contact time. Let air dry. *[When TB is included on the label, add the following or use the TB directions below:]* These directions also apply to Mycobacterium bovis BCG (Tuberculosis) at 68°F (20°C).

[The following sets of directions are optional and for hospital/medical facility use only:

Tuberculosis (TB) [Outbreaks]: This product is effective against Mycobacterium bovis BCG (Tuberculosis) in one minute at 68°F (20°C). Use a wipe to remove heavy soil. Unfold a clean wipe and thoroughly wet surface. Treated surface must remain visibly wet for one (1) full minute. Use additional wipe(s) if needed to assure continuous one (1) minute wet contact time. Let air dry. Although efficacy at a one (1) minute contact time has been shown to be adequate against TB, this time is not sufficient for other organisms on the label. Therefore a 2 minute wet contact time must be used for all listed non TB organisms.

and/or

Ultrasound transducers and probes: For use on external surfaces of ultrasound transducers and probes only. To disinfect: Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected by using a clean washcloth. Then thoroughly wet the probe with a wipe. Treated surface must remain visibly wet for a full two (2) minutes. Use additional wipe(s) if needed to assure continuous two (2) minute wet contact time. Let air dry.

†SPECIAL INSTRUCTIONS FOR CLEANING AGAINST HIV-1, HEPATITIS B VIRUS AND HEPATITIS C VIRUS OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS

Personal protection: When using this product, wear disposable protective gloves, protective gowns, face masks, or eye coverings as appropriate when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

Cleaning procedure: All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold and use first germicidal wipe to remove heavy soil.

Disposal of infectious materials: Dispose of used wipe in accordance with local regulations for infectious waste disposal.

Contact time: Use second germicidal wipe to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[Use storage/disposal statement A if container or canister directions used. Use storage/disposal statement B if travel pack or [refill] packet/package directions used.]

Storage and disposal

Statement A: For [container or] [canisters and] refill dispensers:

Storage: Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss.

Towelette Disposal: Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet.

Dispenser Disposal: Do not reuse empty container. Triple rinse and offer for recycling. If recycling is not available, put in trash collection.


[If refillable container is being used, substitute the following: **Refill Dispenser Disposal:** Triple rinse cracked or broken containers and offer for recycling. If recycling is not available, discard in trash.]

Statement B: For packets, travel packs, and [refill] packets or packages:

Storage: Do not store near heat or open flame.

Towelette Disposal: Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet.

Travel Pack [or Solo® Pak or Refill Packet] [or [bulk] package]] Disposal: Put empty travel pack [or Solo® Pak or refill packet] [or packaging]] in trash collection.



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Optional Claims and Descriptors [Note: Deodorizing directions must be used when deodorizing claims (claims involving odor control) are used]

X" x Y" [wipes]

___% bigger [[quilted] wipes] [vs. X" x Y" wipes]

___# pre-saturated towelettes X" x Y"

© [insert year] [insert company name]

[Absorbs] [Absorbent] [wipe]

After use, just toss wipe(s) away

Always ready to use

An effective method to deodorize and disinfect against odor-causing organisms on hard non-porous surfaces listed on this label.

An effective way to disinfect hard, non-porous environmental surfaces listed on this label.

A non-woven disposable cloth for use in hospitals and other critical care areas where control between treated surfaces is required.

Antibacterial

Antimicrobial

[Antibacterial] [Antimicrobial] health club equipment wipe

A Quaternary [ammonium chloride] solution impregnated in a wiping cloth

Are pre-measured to have the recommended amount of germicidal solution, wipe after wipe.

A Sani-Surface® Brand

Assures proper dilution strength, product effectiveness and standardized technique.

Bleach-free [formula]

Can help reduce the risk of cross contamination on treated surfaces

Cellulosic non-woven disposable cloth

[Cloth-like texture] [Textured] [strong] [for] [with] scrubbing power

Contains [name of active ingredient], the [same] active ingredient used in [Brand-name]

Contains [number] canisters [plus] [and] [number] refills

Convenient

Deodorize hard, non-porous surfaces [Use to]

Deodorizes, and kills common household germs*

Diamond-embossed wipe is low-linting [for sensitive equipment surfaces]

Disinfect external surfaces of ultrasound transducers

Disinfect the precleaned surface

Disinfecting towelettes

Disinfects and deodorizes in 2 minutes

Disinfects grocery cart handles or shopping cart handles or grocery cart child seats or shopping cart child seats

Disinfects hard non-porous environmental surfaces listed on this label [Use to]

Disinfects hard, non-porous surface in 2 minutes

[Disinfects] hard surfaces [that people touch]

Disinfects [health club equipment] [treadmill handles][elliptical handles][exercise bike handles][weight benches]

Disinfects [in [just] two [2] minutes]

Disinfects on the go

Dispenses one wipe at a time

Does not contain bleach

Does not produce mist or overspray

Double-sided scrubbing wipes

Durable, low-lint wipe designed for sensitive equipment surfaces

Easy

Easy [and convenient] way to [touch-up] [bathroom] [office] everyday

[Easy] [convenient] [dispensing] lid [closure] [cap] [package]

Easy drop-in refill[s]

[Easy] pop-up-tub


Easy start lid [package]

Easy to use [easy to refill]

Easy to use package [canister] [container] keeps [your] wipes moist

[Economical] Refill[s] [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] -pack

[Economical] Refills [now] available



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

12/17

Economy

Effective against bathroom bacteria and viruses* [use at least two bacteria and two viruses from this list:]
(Escherichia coli O157:H7, Salmonella enterica, Staphylococcus aureus, Herpes Simplex Type 2
(causative agent of genital herpes), *Influenza A2/Hong Kong (representative of the common flu virus),
Rhinovirus, Rotavirus)

Effective against [the following micro-organisms on hard, non-porous surfaces] [list organisms]

Effective in the presence of 5% blood serum contamination

Eliminates odors

Eliminates potential for overspray

Extra large [disinfecting] [Y" x Z" compared to our A" x B" size] wipes

Fast [Convenient] [Portable] way to disinfect hard non-porous surfaces

For daily use

Gets rid of –or- eliminates odors

[Helps] [kill][s] [eliminate][s] odor-causing bacteria

Helps prevent cross contamination on treated surfaces

Hospital disinfectant [disinfection] [use only when hospital use sites on label]

Includes ___ count refill [pack] [contains] ___ wipes

Just –or- simply wipe exterior toilet surface –or- bathroom germs+ away

Just [pull] wipe and toss

Kill germs+ and bacteria instead of germ paper towels, sponges and rags [Use this product to]

Kills bacteria and viruses* in the bathroom [use at least two bacteria and two viruses from this list:] (Escherichia coli O157:H7, Salmonella enterica, Staphylococcus aureus, Herpes Simplex Type 2 (causative agent of genital herpes), *Influenza A2/Hong Kong (representative of the common flu virus), Rhinovirus, Rotavirus)
–or- on the exterior toilet surfaces

Kills bacteria in 2 minutes

[Kills] [Effective against] [MRSA] [(Methicilin Resistant Staphylococcus aureus)]

Kills household –or- exterior toilet surface –or- bathroom germs+ [use at least two bacteria and two viruses from this list:] (Escherichia coli O157:H7, Salmonella enterica, Staphylococcus aureus, Herpes Simplex Type 2 (causative agent of genital herpes), *Influenza A2/Hong Kong (representative of the common flu virus), Rhinovirus, Rotavirus)

Kills odor causing bacteria

Leaves [bathroom] [office] disinfected [and smelling fresh]

Less [reduced] packaging [less to throw away]

Moisture lock lid

Multi-pack

Multi-purpose disinfectant [Multi [-refill] pack]

New larger, thicker wipe [with diamond-embossed design]

New thicker wipe with diamond-embossed design

Next wipe pops up automatically [everytime] [!]

No canister needed

No measuring

No mixing

No [pinched fingers] [rips] [tears]

No sponges, no rags, no buckets, no mess

No threading [required]

Non-scratching [wipe]

Non-woven disposable cloth

Not for sale in [insert name of state(s)]

Odor [elimination] [neutralization]

One [wipe] at a time dispensing

On-the-go disinfectant

[Package][insert type of package] made from [x%] recycled materials


Perfect for disinfecting exterior toilet surfaces –or- around the exterior toilet surface

Picks up [particles] [dust] [dirt] [hair] [crumbs]

Pop-up [dispensing] [top]

Portable

[Premixed], [premeasured], [premoistened]



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Pull n/ Go [pack]
 Quick
 Quick [and easy] [and convenient] way to disinfect
 Quick and easy way to disinfect [bathroom –or- exterior toilet surfaces] [hard, nonporous surfaces]
 [Quickly] eliminates [lingering] [food] odors on surfaces
 Ready to Use
 Recyclable [package]
 Reduce [risk of cross contamination] [on hard non-porous surfaces]
 Reduces [packaging] waste [by x%] [when compared to regular canister or bucket or pail]
 [Refill] package contains [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] [wet] [pre-moistened] [wipe]
 [refill[s]] [roll[s]]
 Refill [size] [pack] [value pack]
 Refillable [with ___# wipes] [dispenser] refill [pack] [for dispenser]
 Requires no mixing or measuring and are ready-to-use
 Resealable [stay-moist] pouch [package] [container] inside
 Reuse your [empty] canisters [again [and again] [!]
 Sample
 Scrubbing power of a sponge in a disinfecting wipe
 Strong, durable diamond-embossed textured wipe deodorizes
 Strong, durable low-lint wipe [designed for sensitive equipment surfaces]
 Tear at perforation
 Textured diamond-embossed design [deodorizes hard, non-porous surfaces]
 The 2-minute germicidal wipe [disinfectant wipe]
 Thick, heavy duty wipes are bactericidal, tuberculocidal and virucidal* in 2 minutes
 Thick, large wipe disinfects
 Thick, textured wipe with diamond-embossed design
 This [outer] [selling] [package] [container][canister] contains x% post consumer [material][plastic][resin]
 Time to reorder
 Trial
 Tub refill
 Tub with pop-up action
 Twin pack
 Use on hard, non-porous surfaces only: [stainless steel] [Formica®] [glass] [glass tables] [glazed tile] [Test wipe on small inconspicuous area first.]
 Use to preclean or decontaminate critical or semi-critical medical devices prior to sterilization / high-level disinfection
 [Use to] refill[s] [count] [wipe] canisters or tubs
 Value
 Versatile
 Wet wipes
 Where control of the hazards of cross contamination between treated surfaces is of prime importance
 Wipe [material] made from [x] % [sustainable] [ingredients] [fibers] [*May be used only when cellulose wipe material used*]
 [Wipes out] [Kills] [common] [household] [bacteria] [viruses*] [germs*]
 [Wipes out] [Kills] [eliminates] [*list organisms*] in [just] [2 minutes]
 [Wipes out] [Kills] [the] flu virus* (*Influenza A2/Hong Kong)
 [Wipes out] [Kills] viruses* and bacteria [including the flu virus* (*Influenza A2/Hong Kong)]
 [Wiping] [soft] side for disinfecting
 With a bleach-free formula
 X% [more] [wipes] [free] [vs. Z count] [canisters] [packs]
 X% reduction in packaging [material]

The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

14/17

[Note: The following claims can be used only for 6 months]

[Better] [New] [Improved] [Stronger] [Thicker] Wipe [!]

[Better] [new] [&] [improved] [dispensing] lid [closure] [cap] [package] [dispensing] [better] [than before]

New and improved

[New] [Package] [packaging] uses [x%] less material (plastic) [then] [when compared to] [regular] containers [or canisters]

New thicker and larger diamond-embossed wipe cleans more surface area

[Now] [Reinforced] [With] Quilted [strength]

Now with [%] more towel

Now with [#] more wipes

[Scrubs] (&) [cleans] [better than ever]

[Stronger] more durable [wipe]

Thick[er] [wipe] [with] cloth-like texture

X% [Bigger] [larger] [quilted] wipes

[x%] thicker wipes

Master Label Reader Notes

[A footnote can be used to list organisms shown on page 2. Bacteria and virus can be separated into footnotes 1 and 2 if desired.]

[When organisms are listed next to a claim, the organisms may be listed in a footnote instead of directly next to the claim]

[* All pesticidal "germ" claims must be linked to bacteria (minimum: Salmonella and Staph) and at least one virus. For example: Front panel: "Kills germs*" Back panel: "Kills Salmonella enterica, Staphylococcus aureus, Herpes Simplex Type 2, and Influenza A/Hong Kong".]

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

15/17

CLEANING: [cleaning directions and claims are grouped on this page]

DIRECTIONS FOR USE:

[CLEANING] [and] [DEODORIZING]: Wipe surface with towelette until clean. Let air dry.

CLAIMS:

An effective method to clean against odor-causing organisms on hard non-porous surfaces listed on this label.

A One Step▲ Cleaner [and disinfectant] (except precleaning is required for blood and body fluids)

Built-in [tough] scrubbing fibers help clean away tough

Cleans

Cleans and disinfects [on the go]

Cleans grocery carts and shopping carts

Diamond-embossed pattern [texture] provides effective cleaning

Easy One Step▲ Cleaning [and disinfecting] (except precleaning is required for heavy soil and blood and body fluids)

Leaves bathroom fresh and clean in [just] minutes

Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]

Multi-purpose cleaner [Multi [-refill] pack]

Perfect for cleaning [and disinfecting] around the toilet

Premoistened [wet] household cleaning wipe

Premoistened with specially formulated cleaners to handle daily clean-ups without bottles, sponges or buckets

[This product is]

Quick [and easy] [and convenient] way to clean [and disinfect]

Quick and easy way to clean [and disinfect] [bathroom –or- exterior toilet surfaces] [hard, nonporous surfaces]

[So] surfaces are more than [just] [clean] [they're] [disinfected]

[Soft] side for [cleaning] dirt away

Strong, durable diamond-embossed textured wipe cleans [and deodorizes]

Textured diamond-embossed design [cleans hard, non-porous surfaces] [helps remove soil from surfaces]

Textured for better [cleaning] [scrubbing] power

[Textured] [scrubby] [rough] side for [scrubbing] [cleaning] [removing] [wiping] [tough] [stains] [soils] [messes]

Textured wipe for effective cleaning

Thick, large wipe cleans

Virtually [Lint Free] [streak free] [clean] [shine]

Wipe away odors

[Wipes out] [Kills] germs+ [or bacteria] [while it cleans] [that cause odor]

[Wiping] [soft] side for [wiping] dirt away [and disinfecting]

▲ If one step cleaner claims used and HIV-1, HBV, and/or HBC claims and/or directions for ultrasound transducers are included on the label, then the following also must be included on the same panel of the label:

▲ a pre-cleaning step is required for HIV-1 [and] [or] use on ultrasound transducers [or probes].
If HIV-1, HBV, and/or HBC or ultrasound transducers are not listed on the label, any of the one step cleaning/disinfecting claims may be used without the phrase requiring pre-cleaning.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.


[The following graphics and statements are optional. They may be placed anywhere on label/container.]




NOT A BABY WIPE

NOT FOR USE ON SKIN

This graphic shows use of the wipe dispenser.



FINISHED FLOORS



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.


[The following graphics and statements are optional. They may be placed anywhere on label/container.]




Not a Skin Wipe.
For Use on Hard Surfaces Only.



Not a Skin Wipe.
For Use on Hard Surfaces Only.



[One side scrubs]
[One side wipes]



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.