



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY  
WASHINGTON, D.C. 20460

OFFICE OF PREVENTION, PESTICIDES AND TOXIC SUBSTANCES

OCT 30 2007

Ms. Cristina Griffin  
Agent for PDI  
Delta Analytical Corp  
12510 Prosperity Drive, Suite 160  
Silver Spring, MD 20904

Subject: Sani-Cloth Plus Germicidal Disposable Cloth  
EPA Reg. No.: 9480-6  
Amendment Date: June 14, 2007  
EPA Receipt Date: June 15, 2007

Dear Ms. Griffin,

The following amendment submitted in connection with registration under FIFRA section 3(c)(7)(A) is acceptable with the conditions listed below.

- Removal of Isopropanol as an active ingredient
- Add organisms, directions, label claims

**Conditions**

Revise as follows:

- 1) On page 2, delete C. difficile (ATCC 9689), Methicillin Resistant S. aureus (ATCC 33592), and Vancomycin Resistant E. faecalis (ATCC 51299) because the submitted efficacy data does not support claims against the organism.
- 2) On page 2, delete the "Kills HIV-1 (AIDS VIRUS), HBV, and HCV" paragraph, and place it on page 5. It must precede the "Special Instruction for cleaning and decontamination of HIV, HBV, and HCV" directions.
- 3) On page 5, revise the "Disinfecting" heading to read "Deodorizing and Disinfecting."
- 4) On page 7, delete the 34<sup>th</sup> claim "Powerful enough to disinfect" because its heightened efficacy.
- 5) On page 8, revise the 15<sup>th</sup> claim by deleting C. difficile (vegetative cells).

6) On page 8, revise the 20<sup>th</sup> and 21<sup>st</sup> claims by stating "[Wipes out][Kills] the flu virus, (Influenza A2/Hong Kong)" and "[Wipes Out][kills viruses\* and bacteria [including the flu virus, (Influenza A2/Hong Kong)."

### Data Summary

| Data Requirement   | Means of Support    | Status  |
|--|---------------------|---|
| Single Use Towelette for Hard Surface Disinfection – S. aureus (ATCC 6538)                         | MRID No.: 471562-01 | Acceptable, hard nonfood surface in 5% soil for 3 minutes   |
| Single Use Towelette for Hard Surface Disinfection- S. aureus (ATCC 33592) (Methicillin Resistant) | MRID No.: 471562-02 | Unacceptable, need to submit antibiotic resistant data  |
| Single Use Towelette for Hard Surface Disinfection – S. choleraesuis (ATCC 10708)                  | MRID No.: 471562-03 | Acceptable, hard nonfood surface in 5% soil for 2 minutes   |
| Single Use Towelette for Hard Surface Disinfection – P. aeruginosa (ATCC 15442)                    | MRID No.: 471562-04 | Acceptable, hard nonfood surface in 5% soil for 2 minutes   |
| Single Use Towelette for Hard Surface Disinfection- E. coli (ATCC 11229)                           | MRID No.: 471562-05 | Acceptable, hard nonfood surface in 5% soil for 2 minutes   |
| Single Use Towelette for Hard Surface Disinfection – E. coli (ATCC 35150) Serotype O157:H7         | MRID No.: 471562-06 | Acceptable, hard nonfood surface in 5% soil for 2 minutes   |
| Single Use Towelette for Hard Surface Disinfection – E. faecalis (ATCC 51299) VRE                  | MRID No.: 471562-07 | Unacceptable, need antibiotic resistant data. Testing must be completed w/ minimum of $1 \times 10^5$ CFU/ft <sup>2</sup> (per carrier) |
| Hard Surface Disinfection 1 Presaturated Towelette – C. difficile (ATCC 9689) vegetative cells     | MRID No.: 471562-07 | Not reviewed. Agency no longer accepting data on vegetative cells only against C. difficile spores                                      |

### General Comments

The Confidential Statements of Formulas for the Basic and Alternate Formula #1-3 are acceptable. They are in agreement with the label and PR Notice 91-2.

A stamped copy of the accepted labeling is enclosed. Submit one (1) copy of your final printed labeling before distributing or selling the product bearing the revised labeling.


Submit and/or cite all data required for registration/reregistration of your product under FIFRA section 3 (c) (5) and section 4 (a) when the Agency requires all registrants of similar products to submit such data.

If the above conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6 (e). Your release for shipment of the product bearing the amended labeling constitutes acceptance of these conditions.

3 816

Should you have any questions concerning this letter, please contact Drusilla Copeland at (703) 308-6224.

Sincerely,

A handwritten signature in black ink, appearing to read 'Velma Noble', with a stylized, flowing script.

Velma Noble  
Product Manager (31)  
Regulatory Management Branch I  
Antimicrobials Division (7510C)

4 816

# SANI-CLOTH® PLUS

## GERMICIDAL DISPOSABLE CLOTH

### ACTIVE INGREDIENTS:

|  |                 |
|--|-----------------|
| n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides.....                                      | 0.125%          |
| n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> ) dimethyl benzyl ammonium chlorides..... | 0.125%          |
| <u>OTHER INGREDIENTS</u> .....   | <u>99.750%</u>  |
| <b>TOTAL</b> .....   | <b>100.000%</b> |

KEEP OUT OF REACH OF CHILDREN

### CAUTION

See back [side] panel for additional precautionary statements

### NET CONTENTS

\_\_\_\_ # OF CLOTHS [SIZE of each cloth]

EPA Reg. No.: 9480-6

EPA Est. No. : 9480-NY-1

or:

EPA Est. No. A=9480-NY-1, B=9480-IN-1

Alpha character will precede batch code on product

Manufactured by:

PDI, The Healthcare Division of Nice-Pak Products, Inc.  
Two Nice-Pak Park, Orangeburg, NY 10962-1376

### Precautionary Statements

#### Hazards to humans and domestic animals

Caution. Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

#### First aid

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

#### Physical or chemical hazard

Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:

**OCT 30 2007**

Under the Federal Insecticide,  
Fungicide, and Rodenticide Act as  
amended, for the pesticide,  
registered under EPA Reg. No. **9480-6**

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

F:\shared\labels\pdi\9480-6\2007\draft label 061207.doc

June 12, 2007

## Bactericidal

## Virucidal\*

\* Herpes Simplex Type 2 (ATCC VR-734), Influenza A2/Hong Kong (ATCC VR-544), Respiratory Syncytial Virus (RSV)

Virucidal against Herpes Simplex Type 2, Influenza A2/Hong Kong, Respiratory Syncytial Virus (RSV), HIV-1<sup>†</sup>, Hepatitis B Virus<sup>†</sup>, Hepatitis C Virus<sup>†</sup>

## Tuberculocidal

[Product name] kills bacteria listed on this label

*Clostridium difficile* (ATCC 9689: Vegetative cells)

*Escherichia coli* (E. coli) O157:H7 (ATCC 35150)

Methicillin Resistant *Staphylococcus aureus* (MRSA) (ATCC 33592)

*Mycobacterium bovis* BCG (Tuberculosis) (ATCC 35743)

*Pseudomonas aeruginosa* (ATCC 15442)

*Salmonella enterica* (ATCC 10708)

*Staphylococcus aureus* (ATCC 6538)

Vancomycin Resistant *Enterococcus faecalis* (VRE) (ATCC 51299)

[Product name] is virucidal against Herpes Simplex Type 2 (ATCC VR-734), Influenza A2/Hong Kong (ATCC VR-544), Respiratory Syncytial Virus (RSV), HIV-1<sup>†</sup>, Hepatitis B Virus<sup>†</sup>, Hepatitis C Virus<sup>†</sup>

<sup>†</sup>KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68°- 77°F) or (20°- 25°C) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

A quaternary/alcohol solution impregnated in a wiping cloth. A non-woven disposable cloth for use in hospitals and other critical care areas where the control of the hazards of cross-contamination between treated surfaces is required. Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica® and glass. Some organisms are removed from the surface by thoroughly wiping the surface with the wipe. Most remaining organisms are killed within three (3) minutes by exposure to the liquid in the wipe.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

F:\shared\labels\pdi\9480-6\2007\draft label 061207.doc

June 12, 2007

**Areas of Use:**

**Hospital and Healthcare Settings:** Ambulances, Ambulatory Surgical Centers (ASC), Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctors Offices, Donation Centers [blood][plasma][semen][milk][apheresis], Emergency Medical Settings, Emergency Rooms, E.R., Emergency Vehicles, Eye Surgical Centers, Healthcare Settings, Health Care Facilities, Home Health Care, Hospices, Hospitals, ICU, Isolation Areas, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating Rooms (OR), Ophthalmic Offices, Orthopedics, Out-Patient Surgical Centers (OPSC), Patient Care Areas, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory therapy, Surgery, Surgical Centers, Transport vehicles and X-Ray.

and/or

**Critical Care Areas:** CCU, Emergency Rooms, E.R., ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU) and Surgery.

and

**Hospital, Healthcare, and Critical Care Use Sites:** May be used on hard non-porous surfaces of: ambulance equipment surfaces; bed railings; cabinets; carts; chairs; counters; dental chairs; dental countertops; dental unit instrument trays; endodontic equipment such as apex locators; exam tables; gurneys; infant incubators; [interior and exterior surfaces of], infant warmers [interior and exterior surfaces of]; isolettes; IV poles; laboratory equipment and surfaces; loupes; operating room tables and lights; operator light switches; oxygen hoods; physical therapy (PT) equipment surfaces; slit lamps; spine back boards; stethoscopes; stretchers; stools; tables; telephones; toilet seats; and hard non-porous outside surfaces of: amalgamators and dental curing lights; anesthesia machines and respiratory therapy equipment; apheresis machines; diagnostic equipment; dialysis machines; patient monitoring equipment; patient support and delivery equipment; pulp testers and motors; toilets; and ultrasound transducers and probes.

**Grocery Stores or Shopping Centers or Supermarkets or Delicatessens:** Disinfects grocery cart handles and grocery cart child seats. Disinfects shopping cart handles and child seats. Cleans grocery carts or shopping carts.

**Nursery or Daycare Centers:** Bassinets [interior and exterior surfaces of], changing tables, cribs; diaper changing stations; diaper pails; infant/child care equipment surfaces; trash cans; toys and other hard, non-porous surfaces listed on this label.

**Bathroom[s] or Public Restrooms:** Basins, [bath]tubs, bathroom fixtures, cabinets, countertops, doorknobs, faucet[s], handles, hampers, handrails, shower[s] [stalls], showers, sinks, tiled walls, tiles, toilet [handle] [rim] [seats] [tank] [and other exterior surfaces of], urinals [exterior surfaces of] and vanity tops, and other hard non-porous surfaces listed on this label.

**Health club facilities and Gyms:** Gym equipment, weight-lifting surfaces and other hard non-porous surfaces listed on this label.

**Schools or Classrooms or Office Buildings or Public Facilities or Manufacturing Facilities or Libraries or Public Telephones:** Cabinets, chairs, computers, computer monitors, computer mouses, counters, countertops, desks, desktops, doorknobs, fax machines, filing cabinets, filing doors, filing handles, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, telephones, toys, trash cans, work stations, and other hard non-porous surfaces listed on this label.

**Transport vehicles, school busses, public busses, fire trucks, police cars:** Hard non porous surfaces of vehicles including steering wheel, hard plastic or metal flooring, seats, and walls.

**Restaurants:** Hard non porous, non food contact surfaces only including chairs, display cabinets, garbage cans, booster seats, floor, sinks, bathroom surfaces, tiled walls, telephones, cash registers, and doorknobs.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

#### Directions for use

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

#### Dispenser directions:

##### Canister directions:

To start feed: Remove lid and discard inner seal from canister. Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep center cap of lid closed to prevent moisture loss.

##### Travel pack directions:

Peel back label slowly. Pull out wipes as needed and reseal label. Unfold wipe.


##### Box directions:

Open along perforations as indicated. Remove packets as needed.

##### Packet:

Tear open packet. Unfold wipe and use as directed. Discard after single use.

#### SOLO® PAK:



**DO NOT FLUSH • STORE AT ROOM TEMPERATURE  
DISPOSE OF USED WIPES IN TRASH RECEPTACLE • KEEP SOLO®  
LID CLOSED WHEN NOT IN USE**

#### REFILL PACKET

This refill package is made for use only in a [product name] canister, do not use this product in a canister for another product. Do not use the refill wipes in a container that is cracked or broken. Discard the cracked or broken container in accordance with the container disposal instructions on the label.

##### [Option 1]

Remove [entire] lid from canister. Tear open refill [package] at [notch] [mark] [arrow] and place [wipe] roll in empty canister. Follow directions on canister to dispense wipes.

##### [Option 2]

Remove [entire] lid from canister. Tear open refill [package] at [notch] [mark] [arrow]. Place [entire] refill package of wipes in empty canister. Follow directions on canister to dispense wipes.

#### GENERIC:

Open lid. Pull out wipe. Close lid tightly.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

**DISINFECTING:** To disinfect nonfood contact surfaces only:

Use a wipe to remove heavy soil. Unfold a clean wipe and thoroughly wet surface. Treated surface must remain visibly wet for a full three (3) minutes. Use additional wipe(s) if needed to assure continuous three (3) minute wet contact time. Let air dry. *[When TB is included on the label, add the following:]* These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F (20°C).

*[optional directions:]* Ultrasound transducers and probes:

For use on external surfaces of ultrasound transducers and probes only.

To disinfect: Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected by using a clean washcloth. Then thoroughly wet the probe with a towelette. Treated surface must remain visibly wet for a full three (3) minutes. Use additional wipe(s) if needed to assure continuous three (3) minute wet contact time. Let air dry.

**†SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS**

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, face masks, or eye coverings as appropriate when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

**Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold first germicidal cloth to remove heavy soil.

**Disposal of infectious materials:** Dispose of used towelette in accordance with local regulations for infectious waste disposal.

**Contact time:** Use second germicidal cloth to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry. Although efficacy at a two (2) minute contact time has been shown to be adequate against HIV-1, HBV and HCV this time is not sufficient for the other organisms listed on this label. Therefore a three (3) minute wet contact time must be used for other listed organisms.

*[Use storage/disposal statement A if canister directions used. Use storage/disposal statement B if travel pack or packet directions used.]*

**Storage and disposal**

**Statement A: For canisters and refill dispensers:**

**Storage:** Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss.

**Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet.

**Dispenser Disposal:** Do not reuse empty container. Triple rinse and offer for recycling. If recycling is not available, put in trash collection.

*[If refillable container is being used, substitute the following:]* **Refill Dispenser Disposal:** Triple rinse cracked or broken containers and offer for recycling. If recycling is not available, discard in trash.]

**Statement B: For packets, travel packs, and [refill] packets:**

**Storage:** Do not store near heat or open flame.

**Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet.

**Travel Pack [or Solo® Pak or Refill Packet] Disposal:** Put empty travel pack [or Solo® Pak or refill packet] in trash collection.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Optional Claims and Descriptors

X" x Y" [wipes]

\_\_\_% bigger [[quilted] wipes] [vs. X" x Y" wipes]

\_\_\_# pre-saturated towelettes X" x Y"

After use, just toss wipe(s) away

Always ready to use

An effective method to deodorize and disinfect against odor-causing organisms on hard non-porous surfaces listed on this label.

An effective way to disinfect hard, non-porous environmental surfaces listed on this label.

A non-woven disposable cloth for use in hospitals and other critical care areas where control of the hazards of cross contamination between treated surfaces are required.

Antibacterial

Antimicrobial

**A Sani-Surface® Brand**

Are pre-measured to have the recommended amount of germicidal solution, wipe after wipe.

Assures proper dilution strength, product effectiveness and standardized technique.

**Bleach-free [formula]****Bonus****Can help reduce the risk of cross contamination on treated surfaces**

Cellulosic non-woven disposable cloth

**Contains [name of active ingredient], the [same] active ingredient used in [Brand name]****Contains [number] canisters [plus] [and] [number] refills**

Convenient

**Deodorize hard, non-porous surfaces [Use to]**

Deodorizes, and kills common household germs

**Diamond-embossed wipe is low-linting [for sensitive equipment surfaces]**

Disinfect and scrub clean

Disinfect [and for touch-up clean]

Disinfects in 3 minutes

Disinfects [in [just] three [3] minutes]

Disinfect external surfaces of ultrasound transducers

Disinfect the precleaned surface

Disinfects and deodorizes

Disinfects hard, non-porous surface in 3 minutes

Disinfects and deodorizes in 3 minutes

Disinfecting towelettes

Disinfects grocery cart handles or shopping cart handles or grocery cart child seats or shopping cart child seats

Disinfect[s] hard non-porous environmental surfaces listed on this label **[Use to]****Dispenses one wipe at a time****Does not contain bleach**

Does not produce mist or overspray

**Double-sided scrubbing wipes****Durable, low-lint wipe designed for sensitive equipment surfaces****Easy****Easy [and convenient] way to [touch-up] [bathroom] [office] everyday****Easy drop-in refill[s]****[Easy] pop-up-tub****Easy start lid [package]****Easy to use [easy to refill]****[Economical] Refill[s] [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] -pack****[Economical] Refills [now] available****Economy****Effective against bathroom bacteria and viruses****Effective against [the following micro-organisms on hard, non-porous surfaces] [list organisms]****Eliminates odors****Extra large disinfecting wipes**

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Fast [Convenient] [Portable] way to disinfect hard non-porous surfaces  
 For daily use  
 Gets rid of [or] eliminates odors  
 Helps prevent [bacterial] [viral] cross-contamination [between] [hard non-porous] [surfaces]  
 Includes \_\_\_\_\_ count refill [pack] [contains] \_\_\_\_\_ wipes  
 Just [or] simply wipe exterior toilet surface [or] bathroom germs away  
 Just [pull] wipe and toss  
 Kill germs and bacteria instead of germey paper towels, sponges and rags [Use this product to]  
 Kills bacteria in 3 minutes  
 Kills germs and bacteria and viruses in the bathroom [or] on the exterior toilet surfaces  
 Kills household [or] exterior toilet surface [or] bathroom germs  
 Kills [list one or more organisms]  
 Kills odor-causing germs  
 Kills TB in 3 minutes  
 Leaves [bathroom] [office] disinfected [and smelling fresh]  
 Multi-pack  
 Multi-purpose disinfectant [Multi [-refill] pack]  
 New and improved [use this for 6 months after announced improvement]  
 New thicker wipe with diamond-embossed design  
 New larger, thicker wipe [with diamond-embossed design]  
 Next wipe pops up automatically [everytime] [!]  
 No canister needed  
 No measuring  
 No mixing  
 No sponges, no rags, no buckets, no mess  
 Non-woven disposable cloth  
 Now with [%] more towel [use this for 6 months after converting to larger towel]  
 Now with [#] more wipes [use this for 6 months after converting to a larger number of towels]  
 On-the-go disinfectant  
 Perfect for disinfecting exterior toilet surfaces –or– around the exterior toilet surface  
 Pop-up dispensing  
 Pop-up [top]  
 Portable  
 Powerful enough to disinfect  
 Premeasured  
 Premixed  
 Premixed, premeasured, premoistened  
 Premoistened  
 Pull n/ Go [pack]  
 Quick  
 Quick [and easy] [and convenient] way to disinfect  
 Quick and easy way to disinfect bathroom [or] exterior toilet surfaces  
 Quick and easy way to disinfect hard, nonporous surfaces  
 Ready to Use [on the following hard non porous surfaces]  
 [Refill] package contains [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] [wet] [pre-  
 moistened] [wipe] [refill[s]] [roll[s]]  
 Refill [size] [pack] [value pack]  
 Refillable [with \_\_\_\_\_ # wipes] [dispenser] refill [pack] [for dispenser]  
 Requires no mixing or measuring and are ready-to-use  
 Resealable pouch [package] [container] inside  
 Reuse your [empty] canisters [again [and again]] [!]  
 Sample  
 Scrubbing power of a sponge in a disinfecting wipe  
 Starts killing [bacteria] on contact  
 Strong, durable diamond-embossed textured wipe deodorizes  
 Strong, durable low-lint wipe [designed for sensitive equipment surfaces]  
 Tear at perforation  
 Textured diamond-embossed design deodorizes hard, non-porous surfaces

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Textured diamond-embossed design

The 3-minute germicidal wipe [disinfectant wipe]

Thick, heavy duty wipes are bactericidal, tuberculocidal and virucidal in 3 minutes

Thick, large wipe disinfects

Thick, textured wipe with diamond-embossed design

Time to reorder

Trial

Tub with pop-up action

Tub refill

Twin pack

Use to preclean or decontaminate critical or semi-critical medical devices prior to sterilization /  
high-level disinfection

[Use to] refill[s] [count] [wipe] canisters or tubs

Value

Versatile

Virucide\* including HIV-1, HBV and HCV applications and effective against Clostridium difficile  
(vegetative cells) and Staphylococcus aureus

Wet wipes

Where control of the hazards of cross contamination between treated surfaces is of prime importance

[Wipes out] [Kills] [common] [household] [bacteria] [viruses\*] [germs]

[Wipes out] [Kills] [eliminates] [list organisms] in [just] [3 minutes]

[Wipes out] [Kills] [the] flu virus

[Wipes out] [Kills] viruses\* and bacteria [including the flu virus]

[Wiping] [soft] side for disinfecting

With a bleach-free formula

© [insert year] [insert company name]

*1 [A footnote can be used to list organisms shown on page 2. Bacteria and virus can be separated into footnotes 1 and 2 if desired.]*

**CLEANING:** [cleaning directions and claims have been grouped on this page (they were previously incorporated into other parts of the label)]

**DIRECTIONS FOR USE:**

**CLEANING:** Wipe surface with towelette until clean. Let air dry.

**CLAIMS:**

An effective method to clean against odor-causing organisms on hard non-porous surfaces listed on this label.

A One Step<sup>†</sup> Cleaner [and disinfectant] (except precleaning is required for blood and body fluids)

Cleans

Cleans and disinfects

Cleans, deodorizes, and kills common household germs

Cleans, disinfects and deodorizes

Cleans grocery carts and shopping carts

**Diamond-embossed pattern [texture] provides effective cleaning**

Easy One Step<sup>†</sup> Cleaning [and disinfecting] (except precleaning is required for blood and body fluids)

**Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]**

**Leaves bathroom fresh and clean in [just] minutes**

**Multi-purpose cleaner [Multi [-refill] pack]**

**New thicker and larger diamond-embossed wipe/cleans more surface area [use this for 6 months after converting to larger towel]**

Perfect for cleaning [and disinfecting] around the toilet

Premoistened [wet] household cleaning wipe

Premoistened with specially formulated cleaners to handle daily clean-ups without bottles, sponges or buckets [This product is]

Quick [and easy] [and convenient] way to clean [and disinfect]

Quick and easy way to clean [and disinfect] bathroom –or– exterior toilet surfaces

Quick and easy way to clean [and disinfect] hard, nonporous surfaces

[So] surfaces are more than [just] [clean] [they're] [disinfected]

[Soft] side for [cleaning] dirt away

Strong, durable diamond-embossed textured wipe cleans [and deodorizes]

Textured wipe for effective cleaning

Textured diamond-embossed design cleans hard, non-porous surfaces

Textured diamond-embossed design helps remove soil from surfaces

[Textured] [scrubby] [rough] side for [scrubbing] [cleaning] [removing] [wiping] [tough] [stains] [soils] [messes]

Thick, large wipe cleans

Wipe away odors

[Wipes out] [Kills] germs that cause odors

[Wipes out] [Kills] germs [while it cleans]

[Wiping] [soft] side for [wiping] dirt away [and disinfecting]

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

F:\shared\labels\pdi\9480-6\2007\draft label 061207.doc

June 12, 2007

[The following graphics and statements are optional.  
They may be placed anywhere on label/container.]



**NOT A BABY WIPE**

**NOT FOR USE ON SKIN**



*This graphic shows use of the wipe dispenser.*

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

F:\shared\labels\pdi\9480-6\2007\draft label 061207.doc

June 12, 2007

[The following graphics and statements are optional.  
They may be placed anywhere on label/container.]



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[The following graphics and statements are optional.  
They may be placed anywhere on label/container.]



Not a Skin Wipe.  
For Use on Hard Surfaces Only.




Not a Skin Wipe.  
For Use on Hard Surfaces Only.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

F:\shared\labels\pdi\9480-6\2007\draft label 061207.doc

June 12, 2007

[The following graphics and statements are optional.  
They may be placed anywhere on label/container.]



[One side scrubs]

[One side wipes]

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

F:\shared\labels\pdf\9480-6\2007\draft label 061207.doc

June 12, 2007