

9480-6

12/10/2009

1/17

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY  
WASHINGTON, D.C. 20460

OFFICE OF PREVENTION, PESTICIDES AND TOXIC SUBSTANCES

Christina Griffin  
Agent for Professional Disposables International, Inc.  
Delta Analytical Corporation  
12510 Prosperity Drive, Suite 160  
Silver Spring, MD 20904

DEC 10 2009

Subject: Sani-Cloth Plus Germicidal Disposable Cloth  
EPA Registration No.: 9480-6  
Notification Date: November 9, 2009  
EPA Receipt Date: November 10, 2009

Dear Ms. Griffin,

This letter acknowledges receipt of your notification submitted under the provisions of FIFRA section 3(c)9 and PR Notice 98-10.

**Proposed Notification**

Add Pandemic 2009 H1N1 Influenza A Virus claim

**General Comments**

Based on a review of the submitted information, your notification to add 2009 H1N1 Influenza A Virus marketing claims in accordance with the Agency guidance for testing labeling claims against Pandemic 2009 H1N1 Influenza A (Formerly called Swine Flu) is acceptable. A copy has been placed in our records for future reference.

Should you have any questions or comments concerning this letter, please contact Velma Noble at (703) 308-6233.

Sincerely,

Velma Noble  
Product Manager (31)

Regulatory Management Branch I  
Antimicrobial Division (7510P)

CONCURRENCES							
SYMBOL							
SURNAME							
DATE							



United States  
**Environmental Protection Agency**  
 Washington, DC 20460

Registration  
 Amendment  
 Other

OPP Identifier Number  
**NOTIF**

**Application for Pesticide - Section I**

1. Company/Product Number <b>9480-6</b>	2. EPA Product Manager <b>Velma Noble</b>	3. Proposed Classification <input type="checkbox"/> None <input type="checkbox"/> Restricted
4. Company/Product (Name) <b>Professional Disposables International, Inc. (PDI)/ Sani-Cloth Plus Germicidal Disposable Cloth</b>	PM# <b>31</b>	
5. Name and Address of Applicant (Include ZIP Code) <b>Professional Disposables International, Inc. (PDI) c/o Delta Analytical Corp. 12510 Prosperity Drive, Suite 160 Silver Spring, MD 20904</b> <input type="checkbox"/> Check if this is a new address	6. <b>Expedited Review.</b> In accordance with FIFRA Section 3(c)(3)(b)(i), my product is similar or identical in composition and labeling to:  EPA Reg. No. _____ Product Name _____	

**Section - II**

Amendment - Explain below  
 Resubmission in response to Agency letter dated \_\_\_\_\_  
 Notification - Explain below.

Final printed labels in response to Agency letter dated \_\_\_\_\_  
 "Me Too" Application.  
 Other - explain below.

Explanation: Use additional page(s) if necessary. (For section I and Section II.)

**Notification of additional influenza A virus Pandemic 2009 H1N1 statements per PR Notice 95-2.**

This notification is consistent with the provisions of PR Notice 95-2 and EPA regulations at 40 CFR 152.46, and no other changes have been made to the labeling or the confidential statement of formula of this product. I understand that it is a violation of 18 U.S.C. Sec. 1001 to willfully make any false statement to EPA. I further understand that if this notification is not consistent with the terms of PR Notice 95-2 and 40 CFR 152.46, this product may be in violation of FIFRA and I may be subject to enforcement action and penalties under sections 12 and 14 of FIFRA.

**Section - III**

1. Material this Product will be Packaged in:

Child-Resistant Packaging <input type="checkbox"/> Yes* <input type="checkbox"/> No	Unit Packaging <input type="checkbox"/> Yes <input type="checkbox"/> No	Water Soluble Packaging <input type="checkbox"/> Yes <input type="checkbox"/> No	2. Type of Container <input type="checkbox"/> Metal <input type="checkbox"/> Plastic <input type="checkbox"/> Glass <input type="checkbox"/> Paper <input type="checkbox"/> Other (Specify) _____
* Certification must be submitted.		If "Yes," Unit Package wgt.    No. per container	If "Yes," Unit Package wgt.    No. per container

3. Location of Net Contents Information  
 Label  Container

4. Size(s) of Retail Container

5. Location of Label Direction  
 On Label  
 On Labeling accompanying product

6. Manner In Which Label Is Affixed to Product  
 Lithograph  Other \_\_\_\_\_  
 Paper glued  
 Stenciled

**Section - IV**

1. Contact Point (Complete items directly below for identification of individual to be contacted, if necessary, to process this application.)

Name <b>Cristina Griffin</b>	Title <b>Agent for Professional Disposables International, Inc. (PDI)</b>	Telephone No. (Include Area Code) <b>301-680-7971</b>
---------------------------------	--	--

**Certification**  
 I certify that the statements I have made on this form and all attachments thereto are true, accurate and complete. I acknowledge that any kind of knowingly false or misleading statement may be punishable by fine or imprisonment or both under applicable law.

2. Signature 	3. Title <b>Agent for Professional Disposables International, Inc. (PDI)</b>	6. Date Application Received  <b>(Stamped)</b>
4. Typed Name <b>Cristina Griffin</b>	5. Date <b>November 9, 2009</b>	



November 9, 2009

Document Processing Desk (NOTIF)  
Office of Pesticide Programs (7504P)  
U.S. Environmental Protection Agency  
One Potomac Yard  
2777 S. Crystal Drive, Room S-4900  
Arlington VA 22202

Attn: Velma Noble, PM 31

**RE: Notification of additional influenza A virus Pandemic 2009 H1N1 statements per PR Notice 95-2 and EPA Guidance of 10/21/09**  
**Product: Sani-Cloth Plus Germicidal Disposable Cloth, EPA Reg. No. 9480-6**  
**Company: Professional Disposables International, Inc. (PDI)**

Dear Ms. Noble:

On behalf of PDI, I am submitting a notification per PR Notice 95-2 and EPA's "Guidance for Testing and Labeling Claims against Pandemic 2009 H1N1 Influenza A Virus (Formerly called Swine Flu)" for the product **Sani-Cloth Plus Germicidal Disposable Cloth**, EPA Reg. No. 9480-6.

This is a notification to add additional Influenza A virus Pandemic 2009 H1N1 statements as outlined in EPA's Guidance. The product qualifies for addition of the statements because the product has been tested against Influenza A and EPA previously approved Influenza A as one of the organism claims on the label. No other changes have been made to the label. The change is highlighted.

Please provide electronic approval to [cgriffin@delta-ac.com](mailto:cgriffin@delta-ac.com) .

**Enclosures**


- EPA form 8570-1
- 1 copy of the revised label (additions are on page 2 of the label and highlighted)

If you have any questions regarding this submission, please contact me at 301-680-7971 or [cgriffin@delta-ac.com](mailto:cgriffin@delta-ac.com).

Sincerely,

Cristina Griffin  
Agent for PDI

Enclosures  
cc: David Jones, PDI



4/17

# SANI-CLOTH® PLUS

## GERMICIDAL DISPOSABLE CLOTH

**ACTIVE INGREDIENTS:**

n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides.....	0.125%
n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> ) dimethyl benzyl ammonium chlorides.....	0.125%
<b>OTHER INGREDIENTS</b> .....	<b>99.750%</b>
<b>TOTAL</b>	<b>100.000%</b>

[Does not include the weight of the [cloth] [towelette] [or] [wipe]]

**KEEP OUT OF REACH OF CHILDREN**  
**CAUTION**

See back [side] panel for additional precautionary statements

NET CONTENTS [and/or Net Wt.] \_\_\_\_\_ oz (g) [or lbs. (kg)]  
\_\_\_\_ # OF CLOTHS [SIZE of each cloth]

EPA Reg. No.: 9480-6

EPA Est. No. : 9480-NY-1

or

EPA Est. No. A=9480-NY-1, B=72956-WI-1, C=72956-AR-1

Alpha character will precede batch code on product

**NOTIFICATION**  
Date Reviewed: 12/17/09  
Reviewed By: [Signature]

[Manufactured by [or for]:]  
Professional Disposables International, Inc.  
Two Nice-Pak Park  
Orangeburg, NY 10962-1376  
[For information call: [phone number to be inserted]]

[Made in USA]

**PRECAUTIONARY STATEMENTS:**

**Hazards to humans and domestic animals.**

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes and clothing. Wash thoroughly with soap and water after handling, and before eating, drinking, chewing gum, using tobacco, or using toilet.

**FIRST AID:**

**If in eye:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.  
[For additional information in case of emergency call toll free [phone number to be inserted]].  
Have the product container or label with you when calling a poison control center or doctor for going for treatment.

**PHYSICAL OR CHEMICAL HAZARDS:**

Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably, except for within the product primary brand name.

[Use one or more organisms listed below]  
[ATCC numbers are not required to be listed on production label.]

- [Bactericidal] [in 3 minutes] [Tuberculocidal][in 5 minutes]
- [[Product name] kills bacteria listed on this label:]
- Campylobacter jejuni [ATCC 29428]
- Escherichia coli (E. coli) O157:H7 [ATCC 35150]
- Escherichia coli (E. coli) [ATCC 11229]
- Methicillin Resistant Staphylococcus aureus (MRSA) [ATCC 33592]
- Mycobacterium bovis BCG (Tuberculosis) [ATCC 35743]^
- Pseudomonas aeruginosa [ATCC 15442]
- Salmonella enterica [ATCC 10708]
- Staphylococcus aureus [ATCC 6538]
- Vancomycin Resistant Enterococcus faecalis (VRE) [ATCC 51299] [in 3 minutes]

- [Virucidal\*] [ATCC # not available unless shown]
- [[Product name] is virucidal against:]
- \*Hepatitis B Virus<sup>†</sup>
- \*Hepatitis C Virus<sup>†</sup>
- \*Herpes Simplex Type 2 [ATCC VR-734]
- \*HIV-1<sup>†</sup>
- \*Influenza A (H1N1) virus [ATCC VR-98] (Strain A/Malaya/302/54)
- \*Influenza A/Hong Kong [ATCC VR-544]
- \*Respiratory Syncytial Virus (RSV) (Strain Long) [ATCC VR-26]

A quaternary/alcohol solution impregnated in a wiping cloth. A non woven disposable cloth for use in hospitals and other critical care areas where the control of the hazards of cross-contamination between treated surfaces is required. Suitable for use on equipment requiring alcohol-based products. Use on hard non porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica® and glass. Some organisms are removed from the surface by thoroughly wiping the surface with the wipe. Most remaining organisms are killed within three (3) minutes by exposure to the liquid in the wipe^.

^ 5 minute contact time for TB

[Optional: Choose one of the following statements:]

[Respiratory illnesses attributable to Pandemic 2009 H1N1 are caused by influenza A virus. This product is a broad-spectrum hard surface disinfectant that has been shown to be effective against (influenza A virus tested and listed on the label) and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 (formerly called swine flu).

This product has demonstrated effectiveness against influenza A virus and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 influenza A virus.

This product has demonstrated effectiveness against (influenza A virus tested and listed on the label) and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 (formerly called swine flu).

Kills Pandemic 2009 H1N1 influenza A virus (formerly called swine flu).

Kills Pandemic 2009 H1N1 influenza A virus.]

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

**Areas of Use:**

**Hospital and Healthcare Settings:** Ambulances, Ambulatory Surgical Centers (ASC), Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctors Offices, Donation Centers [blood][plasma][semen][milk][apheresis], Emergency Medical Settings, Emergency Rooms, E.R., Emergency Vehicles, Eye Surgical Centers, Healthcare Settings, Health Care Facilities, Home Health Care, Hospices, Hospitals, ICU, Isolation Areas, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating Rooms (OR), Ophthalmic Offices, Orthopedics, Out-Patient Surgical Centers (OPSC), Patient Care Areas, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory therapy, Surgery, Surgical Centers, Transport vehicles and X-Ray.

and/or

**Critical Care Areas:** CCU, Emergency Rooms, E.R., ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU) and Surgery.

and

**Hospital, Healthcare, and Critical Care Use Sites:** May be used on hard non porous surfaces of: ambulance equipment surfaces; bed railings; cabinets; carts; chairs; counters; dental chairs; dental countertops; dental unit instrument trays; endodontic equipment such as apex locators; exam tables; gurneys; infant incubators [interior and exterior surfaces of]; infant warmers [interior and exterior surfaces of]; isolettes; IV poles; laboratory equipment and surfaces; loupes; operating room tables and lights; operatory light switches; oxygen hoods; physical therapy (PT) equipment surfaces; slit lamps; spine back boards; stethoscopes; stretchers; stools; tables; telephones; toilet seats; and hard non porous outside surfaces of: amalgamators and dental curing lights; anesthesia machines and respiratory therapy equipment; apheresis machines, diagnostic equipment, dialysis machines; patient monitoring equipment, patient support and delivery equipment, pulp testers and motors; toilets; and ultrasound transducers and probes.

**Grocery Stores or Shopping Centers or Supermarkets or Delicatessens:** Disinfects grocery cart handles and grocery cart child seats. Disinfects shopping cart handles and child seats. Cleans grocery carts or shopping carts.

**Nursery or Daycare Centers:** Bassinets [interior and exterior surfaces of], changing tables, cribs, diaper changing stations, diaper pails, infant/child care equipment surfaces, trash cans, toys and other hard, non porous surfaces listed on this label.

**Bathroom[s] or Public Restrooms:** Basins, [bath]tubs, bathroom fixtures, cabinets, countertops, doorknobs, faucet[s], handles, hampers, handrails, shower[s] [stalls], showers, sinks, glazed tiled walls, glazed tiles, toilet [handle] [rim] [seats] [tank] [and other exterior surfaces of], urinals [exterior surfaces of] and vanity tops, and other hard non porous surfaces listed on this label.

**Health club facilities and Gyms:** Gym equipment, weight-lifting surfaces and other hard non porous surfaces listed on this label.

**Schools or Classrooms or Office Buildings or Public Facilities or Manufacturing Facilities or Libraries or Public Telephones:** Cabinets, chairs, computers, computer monitors, computer mice, counters, countertops, desks, desktops, doorknobs, fax machines, filing cabinets, filing doors, filing handles, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, telephones, toys, trash cans, work stations, and other hard non porous surfaces listed on this label.

**Transport vehicles, school busses, public busses, fire trucks, police cars:** Hard non porous surfaces of vehicles including steering wheel, hard plastic or metal flooring, seats, and walls.

**Restaurants:** Hard non porous, non food contact surfaces only including chairs, display cabinets, garbage cans, booster seats, floor, sinks, bathroom surfaces, tiled walls, telephones, cash registers, and doorknobs.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[Use one of the following as applicable to type of container:]

**[To Dispense TOWELETTES] [or WIPES] [or TOWELS] [or CLOTHS]:**

Remove lid [and discard inner seal] [from container or canister] Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the [container or] canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep [center cap of] lid closed to prevent moisture loss.

**TRAVEL PACK:**

Peel back label slowly. Pull out wipes as needed and reseal label. Unfold wipe.

**BOX DIRECTIONS:**

Open along perforations as indicated. Remove packets as needed.


**PACKET/PACKAGE:**

Tear open packet. Unfold wipe and use as directed. Discard after single use.

**To Dispense [or To use] [SOLO® PAK] [or soft pack] [or product]:**

Open package lid. Completely remove inner seal. Pull out wipe [at an angle] [towards you through orifice or opening]. [The next wipe pops up automatically.] Snap lid closed after use to retain moisture.

**SOLO® PAK:**



DO NOT FLUSH • STORE AT ROOM TEMPERATURE  
DISPOSE OF USED WIPES IN TRASH RECEPTACLE • KEEP SOLO® LID CLOSED WHEN NOT IN USE

SOLO® is a registered trademark of Nice-Pak Products, Inc.

**REFILL PACKET [or PACKAGE] [or POUCH]:**

This refill package is made for use only in a [product name] container; do not use this product in a container for another product. Do not use the refill wipes in a container that is cracked or broken. Discard the cracked or broken container in accordance with the container disposal instructions on the label.

**[Option 1]**

Remove [entire] lid from container. Tear open refill [package] at [notch] [mark] [arrow] and place [wipe] roll in empty container. Follow information on container to dispense wipes.

**[Option 2]**


Remove [entire] lid from container. Tear open refill [package] at [notch] [mark] [arrow]. Place [entire] refill package of wipes in empty container. Follow information on container to dispense wipes.

**GENERIC:**

Open lid. Pull out wipe. Close lid tightly.

**Bucket [or pail or container]**

1. Remove the lid from the [pail or bucket or container.]
2. Press open the dispenser at center of the lid.
3. Open the refill pack.
4. Locate the first wipe at the center of the roll.
5. Feed the first wipe through the dispensing orifice.
6. Reapply the lid to the [pail or bucket or container.]
7. Dispense the wipe.



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably, except for within the product primary brand name.

- For the refills:
1. Remove the lid from the [pail or bucket or container.]
  2. Remove the empty refill pack material from the [pail or bucket or container.]
  3. Press open the dispenser at center of the lid.
  4. Open the refill pack.
  5. Locate the first wipe at the center of the roll.
  6. Feed the first wipe though the dispensing orifice.
  7. Reapply the lid to the pail.
  8. Dispense the wipe.

**DIRECTIONS FOR USE**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**DEODORIZING AND DISINFECTING:** To disinfect non food contact surfaces only:  
Use a wipe to remove heavy soil. Unfold a clean wipe and thoroughly wet surface. Treated surface must remain visibly wet for a full three (3) minutes. Use additional wipe(s) if needed to assure continuous three (3) minute wet contact time. Let air dry. *[When TB is included on the label, add the following:]* These directions also apply to Mycobacterium bovis BCG (Tuberculosis) at a five (5) minute contact time at 68°F (20°C).

*[optional directions:]* Ultrasound transducers and probes:  
For use on external surfaces of ultrasound transducers and probes only.  
To disinfect: Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected by using a clean washcloth. Then thoroughly wet the probe with a towelette. Treated surface must remain visibly wet for a full three (3) minutes. Use additional wipe(s) if needed to assure continuous three (3) minute wet contact time. Let air dry. *[When TB is included on the label, add the following:]* These directions also apply to Mycobacterium bovis BCG (Tuberculosis) at a five (5) minute contact time at 68°F (20°C).

†KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68°- 77°F) or (20°- 25°C) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).


†SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, face masks, or eye coverings as appropriate when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

**Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold first germicidal cloth to remove heavy soil.

**Disposal of infectious materials:** Dispose of used towelette in accordance with local regulations for infectious waste disposal.

**Contact time:** Use second germicidal cloth to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry. Although efficacy at a two (2) minute contact time has been shown to be adequate against HIV-1, HBV and HCV this time is not sufficient for the other organisms listed on this label. Therefore a three (3) minute wet contact time must be used for other listed organisms. *[When TB is included on the label, add the following:]* five (5) minutes for TB.



The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.



[Use storage/disposal statement A if canister, tub, bucket, pail or refill dispenser directions used. Use storage/disposal statement B if travel pack, Solo® Pak, soft pack, pouch, or [refill or [bulk]] packet/package] directions used.]

**STORAGE AND DISPOSAL**


Do not contaminate water, food, or feed by storage and disposal. [Statement not required for residential or household use products.]

[Statement A:]

**Storage:** Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Dispenser [or Container] Disposal:** Nonrefillable container. Do not reuse or refill this container. Triple rinse and offer for recycling. If recycling is not available, put in trash collection. [If refillable dispenser is being used, substitute the following: **Refill Dispenser [or Container] Disposal:** Triple rinse cracked or broken containers and offer for recycling. If recycling is not available, discard in trash.]

[Statement B:]


**Storage:** Do not store near heat or open flame. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Dispenser [or Container] Disposal:** Put empty travel pack [or Solo® Pak or [refill] packet [or packaging]] in trash collection.




The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

**Optional Claims and Descriptors**

- X" x Y" [wipes]
- \_\_\_% bigger [[quilted] wipes] [vs. X" x Y" wipes]
- \_\_\_# pre-saturated towelettes X" x Y"
- © [insert year] [insert company name]
- A bleach-free (formula) [a phenolic free (formula)]
- After use, just toss wipe(s) away
- Always ready to use
- An effective method to deodorize and disinfect against odor-causing organisms on hard non porous surfaces listed on this label.
- An effective way to disinfect hard, non porous environmental surfaces listed on this label.
- A non woven disposable cloth for use in hospitals and other critical care areas where control of the hazards of cross contamination between treated surfaces is required.
- Antibacterial
- Antimicrobial
- A Sani-Surface® Brand
- Are pre-measured to have the recommended amount of germicidal solution, wipe after wipe.
- Assures proper dilution strength, product effectiveness and standardized technique.
- Bleach-free [formula]
- Bonus
- Can help reduce the risk of cross contamination on treated surfaces
- Cellulosic non woven disposable cloth
- Contains [name of active ingredient], the [same] active ingredient used in [Brand-name]
- Contains [number] canisters [plus] [and] [number] refills
- Convenient
- Deodorize hard, non porous surfaces [Use to]
- Deodorizes, and kills common household germs\*
- Diamond-embossed wipe is low-linting
- Disinfect and scrub clean
- Disinfect [and for touch-up clean]
- Disinfect external surfaces of ultrasound transducers
- Disinfect the precleaned surface
- Disinfects and deodorizes
- Disinfecting towelettes
- Disinfects grocery cart handles or shopping cart handles or grocery cart child seats or shopping cart child seats
- Disinfect[s] hard non porous environmental surfaces listed on this label [Use to]
- Dispenses one wipe at a time
- Does not contain bleach
- Does not produce mist or overspray
- Double-sided scrubbing wipes
- Durable, low-lint wipe
- Easy
- Easy [and convenient] way to [touch-up] [bathroom] [office] everyday
- Easy drop-in refill[s]
- [Easy] pop-up-tub
- Easy start lid [package]
- Easy to use [easy to refill]
- [Economical] Refill[s] [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] –pack
- [Economical] Refills [now] available
- Economy
- Effective against bathroom bacteria and viruses
- Effective against [the following micro-organisms on hard, non porous surfaces] *[list organisms]*
- Eliminates odors
- Extra large disinfecting wipes
- [Convenient] [Portable] way to disinfect hard non porous surfaces
- For daily use
- Gets rid of [or] eliminates odors
- The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.



Helps prevent [bacterial] [viral] cross contamination [between] [hard non porous] [surfaces]  
 Includes \_\_\_ count refill [pack] [contains] \_\_\_ wipes  
 Just [or] simply wipe exterior toilet surface [or] bathroom germs\* away  
 Just [pull] wipe and toss  
 Kill germs\* and bacteria instead of germ paper towels, sponges and rags [Use this product to]  
 Kills germs\* and bacteria and viruses in the bathroom [or] on the exterior toilet surfaces  
 Kills household [or] exterior toilet surface [or] bathroom germs\*  
 Kills odor-causing germs\*  
 Leaves [bathroom] [office] disinfected [and smelling fresh]  
 Multi-pack  
 Multi-purpose disinfectant [Multi [-refill] pack]  
 New and improved [use this for 6 months after announced improvement]  
 New thicker wipe with diamond-embossed design  
 New larger, thicker wipe [with diamond-embossed design]  
 Next wipe pops up automatically [everytime] [!]  
 No canister needed  
 No measuring  
 No mixing  
 Non woven disposable cloth  
 No sponges, no rags, no buckets, no mess  
 Not for sale in [insert name of state(s)]  
 Now with [%] more towel [use this for 6 months after converting to larger towel]  
 Now with [#] more wipes [use this for 6 months after converting to a larger number of towels]  
 On-the-go disinfectant  
 [A][An] [Easy] One-step disinfectant [or disinfecting] (except precleaning is required for heavy soil, blood and body fluids)  
 Perfect for disinfecting exterior toilet surfaces –or- around the exterior toilet surface  
 Pop-up dispensing  
 Pop-up [top]  
 Portable  
 Premeasured  
 Premixed  
 Premixed, premeasured, premoistened  
 Premoistened  
 Pull n/ Go [pack]  
 Easy [and convenient] way to disinfect  
 Easy way to disinfect bathroom [or] exterior toilet surfaces  
 Easy way to disinfect hard, non porous surfaces  
 Ready to Use [on the following hard non porous surfaces]  
 [Refill] package contains [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] [wet] [pre-moistened] [wipe] [refill[s]] [roll[s]]  
 Refill [size] [pack] [value pack]  
 Refillable [with \_\_\_# wipes] [dispenser] refill [pack] [for dispenser]  
 Requires no mixing or measuring and are ready-to-use  
 Resealable pouch [package] [container] inside  
 Reuse your [empty] canisters [again [and again] [!]  
 Sample  
 Scrubbing power of a sponge in a disinfecting wipe  
 Strong, durable diamond-embossed textured wipe deodorizes  
 Strong, durable low-lint wipe  
 Tear at perforation  
 Textured diamond-embossed design deodorizes hard, non porous surfaces  
 Textured diamond-embossed design  
 Thick, large wipe disinfects  
 Thick, textured wipe with diamond-embossed design  
 Time to reorder  
 Trial



The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

Tub with pop-up action  
 Tub refill  
 Twin pack  
 Use to preclean or decontaminate critical or semi-critical medical devices prior to sterilization / high-level disinfection  
 [Use to] refill[s] [count] [wipe] canisters or tubs  
 Value  
 Versatile  
 Virucide\* including HIV-1, HBV and HCV applications and effective against Staphylococcus aureus  
 Wet wipes  
 Where control of the hazards of cross contamination between treated surfaces is of prime importance  
 [Wipes out] [Kills] [common] [household] [bacteria] [viruses] [germs\*]  
 [Wipes out] [Kills] [the] flu virus, (Influenza A/Hong Kong)  
 [Wipes out] [Kills] viruses\* and bacteria [including the flu virus, (Influenza A/Hong Kong)]  
 [Wiping] [soft] side for disinfecting

**3 minute claims:** *[if TB is included on the label add footnote for TB as shown, otherwise leave off footnote and mark]*

Disinfects [in [just] three [3] minutes^]  
 Disinfects hard, non porous surfaces [in 3 minutes^]  
 Disinfects and deodorizes [in 3 minutes^]  
 Kills bacteria in 3 minutes ^  
 Kills *[list one or more organisms]*[in 3 minutes^]  
 The 3-minute^ germicidal wipe [disinfectant wipe]  
 Thick, heavy duty wipes are bactericidal and virucidal in 3 minutes^  
 [Wipes out] [Kills] [eliminates] *[list organisms]* in [just] [3 minutes ^]


^ 5 minutes for TB

**5 minute TB claim:**

Kills TB in 5 minutes

1 *[A footnote can be used to list organisms shown on page 2. Bacteria and virus can be separated into footnotes 1 and 2 if desired.]*

\* All "germ" claims must be linked to bacteria (minimum: Salmonella and Staph) and at least one virus with an asterisk or footnote. For example: Front panel: "Effective against household germs\*" Back panel: "\*\*Kills Salmonella enterica, Staphylococcus aureus, Herpes Simplex Type 2, and Influenza A/Hong Kong.



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

**CLEANING:** [cleaning directions and claims have been grouped on this page (they were previously incorporated into other parts of the label)]

**DIRECTIONS FOR USE:**

**CLEANING:** Wipe surface with towelette until clean. Let air dry.

**CLAIMS:**

An effective method to clean against odor-causing organisms on hard non porous surfaces listed on this label.

Cleans

Cleans and disinfects

Cleans, deodorizes, and kills common household germs\*

Cleans, disinfects and deodorizes

Cleans grocery carts and shopping carts

Diamond-embossed pattern [texture] provides effective cleaning

Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]

Leaves bathroom fresh and clean in [just] minutes

Multi-purpose cleaner [Multi [-refill] pack]

New thicker and larger diamond-embossed wipe cleans more surface area [use this for 6 months after converting to larger towel]

Perfect for cleaning [and disinfecting] around the toilet

Premoistened [wet] household cleaning wipe

Premoistened with specially formulated cleaners to handle daily clean-ups without bottles, sponges or buckets

[This product is]

Easy [and convenient] way to clean [and disinfect]

Easy way to clean [and disinfect] bathroom –or- exterior toilet surfaces

Easy way to clean [and disinfect] hard, non porous surfaces

[So] surfaces are more than [just] [clean] [they're] [disinfected]

[Soft] side for [cleaning] dirt away

Strong, durable diamond-embossed textured wipe cleans [and deodorizes]

Textured wipe for effective cleaning

Textured diamond-embossed design cleans hard, non porous surfaces

Textured diamond-embossed design helps remove soil from surfaces

[Textured] [scrubby] [rough] side for [scrubbing] [cleaning] [removing] [wiping] [tough] [stains] [soils] [messes]


Thick, large wipe cleans

Wipe away odors

[Wipes out] [Kills] germs\* that cause odors


[Wipes out] [Kills] germs\* [while it cleans]

[Wiping] [soft] side for [wiping] dirt away [and disinfecting]



The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.


[The following graphics and statements are optional. They may be placed anywhere on label/container.]



**NOT A BABY WIPE**



**NOT FOR USE ON SKIN**

*This graphic shows use of the wipe dispenser.*



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[The following graphics and statements are optional.  
They may be placed anywhere on label/container.]



The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.


[The following graphics and statements are optional.  
They may be placed anywhere on label/container.]



Not a Skin Wipe.  
For Use on Hard Surfaces Only.




Not a Skin Wipe.  
For Use on Hard Surfaces Only.




The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.



[The following graphics and statements are optional.  
They may be placed anywhere on label/container.]



[One side scrubs]  
[One side wipes]



Terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for the product primary brand name.