



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

November 15, 2016

Melissa Louis-Juste
Regulatory Affairs Specialist
Professional Disposables International (PDI) Inc.
100 Philips Parkway
Montvale, NJ 07645

Subject: Label Amendment – Additional alternate brand name and minor label
modifications
Product Name: Sani-Cloth Plus Germicidal Disposable Cloth
EPA Registration Number: 9480-6
Application Date: October 18, 2016
Decision Number: 522942

Dear Ms. Louis-Juste:

The amended label referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act, as amended, is acceptable. This approval does not affect any conditions that were previously imposed on this registration. You continue to be subject to existing conditions on your registration and any deadlines connected with them.

A stamped copy of your labeling is enclosed for your records. This labeling supersedes all previously accepted labeling. You must submit one copy of the final printed labeling before you release the product for shipment with the new labeling. In accordance with 40 CFR 152.130(c), you may distribute or sell this product under the previously approved labeling for 18 months from the date of this letter. After 18 months, you may only distribute or sell this product if it bears this new revised labeling or subsequently approved labeling. "To distribute or sell" is defined under FIFRA section 2(gg) and its implementing regulation at 40 CFR 152.3.

The alternate brand name "Sani-Cloth Plus Germicidal Disposable Wipe" has also been added to the product record.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) list examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Compliance.

Your release for shipment of the product constitutes acceptance of these conditions. If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6. If you have any questions, please contact Tara Flint at flint.tara@epa.gov or Eric Miederhoff at Miederhoff.eric@epa.gov.

Sincerely,

A handwritten signature in blue ink that reads "E. Miederhoff". The signature is written in a cursive, flowing style.

Eric Miederhoff
Product Manager 31
Regulatory Management Branch I
Antimicrobials Division (7510P)
Office of Pesticide Programs

Enclosure

ACCEPTED

11/15/2016

Under the Federal Insecticide, Fungicide
and Rodenticide Act as amended, for the
pesticide registered under
EPA Reg. No. 9480-6

SANI-CLOTH® PLUS

GERMICIDAL DISPOSABLE CLOTH

Alternate brand name: SANI-CLOTH® PLUS GERMICIDAL DISPOSABLE WIPE

ACTIVE INGREDIENTS:

n-Alkyl (68% C₁₂, 32% C₁₄) dimethyl ethylbenzyl ammonium chlorides..... 0.125%

n-Alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈) dimethyl benzyl ammonium chlorides..... 0.125%

OTHER INGREDIENTS..... 99.750%

TOTAL..... 100.000%

[Does not include the weight of the [cloth] [towelette] [or] [wipe]]

KEEP OUT OF REACH OF CHILDREN

CAUTION

See [back] [outer carton] [side] [container inside] [panel] for additional precautionary statements

[Contains:] NET CONTENTS [and/or Net Wt.] _____ oz (g) [or lbs. (kg)]
____ [# OF CLOTHS] [SIZE of each cloth]

EPA Reg. No.: 9480-6

EPA Est. No. : 9480-NY-1

or

EPA Est. No. A=9480-NY-1, B=72956-WI-1, C=72956-AR-1 [or others]

Alpha character will precede batch code on product

[Manufactured by [or for]:]

Professional Disposables International, Inc.

Two Nice-Pak Park

Orangeburg, NY 10962-1376 [USA]

[For information call: [phone number to be inserted]]

[Made in [USA] [with domestic and [foreign] [imported] materials] [insert country]

PRECAUTIONARY STATEMENTS:

Hazards to humans and domestic animals.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes and clothing. Wash hands thoroughly with soap and water after handling, and before eating, drinking, chewing gum, using tobacco, or using the toilet.

FIRST AID:

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment. [For additional information in case of emergency call toll free [phone number to be inserted]].

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

PHYSICAL OR CHEMICAL HAZARDS:

Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

[Symbols including asterisks, bullets or footnotes and similar may be used to designate categories of organisms.]
[ATCC numbers are not required to be listed on production label. Organisms may be listed in a table and in any order.] [ATCC # not available unless shown] [Use one or more organisms listed below]

Bacteria [(3 Minute Contact)] [Time)]

Campylobacter jejuni [ATCC 29428]
Escherichia coli (*E. coli*) [ATCC 11229]
Escherichia coli O157:H7 [ATCC 35150]
Pseudomonas aeruginosa [ATCC 15442]
Salmonella enterica [ATCC 10708]
Staphylococcus aureus [ATCC 6538]

Multi-Drug Resistant Bacteria [(3 Minute Contact)] [Time)]

Methicillin Resistant *Staphylococcus aureus* (MRSA) [ATCC 33592]
Vancomycin Resistant *Enterococcus faecalis* (VRE) [ATCC 51299] [in 3 minutes]

Viruses* [(3 Minute Contact)] [Time)]

*Herpes Simplex virus Type 2 [ATCC VR-734]
*Influenza A virus /Hong Kong [ATCC VR-544]
*Influenza A (H1N1) virus [ATCC VR-98] [Strain A/Malaya/302/54]
*Respiratory Syncytial Virus [RSV] (Strain Long) [ATCC VR-26]

Bloodborne Pathogens [(2 Minute Contact)] [Time)]

*Hepatitis B virus (HBV) -Duck HBV[†]
*Hepatitis C virus (Human) (HCV) -Bovine Diarrhea Virus[†]
*HIV-1 (AIDS virus)[†]

[A quaternary/alcohol solution impregnated in a wiping cloth. A nonwoven disposable cloth for use in hospitals and other critical care areas where the control of the hazards of cross-contamination between treated surfaces is required. Suitable for use on equipment requiring alcohol-based products. Use on hard nonporous surfaces and equipment listed on this label made of stainless steel, plastic, Formica® and glass.[Formica® is a registered trademark of The Diller Corporation.] Some organisms are removed from the surface by thoroughly wiping the surface with the wipe. Most remaining organisms are killed within [three (3) minutes by exposure to the liquid in the wipe.]

[Optional: Choose one of the following statements:]

- [Respiratory illnesses attributable to Pandemic 2009 H1N1 are caused by influenza A virus. This product is a broad-spectrum hard surface disinfectant that has been shown to be effective against (influenza A virus tested and listed on the label) and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 (formerly called swine flu).
- This product has demonstrated effectiveness against influenza A virus and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 influenza A virus.
- This product has demonstrated effectiveness against (influenza A virus tested and listed on the label) and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 (formerly called swine flu).
- Kills Pandemic 2009 H1N1 influenza A virus (formerly called swine flu).
- Kills Pandemic 2009 H1N1 influenza A virus.]

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Areas of Use:

Hospital and Healthcare Settings: Ambulances, Ambulatory Surgical Centers (ASC), Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctors Offices, Donation Centers [blood][plasma][semen][milk][apheresis], Emergency Medical Settings, Emergency Rooms, E.R., Emergency Vehicles, Eye Surgical Centers, Healthcare Settings, Healthcare Facilities, Home Healthcare, Hospices, Hospitals, ICU, Isolation Areas, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating Rooms (OR), Ophthalmic Offices, Orthopedics, Out-Patient Surgical Centers (OPSC), Patient Care Areas, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory therapy, Surgery, Surgical Centers, Transport vehicles and X-Ray.

and/or

Critical Care Areas: CCU, Emergency Rooms, E.R., ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU) and Surgery.

and

Hospital, Healthcare, and Critical Care Use Sites: May be used on hard nonporous surfaces of: ambulance equipment surfaces; bed railings; cabinets; carts; chairs; counters; dental chairs; dental countertops; dental unit instrument trays; endodontic equipment such as apex locators; exam tables; gurneys; infant incubators [interior and exterior surfaces of]; infant warmers [interior and exterior surfaces of]; isolettes; IV poles; laboratory equipment and surfaces; loupes; operating room tables and lights; operatory light switches; oxygen hoods; physical therapy (PT) equipment surfaces; slit lamps; spine back boards; stethoscopes; stretchers; stools; tables; telephones; toilet seats; and hard nonporous outside surfaces of: amalgamators and dental curing lights; anesthesia machines and respiratory therapy equipment; apheresis machines, diagnostic equipment, dialysis machines; patient monitoring equipment, patient support and delivery equipment, pulp testers and motors; toilets; and ultrasound transducers and probes.

Grocery Stores or Shopping Centers or Supermarkets or Delicatessens: Disinfects grocery cart handles and grocery cart child seats. Disinfects shopping cart handles and child seats. Cleans grocery carts or shopping carts.

Nursery or Daycare Centers: Bassinets [interior and exterior surfaces of], changing tables, cribs, diaper changing stations, diaper pails, infant/child care equipment surfaces, trash cans, toys and other hard, nonporous surfaces listed on this label.

Bathroom[s] or Public Restrooms: Basins, [bath]tubs, bathroom fixtures, cabinets, countertops, doorknobs, faucet[s], handles, hampers, handrails, shower[s] [stalls], showers, sinks, glazed tiled walls, glazed tiles, toilet [handle] [rim] [seats] [tank] [and other exterior surfaces of], urinals [exterior surfaces of] and vanity tops, and other hard nonporous surfaces listed on this label.

Health club facilities and Gyms: Gym equipment, weight-lifting surfaces and other hard nonporous surfaces listed on this label.

Schools or Classrooms or Office Buildings or Public Facilities or Manufacturing Facilities or Libraries or Public Telephones: Cabinets, chairs, computers, computer monitors, computer mouses, counters, countertops, desks, desktops, doorknobs, fax machines, filing cabinets, filing doors, filing handles, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, telephones, toys, trash cans, work stations, and other hard nonporous surfaces listed on this label.

Transport vehicles, school busses, public busses, fire trucks, police cars: Hard nonporous surfaces of vehicles including steering wheel, hard plastic or metal flooring, seats, and walls.

Restaurants: Hard nonporous, non food contact surfaces only including chairs, display cabinets, garbage cans, booster seats, floor, sinks, bathroom surfaces, tiled walls, telephones, cash registers, and doorknobs.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to [sterilization] [or] [high level disinfection].

[Use one of the following as applicable to type of container:]

[To Dispense TOWELETTES] [or WIPES] [or TOWELS] [or CLOTHS]:

Remove lid [and discard inner seal] [from container or canister] Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the [container or] canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep [center cap of] lid closed to prevent moisture loss.

TRAVEL PACK:

Peel back label slowly. Pull out wipes as needed and reseal label. Unfold wipe.

BOX DIRECTIONS:

Open along perforations as indicated. Remove packets as needed.


PACKET/PACKAGE:

Tear open packet [or package]. Unfold wipe and use as directed. Discard after single use.

To Dispense [or To use] [SOLO® PAK] [or soft pack] [or product]:

Open package lid. Completely remove inner seal. Pull out wipe [at an angle] [towards you through orifice or opening]. [The next wipe pops up automatically.] Snap lid closed after use to retain moisture.

SOLO® PAK:



DO NOT FLUSH • STORE AT ROOM TEMPERATURE
DISPOSE OF USED WIPES IN TRASH RECEPTACLE • KEEP SOLO® LID CLOSED WHEN NOT IN USE
SOLO® is a registered trademark of Nice-Pak Products, Inc.

REFILL [PACKET] [or PACKAGE] [or POUCH]:

This refill package is made for use only in a [product name] [container or] canister; do not use this product in a [container or] canister for another product. Do not use the refill wipes in a [container or] canister that is cracked or broken. Discard the cracked or broken [container or] canister in accordance with the [container or] canister disposal instructions on the label.

[Option 1]

Remove [entire] lid from [container or] canister. Tear open refill [package] at [notch] [mark] [arrow] and place [wipe] roll in empty [container or] canister. Follow information on [container or] canister to dispense wipes.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[Option 2]

Remove [entire] lid from [container or] canister. Tear open refill [package] at [notch] [mark] [arrow]. Place [entire] refill package of wipes in empty [container or] canister. Follow information on [container or] canister to dispense wipes.

GENERIC:

Open lid. Pull out wipe. Close lid tightly.

Bucket [or pail or container] *[[numbers optional]*

1. Remove [the] lid from [the] [pail or bucket or container.]
2. Press open [the] dispenser at center of [the] lid.
3. Open [the] refill pack.
4. Locate [the] first wipe at [the] center of [the] roll.
5. Feed [the] first wipe through [the] dispensing orifice.
6. Reapply [the] lid to [the] [pail or bucket or container.]
7. Dispense [the] wipe.

[For] [the] refill[s] [pouch]:

1. Remove [the] lid from [the] [pail or bucket or container.]
2. Remove [the] empty refill pack material from [the] [pail or bucket or container.]
3. Press open [the] dispenser at center of [the] lid.
4. Open [the] refill pack.
5. Locate [the] first wipe at [the] center of [the] roll.
6. Feed [the] first wipe through [the] dispensing orifice.
7. Reapply [the] lid to [the] pail.
8. Dispense [the] wipe.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

[See directions for use on the container inside]

TO [CLEAN] [AND] [,] [DISINFECT [AND DEODORIZE]: Use a wipe to remove [gross filth] [and] [or] [heavy soil] Unfold a clean wipe and thoroughly [wet] [treat] surface. Allow surface to remain [wet] [treated] for (3) minutes. [Use additional wipe(s) if needed to assure continuous three (3) minute [wet] [treated] contact time]. Let air dry.

Or

TO DISINFECT [CLEAN] [AND DEODORIZE]: Unfold a clean wipe and thoroughly [wet] [treat] surface. Allow surface to remain [treated] [wet] for three (3) minutes. Let air dry. [A potable water rinse is required for food contact surfaces.] If present, remove [gross filth] [and] [heavy] [soil loads] prior to disinfecting.

Or

TO DISINFECT [AND DEODORIZE]: To disinfect nonfood contact surfaces only: Use a wipe to remove [heavy] soil as needed. Unfold a clean wipe and thoroughly [wet] [treat] surface. Allow surface to remain [wet] [treated] for three (3) minutes. [Use additional wipe(s) if needed to assure continuous three (3) minute [wet] [treated] contact time.] Let air dry.

Or

TO DISINFECT [AND DEODORIZE]: To disinfect nonfood contact surfaces only: Unfold a clean wipe and thoroughly [wet]] [treat] surface. Allow surface to remain [wet]] [treated] for three (3) minutes. Let air dry. For [heavily] soiled surfaces, use a wipe to [pre-] clean prior to disinfecting.

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[Optional directions:] Ultrasound transducers and probes:

For use on external surfaces of ultrasound transducers and probes only.

To disinfect: Preclean gel from external surface of ultrasound transducer or probe to be disinfected by using a clean washcloth. Then thoroughly [wet] [treat] the probe with a towelette. Allow surface to remain [wet] [treated] for three (3) minutes. [Use additional wipe(s) if needed to assure continuous three (3) minute [wet] [treated] contact time]. Let air dry.

†KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRECLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68°- 77°F) or (20°- 25°C) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

†SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS

Personal protection: When using this product, wear disposable protective gloves, protective gowns, face masks, and eye coverings [per facility protocol] when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

Cleaning procedure: All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold first germicidal cloth to remove heavy soil.

Disposal of infectious materials: Dispose of used towelette in accordance with local regulations for infectious waste disposal.

Contact time: Use second germicidal cloth to thoroughly [wet] [treat] surface. Allow to remain [wet] [treated] [for] two (2) minutes, let air dry. [Although efficacy at a two (2) minute contact time has been shown to be adequate against HIV-1, HBV and HCV this time is not sufficient for the other organisms listed on this label. Therefore a three (3) minute [wet] [treated] contact time must be used for other listed organisms] or [For all other organisms see directions for contact times.]

[Use storage/disposal statement A if canister, tub, bucket, pail or refill dispenser directions used. Use storage/disposal statement B if travel pack, Solo® Pak, soft pack, pouch, or [refill or [bulk]] packet/package] directions used.]

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal. [Statement not required for residential or household use products.]

[Statement A:]

Storage: Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Dispenser [or Container] [or Package] Disposal:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling. If recycling is not available, put in trash collection.

[If refillable dispenser is being used, substitute the following: **Refill Dispenser [or Container] [or Package] Disposal:** Offer for recycling. If recycling is not available, discard in trash.]

[Statement B:]

Storage: Do not store near heat or open flame. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Dispenser [or Container] [or Package] Disposal:** Put empty [travel pack] [or Solo® Pak or [refill] packet [or packaging]] in trash collection.

The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

Optional Claims and Descriptors [Note: Deodorizing directions must be used when deodorizing claims (claims involving odor control) are used]

X" x Y" [wipes]

___% bigger [[quilted] wipes] [vs. X" x Y" wipes]

___# pre-saturated towelettes X" x Y"

© [insert year] [insert company name]

[A][An] [Easy] One-step disinfectant [or disinfecting] (except precleaning is required for heavy soil, blood and body fluids)

A bleach-free (formula) [a phenolic free (formula)]

[Absorbs] [Absorbent] [wipe]

After use, just toss wipe(s) away

Always ready to use

An effective method to deodorize and disinfect against odor-causing organisms on hard nonporous surfaces listed on this label.

An effective way to disinfect hard, nonporous environmental surfaces listed on this label.

A nonwoven disposable cloth for use in hospitals and other critical care areas where control of the hazards of cross contamination between treated surfaces is required.

[Antibacterial] [Antimicrobial] health club equipment wipe

A Quaternary [ammonium chloride] solution impregnated in a wiping cloth

A Sani-Surface® Brand

Are pre-measured to have the recommended amount of germicidal solution, wipe after wipe.

Assures proper dilution strength, product effectiveness and standardized technique.

[Bactericidal]

[Bactericidal] [and] [Virucidal*] in [only] [just] 3 minutes

Bleach-free [formula]

Bonus

Can help reduce the risk of cross contamination on treated surfaces

Cellulosic nonwoven disposable cloth

[Cloth-like texture] [Textured] [strong] [for] [with] scrubbing power

Compatible for use on hard, nonporous surfaces commonly found in healthcare settings listed on this label

Contains [name of active ingredient], the [same] active ingredient used in [Brand-name]

Contains [number] canisters [plus] [and] [number] refills

Convenient

[Convenient] [Portable] [way to disinfect hard nonporous surfaces]

Deodorize hard, nonporous surfaces [Use to]

Deodorizes, and kills common household germs*

[Diamond][_] [embossed] wipe is low-linting

Disinfect [and for touch-up clean]

Disinfect and scrub clean

Disinfect external surfaces of ultrasound transducers

Disinfect the precleaned surface

Disinfecting towelettes

Disinfects and deodorizes

Disinfects grocery cart handles or shopping cart handles or grocery cart child seats or shopping cart child seats

Disinfect[s] hard nonporous environmental surfaces listed on this label [Use to]

Disinfects hard nonporous surfaces [that people touch]

Disinfects hard nonporous [surfaces:] [health club equipment] [treadmill handles] [elliptical handles] [exercise bike handles] [weight benches]

Disinfects on the go

Dispenses one wipe at a time

Does not contain bleach

Does not produce mist or overspray

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

Double-sided scrubbing wipes
 [Durable] [,] [low-lint] wipe
 Easy
 Easy [and convenient] way to disinfect bathroom [or] exterior toilet surfaces
 Easy [and convenient] way to [touch-up] [bathroom] [office] everyday
 [Easy] [convenient] [dispensing] lid [closure] [cap] [package]
 Easy drop-in refill[s]
 [Easy] pop-up-tub
 Easy start lid [package]
 Easy to use [easy to refill]
 Easy to use package [canister] [container] keeps [your] wipes moist
 Easy way to disinfect hard, nonporous surfaces
 [Economical] Refill[s] [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] –pack
 [Economical] Refills [now] available
 Economy
 Effective against bathroom bacteria and viruses* *[use at least two bacteria and two viruses from the organism list]*
 Effective against [the following micro-organisms on hard, nonporous surfaces] *[list organisms]*
 Eliminates odors
 Eliminates potential for overspray
 [Extra] large [disinfecting] [Y" x Z" compared to our A" x B" size] wipe[s]
 For daily use
 Gets rid of [or] eliminates odors
 [Helps] [kill][s] [eliminate][s] odor-causing bacteria
 Helps prevent [bacterial] [viral] cross contamination [between] [hard nonporous] [on treated] [surfaces]
 Hospital disinfectant [disinfection] *[use only when hospital use sites on label]*
 Includes ____ count refill [pack] [contains] ____ wipes
 [Just] [or] simply wipe exterior toilet surface [or] bathroom germs* away
 [Just] [pull] wipe and toss
 Kills bacteria and viruses* in the bathroom *[use at least two bacteria and two viruses from the organism list:]* –or-
 on the exterior toilet surfaces
 [Kills] [Effective against] [MRSA] [(Methicillin Resistant Staphylococcus aureus)]
 Kill germs* and bacteria instead of germ paper towels, sponges and rags [Use this product to]
 Kills germs* and bacteria and viruses in the bathroom [or] on the exterior toilet surfaces
 Kills household [or] exterior toilet surface [or] bathroom germs* + *[use at least two bacteria and two viruses from the organism list]*
 Kills odor causing [bacteria] [germs*]
 Leaves [bathroom] [office] disinfected [and smelling fresh]
 Less [reduced] packaging [less to throw away]
 Moisture lock lid
 Multi-pack
 Multi-purpose disinfectant [Multi [-refill] pack]
 Next wipe pops up automatically [everytime] [!]
 No canister needed
 No measuring
 No mixing
 Nonwoven disposable cloth
 No [pinched fingers] [rips] [tears]
 No sponges, no rags, no buckets, no mess
 No threading [required]
 Non-scratching [wipe]
 Not for sale in *[insert name of state(s)]*
 Not for in-home use
 Not for sale [or use] after *[insert date]*
 Odor [elimination] [neutralization]

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

One [wipe] at a time dispensing
 On-the-go disinfectant
 [Package] *[insert type of package]* made from [x%] recycled materials
 Perfect for disinfecting exterior toilet surfaces –or- around the exterior toilet surface
 Picks up [particles] [dust] [dirt] [hair] [crumbs]
 Pop-up dispensing
 Pop-up [top]
 Portable
 Premeasured
 Premixed
 Premixed, premeasured, premoistened
 Premoistened
 [[Product name] is virucidal against:] *[use at least two bacteria and two viruses from the organism list:]*
 [[Product name] kills bacteria listed on this label:] *[use at least two bacteria and two viruses from the organism list:]*
 Pull n/ Go [pack]
 Ready to Use [on the following hard nonporous surfaces]
 Recyclable [package]
 Reduce [risk of cross contamination] [on hard nonporous surfaces]
 Reduces [packaging] waste [by x%] [when compared to regular canister or bucket or pail]
 [Refill] package contains [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] [wet] [pre-moistened] [wipe] [refill[s]] [roll[s]]
 Refill [size] [pack] [value pack]
 Refillable [with ____# wipes] [dispenser] refill [pack] [for dispenser]
 Requires no mixing or measuring and are ready-to-use
 Resealable [stay-moist] pouch [package] [container] inside
 Reuse your [empty] canisters [again [and again] [!]
 Sample
 Scrubbing power of a sponge in a disinfecting wipe
 Strong, durable [diamond] [-] [embossed] textured wipe deodorizes
 Strong, durable low-lint wipe [designed for sensitive nonporous equipment surfaces]
 Suitable for use on equipment requiring alcohol-based products
 Tear at perforation
 Textured [diamond] [-] [embossed] design [deodorizes hard nonporous surfaces]
 Thick, large wipe disinfects
 Thick, textured wipe with [diamond] [-] [embossed] design
 This [outer] [selling] [package] [container][canister] contains x% post consumer [material][plastic][resin]
 Time to reorder
 Trial
 Tub refill
 Tub with pop-up action
 Twin pack
 Use by *[insert date]*
 Use on hard nonporous surfaces only: [stainless steel] [Formica®] [glass] [glass tables] [glazed tile]
 [Test wipe on small inconspicuous area first.] [Formica® is a registered trademark of The Diller Corporation.]
 Use to preclean or decontaminate critical or semi-critical medical devices prior to sterilization / high-level disinfection
 [Use to] refill[s] [count] [wipe] canisters or tubs
 Value
 Versatile
 [Virucidal*]
 Virucide* including *HIV-1, HBV and HCV* applications and effective against *Staphylococcus aureus*
 Wet wipes
 Where control of the hazards of cross contamination between treated surfaces is of prime importance

The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

Wipe [material] made from [x] % [sustainable] [ingredients] [fibers] *[May be used only when cellulose wipe material used]*

[Wipes out] [Kills] [common] [household] [bacteria] [viruses] [germs*]

[Wipes out] [Kills] [the] flu virus, (Influenza A/Hong Kong)

[Wipes out] [Kills] viruses* and bacteria [including the flu virus, (Influenza A/Hong Kong)]

[Wiping] [soft] side for disinfecting

X% [more] [wipes] [free] [vs. Z count] [canisters] [packs]

X% reduction in packaging [material]

[Note: The following claims can be used only for 6 months]

[Better] [New] [Improved] [Stronger] [Thicker] Wipe [!]

[Better] [new] [&] [improved] [dispensing] lid [closure] [cap] [package] [dispensing] [better] [than before]

[New] [Package] [packaging] uses [x%] less material (plastic) [then] [when compared to] [regular] containers [or canisters]

New and improved

New larger, thicker wipe [with [diamond] [-] [embossed] design]

New thicker and [or with] larger [diamond] [-] [embossed] [wipe cleans more surface area]

[Now] [Reinforced] [With] Quilted [strength]

Now with [%] more towel *[use this for 6 months after converting to larger towel]*

Now with [#] more wipes *[use this for 6 months after converting to a larger number of towels]*

[Scrubs] (&) [cleans] [better than ever]

[Stronger] more durable [wipe]

Thick[er] [wipe] [with] cloth-like texture

X% [Bigger] [larger] [quilted] wipes

[x%] thicker wipes

3 minute claims:

Disinfects and deodorizes [in 3 minutes]

Disinfects hard, nonporous surfaces [in 3 minutes]

Disinfects [in [just] three [3] minutes]

Kills bacteria in 3 minutes

Kills *[list one or more organisms]* [in 3 minutes]

The 3-minute germicidal wipe [disinfectant wipe]

Thick, heavy duty wipes are bactericidal and Virucidal* in 3 minutes

[Wipes out] [Kills] [eliminates] *[list organisms]* in [just] [3 minutes]

1 *[A footnote can be used to list organisms shown on page*

2. *Bacteria and virus can be separated into footnotes 1 and 2 if desired.]*

* All “germ” claims must be linked to bacteria (minimum: *Salmonella* and *Staph*) and at least one virus with an asterisk or footnote. For example: Front panel: “Effective against household germs*” Back panel: “*Kills *Salmonella enterica*, *Staphylococcus aureus*, *Herpes Simplex Type 2*, and *Influenza A/Hong Kong*.”

The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

CLEANING: *[cleaning directions and claims have been grouped on this page (they were previously incorporated into other parts of the label)]*

DIRECTIONS FOR USE:

CLEANING: Wipe surface with towelette until clean. Let air dry.

CLAIMS:

An effective method to clean against odor-causing organisms on hard nonporous surfaces listed on this label.

Cleans

Cleans and disinfects

Cleans, deodorizes, and kills common household germs*

Cleans, disinfects and deodorizes

Cleans grocery carts and shopping carts

[Diamond][-][embossed] pattern [texture] provides effective cleaning

Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]

Leaves bathroom fresh and clean in [just] minutes

Multi-purpose cleaner [Multi [-refill] pack]

New thicker and larger [diamond] [-] [embossed] wipe cleans more surface area *[use this for 6 months after converting to larger towel]*

Perfect for cleaning [and disinfecting] around the toilet

Premoistened [wet] household cleaning wipe

Premoistened with specially formulated cleaners to handle daily clean-ups without bottles, sponges or buckets

[This product is]

Easy [and convenient] way to clean [and disinfect]

Easy way to clean [and disinfect] bathroom –or- exterior toilet surfaces

Easy way to clean [and disinfect] hard, nonporous surfaces

[So] surfaces are more than [just] [clean] [they're] [disinfected]

[Soft] side for [cleaning] dirt away

Strong, durable [diamond] [-] [embossed] textured wipe cleans [and deodorizes]

Textured wipe for effective cleaning

Textured [diamond] [-] [embossed] design cleans hard nonporous surfaces

Textured [diamond] [-] [embossed] design helps remove soil from surfaces

[Textured] [scrubby] [rough] side for [scrubbing] [cleaning] [removing] [wiping] [tough] [stains] [soils] [messes]

Thick, large wipe cleans

Wipe away odors

[Wipes out] [Kills] germs* that cause odors

[Wipes out] [Kills] germs* [while it cleans]

[Wiping] [soft] side for [wiping] dirt away [and disinfecting]

The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[The following graphics and statements are optional.
They may be placed anywhere on label/container.]





NOT A BABY WIPE

NOT FOR USE ON SKIN



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.

[The following graphics and statements are optional. They may be placed anywhere on label/container.]



The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.



Not a Skin Wipe.
For Use on Hard Surfaces Only.



Not a Skin Wipe.
For Use on Hard Surfaces Only.

The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

[The following graphics and statements are optional.
They may be placed anywhere on label/container.]



[One side scrubs]
[One side wipes]



The terms “wipe(s)”, “towel(s)”, “towelette(s)”, and “cloth(s)” and the product name may be used interchangeably except for within the product primary brand name.

[The following are company logos for use on PDI product labels. The purpose of the logo is for easy recognition of PDI products in the marketplace.]



The terms "wipe(s)", "towel(s)", "towelette(s)", and "cloth(s)" and the product name may be used interchangeably except for within the product primary brand name.