SAVING WITTED STATED. ID NOON	U.S. ENVIRONMENTAL PROTECTION AGENCY Office of Pesticide Programs Antimicrobials Division (7510P) 1200 Pennsylvania Ave., N.W. Washington, D.C. 20460 NOTICE OF PESTICIDE: <u>X</u> Registration <u>Reregistration</u> (under FIFRA, as amended)	EPA Reg. Number: 9480-14	Date of Issuance:
		Term of Issuance: Conditional	
		Name of Pesticide Product: Project Flash Spray	
Jean Claude Mar Manager, Regula Nice-Pak Produc	atory Affairs EPA Services ets, Inc. sposables International (PDI) Inc. sway		
Note: Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Antimicrobials Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.			
On the basis of information furnished by the registrant, the above named pesticide is hereby registered under the Federal Insecticide, Fungicide and Rodenticide Act. Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.			
This product is conditionally registered in accordance with FIFRA section $3(c)(7)(A)$. You must comply with the following conditions:			
1. Submit and/or cite all data required for registration/reregistration/registration review of your product under FIFRA when the Agency requires all registrants of similar products to submit such data.			
Signature of Approving Official:		Date:	
Zeno Bain, Actir Regulatory Mana Antimicrobials D Office of Pesticid	ivision (7510P)		

Registration Notice Conditional v.20150320

Page 2 of 2 EPA Reg. No. 9480-14 Decision No.529432

EPA Form 8570-6

- 2. You are required to comply with the data requirements described in the DCI identified below:
 - a. Hydrogen Peroxide GDCI-000595-1127

You must comply with all of the data requirements within the established deadlines. If you have questions about the Generic DCI or EDSP Order listed above, you may contact the Reevaluation Team Leader (Team 36): <u>http://www2.epa.gov/pesticide-contacts/contacts-office-pesticide-programs-antimicrobial-division</u>

- 3. The data requirements for storage stability and corrosion characteristics (Guidelines 830.6317 and 830.6320) are not satisfied. A one year study is required to satisfy these data requirements. You have 18 months from the date of registration to provide these data.
- 4. Make the following label changes before you release the product for shipment:
 - Revise the EPA Registration Number to read, "EPA Reg. No. 9480-14."
- 5. Submit one copy of the final printed label for the record before you release the product for shipment.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) list examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Compliance.

If you fail to satisfy these data requirements, EPA will consider appropriate regulatory action including, among other things, cancellation under FIFRA section 6(e). Your release for shipment of the product constitutes acceptance of these conditions. A stamped copy of the label is enclosed for your records. Please also note that the record for this product currently contains the following CSFs:

- Basic CSF dated July 18, 2017
- Alternate CSF 1 dated July 18, 2017

If you have any questions, please contact Terria Northern by phone at 703-347-0265, or via email at northern.terria@epa.gov.

Enclosure: Approved label

ACCEPTED

Dec 05, 2017 Under the Federal Insecticide, Fungicide and Rodenticide Act as amended, for the

100.00%

EPA Reg. No. 9480-14

PROJECT FLASH SPRAY

ACTIVE INGREDIENTS:

Hydrogen Peroxide
OTHER INGREDIENTS
TOTAL

KEEP OUT OF REACH OF CHILDREN CAUTION

See back [side] panel for additional precautionary statements

NET CONTENTS _____fl. oz. [or qt. or mL or L]

EPA Reg. No.: 9480-EPA Est. No. : 9480-NY-1 *Or* EPA Est. No. A=9480-NY-1, B=72956-AR-1 [*or others]* Alpha character will precede batch code on product

> [Manufactured by [*or* for]]: Professional Disposables International, Inc. Two Nice-Pak Park, Orangeburg, NY 10962-1376 [USA] [For information call: [*phone number to be inserted*]] [Made in [USA] [with domestic and imported materials] [*insert country*]

[On any panel]

PRECAUTIONARY STATEMENTS: Hazards to humans and domestic animals

CAUTION: Causes moderate eye irritation. Avoid contact with eyes and clothing. Wash hands thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet [*or* restroom].

FIRST AID:

Call a poison control center [or healthcare professional] [or doctor] for treatment advice. Have the product container or label with you when calling a poison control center [or healthcare professional] [or doctor] or going for treatment. [For additional information in case of emergency call the Poison Control Center toll free at 1-800-222-1222.]

If in eye: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

Physical or chemical hazard:

This product contains an oxidizing agent. Avoid contact with other sanitizers, cleaners or organic substances.

DIRECTIONS FOR USE:

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

[[[This product] [*Insert product name*] is a "ready-to-use" liquid, comprised of 4.04% concentration of hydrogen peroxide.] [[It] should be used at full strength and should not be diluted before use.]]

TO DISINFECT [CLEAN] [AND DEODORIZE]: Spray [*product name*] 6-8 inches from surface until thoroughly wet. Allow surface to remain visibly wet for one (1) minute. For [*Clostridium difficile*] [*C.difficile*] spores allow surface to remain visibly wet for five (5) minutes. Let air dry [or wipe excess liquid using a towel]. [A potable water rinse is required for food contact surfaces.] If present, remove [gross filth] [and] [heavy] [soil loads] prior to disinfecting. A precleaning step is required to kill [*Clostridium difficile*] [or *C.difficile*] spores and *Mycobacterium bovis* BCG (TB).

OR

TO DISINFECT [AND DEODORIZE]: Spray [*product name*] 6-8 inches from surface until thoroughly wet. Allow surface to remain visibly wet for one (1) minute. For [*Clostridium difficile*] [*C.difficile*] spores allow surface to remain visibly wet for five (5) minutes. Let air dry [or wipe excess liquid using a towel]. If present, remove [gross filth] [and] [or] [heavy soil loads] prior to disinfecting. A precleaning step is required to kill [*Clostridium difficile*] [*C.difficile*] spores and *Mycobacterium bovis* BCG (TB). [A thorough rinse with potable water is required for surfaces in direct contact with food.] [If label includes cleaning claims for grease, soap scum, grime, dirt, messes, soil, add the following as the first sentence : A precleaning step is required for [grease], [soap scum], [grime], [dirt], [messes] [soil].]

OR

TO [CLEAN] [AND] [,] [DISINFECT] [AND DEODORIZE]: Remove [gross filth] [and] [or] [heavy soil].Spray [*product name*] 6-8 inches from surface until thoroughly wet. Allow surface to remain visibly wet for one (1) minute. For [*Clostridium difficile*] [*C.difficile*] spores allow surface to remain visibly wet for five (5) minutes. Let air dry [or wipe excess liquid using a towel]. [A potable water rinse is required for food contact surfaces.] A precleaning step is required to kill [*Clostridium difficile*] [*or C.difficile*] spores and *Mycobacterium bovis* BCG (TB).

OR

TO [CLEAN] [AND] [,] [DISINFECT] [AND DEODORIZE]: Spray [*product name*] 6-8 inches from surface until thoroughly wet. Allow surface to remain visibly wet for one (1) minute. For [*Clostridium difficile*] [*C.difficile*] spores allow surface to remain visibly wet for five (5) minutes. Let air dry [or wipe excess liquid using a towel]. [A potable water rinse is required for food contact surfaces.] For heavily soiled surfaces, preclean prior to disinfecting. A precleaning step is required to kill [*Clostridium difficile*] [*or C.difficile*] spores and *Mycobacterium bovis* BCG (TB).

OR

TO DISINFECT [AND DEODORIZE]: To disinfect nonfood contact surfaces only: Remove heavy soil from surface. Spray [*product name*] 6-8 inches from surface until thoroughly wet. Allow surface to remain visibly wet for one (1) minute. For [*Clostridium difficile*] [*C.difficile*] spores allow surface to remain visibly wet for five (5) minutes. Let air dry [or wipe excess liquid using a towel]. A precleaning step is required to kill [*Clostridium difficile*] [*or C.difficile*] spores and *Mycobacterium bovis* BCG (TB).

OR

TO DISINFECT [AND DEODORIZE]: To disinfect nonfood contact surfaces only: Spray [*product name*] 6-8 inches from surface until thoroughly wet. Allow surface to remain visibly wet for one (1) minute. For [*Clostridium difficile*] [*C.difficile*] spores allow surface to remain visibly wet for five (5) minutes. Let air dry [or wipe excess liquid using a towel]. For [heavily] soiled surfaces, [pre] clean prior to disinfecting. A precleaning step is required to kill [*Clostridium difficile*] [*or C.difficile*] spores and *Mycobacterium bovis* BCG (TB).

[The following sets of directions are optional and for hospital/medical facility use only:]

[Norovirus Disinfection] [*or* Tuberculosis (TB) [Outbreaks]]: This product is effective against [Norovirus (as Feline calicivirus)] [*or Mycobacterium bovis* BCG (Tuberculosis)] in one minute at 68°F (20°C). Spray [*product name*] 6-8 inches from surface until thoroughly wet. Allow surface to remain visibly wet for one (1) minute [*or* 60 seconds]. Let air dry [or wipe excess liquid using a towel]. [A potable water rinse is required for food contact surfaces.] If present, remove [gross filth] [and] [heavy] [soil loads] prior to disinfecting. Although efficacy at a 1 (one) minute contact time has been shown to be adequate against [Norovirus] [*or* TB], this time is not sufficient for [*Clostridium difficile*] [*C.difficile*] spores listed on the label. Therefore a 5 minute wet contact time must be used if [*Clostridium difficile*] [*c.difficile*] spores disinfection is required. A precleaning step is required to kill [*Clostridium difficile*] [*or C.difficile*] spores and] *Mycobacterium bovis* BCG (TB).

and/or

Ultrasound transducers and probes: For use on external surfaces of ultrasound transducers, scopes [and] [or] probes only prior to high-level disinfection [(HLD)].

To clean [and decontaminate] prior to high-level disinfection [HLD]: Preclean gel from external surface of ultrasound transducer or probe by using a clean washcloth. Then, follow instructions to disinfect.

[†]KILLS HIV-1 (AIDS VIRUS) AND HEPATITIS B VIRUS (HBV) IN 1 MINUTE ON THE PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRECLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in 1 minute at room temperature (68°- 77°F) *or* (20°- 25°C) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) and Hepatitis B Virus (HBV).

[†]SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 and HEPATITIS B VIRUS (HBV) ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal protection: When handling HIV-1 (AIDS Virus) and HBV infected blood or body fluids; wear disposable protective gloves, protective gowns, masks and eye coverings, per facility protocol.

Cleaning procedure: All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the [germicidal] spray.

Disposal of infectious materials: Used towels that may be contaminated with blood and other body fluids should be disposed of according to local regulation for infectious waste disposal.

Contact time: Use [germicidal] spray to thoroughly wet surface. Allow surface to remain visibly wet for one (1) minute. Let air dry.

Special Instructions for cleaning prior to disinfection against *Clostridium difficile* spores.

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the sporicidal product. Cleaning is to include [vigorous] wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Waste Disposal: Materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for disposal of infectious materials.

To Disinfect: Spray 6-8 inches from surface until surface is thoroughly wet. Allow surface to remain visibly wet for five (5) minutes, let air dry.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean and decontaminate critical or semi-critical medical devices prior to sterilization or [high level disinfection] [or HLD].

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

Storage: Store this product in a cool, dry area, away from direct sunlight.

For a non-refillable container

Spray Bottle [or Container] [or Jug] Disposal: Nonrefillable container. Do not reuse or refill this container. Recycle empty container or discard in trash.

Or, if refillable dispenser is being used, use the following:

Refill Bottle [or Container] [or Jug] Disposal: Rinse cracked or broken containers and offer for recycling. If recycling is not available, discard in trash.

Areas of Use:

Hospital and Healthcare Settings: May be used on hard, nonporous surfaces in: Ambulances, Ambulatory Surgical Centers (ASC), Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctors' Offices, Donation Centers [blood][plasma][semen][milk][apheresis], Emergency Medical Settings, Emergency Rooms, E.R., Emergency Vehicles, Eye Surgical Centers, Healthcare Settings, Health Care Facilities, Home Health Care, Hospices, Hospitals, ICU, Isolation Areas, Laboratory, Laundry rooms, Long Term Care Centers, Medical Intensive Care Units (MICU), Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating Rooms (OR), Oncology Facilities, Ophthalmic Offices, Orthopedics, Out-Patient Surgical Centers (OPSC), Patient Care Areas, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physicians' offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory therapy, Skilled Nursing Facility, Surgery, Surgery Intensive Care Unit (SICU), Surgical Centers, Transport vehicles, Urgent Care Centers, and

X-Ray.

and/or

Critical Care Areas: CCU, Emergency Rooms, E.R., ICU, Medical Intensive Care Units (MICU), Neonatal Intensive Care Units (NICU), Operating Rooms, Post-Operative Care Unit (PACU), Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Unit (SICU).

And

Hospital, Healthcare, and Critical Care Use Sites: May be used on hard nonporous surfaces of: ambulance equipment surfaces; bed railings; bladder scan; blood glucose meters; blood pressure cuffs; cabinets; cardiac gym equipment; carts; chairs; cardiac monitors; workstations on wheels (WOW); counters; defibrillators; dental chairs; dental countertops; dental unit instrument trays; dopplers; electronic thermometers; endodontic equipment such as apex locators; exam tables; gurneys; infant incubators [interior and exterior surfaces of]; isolettes; IV poles; IV pumps; laboratory equipment and surfaces; loupes; operating room tables and lights; operatory light switches; Ophthalmoscopes; Otoscopes; oxygen hoods; physical therapy (PT) equipment surfaces; patient lift equipment; pulse oximeter; radiology equipment; slit lamps; spine back boards; stethoscopes; stretchers; stools; tables; telemetry equipment; telephones; toilet seats; anesthesia machines and respiratory therapy equipment; apheresis machines, diagnostic equipment, dialysis machines; patient monitoring equipment, patient support and delivery equipment, pulp testers and motors; toilets; ultrasound monitors; and ultrasound transducers and probes.

Grocery Stores *or* **Shopping Centers** *or* **Supermarkets** *or* **Delicatessens:** Disinfects grocery cart handles and grocery cart child seats. Disinfects shopping cart handles and child seats.

Nursery *or* **Daycare Centers:** Bassinets [interior and exterior surfaces of], changing tables, cribs, diaper changing stations, diaper pails, infant/child care equipment surfaces, trash cans, toys and other hard, nonporous surfaces listed on this label.

Bathroom[s] or Public Restrooms: Basins, [bath]tubs, bathroom fixtures, cabinets, countertops, doorknobs, faucet[s], handles, hampers, handrails, shower[s] [stalls], showers, sinks, glazed tiled walls, glazed tiles, toilet [handle] [rim] [seats] [tank] [and other exterior surfaces of], urinals [exterior surfaces of] and vanity tops, and other hard nonporous surfaces listed on this label.

Health club facilities, Fitness Centers, and Gyms: Gym equipment, weight-lifting surfaces and other hard nonporous surfaces listed on this label.

Schools or Classrooms or Churches or Office Buildings or Public Facilities or Manufacturing Facilities or Libraries or Public Telephones: Cabinets, chairs, computers, computer monitors, computer mouses, counters, countertops, desks, desktops, doorknobs, fax machines, filing cabinets, filing doors, filing handles, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, telephones, toys, trash cans, work stations, and other hard nonporous surfaces listed on this label.

Transport vehicles, school buses, public buses, fire trucks, police cars, airplanes, cruise ships, cargo vessels: Hard nonporous surfaces of vehicles including steering wheel, hard plastic or metal flooring, seats, and walls.

Restaurants: Hard nonporous, nonfood contact surfaces only including chairs, display cabinets, garbage cans, booster seats, floor, sinks, bathroom surfaces, glazed tiled walls, telephones, cash registers, and doorknobs.

[Symbols including asterisks, bullets or footnotes and similar may be used to designate categories of organisms.] [ATCC numbers are not required to be listed on production label. Organisms may be listed in a table and in any order.]

Antibiotic Desistant Destaria [/d Minute Contest)]

Antibiotic Resistai	nt Bacteria [(1 Minute Contact)]
Acinetobacter baumannii - Multidrug Resistant (MI	DR) [ATCC 19606] [Effective against organism resistant to
Gentamicin, Trimethoprim/Sulfa and Intermedi	ate resistance to Cefepime, Ceftazidime, Ceftriaxone]
Carbapenem Resistant - Klebsiella pneumoniae [(
ESBL Positive Escherichia coli [ATCC BAA-196]	
Methicillin Resistant Staphylococcus aureus (MRS	A) [ATCC 33592]
NDM-1- Positive Enterobacter cloacae [CDC 1000	
Vancomycin Resistant Enterococcus faecalis (VRE	
Bacteria [(1 Minute Contact)]	Non-Enveloped Viruses* [(1 Minute Contact)]
Pseudomonas aeruginosa [ATCC 15442]	*Adenovirus type 5 [ATCC VR-5] [Strain Adenoid 75]
Salmonella enterica [ATCC 10708]	*Rhinovirus type 1a [ATCC VR-1559] [Strain 2060]
Staphylococcus aureus [ATCC 6538]	*Rotavirus [Strain WA]
	* Norovirus (Feline Calicivirus) [ATCC VR-782]
Bloodborne Pathogens [(1 Minute Contact)]	Enveloped Viruses* [(1 Minute Contact)]
*Hepatitis B virus (HBV) - Duck HBV [†]	*Herpes Simplex virus type 2 [ATCC VR-734] [Strain G]
[Strain 11/04/12]	*Influenza A virus (H3N2) / Strain Hong Kong
*HIV-1+ (AIDS virus)[Strain HTLV-III _B]	[ATCC VR-544]
	*Respiratory Syncytial virus (RSV) [ATCC VR-26],
	[Strain Long]
	Pathogenic Fungi [(1 Minute Contact)]
	Candida albicans [ATCC 10231]
	Trichophyton intergitale [[(Formerly known as] [(Tested as]
	Trichophyton mentagrophytes)] [ATCC 9533]
	Mycobacterium (TB) [Tuberculocidal] [(1 Minute Contact)]
	Mycobacterium bovis BCG (Tuberculosis) (TB) [Organon
	Teknika]
	C.difficile Spores [(5 Minute Contact)]
	Clostridium difficile spores [ATCC 43598]

[Effective against all listed organisms in the presence of organic soil (5% blood serum)]

[Optional organism listings and claims:]

[*][Product name] is virucidal against Adenovirus type 5, Hepatitis B virus (HBV),), Herpes Simplex Virus type 2 (causative agent of genital herpes), HIV-1, Influenza A virus (H3N2)/Hong Kong strain (representative of the common flu virus), Respiratory Syncytial Virus [RSV], Rhinovirus type 1a, Rotavirus (Strain WA), in the presence of organic soil (5% blood serum).

[Product name] is a germicidal liquid containing a Hydrogen peroxide based solution. Recommended for use in hospitals and other critical care areas where the control of the hazards of cross-contamination between treated surfaces is required. Use on hard nonporous surfaces and equipment listed on this label made of stainless steel, plastic, Formica[®] and glass. [Formica[®] is a registered trademark of The Diller Corporation.]

The Following Emerging Viral Pathogen Claims Are Optional:

Emerging Viral Pathogen Claims

This product qualifies for emerging viral pathogen claims per the EPA's 'Guidance to Registrants: Process for Making Claims Against Emerging Viral Pathogens not on EPA-Registered Disinfectant Labels' when used in accordance with the appropriate use directions indicated below.

This product meets the criteria to make claims against certain emerging viral pathogens from the following viral category[ies]:

- Enveloped Viruses
- Large, Non-Enveloped Viruses
- Small, Non-Enveloped Viruses

For an emerging viral pathogen	follow the directions for use for the
that is a/an	following organisms on the label:
Enveloped virus	One of the following:
	Adenovirus Type 5, Strain Adenoid 75 (ATCC
	VR-5)
	Rotavirus, Strain WA (ATCC VR-2018)
	Rhinovirus Type 1A, Strain 2060 (ATCC VR-
	1559)
	Feline Calicivirus, Strain F-9 (ATCC VR-782)
Large, non-enveloped virus	One of the following:
	Rhinovirus Type 1A, Strain 2060 (ATCC VR- 1559)
	Feline Calicivirus, Strain F-9 (ATCC VR-782)
Small, non-enveloped virus	Two of the following:
	Rhinovirus Type 1A, Strain 2060 (ATCC VR-
	1559)
	Feline Calicivirus, Strain F-9 (ATCC VR-782)

[Product name] has demonstrated effectiveness against viruses similar to [name of emerging virus] on hard, [porous and/or nonporous surfaces]. Therefore, [product name] can be used against [name of emerging virus] when used in accordance with the directions for use against [name of supporting virus(es)] on [hard, porous/nonporous surfaces]. Refer to the [CDC or OIE] website at [pathogen specific website address] for additional information.

[*Name of illness/outbreak*] is caused by [*name of emerging virus*]. [*Product Name*] kills similar viruses and therefore can be used against [*name of emerging virus*] when used in accordance with the directions for use against [*name of supporting virus*(es)] on [hard, porous/nonporous surfaces]. Refer to the [*CDC or OIE*] website at [*website address*] for additional information.

<u>Optional Claims and Descriptors</u> [Note: Deodorizing directions must be used when deodorizing claims (claims involving odor control) are used]

© [insert year] [insert company name]

Always ready to use

An effective method to deodorize and disinfect against odor-causing organisms on hard, nonporous surfaces [listed on this label]

An effective way to disinfect hard, nonporous environmental surfaces listed on this label.

Antibacterial

Antimicrobial

[Antibacterial] [Antimicrobial] health club equipment spray

A Sani-Surface[®] Brand

[Bactericidal], [Tuberculocidal], [Fungicidal] [Virucidal*] [in 1 [one] minute]

Bactericidal

Bleach-free [formula] [this is included as useful info for users needing a bleach-free product]

Broad spectrum [microorganism coverage]

Broad Spectrum Disinfectant

Broad spectrum efficacy

Brought to you by [Sani-Cloth[®] Brand] [PDI] [Sani Professional]

Can help reduce the risk of cross contamination on treated hard, nonporous surfaces

Cleans [as it disinfects] [and disinfects]

Cleans [,] [Deodorizes,] and kills common household germs⁺ [use at least two bacteria and two viruses from the organism list]

Compatible for use on [hard, nonporous] surfaces commonly found in healthcare settings [listed on this label]

Completely inactivates [viruses*] [*insert virus(es)**]

Contains [number] bottles [plus] [and] [number] [refills]

Contains no VOC-emitting ingredients

Convenient

Demonstrates 6 [six] log reduction against [*Clostridium difficile*] [*C. difficile*] [C.diff] spores [in [5] [five] minutes] Demonstrates 99.9% reduction of [bacteria] [*insert bacteria*]

Deodorize hard, nonporous surfaces [Use to]

Deodorizes [,] [cleans,] and kills common household germs⁺ [use at least two bacteria and two viruses from the organism list]

Designed for killing [*Clostridium difficile*] [*C. difficile*] [C.diff] spores on hard, nonporous surfaces in hospitals Destroys [Deactivates] [[Helps] Eliminate] [*Clostridium difficile* spores] [*C.difficile* spores] [C.diff spores] [on hard, nonporous surfaces in hospital]

Disinfect external hard, nonporous surfaces of ultrasound transducers

Disinfects precleaned hard, nonporous surfaces

Disinfectant spray

Disinfects [,] [cleans] [and deodorizes] [hard] [nonporous surface(s)] [that people touch] [in [just] 1 minute]

Disinfects [grocery cart handles] [shopping cart handles] [grocery cart child seats] [shopping cart child seats]

[Disinfects] hard, nonporous surfaces

Disinfects [health club equipment] [treadmill handles] [elliptical handles] [exercise bike handles] [weight benches]

Disinfects on the go

Do not use on natural marble, windows, unpainted wood, brass, clear plastic, or colored grout.

Does not contain bleach [this is included as useful info for users needing a bleach-free product]

[Easy] [Convenient] [Portable] way to disinfect hard nonporous surfaces

Easy [and convenient] way to [touch-up] [bathroom] [office] [and other sites listed on this label] everyday Easy [to use] [to refill]

Easy to use trigger spray

[Economical] Refills [now] available

Effective against bacteria and viruses* [use at least two bacteria and two viruses from the organism list] Effective against [Clostridium difficile] [C.diff] spores [[in [5] [five] minutes]

Effective against [:] [the following micro-organisms on hard, nonporous surfaces] [*list organisms*]

Effective in the presence of 5% blood serum contamination Eliminates odors Eliminates daily mixing and measuring Eliminates [lingering] [food] odors on surfaces Environmental services use Exclusive to PDI. Inc. Expiration date [or EXP] [insert date] [For] daily use For [hospital] [healthcare] use [For] nonfood contact surfaces [only] From the makers of Sani-Cloth® Brand Fungicidal [Gets rid of] [Eliminates] odors Germicidal* [Spray] [use at least two bacteria and two viruses from the organism list] Hard, nonporous surface disinfectant [Helps] [kill(s)] [eliminate(s)] odor-causing bacteria Helps prevent cross contamination on hard nonporous treated surfaces Hospital [disinfectant] [disinfection] [disinfectant-cleaner] [use only when hospital use sites on label] [Ideal] for use on high-touch, hard, nonporous surfaces Kills bacteria and viruses* [in the bathroom] [on the exterior toilet surfaces] [use at least two bacteria and two viruses from the organism list:] Kills bacteria [in 1 minute] Kills [Clostridium difficile] [C. difficile] [C.diff] spores [in [5] [five] minutes] [Kills] [Disinfects] Norovirus [Feline Calicivirus] [in] 1 [one] minute [Kills] [Effective against] 99.9999% of [Clostridium difficile] [C. difficile] [C.diff] spores [in [5] [five] minutes] [Kills] [Effective against] MRSA [(Methicillin Resistant Staphylococcus aureus)] [in 1 (one) minute] [Kills] [Effective against] CRKP [(Carbapenem Resistant Klebsiella pneumoniae)] [in 1(one) minute] Kills odor causing bacteria [Kills] [Effective against] TB [in [1] [one] minute] [Kills] [Effective against] the spore form of Clostridium difficile [or C. difficile] [Kills] [Effective against] viruses* [in [only] 1 [one] minute] Leaves [bathroom] [office] [and other sites listed on this label] disinfected [and smelling fresh] [Made in USA] [Of Domestic and [Foreign] [Imported] Materials] [Made in [insert country]] Meets OSHA Bloodborne Pathogen Standard [or AHE Practice Guidance for environmental cleaning] [CDC Guidelines for disinfection and sterilization in healthcare settings] Multi-pack Multi-purpose disinfectant [Multi [-refill] pack] New (may be used for 6 months on marketed label) No diluting or mixing required Non-Bleach [this is included as useful info for users needing a bleach-free product] Non-Bleach Sporicidal[^] ([^] [Clostridium difficile] or [C. difficile] spores) [hospital] [disinfectant] [this is included as useful info for users needing a bleach-free product] Non-corrosive Not available in private label Not for sale in [insert name of state(s)] Not for sale [or use] after [insert date] Odor [eliminator] [neutralizer] One-Step Cleaner/Disinfectant [Grease and soap scum claims are not to be used on labeling with onestep cleaning/disinfection claims. One step cleaner/disinfectant claim should not be used for C.difficile.] Overall [5] [five] minute contact time [Package] [insert type of package] made from [x%] recycled materials Perfect for disinfecting [exterior toilet surfaces] [around the exterior toilet surface] Powers through tough messes Proprietary Formulation [product]

Ready to Use [disinfectant] [spray] Recyclable [package] Reduce [risk of] cross contamination [on hard, nonporous surfaces] Reduces exposure to microorganism(s) from treated surface [Refill] package contains [various numbers] [refill(s)] Refill [size] [pack] [value pack] Removes [odors] [dirt] Shelf-Stable formula Sporicidal[^] ([^] [*Clostridium difficile*] *or* [*C. difficile*] spores) The 5-minute [germicidal* spray] [disinfectant spray] Trial [pack] Tuberculocidal Twin pack Use by [insert date] Use on hard, nonporous surfaces only: [stainless steel] [Formica®] [glass [glass tables] [glazed tile] [Test solution on a small inconspicuous area first.] [Formica® is a registered trademark of The Diller Corporation.1 [Use to] disinfect(s) hard, nonporous environmental surfaces listed on this label Use to preclean or decontaminate critical or semi-critical medical devices prior to sterilization / high-level disinfection Value Versatile Virucidal* Where control of the hazards of cross contamination between treated surfaces is of prime importance With a bleach-free formula *Ithis is included as useful info for users needing a bleach-free product*

Master Label Reader Notes

Statements in brackets [] are optional or instructional.

[* All pesticidal [•]germ" claims must be linked to bacteria (minimum: Salmonella and Staph) and at least one virus. For example: Front panel: "Kills germs*" Back panel: "*Kills Salmonella enterica, Staphylococcus aureus, Herpes Simplex Virus Type 2, and Influenza A/Hong Kong"].

[[]A footnote can be used to list organisms. Bacteria and virus can be separated into footnotes 1 and 2 if desired.]

[[]When organisms are listed next to a claim, the organisms may be listed in a footnote instead of directly next to the claim]

[Graphics depicting use sites listed on this label can appear on product label and labeling.] [The following graphics and statements are optional. They may be placed anywhere on label/container.]

FOR USE ON HARD, NONPOROUS SURFACES ONLY

Do not freeze

Do not allow to freeze





[The following are company logos for use on PDI product labels. The purpose of the logo is for easy recognition of PDI products in the marketplace.]

