

BOTTLE LABEL (Approx. 4½" high x 3¼" wide)

CARTON LABEL (Approx. 5" high x 3¼" wide)

FRONT PANEL

1 Quart

LITTON
GERMICIDAL
CONCENTRATE

Brand of 12.9% Benzalkonium Chloride

MUST BE DILUTED

FOR INSTRUMENT
DISINFECTION

ACTIVE INGREDIENT:

n-Alkyl (C₁₂ 50%, C₁₄ 30%, C₁₆ 17%, C₈-C₁₈ 3%)
Dimethyl Benzyl Ammonium Chloride
(Benzalkonium Chloride)

12.9%

INERT INGREDIENTS:

(including 17% Sodium Nitrite as
rust inhibitor)

87.1%

CAUTION: KEEP OUT OF THE REACH OF CHILDREN

(See back panel for other cautions.)

Manufactured for

LITTON DENTAL PRODUCTS
DES PLAINES, ILLINOIS

ACCEPTED

JUL 26 1968

UNDER THE FEDERAL INSECTICIDE
FUNGICIDE AND ROENTGENIC ACT
FOR LISTING PURPOSES
EPA REG. NO. 10223 / 1
TO ATTACHED COMMENTS

BOTTLE LABEL (Approx. $4\frac{1}{2}$ " high x 4" wide)

CARTON LABEL (Approx. 5" high x $3\frac{1}{4}$ " wide)

BACK PANEL

DIRECTIONS:

To prepare a one-quart solution pour 10 ml (a 10-50 ml measuring cup is included) into 1 quart of tap or distilled water. To prepare a one-gallon solution pour 40 ml into a gallon of tap or distilled water. In known hard water areas always use distilled water. Be sure to wash instruments removing all blood, debris and soap. Immerse instruments for at least 15 minutes and remove ready to use. Chemical disinfection should not be relied upon to disinfect hinged instruments, those with deep narrow crevices, hypodermic needles and syringes. However, to store such sterile instruments a 1:750 solution which has not been previously used may be employed. Precautions should be observed to prevent contamination of these solutions with resistant organisms. Quaternary ammonium type disinfectants will not eliminate spore bearing organisms, Mycobacterium tuberculosis or viral hepatitis and any instruments thought to be contaminated by such organisms should be sterilized by autoclaving.

Do not use LITTON GERMICIDAL CONCENTRATE for disinfection of optical instruments with cemented parts.

LITTON GERMICIDAL CONCENTRATE is a formula consisting of a very efficient germicidal concentrate with an anti-rust factor included and has been shown to be effective in the chemical disinfection of professional and hospital instruments and equipment. In the recommended dilutions it is odorless, non-staining and non-irritating to normal skin and will inhibit rust formation. Its extreme wetting action and low surface tension give it a distinct advantage for disinfection of hospital and dental equipment.

CAUTION: Concentrate is irritating. Avoid contamination of foodstuff. Harmful if swallowed. Avoid contact with skin and eyes. In case of contact, wash thoroughly with water. If eye irritation persists, consult physician promptly.

USDA Reg. No. XXXX-X

* * * * *

NOTE: The established name "12.9% Benzalkonium Chloride" will be at least one-half the type size of the trade name.

The signal word "CAUTION" will be 14 point or larger.

"Keep out of the reach of children" will be 10 point or larger.