

For Effective Disinfection of Surfaces and Equipment

BETASEPTIC

DIRECTIONS FOR USE

• For Routine Environmental Disinfection

1. Add 2 fluid ounces (approximately 4 tablespoonfuls) of BETASEPTIC Germicidal Concentrate per each gallon of water to obtain a 150 ppm solution of available iodine.
2. Remove gross filth and heavy soil deposits from surfaces to be cleaned prior to disinfecting with the diluted solution of BETASEPTIC Germicidal Concentrate.
3. Scrub the surfaces with the diluted solution of BETASEPTIC Germicidal Concentrate. If surface treated will come in contact with food, rinse thoroughly with potable water after disinfection.
4. Rinse all cleaning equipment with clear water after each use.

• For Cold Disinfection of Surgical and Hospital Instruments

1. Remove all adhering blood and serous exudates from surgical and hospital instruments prior to immersion in the diluted solution of BETASEPTIC Germicidal Concentrate.
2. Dilute BETASEPTIC Germicidal Concentrate as previously described or in accordance with dilution table listed below. Immerse instruments in the diluted solution for at least 10 minutes and then rinse thoroughly with sterile water.

Dilution Table

Amount of BETASEPTIC Germicidal Solution desired	1 Quart	1/2 Gallon	1 Gallon	5 Gallons
Amount of BETASEPTIC Germicidal Concentrate to yield 150 parts per million iodine	1/2 fl. oz.	1 fl. oz.	2 fl. oz.	10 fl. oz.

Note: 1 fluid oz. = 2 tablespoonfuls or 6 teaspoonfuls or 30 cc.

AVOID STORING AT EXCESSIVE HEAT

dist. by

THE SUPPLY COMPANY

For Environmental Use

Professional Detergent Concentrate

Germicide / Disinfectant / Detergent

For Effective Use In:

- Hospital Disinfection
- Food Service Sanitizing
- Laboratory Disinfection

Active Ingredient: Polyvinylpyrrolidone-iodine	7.5%
(providing .75% of titratable iodine)	
Inert Ingredients:	92.5%

Caution: Keep out of reach of children. Avoid contact with eyes. In case of contact, flush thoroughly with plenty of water. If irritation persists, get medical attention. Harmful if swallowed.

Directions for use and dilution instructions appear on the back panel.

Net Contents: One Gallon

Professional Uses:

- Operating Room Disinfection
- Surgical and Hospital Instruments
- Isolation Techniques and Nursing Procedures
- Thermometer Disinfection, Stethoscopes
- Physical Therapy Departments
- Autopsy Facilities
- Food Service
- Laboratories

Housekeeping Uses:

- Walls, Floors, Tables
- Woodwork, Beds and Other Inanimate Surfaces
- Telephones
- Combs

Other Uses:

- Diet Kitchens and Cafeterias