



### **CONCENTRATED GERMICIDAL SOLUTION**

A most effective, stable chemical solution for the cold disinfection of Dental and Medical instruments, prosthetic and orthodontic appliances, plastic and rubber apparatus, porcelain bowls and non-porous surfaces such as stainless steel table tops.

### **DIRECTIONS FOR USE**

Wash instruments to be disinfected with warm water and soap to remove all debris, pus, blood or exudates. Thoroughly rinse away all soap with clear water, immerse cleaned instruments in properly diluted Germicidal Solution (1:100) for from 10 to 15 minutes. No rust inhibitor need be added. Leave instruments in solution until ready to be used. Replace solution frequently.

### **WARNING**

This solution **MUST** come in direct contact with all surfaces of instruments to insure disinfection, therefore, it should **NOT** be relied upon to disinfect hinged instruments, or any instruments with deep, narrow crevices, hypodermic needles or syringes. Nor should it be relied upon to disinfect instruments or surfaces which have been exposed to Mycobacterium tuberculosis, to spores or to the etiologic agent for viral hepatitis.

Manufactured for:

**THE J. BIRD MOYER COMPANY, INC.**  
Philadelphia, Pa. 19132  
Made in U.S.A.



These MOYCO products bear the Seal of Acceptance of the Council on Dental Therapeutics of the American Dental Association:

Camphor-Paramonochlorophenol  
Clove Oil  
Creosote  
Eucalyptol  
Eucalyptus Oil  
Eugenol  
Firmdent Denture Adhesive  
Formalin  
Iodine Tincture  
Paramonochlorophenol crystals  
Paramonochlorophenol Liquified 98%  
Paraformaldehyde  
Phenol  
Phenol Liquified  
Sulfuric Acid 50% Solution  
Tannic Acid  
Tannic Acid Solution 25% w/w  
Thymol  
Trichloroacetic Acid Solutions  
Zinc Chloride  
Zinc Chloride Solution 50% w/w  
Zinc Oxide  
  
Zinc-Oxide-Rosin-Eugenol Cement  
with Menthol and Thymol Iodide  
"ZRE CEMENT"



## CONCENTRATED GERMICIDAL SOLUTION

H. I. 4624 - 1835

A stable chemical solution for the cold disinfection of Dental and Medical instruments and appliances.

Requires No rust inhibitor.

Mixes with ordinary tap water.

Keeps instruments looking new and bright.

Non-flammable — Non-Staining.

---

When one ounce of this concentrate is mixed with 99 ounces of water (approximately 3 quarts) the result is a dilution of 1:1000 of the active ingredient.

U.S.D.A. REG. NO. 10000

Manufactured for  
**THE J. BIRD MOYER COMPANY, INC.**  
Philadelphia, Pa. 19132  
Made in USA