

DisCide Ultra Disinfecting Spray

Kills A Wide Range of Pathogens In One Minute*

Bactericidal ♦ Tuberculocidal ♦ Virucidal ♦ Fungicidal ♦ Herpecidal**

DisCide Ultra Disinfecting Spray is an herbal scented, dual chain quaternary ammonium/high level alcohol solution.

Use in Health Care & Alternate Care Settings such as hospitals, intensive care units, emergency medical settings, central supply, laboratories, clinics, nursing homes and dental suites.

DisCide Ultra Disinfecting Spray is effective against the following pathogens on hard, non-porous surfaces when used as directed:

- ♦ HIV-1** (AIDS virus)
- ♦ Methicillin Resistant Staphylococcus aureus (MRSA)
- ♦ Vancomycin Resistant Enterococcus faecalis (VRE)
- ♦ Staphylococcus aureus (Staph)
- ♦ Pseudomonas aeruginosa
- ♦ Salmonella choleraesuis
- ♦ Escherichia coli (E. coli)
- ♦ Adenovirus type 2
- ♦ Herpes Simplex Virus type 2** (genital herpes virus)
- ♦ Mycobacterium bovis (Tuberculosis bacteria) (5-minute contact time required)*
- ♦ Trichophyton mentagrophytes (Athlete's Foot fungus)

ACTIVE INGREDIENTS

n-alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride	0.12%
n-alkyl (68% C12, 32% C14) dimethyl ethyl benzyl ammonium chloride	0.12%
Isopropyl alcohol	63.25%
INERT INGREDIENTS	36.51%
TOTAL INGREDIENTS	100.00%

WARNING
KEEP OUT OF REACH OF CHILDREN
 See back panel for additional precautionary statements

EPA Reg. No. 10492-5
EPA Est. No. 8325-PA-001

DisCide Ultra Disinfecting Spray meets the disinfection requirements of OSHA's Bloodborne Pathogen Standard (29 CFR §1910.1030).

NET CONTENTS:

Rev. 10/29/1998

ACCEPTED
 with COMMENTS
 in EPA Letter Dated:
 NOV 19 1998
 Under the Federal Insecticide,
 Fungicide, and Rodenticide Act as
 amended, for the pesticide,
 registered under EPA Reg. No.
 10492-5

20/2

DIRECTIONS FOR USE:

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

TO DISINFECT: Completely preclean surfaces to be disinfected. DisCide Ultra Disinfecting Spray can be used for this purpose. Then, hold the DisCide bottle approximately 12 inches from the surface and thoroughly wet. Allow to remain wet for 1 minute. When used as a tuberculocidal disinfectant, allow the surface to remain wet for 5 minutes at 20° C.

For Disinfection of: inanimate, precleaned, hard, non-porous surfaces such as chairs, counter tops, drawer pulls, carts, baskets, tables, cabinets, telephones and other related hard surfaces in hospitals, operating rooms, intensive care and emergency units, dental or veterinary offices, surgery, recovery, anesthesia, X-Ray, Cat. Lab, orthopedics, newborn nursery, respiratory therapy, radiology, central supply, emergency medical vehicles, police and fire vehicles, correctional facilities, schools and school medical facilities, school buses, day care and non-food contact surfaces. Use on stainless steel, glass, formica and other similar surfaces. Use care when using on varnished surfaces, plexi-glass or leather.

KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings (hospitals, nursing homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (AIDS VIRUS) FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS

PERSONAL PROTECTION: Wear disposable latex gloves, protective gown, face mask and eye covering, as appropriate, when handling HIV-1 infected blood or body fluids.

CLEANING PROCEDURE: All blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this disinfectant. DisCide Ultra Disinfecting Spray can be used for this purpose.

CONTACT TIME: Thoroughly wet surface with DisCide Ultra Disinfecting Spray. Allow to air dry. The efficacy of a 1 minute contact time has been shown to be adequate against HIV-1 (AIDS virus).

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

PRECAUTIONARY STATEMENTS:

Hazards to humans and domestic animals: **WARNING.** Causes substantial but temporary eye injury. Do not get in eyes or on clothing. Wear protective eyewear. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse. **IF IN EYES:** Hold eyelids open and flush with a steady, gentle stream of water for 15 minutes. Get medical attention.

Physical or chemical hazards: Flammable. Keep away from heat or open flame.

STORAGE AND DISPOSAL

DISPOSAL: Do not reuse empty container. Discard in trash.

For product emergency, contact Chem•Tel Inc. at 1-800-255-3924, 24 hours.

For technical information on this and other fine PHC spray disinfectants, antimicrobial hand asepsis products and other infection control products, call 1-800-344-6424 or fax 203-377-8988.

PALMERO HEALTH CARE, 120 GOODWIN PLACE, STRATFORD, CT 06497

Rev. 10/29/1998

ACCEPTED
with COMMENTS
in EPA Letter Dated:

NOV 19 1998

Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No.

10492-5