ETHYLENE OXIDE & HALOCARBON 12 STERILIZING GAS

12% ETHYLENE OXIDE & 88% HALOCARBON 12 STERILIZING GAS

(DICHLORODIFLUOROMETHANE AND ETHYLENE OXIDE)
UN 3070
RQ 10 LBS./4.54 KG. (FOR ETHYLENE OXIDE)

NON-FLAMMABLE GAS

ACCEPTED APR 26 1994

Under the Federal Insecticid, Fungicide, and Rodenticide of, as amended, for the posticid registered under 1491–2.

ACTIVE INGREDIENT: ETHYLENE OXIDE (CAS NO. 75-21-8)

INERT INGREDIENT: DICHLORODIFUOROMETHANE (CAS NO. 75-71-8)_

. 75-71-8)<u>88%</u>



KEEP OUT OF REACH OF CHILDREN DANGER! PELIGRO!



PRECAUCION AL USARIO: SI USTED NO LEE INGLES, NO USE ESTÉ PRODUCTO HASTA QUE LA ETIQUETA LE HAYA SIDO EXPLICADO AMPLIAMENTE.

STATEMENT OF PRACTICAL TREATMENT

IF INHALED OR SUSPICION OF INHALATION: Remove person to fresh as fleet warm. If breathing has stopped, give arbificial respiration. If breathing is eifficial, give Oxygen, Call a physician even if no symptome are present. Keep under medical observation. Symptoms may be delayed.

IN CASE OF CONTACT WITH EYES, SKIN OR CLOTHING: immediately flush eyes or size with plonty of water for at least 15 minutes white removing contentinated clothing and shoes. Call a physician, Wash clothing before re-use. Discard contaminated shoes

NOTE TO PHYSICIAN: See attached cylinder tag

PRECAUTIONARY STATEMENTS:

HAZARDS TO HUMANS AND DOMESTIC ANIMALS:

DANGER I Non-flommable liquid and gas under pressure. Can cause eye and side borne. Be not get liquid in eyes, on side, or on alothing. May be borneled or fatel it linksied. Do not breathe gas. Can excee irritation of the respiratory trant. Symptoms of amount may be delayed. Cancer and reproductive beauty on the basis of Ethylene Guide casteet.

ENVIRONMENTAL HAZAROS:

Do not discharge effluent containing this product into lates, etwame, ponde, estuaries oceans or public waters inless the is specifically identified and addressed in an N P D.E.S. permit. Do not discharge effluent containing this product to sewer systems without prinviously notifying the sewage treatment plant authority. For guidance, contact your state water board or regional office of this E.P.A.

PHYSICAL AND CHEMICAL HAZARDS:

Contents under pressure. Back, flow into container could cause cylinder raphyre. Alway use a back flow preventing device in the container decharge, pripring. Container temperature should not exceed 125°F (52°C). Secure the cylinder in an upright position during storage and use. Keep protective valve cap in place when not in use. Use only in close systems which are raised for cylinder pressure (67 PSIG). Close valve after each use and when empty. Cylinder must not be refitted except by or with the consent of Legad Carbonic.

CONTAINS CHLOROFLUDROCARSON 12

STORAGE AND DISPOSAL:

STORAGE Do not etore near equirces of heat, open flames, or in direct sun. Cylunder storage temperature must not exceed 125°F (52°C). To minimize polymerization of Ethylene Oxide, storage of less than lambsent temperature is recommended. Do not puncture of drop cylinders.

PESTICIDE DISPOSAL. The cumterts of this cylinder World wither be used for the label purpose or must be returned in the original cylinder to. Liquid Carbonic for handling. Do not attempt disposal of this product without the express direction of Liquid Carbonic.

CONTAINER DISPOSAL. Return empty and out-dated cylinders for re-use. Prior to return, check valves for leaks. If leaks are detected, contact supplier.

BEST AVAILABLE CODY

WARNING: Contains Dichlorodifluoromethane,

a substance which harms public health and the environment by destroying ozone in the upper atmosphere

DIRECTION FOR USE: IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PROBUCT IN A MANNER INCONSISTENT WITH ITS LABELING. USE IN ACCORDANCE WITH STEINLIZER MANUFACTURERS OF PERATURE INSTRUCTIONS.
THIS PROBLICT SHOULD ONLY BE USED BY PERSONNE EXPERIENCED IN THE USE OF ETHYLLING SETENLIANS OAS INSTRUCES OR UNDER THE SUPERVISION OF SUCH AN EXPERIENCED PERSON. SEE ATTACHES CYLINDER TAIL FOR FURTHERS INFECTIONS FOR USE. USE IN ACCORDANCE WITH THE MATERIAL SAFETY DATA SHEET FOR THIS PRODUCT.

EPA REG. NO. 11491-2 EPA ESTABLISHMENT NUMBER :

11491-CA-1 11491-CA-2 11491-FL-1

11491-PA-2 11491-TX-1

11941-IL -2

NET CONTENT

FOR INDUSTRIAL AND MEDICAL USE
IN COMMERCIAL GAS STERILIZERS
DESIGNED FOR USE WITH ETHYLENE OXIDE

EMERGENCY RESPONSE PHONE# (504) 673-8831 CHEMTREC (800) 424-9300

DO NOT REMOVE THIS LABEL



Made in U.S.A.

LIQUID CARBONIC

CYLINDER GAS PRODUCTS
#19 JOINE BOULEVARD OAK \$5000K, IL 40021-5219

Form 03120

Rev. 7/93

LCCO

INTEGRAL PART OF THE LABEL FOR THIS PRODUCT.

NOTE TO PHYSICIAN: THIS PRODUCT CONTAINS ETHYLENE OXIDE AS THE ACTIVE INGREDIENT. EXPOSURE OF EYES OR SKIN MAY RESULT IN IMMEDIATE OR DELAYED SKIN IRRITATION AND BLISTERING WITH POSSIBLE ALLERGIC SKIN REACTION. THE THOROUGHLY AFFECTED AREAS SHOULD BE CLEANED AND TREATED SYMPTOMATICALLY.INHALATION MAY PRODUCE A NUMBER OF EFFECTS. THE MOST IMMEDIATELY SERIOUS OF WHICH IS PULMONARY EDEMA. A RESULT OF ETHYLENE OXIDE'S RESPIRATORY TRACT IRRITANCY AND WHICH MAY BE DELAYED IN ONSET. IF THIS OCCURS TO A SIGNIFICANT DEGREE IT CAN BE FATAL MANAGEMENT OF INHALATION EXPOSURES SHOULD INCLUDE EXTENDED OBSERVATION AND RESPIRATORY SUPPORT AS INDICATED BY THE PATIENTS PROGRESS AND PRESENTATION. DUE TO THE UNDERLYING RESPIRATORY IRRITATION IN CASES OF PULMONARY EDEMA THE USE OF HYPERBARIC OXYGEN MAY, HOWEVER BE CONTRAINDICATED AND ITS USE MUST BE CAREFULLY CONSIDERED.OTHER SEQUELAE OF INHALATION MAY INCLUDE TIGHTNESS OF THE CHEST, HEADACHE, NAUSEA, VOMITING, DIARRHEA, VERTIGO, WEAKNESS, DROWSINESS CYANOSIS, INCOORDINATION, CONVULSION AND/OR COMA DEPENDANT ON THE DEGREE OF EXPOSURE.ANY OF THESE MAY BE DELAYED IN ONSET.MANAGEMENT IS SUPPORTIVE AND SYMPTOMATIC.

BEST AVAILABLE COPY

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

THIS PRODUCT IS FOR INDUSTRIAL USE ONLY IN COMMERCIAL GAS STERILIZERS DESIGNED FOR USE WITH ETHYLENE OXIDE. THIS PRODUCT SHOULD BE USED ONLY BY PERSONS EXPERIENCED IN STERILIZATION OR FUMIGATION USING THIS MIXTURE, OR ARE UNDER THE DIRECT SUPERVISION OF PERSONS WHO ARE EXPERIENCED IN STERILIZATION OR FUMIGATION USING THIS MIXTURE.USE THIS PRODUCT IN ACCORDANCE WITH THE MANUFACTURER'S OPERATING INSTRUCTIONS FOR THE PARTICULAR STERILIZER AND ONLY AFTER READING LIQUID CARBONIC'S MATERIAL SAFETY DATA SHEET FOR 12% ETHYLENE OXIDE/88% HALOCARBON 12.ALWAYS CHECK FOR LEAKS OF THE CYLINDER VALVE BEFORE MOVING THE CYLINDER INTO THE STERILIZING FACILITY.DO NOT USE ANY OIL OR GREASE ON VALVES , REGULATOR OR FITTINGS FOR DELIVERY OF THIS GAS INTO THE STERILIZER, BACK FLOW INTO THE CONTAINER COULD CAUSE CONTAINER RUPTURE. ALWAYS USE A BACK FLOW PREVENTING DEVICE IN THE CONTAINER DISCHARGE PIPING. CHECK FOR LEAKS AFTER CONNECTING CYLINDER TO THE STERILIZER INLET PIPING. CONTAINER TEMPERATURE SHOULD NOT EXCEED

125°F (52°C) DURING STORAGE OR USE.KEEP PROTECTIVE CAP ON WHEN NOT IN USE. SECURE THE CYLINDER IN AN UPRIGHT POSITION DURING STORAGE AND USE.USE ONLY IN CLOSED SYSTEMS WHICH ARE RATED FOR THE CYLINDER PRESSURE (67 PSIG).CLOSE VALVE AFTER EACH USE AND WHEN EMPTY.

ITEMS TO BE STERILIZED SHOULD BE THOROUGHLY CLEANED OF SOIL BEFORE BEING STERILIZED. THIS PRODUCT , WHEN USED IN ACCORDANCE WITH THE STERILIZER MANUFACTURER'S OPERATING PROCEDURES, WILL STERILIZE AND/OR REDUCE MICROBIAL LOAD ON MEDICAL AND LABORATORY ITEMS, PHARMACEUTICALS, ASEPTIC PACKAGING, WHOLE AND GROUND SPICES, AND OTHER PROCESSED NATURAL SEASONING MATERIALS EXCEPT MIXTURES TO WHICH SALT HAS BEEN ADDED (SEE 40 CFR185.2850) and COSMETICS.

ACCEPTED APR 26 1994

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the postlaide registered under III G - C EPA Reg. No.

ST AVAILABLE COPY