ACCEPTED

OCT 2 3 2000

Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide, registered under EPA Reg. No. 116*33-3*4

ENVICIDE

Hospital Spray Disinfectant

Kills HIV-1 (AIDS Virus) on Pre-Cleaned Environmental Surfaces Hospital Use and Broad-Spectrum Disinfectant & Anti-Microbial Agent

An Effective Surface Disinfectant Against:

- . HIV-1 (AIDS Virus)
- Herpes Simplex II
- Streptococci

- Staphylococci
- Pseudomonas
- Salmonella
- · Mold and Mildew Influenza A, Virus (Hong Kong)

ACTIVE INGREDIENTS:

п-Alkyl (60% С , 30% С , 5% С , 5% С	
dimethyl benzyl ammonium chlorides	0.1°8
n-Alkyl (68% C.,., 32% C.,.) dimethyl ethylbenzyl	
ammonium chlorides	
Isopropanol	
INERT INGREDIENTS	46 8%
TOTAL	100 0%
EPA EST. NO. 11623-GA-1	EPA REG. NO. 11623-24

CAUTION KEEP OUT OF REACH OF CHILDREN PHYSICAL AND CHEMICAL HAZARDS

Flammable! Contents under pressure. Do not use near fire, sparks, or flame. Do not puncture or incinerate container. Exposure to temperatures above 120° F may cause bursting. See back panel for additional cautions and health hazards.

NET WT. 15 OZ. (425 grams)

MANUFACTURED BY:

APOLLO INDUSTRIES, INC.

1850 S. Cobb Industrial Blvd. Smyrna, Georgia 30082 (770) 433-0210 www.apolloind.com





EPA REV Envicide II 2660 (B)

(single

7/25.00. 12.53 AM

ambulance, then give artificial respiration, preferably mouth-to-mouth, if possible. Call a poison control center or doctor for further treatment advice

Do not contaminate water, food, or feed by storage and disposal, Pesticide Storage: Store in a cool area away from heat or open flame.

before eating, drinking, chewing gum, using tobacco or using the foilet

Pesticide Disposal. Wastes resulting from the use of this product may be disposed of according.

to applicable Federal, State, or local procedures.

Container Disposal. This container may be recycled in aerosol recycling centers. At present, There are only a few such centers in the country. Before offering for recycling, empty the can by using the product according to the label (DO NOT PUNCTURE!). If recycling is not available, wrap the container in several layers of newspaper and discard in the trash

. Use in the following settings industrial, commercial, residential, medical and health care

and motels, ambulances, equipment, and vehicles.

aeruginosa, and Trichophylon mentagrophyles.

surfaces commonly contacted by the bare loot.

to remain wet for 10 minutes, do not wipe, allow to air dry

devices prior to sterifization or high level disinfection.

soiled with blood or body fluids, in medical or healthcare settings.

Prevents and controls mold and mildew and the odors they cause

bed pans, and hampers.

. Designed for use in homes, hospitals, institutions, schools, nursing homes, offices, hotels

· Used as directed, disinfects manimate hard surfaces such as sinks, toilet seats, shower stalls, garbage cans, walls, doors, doorknobs, desktops, lelephone receivers, diaper pails,

· Deodorizes surfaces where odors are caused by microorganisms. Deodorizes shoes. . Kills HIV-1 (AIDS virus) on pre-cleaned environmental surfaces or objects previously

· Disinfects against Staphylococcus aureus, Salmonella choleraesuis, Pseudomonas

· Fungicidal against Pathogenic lungi (Athletes Foot Fungus) when used as directed on hard

DIRECTIONS FOR USE It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

SURFACE DISINFECTION AND DEODORIZATION: First remove soil deposits from surface to be treated. Hold can six inches from surface and spray until completely well. Allow surface

MOLD AND MILDEWSTAT: Spray the surface to be treated, making sure to wel completely. Let

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST

HIV-1 OF SURFACES / OBJECTS SOILED WITH BLOOD / BODY FLUIDS, PERSONAL PROTECTION: Use disposable gloves, gowns, masks, eye covering, etc. as needed to prevent contact with blood or body fluids. CLEANING PROCEDURE: Blood and body fluids

must be thoroughly cleaned from surfaces and objects before application. Then apply this

product as directed above for surface disinfection. DISPOSAL OF INFECTIOUS

MATERIALS: Blood and body fluids should be autoclaved and disposed of according to

federal, state, and local regulations for infectious waste disposal. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally

sterile areas of the body or (2) contact intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas or the body.

This product may be used to preclean or decontaminate critical or semi-critical medical

PRECAUTIONARY STATEMENTS HAZARD TO HUMANS AND DOMESTIC ANIMALS CAUTION

KEEP OUT OF REACH OF CHILDREN Caution. Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash hands

FIRST AID

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove

contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control

center or doctor for treatment advices. IF SWALLOWED: Call a poison control center or doctor immediately for treatment advice. Have a person sip a glass of water if able to swallow. Do not

induce vomiting unless told to by a poison control center or doctor. Do not give anything to an

unconscious person, IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin

immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment. IF INHALED: Move person to fresh air. If person is not breathing, call 911 or an

STORAGE AND DISPOSAL

surface air dry. Repeat at weekly intervals or when mildew growth reappears.