

SONACIDE®

STERILIZING AND DISINFECTING SOLUTION Rapid acting and stable • Use undiluted

Recommended for sterilizing and disinfecting of rubber and plastic objects, lenses and stainless steel surgical and medical instruments.

CAN BE USED IN ULTRASONIC CLEANERS

Active ingredients
Glutaraldehyde (1.5 Pentanedial) 2
Inert ingredients* 98

*Contain Non ionic ethoxylates of isomeric linear alcohols
(CH₂)_n(CH₂)_mO(CH₂-CH₂O)_nH with 11, n = 15² as synergistic agents

WARNING: HARMFUL IF SWALLOWED

E.P.A. Registration No. 15136-1 See other precautions on back label
Net contents: 1 U.S. gallon (128 Fl. Ozs.)

WAVE ENERGY SYSTEMS INC.

NEWTOWN Pa. U.S.A.

SONACIDE

DIRECTIONS FOR USE PREPARATION FOR STERILIZATION OR DISINFECTION

Recommended for use in Hospitals, Dental and Medical Offices and Veterinary Hospitals. Instruments should be thoroughly cleansed, rinsed and drained before immersing in undiluted SONACIDE. Hollow objects and needles should be flushed and filled.

STERILIZATION Immerse objects to be sterilized in SONACIDE solution maintained at 60°C for **60 minutes**. Remove and rinse with sterile water. This is a complete sterilization procedure including the destruction of viruses and tubercule bacilli.

DISINFECTION Immerse objects in SONACIDE solution maintained at 60°C for **5 minutes** to destroy vegetative organisms including Pseudomonas species and pathogenic fungi.

Any standard thermostatically controlled heating bath may be used.


OR immerse objects in the solution for **10 minutes** at room temperature to destroy vegetative organisms including Pseudomonas species and pathogenic fungi. *When used as directed this product will disinfect against Influenza A- (Hong Kong) virus, Herpes Simplex virus, Newcastle Disease Virus, and Echo virus Type 25 on inanimate surfaces. ** To destroy tubercule bacilli immerse objects for **20 minutes**. Instruments should be removed and rinsed. Where necessary store under sterile conditions.

PRECAUTIONS: For Professional Use Only

WARNING: Causes eye irritation. In case of contact with the eye, flush well with water and seek medical aid. Avoid prolonged and repeated contact with skin as possibility of sensitization exists. Avoid contamination of food.

Use in covered containers to avoid evaporation and to minimize odor. For container disposal, rinse thoroughly with water and discard.

*Aqueous solutions of glutaraldehyde are inherently acidic and generally the pH lies between 2.7 to 3.9



Indications: For use in the treatment of...
Availability: ...

No
Postage Stamp
Necessary
If Mailed in the
United States

Wave Energy Systems Inc.
Newtown Commons
Newtown, Pa. 18940
(215) 968-3833



Postage will be paid by Addressee

WAVE ENERGY SYSTEMS INC.
NEWTOWN COMMONS
NEWTOWN, PA 18940

- Please send me more information about Sonacide
- Please have a Wave Energy Systems representative call on me.

Name _____
Position _____
Company _____
Address _____
City _____ State _____ Zip Code _____
Phone _____

Directions for Use: Quaternary ammonium compounds are used to treat a
mammalian animal. For use with Sonacide, the water must be treated
with appropriate disinfectant.

Sterilization: Sonacide is not a sterilant. It is a disinfectant.
It is not a sterilant.

Disinfection: Sonacide is a quaternary ammonium compound. It is used
to disinfect surfaces and equipment.

Read the instructions on the label for the Wave Energy Systems disinfectant.
It is not a sterilant.

SONACIDE

ACCEPTED

19

15136-1

Can be used at room temperature, with moderate heat, and even in autoclaves. Can withstand ultrasonic irradiation. Leaves a smaller residue than buffered alkaline solutions. Kills M. tuberculosis in 20 minutes at room temperature.


Can be time-limited for the broad identity, can be practically used and required 4 to 5 weeks in open container at room temperature or in closed container at elevated temperature. Does not affect elasticity, the sharpness of cutting instruments, or the markings on thermometers. Does not coagulate blood.

The Story of Sonacide

The story of Sonacide begins with the discovery of a new disinfectant... (The text is very faint and difficult to read, but it appears to be a historical account of the product's development.)

Recent experiments have shown that Sonacide is a highly effective disinfectant... (This section discusses the scientific basis and efficacy of the product.)

Recent experiments recently showed that it is possible to reduce sterilization time with Sonacide to 20 minutes when combined with ultrasonics as a synergistic agent.



Additional information and contact details for the manufacturer.

A HIGH SPEED

~~AVAILABLE~~ ~~AVAILABLE TODAY~~
SONACIDE

Kills all microorganisms, including vacuum-dried spores, ~~10 times faster than any other sterilizing solution registered at the Environmental Protection Agency.~~


↓
in 60 minutes. This is a complete sterilization procedure including the destruction of viruses and tubercle bacilli.



Provides complete disinfection, ~~according to the latest Environmental Protection Agency requirements.~~

↓
this product is registered under EPA Reg. No. 15156-1.

Bactericidal action assures death of highly resistant Pseudomonas strains... kills pathogenic fungi and viruses.



Kills all microorganisms, including vacuum dried spores, ~~10 times faster than any other sterilizing solution registered at the Environmental Protection Agency.~~



in 60 minutes. This is a complete sterilization procedure including the destruction of viruses and tubercle bacilli.

provides complete disinfection, according to latest Environmental Protection Agency requirements.

Bactericidal action assures death of highly resistant Pseudomonas strains... kills pathogenic fungi and viruses.

product is registered for EPA Reg. No. 15156-1.




as a disinfectant. Temperature will kill Influenza A₂ Hong Kong, Herpes simplex, Newcastle Disease Virus, and other viruses on surface. It is

The first new sterilizing solution in over a decade

Sets new standards in speed, safety, stability and biocidal capability

COMPARES WITH OTHERS

STERILIZING AGENT	HOW USED	TEMP.	TIME FOR STERILIZATION	Main limitations of today's sterilizing methods
SONACIDE				vapour could be aggressive when heated above 37°C.
STEAM				cannot sterilize heat sensitive materials
FORMALDEHYDE				long exposure time hours for sporicidal action, aggressive vapour, pungent odor and corrosive.
IODINE				poor sporicidal activity, appreciable inactivation by organic matter.
ALKALINE GLUTARALDEHYDE etc.				requires premixing, not stable in presence of heat or ultrasonics, short life after mixing, leaves an alkaline salt residue.
ETHYLENE OXIDE				explosive gas, aggressive chemical, needs long deaeration time after processing, cannot be reused (one time use).



Industrial Applications

For inquiry, contact:
Wave Energy Systems Inc.
Newtown Commons
Newtown, Pa. 18940
(215) 968-3833