

15136-10

12/08/2012

1/2



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, D.C. 20460

DEC 8 2012

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

Mr. Kris Kontis
Medical Chemical Corporation
19430 Van Ness Ave
Torrance, CA 90501

Subject: Med-Chem Pink Germicidal Solution
EPA Registration Number: 15136-10
Amendment Date: July 13, 2012
EPA Receipt Date: July 13, 2012

Dear Mr. Kontis:

The following amendment submitted in connection with registration under FIFRA, as amended, is acceptable.

Proposed Amendment:

Response to ATP Failure by removing Pseudomonas aeruginosa and health care settings from the label

General Comments:

A stamped label with conditions is enclosed for your records. Submit one final printed label to the Agency before selling or distributing the label bearing the revised labeling.

If you have questions concerning this letter, then please contact me at 703-308-6416 or by email at campbell-mcfarlane.jacqueline@epa.gov.

Sincerely,

Jacqueline Campbell
Product Manager (34)
Regulatory Management Branch II
Antimicrobials Division (7510P)

Enclosure: Stamped label

ACTIVE INGREDIENTS:

ALKYL (40% C12, 50% C14, 10% C16) DIMETHYL	
BENZYL AMMONIUM CHLORIDES	0.133%
SODIUM CARBONATE	0.416%
INERT INGREDIENTS	99.451%
Total	100.000%

**CAUTION
KEEP OUT OF
REACH OF
CHILDREN**

Effective against:

Staphylococcus aureus, ATCC 6538

Salmonella enterica, ATCC 10708

Trichophyton mentagrophytes, ATCC 9533

MED-CHEM PINK

Germinicidal Solution

REVISED JANUARY 2012

EPA Reg. No. 15136-10
EPA Est. No. 15136-CA-01
Catalog # 142A-1g1

Net Contents: 1 Gallon

Manufactured By
MCC

Medical Chemical Corporation
19430 Van Ness Ave.
Torrance, CA 90501
1-800-424-9394

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Clean instruments thoroughly and rinse thoroughly before application of Med-Chem Pink. Use Med-Chem Pink at full strength to disinfect hard non-porous surfaces or objects made of metal, or glass and apparatus including equipment and tubing. Expose items to be disinfected for at least 10 minutes at room temperature. Do not use as a mixture with soap or other cleaning agents.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

NOT FOR STERILIZATION of instruments contaminated with spore-bearing pathogens, viral hepatitis, tuberculosis, or rusty instruments. Contaminated instruments and hypodermic needles should always be sterilized by autoclaving. For professional use only.

**PRECAUTIONARY STATEMENTS:
HAZARDOUS TO HUMANS AND DOMESTIC
ANIMALS**

CAUTION:

Causes eye and skin irritation. Do not get in eyes, on skin or on clothing. Harmful if swallowed. Avoid contamination of food. Remove contaminated clothing and wash before reuse. Wash thoroughly with soap and water after handling.

FIRST AID:

IF SWALLOWED:

- Call a poison control center of doctor immediately for treatment advice
- Have person sip a glass of water if able to swallow

FOR LABORATORY USE ONLY

- Do not induce vomiting unless told to do so by a poison control center or doctor
- Do not give anything to an unconscious person

IF ON SKIN:

- Take off contaminated clothing
- Rinse skin immediately with plenty of water for 15-20 minutes
- Call a poison control center or doctor for treatment advice

IF IN EYES:

- Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing
- Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center, doctor, or going for treatment. For emergency information concerning this product, call the National Pesticides Information Center (NPIC) at 1-800-858-7378 seven days a week, 6:30 am to 4:30 pm Pacific Time (NPIC Web site www.npic.orst.edu)

STORAGE AND DISPOSAL:

DO NOT CONTAMINATE WATER, FOOD, OR FEED BY STORAGE DISPOSAL

PESTICIDE STORAGE:

Store in dry place at a temperature between 50°F and 120°F.


PESTICIDE DISPOSAL:

Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.

CONTAINER HANDLING:

Non-refillable container. Do not refill or reuse this container. Triple rinse or equivalent. Offer for recycling or puncture and dispose in a sanitary landfill.

NOT FOR THERAPEUTIC USE



2/2