PROFESSIONAL GERMICIDE CONCENTRATE

Brand of 12.9% Benzalkonium Chloride Must Be Diluted for Instrument Disinfection

ACTIVE INGREDIENTS:

Alkyl (50% C12, 30% C14, 17% C1 3% C18) dimethyl benzyl	6 ·
Ammonium Chloride*	12.9%
Isopropyl Alcohol	13.0%
INERT INGREDIENTS	74.1%
(Formulation Contains 1756) Sodium Nitrite)	
Total Ingredients .	100.0%
*Benzalkonium Chloride	

DANGER:

Keep out of reach of Children. Do not use, pour, spill or store near heat or open flame.

See back panel for antidote statement and other precautions.

NET CONTENTS: 32 FL. OZ. (1 QUART)

Manufactured in the U.S.A by

SULTAN CHEMISTS, INC.

Hackensack, N J. 07601

To prepare a one-quart solution pour 10 ml (a 10-50 ml measuring cup is included) into 1 quart of tap or distilled water. To prepare a onegallon solution pour 40 ml into a gallon of distilled water. In known hard water areas always use distilled water. Be sure to wash instruments, removing all blood, debris and soap. Immerse instruments forat least 15 minutes. Germicides containing a Quaternary Ammonium derivative should not be relied on to destroy spore bearing organisms or Mycobacterium Tuberculosis, or the etiologic agent of viral hepatitis. Thus instruments suspected of such contamination should be sterilized by heat. Needles, hypodermic syringes, corroded instruments or those with deep narrow crevices, as well as hinged or defective plated instruments should be sterilized by autoclaving. To store such sterile instruments and maintain disinfection they may be kept immersed in Professional Germicide solution. Change solution periodically. Do not use Professional Germicidal Concentrate for disinfection of optical instruments with cemented parts. Rinse empty container thoroughly with water and discard it. Professional Germicidal Concentrate when diluted as directed, is a potent germicidal solution combining high bactericidal and disinfecting action against most kinds of pathogenic bacteria. In both concentrated and diluted form, Professional Germicide is stable. In the recommended dilutions it is odorless, non-staining and non-irritating to normal skin and will inhibit rust formation.

Keep Out of Reach of irritation. Do not get in skin when handling co contamination of food. FIRST AID: In case of plenty of water for at move and wash contar If swallowed, drink pro solution or if these ar Avoid alcohol. Call ph NOTE TO PHYSICIAN the use of gastric lava tory depression and co

1975 L'ENSONSCION NERMINENT ACT NER APORTALES 844-2

13.61

DIRECTIONS

DANGER:

Keep Out of Reach of Children. Causes eye damage and severe skin irritation. Do not get in eyes, on skin or on clothing. Protect eyes and skin when handling concentrate. Harmful or fatal if swallowed. Avoid contamination of food.

FIRST AID: In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes call a physician. Remove and wash contaminated clothing before re-use.

If swallowed, drink promptly a large quantity of milk, egg whites gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol, Call physician immediately.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

> E.P.A. No. 18184-2 E.P.A. Est. No. 18184-NJ-1