



WEX-CEL

CONCENTRATED GERMICIDAL DETERGENT

CLEANS ——— DISINFECTS ——— DEODORIZES

**GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL,
FUNGICIDAL and VIRUCIDAL***

in soft or hard water. (A.O.A.C. 400 ppm Hardness as Calcium Carbonate)

ACTIVE INGREDIENTS:

Sodium ortho-Phenylphenate 6.49%, Sodium ortho-Benzyl-para-Chlorophenate 6.50%, Sodium Dodecyl Benzene Sulfonate 1.53%, Tetrasodium Ethylenediamine Tetraacetate 1.21%.

INERT INGREDIENTS:

84.27% includes biodegradable detergents and no phosphates

EPA REGISTRATION NO. 34810-6

DANGER: KEEP OUT OF REACH OF CHILDREN.

See rear panel for additional precautionary and first aid statements.

Wexford Labs, Inc.

1035 South Vandeventer

St. Louis, Missouri 63110

**FOR INDUSTRIAL AND INSTITUTIONAL USE
CONTENTS ONE U. S. GALLON**

ACCEPTED

29 AUG 1975

WEX-CEL

CONCENTRATED GERMICIDAL DETERGENT

WEX-CEL is used for cleaning, disinfecting and deodorizing floors, walls and other similar hard surfaces. Especially suited for use in Operating Rooms, Intensive Care, Nurseries, Emergency Areas, Patient Rooms and Rest Rooms.

When used as directed **WEX-CEL** is effective against *Staphylococcus aureus*, *Salmonella choleraesuis*, *Pseudomonas aeruginosa*, pathogenic fungi, and Lipophilic viruses.

DIRECTIONS

For all disinfection, including hospital and health care facilities, use **WEX-CEL** at a dilution of 1 ounce per gallon of water. *When used at a dilution of 1 ounce per gallon of water on environmental inanimate surfaces. **WEX-CEL** is effective against Influenza A2, Herpes Simplex, Adenovirus Type 5, and Vaccinia viruses.

WEX-CEL solution may be applied with cloth, mop, or brush, or by immersing equipment to be disinfected.

When disinfecting surgical instruments with **WEX-CEL**, these instruments must be pre-cleaned to remove blood and serious exudates prior to immersion in, or treatment with **WEX-CEL**.

Rinse empty container thoroughly with water and detergent before discarding or returning drum to drum conditioner.

WEX-CEL is a synthetic germicidal detergent with broad spectrum kill of both gram positive and gram negative micro-organisms. It is an excellent cleaner which dilutes clear (1:128) in soft or hard water up to 400 ppm water hardness. Contains biodegradable detergents and

is phosphate free. Do not mix with other cleaners.

GERMICIDAL QUALIFICATIONS

(AOAC Use Dilution Confirmation) 1:128

TUBERCULOCIDAL ACTIVITY TEST

(AOAC M. Tuberculosis) 1:128

FUNGICIDAL TEST

(AOAC, T Interdigitale) 1:128

DANGER

KEEP OUT OF REACH OF CHILDREN. Corrosive. Causes eye damage and severe skin irritation. Avoid breathing spray mist. Do not get in eyes, on skin, or on clothing. Protect eyes and skin when handling concentrate. Harmful or fatal if swallowed. Avoid contamination of food.

FIRST AID

In case of contact, immediately flush eyes or skin with plenty of water for a least 15 minutes. For eyes call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution, or if these are not available, drink large quantities of water. Avoid alcohol and call a physician.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

Wexford Labs, Inc.

EPA ESTABLISHMENT NUMBER 34810-MO-1