

WEX-CIDE

**CONCENTRATED GERMICIDAL DETERGENT
CLEANS — DISINFECTS — DEODORIZES
GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL
FUNGICIDAL AND VIRUCIDAL***

in soft or hard water (A.O.A.C. 400 ppm Hardness as Calcium Carbonate) and organic soil (5% serum).

ACTIVE INGREDIENTS:	
ortho-Phenylphenol	8.73%
ortho-Benzyl-para-Chlorophenol	5.75%
INERT INGREDIENTS	87.51%
	100.00%

EPA REGISTRATION NO. 34810-8

DANGER: KEEP OUT OF REACH OF CHILDREN.

See rear panel for additional precautionary and first aid statements.

Wexford Labs, Inc.

325 Lettingwell Avenue

Kirkwood, Missouri 63122

FOR USE AND STORAGE BY MAINTENANCE PERSONNEL ONLY.

Contents: One U.S. Gallon (128 FL. Ozs.)

GERMICIDAL QUALIFICATIONS:

(AOAC Use Dilution Confirmation) in the presence of 400 ppm hardness and 5% serum at 1:256
(A.O.A.C. Ps. aeruginosa, S. aureus, S. choleraesuis).

TUBERCULOCIDAL ACTIVITY TEST:

(AOAC M. Tuberculosis) in the presence of 400 ppm hardness and 5% serum at (20°C for 10 minutes) 1:256

FUNGICIDAL TEST:

(AOAC, T. mentagrophytes) in the presence of 400 ppm hardness and 5% serum at 1:256

VIRUCIDAL ACTIVITY TEST:

(*Influenza A₁ (Hong Kong), Adenovirus Type 2, Herpes simplex Types 1 and 2, and Vaccinia virus) in the presence of 400 ppm hardness and 5% serum at 1:256

STORAGE AND DISPOSAL: Store at room temperature. Do not leave containers of this product open during storage. Rinse thoroughly. Do not reuse empty containers; wrap and put in trash.

ACCEPTED
WITH COMMENTS
BY EPA LABORATORY

SEP 06 1989

Under the authority of the
Federal Insecticide, Fungicide,
and Rodenticide Act, this product is
registered with EPA Reg. No.

34810-8

09.08
12.13.89

WEX-CIDE

CONCENTRATED GERMICIDAL DETERGENT

WEX-CIDE is used for one-step cleaning, disinfecting and deodorizing floors, walls, and other similar hard surfaces. Especially suited for use in Operating Rooms, Intensive Care, Nurseries, Emergency Areas, Patient Rooms and Rest Rooms.

When used as directed WEX-CIDE is effective against Staphylococcus aureus, Salmonella choleraesuis, Pseudomonas aeruginosa, pathogenic fungi, *Influenza A₁ (Hong Kong), Adenovirus Type 2, Herpes simplex Types 1 and 2, and Vaccinia virus.

WEX-CIDE is a synthetic germicidal detergent with broad spectrum kill of both gram positive and gram negative micro-organisms. It is an excellent cleaner. WEX-CIDE contains biodegradable detergents and is phosphate free. Do not mix with other cleaners.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

For all disinfection of hard, inanimate surfaces, including hospital and health care facilities, use WEX-CIDE at a dilution of 1/2 ounce per gallon of water, for walls, floors, bathroom fixtures, surgical instruments, stainless steel, chrome, glass and vinyl. WEX-CIDE solution may be applied with a cloth, mop or brush, or by immersing equipment to be disinfected. Wet all surfaces thoroughly. The surface must remain wet for at least 10 minutes before drying the hard surface with a cloth, or a mop in the case of floors. Gross filth and

heavy soil must be removed prior to use. When disinfecting surgical instruments must be predated prior to immersion. WEX-CIDE is suitable for use on all types of food processing equipment.

PRECAUTIONS

HAZARD TO HUMANS: Keep out of reach of children. Avoid breathing dust, or on clothing. Wear gloves when handling the concentrate.

STATEMENT OF PRACTICE: Immediately flush eyes with water for 15 minutes. For eyes call a physician. If clothing becomes wet, remove and wash before reuse. If of milk, egg whites, gelatin, etc., wash immediately with large quantities of water.

NOTE TO PHYSICIAN: Indicate the use of gas shock, respiratory depression.

Wexford Labs, Inc.

EPA ESTABLISHMENT NUMBER 34810-MO.1

09.08
12.13.89