

PM 32 34810-16

10/2

ACCEPTED
JUN 14 1996
Under the Federal Insecticide, Fungicide, and
Rodenticide Act, as amended, for the
registered under
EPA Reg. No. 34810-16



PHENEX-CEL

CONCENTRATED GERMICIDAL DETERGENT
CLEANS • DISINFECTS • DEODORIZES
GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL,
FUNGICIDAL AND VIRUCIDAL*

In soft or hard water (A.O.A.C. 400 ppm Hardness as Calcium Carbonate) and
organic soil (5% serum)

ACTIVE INGREDIENTS:

p-Phenylphenol, sodium salt	8.45%
o-Benzyl-p-Chlorophenol, sodium salt.....	7.15%
INERT INGREDIENTS:	84.40%
	100.00%

EPA REGISTRATION NO. 34810-16

**KEEP OUT OF REACH OF CHILDREN
DANGER**

See side panel for additional precautionary statements.

Wexford Labs, Inc.

325 Leffingwell Avenue

Kirkwood, Missouri 63122

FOR INDUSTRIAL AND INSTITUTIONAL USE
Contents: One U.S. Gallon (128 Fl. oz.) or 3.785 Liters

GERMICIDAL QUALIFICATIONS:

(A.O.A.C. Ps. aeruginosa, S. aureus, S. choleraesuis) in the presence of 400 ppm hardness as Calcium Carbonate and 5% serum at 1:128

TUBERCULOCIDAL ACTIVITY TEST:

(A.O.A.C. M. Tuberculosis) In the presence of 400 ppm hardness as Calcium Carbonate and 5% serum at (20° C for 10 minutes) 1:128

FUNGICIDAL TEST:

(A.O.A.C. Trichophyton mentagrophytes) in the presence of 400 ppm hardness as Calcium Carbonate and 5% serum at 1:128

VIRUCIDAL ACTIVITY TEST:

*Influenza A3/Hong Kong, Herpes Simplex Type 2, Adenovirus Type 5, and Vaccinia virus (A.O.A.C. Use Dilution Confirmation), and HIV-1 [AIDS virus] in the presence of 400 ppm hardness as Calcium Carbonate and 5% serum at 1:128

PRECAUTIONARY STATEMENTS

HAZARD TO HUMANS AND DOMESTIC ANIMALS, DANGER. Keep out of reach of children. Corrosive. Causes eye damage and skin irritation. Avoid breathing spray mist. Do not get in eyes, on skin, or on clothing. Protect eyes and skin when handling concentrate. Harmful if swallowed. Avoid contamination of food.

STATEMENT OF PRACTICAL TREATMENT: In case of contact immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of water. Avoid alcohol and call a physician.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container. Triple rinse (or equivalent), then offer for recycling or reconditioning, or dispose of in a sanitary landfill or by incineration, if allowed by state and local authorities. Consult federal, state or local disposal authorities for approved alternative procedure.

20/2

PHENEX-CEL CONCENTRATED GERMICIDAL DETERGENT

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Do not mix with other cleaners.

For all disinfection of hard surfaces in hospitals, health-care facilities, schools and nursing homes, such as floors, walls, bathroom fixtures, and shower stalls, use PHENEX-CEL at 1 ounce per gallon of water. PHENEX-CEL use-dilution may be applied with a cloth, mop or brush on surfaces such as stainless steel, chrome, glass, vinyl, ceramic and terrazzo. Wet all surfaces thoroughly. The surface must remain wet for at least ten minutes before drying with a cloth, or a mop in the case of floors. Use-dilution of PHENEX-CEL should be changed daily, or sooner if it becomes over diluted or soiled. For gross filth or heavy soil a precleaning step is required.

DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION: This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

***PHENEX-CEL KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health-care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus). PERSONAL PROTECTION: Wear disposable rubber or latex gloves, gowns, face masks, or eye coverings as appropriate, when performing manual cleaning of these surfaces. **CONTACT TIME:** Leave surfaces wet for 1 minute. This contact time will not control other common types of viruses and bacteria.

CLEANING PROCEDURE: PHENEX-CEL solution may be used to remove and clean blood and other body fluids from surfaces and objects before disinfection. For lightly soiled areas wipe off visual soil/clean with PHENEX-CEL solution, and disinfect surface with a second application of the PHENEX-CEL solution. For large spills the contaminated area should be flooded with the PHENEX-CEL solution before cleaning, then decontaminated with fresh PHENEX-CEL solution.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other body fluids should be autoclaved and disposed of according to Federal, State and local regulations for infectious waste disposal.

Wexford Labs, Inc.

EPA ESTABLISHMENT NUMBER 34810-MO.1

P-Cel-020895