

RECEIVED
MAY 1 1984
34810-17

Loose powder in buckets 10#

184



DRI-CIDE

CONCENTRATED DRY GERMICIDAL DETERGENT
CLEANS • DISINFECTS • DEODORIZES

GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL, FUNGICIDAL, and VIRUCIDAL*
In soft or hard water (A.O.A.C. 400 ppm hardness as Calcium Carbonate) with 5% blood serum.

ACTIVE INGREDIENTS:

o-Phenylphenol 7.33%
o-Benzyl-p-Chlorophenol 6.09%

INERT INGREDIENTS: 85.58%
100.00%

EPA REG. NO. 34810-17 • EST. NO. 34810-MO-1

KEEP OUT OF REACH OF CHILDREN

DANGER

See back panel for additional precautionary statements.

DRI-CIDE is a one-step cleaner-disinfectant, used for cleaning, disinfecting and deodorizing. When used as directed on non-porous, inanimate surfaces it is effective against the following organisms:

GERMICIDAL QUALIFICATIONS:

(AOAC Use Dilution Confirmation) in the presence of 400 ppm hardness and 5% serum at 1:128

(AOAC Pseudomonas aeruginosa, Salmonella choleraesuis and Staphylococcus aureus)

TUBERCULOCIDAL ACTIVITY TEST:

(AOAC M. Tuberculosis) in the presence of 400 ppm hardness and 5% serum (20°C for 10 minutes) at 1:128

FUNGICIDAL TEST:

(AOAC T. Mentagrophytes) in the presence of 400 ppm hardness and 5% blood serum at 1:128

***VIRUCIDAL ACTIVITY TEST:**

(Influenza A₂/Hong Kong, Herpes simplex Type 2, and Vaccinia virus) in the presence of 400 ppm hardness and 5% blood serum, and HIV-1 (AIDS virus) in the presence of 400 ppm hardness and 10% serum at 1:128

AREAS OF USE: Hospitals and nursing homes, veterinary and animal care facilities, industrial plants, office buildings, schools and recreational areas.

Wexford Labs, Inc.

325 LEFFINGWELL • KIRKWOOD, MISSOURI 63122

FOR INDUSTRIAL AND INSTITUTIONAL USE

NET WEIGHT: 10 Pounds (4.53592 kilograms)

PM 32 34810-17

10/2

274



DRI-CIDE
CONCENTRATED DRY GERMICIDAL DETERGENT
DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

DISINFECTION: For disinfection of non-porous hard surfaces such as floors, walls, bathroom fixtures, shower stalls and surgical instruments, stainless steel, chrome, glass, vinyl, vinyl asbestos, ceramic and terrazzo, use 1/2 ounce (1 heaping blue scoop) to 1/2 gallon of warm water; use 1 ounce (1 heaping green scoop) to 1 gallon of warm water; use 1-1/2 ounces (1 heaping clear scoop) to 1 1/2 gallons of warm water. Stir gently for a uniform solution. Solution may be applied with mop, cloth, scrubber or other suitable applicator, or by immersion. Wet surfaces thoroughly for at least ten minutes before drying. For heavily soiled areas, a pre-cleaning step is required. A fresh solution should be prepared at least daily or more often if the solution becomes diluted or soiled. When disinfecting surgical instruments these instruments must be pre-cleaned to remove blood and serous exudates prior to immersion in, or treatment with DRI-CIDE solution. DRI-CIDE is a synthetic germicidal detergent with broad spectrum kill of many gram positive and gram negative microorganisms. It contains biodegradable detergents and is phosphate free. DO NOT MIX WITH OTHER CLEANERS.

*DRI-CIDE KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus). **PERSONAL PROTECTION:** Wear disposable rubber or latex gloves, gowns, face masks, or eye coverings as appropriate, when performing manual cleaning of these surfaces. **CONTACT TIME:** Leave surfaces wet for 1 minute. This contact time will not control other common types of viruses and bacteria. **CLEANING PROCEDURE:** DRI-CIDE solution may be used to remove and clean blood and other body fluids from surfaces and objects before disinfection. For lightly soiled areas wipe off visible soil/clean with DRI-CIDE solution, and disinfect surfaces with a second application of the DRI-CIDE solution. For large spills the contaminated area should be flooded with the DRI-CIDE solution before cleaning, then decontaminated with fresh DRI-CIDE solution. **DISPOSAL OF INFECTIOUS MATERIAL:** Blood and other body fluids should be autoclaved and disposed of according to Federal, State and local regulations for infectious waste disposal.

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS. DANGER: KEEP OUT OF REACH OF CHILDREN. Corrosive. Powder causes eye damage and skin irritation. Harmful if swallowed. Do not get in eyes, on skin or clothing. Wear goggles or face shield and rubber gloves when handling concentrate. Avoid contamination of food.

STORAGE AND DISPOSAL: Store in cool dry place. Do not reuse container. Triple rinse (or equivalent), wrap and place in a trash collection.

STATEMENT OF PRACTICAL TREATMENT: In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of water. Avoid alcohol and call a physician.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

Wexford Labs, Inc.

EPA ESTABLISHMENT NUMBER 34810-MO. 1

374



DRI-CIDE

One-Half Ounce Premeasured Packet

CONCENTRATED DRY GERMICIDAL DETERGENT

CLEANS • DISINFECTS • DEODORIZES

GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL, FUNGICIDAL, and VIRUCIDAL*
in soft or hard water (A.O.A.C. 400 ppm hardness as Calcium Carbonate) with 5% blood serum.

ACTIVE INGREDIENTS:

o-Phenylphenol 7.33%

o-Benzyl-p-Chlorophenol 6.09%

INERT INGREDIENTS: 85.58%
100.00%

EPA REG. NO. 34810-17 • EST. NO. 34810-MO-1

KEEP OUT OF REACH OF CHILDREN

DANGER

See back panel for additional precautionary statements.

DRI-CIDE is a one-step cleaner-disinfectant, used for cleaning, disinfecting and deodorizing. When used as directed on non-porous, inanimate surfaces it is effective against the following organisms:

GERMICIDAL QUALIFICATIONS:

(AOAC Use Dilution Confirmation) in the presence of 400 ppm hardness and 5% serum at 1:128

(AOAC Pseudomonas aeruginosa, Salmonella choleraesuis and Staphylococcus aureus)

TUBERCULOCIDAL ACTIVITY TEST:

(AOAC M. Tuberculosis) in the presence of 400 ppm hardness and 5% serum (20°C for 10 minutes) at 1:128

FUNGICIDAL TEST:

(AOAC T. Mentagrophytes) in the presence of 400 ppm hardness and 5% blood serum at 1:128

***VIRUCIDAL ACTIVITY TEST:**

(Influenza A₂/Hong Kong, Herpes simplex Type 2, and Vaccinia virus) in the presence of 400 ppm hardness and 5% blood serum, and HIV-1(AIDS virus) in the presence of 400 ppm hardness and 10% serum at 1:128


AREAS OF USE: Hospitals and nursing homes, veterinary and animal care facilities, industrial plants, office buildings, schools and recreational areas.

Wexford Labs, Inc.

325 LEFFINGWELL • KIRKWOOD, MISSOURI 63122

FOR INDUSTRIAL AND INSTITUTIONAL USE

NET CONTENTS: 250 HALF OZ. WATER SOLUBLE PACKETS



ts

2/2
484

DRI-CIDE

CONCENTRATED DRY GERMICIDAL DETERGENT

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

DISINFECTION: For disinfection of non-porous hard surfaces such as floors, walls, bathroom fixtures, shower stalls and surgical instruments, stainless steel, chrome, glass, vinyl, vinyl asbestos, ceramic and terrazzo, use 1 ounce (two packets) per gallon of warm water. Stir gently for a uniform solution. Solution may be applied with mop, cloth, scrubber or other suitable applicator, or by immersion. Wet surfaces thoroughly for at least ten minutes before drying. For heavily soiled areas, a pre-cleaning step is required. A fresh solution should be prepared at least daily or more often if the solution becomes diluted or soiled. When disinfecting surgical instruments these instruments must be pre-cleaned to remove blood and serous exudates prior to immersion in, or treatment with DRI-CIDE solution. DRI-CIDE is a synthetic germicidal detergent with broad spectrum kill of many gram positive and gram negative microorganisms. It contains biodegradable detergents and is phosphate free. **DO NOT MIX WITH OTHER CLEANERS.**

***DRI-CIDE KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus). **PERSONAL PROTECTION:** Wear disposable rubber or latex gloves, gowns, face masks, or eye coverings as appropriate, when performing manual cleaning of these surfaces. **CONTACT TIME:** Leave surfaces wet for 1 minute. This contact time will not control other common types of viruses and bacteria. **CLEANING PROCEDURE:** DRI-CIDE solution may be used to remove and clean blood and other body fluids from surfaces and objects before disinfection. For lightly soiled areas wipe off visual soil/clean with DRI-CIDE solution, and disinfect surfaces with a second application of the DRI-CIDE solution. For large spills the contaminated area should be flooded with the DRI-CIDE solution before cleaning, then decontaminated with fresh DRI-CIDE solution. **DISPOSAL OF INFECTIOUS MATERIAL:** Blood and other body fluids should be autoclaved and disposed of according to Federal, State and local regulations for infectious waste disposal.

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS. DANGER: KEEP OUT OF REACH OF CHILDREN. Corrosive. Powder causes eye damage and skin irritation. Harmful if swallowed. Do not get in eyes, on skin or clothing. Wear goggles or face shield and rubber gloves when handling concentrate. Avoid contamination of food.

STORAGE AND DISPOSAL: Store in cool dry place. Do not reuse container. Triple rinse (or equivalent), wrap and place in a trash collection.

STATEMENT OF PRACTICAL TREATMENT: In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of water. Avoid alcohol and call a physician.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

Wexford Labs, Inc.

EPA ESTABLISHMENT NUMBER 34810-MO. 1

