

19 APR 1994

Wexford Labs, Inc.
325 Leffingwell
Kirkwood, MO 63122

Attn: Said Raziq,
Research Director

Subject: Ready To Use Wex-Cide
EPA Registration No. 34810-21
Your Resubmission of November 2, 1993

The proposed, amended product labeling, revised in response to EPA's deficiency letter of September 30, 1993 has been reviewed. Please note the following.

The labeling referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, is acceptable subject to the comments listed below. A stamped copy is enclosed for your records. Submit five copies of the finished printed labeling before releasing the product for shipment bearing the new labels.

- Within 15 days of the date of this letter, you must submit to EPA an updated, accurate Confidential Statement of Formula (CSF) which describes the composition of your product. Note that in your letter of November 2, 1993, it was stated that the un-indexed, proposed dye and corrosion inhibitor are no longer being considered for inclusion in the formulation.

If you should have any questions concerning this letter, you may call Robert Travaglini on (703) 305-6909.

Sincerely,

Ruth G. Douglas
Product Manager (32)
Antimicrobial Program Branch
Registration Division (7505C)

CONCURRENCES

SYMBOL	7505C						
SURNAME	Travaglini						
DATE	19 APR 1994						

with COMMENTS
in EPA Letter Dated:
19 APR 1994

READY TO USE WEX-CIDE

GERMICIDAL DETERGENT

CLEANS • DISINFECTS • DEODORIZES

GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL,
FUNGICIDAL AND VIRUCIDAL*

ACTIVE INGREDIENTS:	
o-Phenylphenol	0.028%
o-Benzyl-p-Chlorophenol	0.023%
INERT INGREDIENTS:	99.951%
Includes biodegradable detergents and no phosphates	100.000%

EPA REG. NO. 34810-21 © EPA EST. NO. 34810-MO-1

KEEP OUT OF REACH OF CHILDREN
CAUTION

GERMICIDAL ACTIVITY: (AOAC Use Dilution Confirmation in the presence of 5% serum.) Ps. aeruginosa, S. aureus, S. choleraesuis; including S. aureus (MRSA) (TNSP) (GRM) and (PF) and S. Epidermidis (MRSE). TUBERCULOCIDAL ACTIVITY: (AOAC M. Tuberculosis in the presence of 5% serum at 20°C for 10 minutes) FUNGICIDAL ACTIVITY: (AOAC T. Mentagrophytes in the presence of 5% serum.) VIRUCIDAL ACTIVITY: (Influenza A(Hong Kong), Adenovirus Type 2, Herpes simplex Types 1 & 2, Vaccinia virus in the presence of 5% serum, and HIV-1 (virus associated with AIDS) in the presence of 10% serum on inanimate surfaces.)

STORAGE and DISPOSAL: Store at room temperature. Do not leave containers of this product open during storage. Do not reuse empty container. Rinse empty container thoroughly and discard in trash.



Wexford Labs, Inc.

325 Lettingwell Avenue

Kinwood, IA 50122

Contents: One Quart

Under the Federal Insecticide, Fungicide, and Rodenticide Act, for registered under EPA Reg. 34810-21

READY TO USE WEX-CIDE is used for one-step cleaning, disinfecting and deodorizing hard, non-porous inanimate surfaces in Operating Rooms, Intensive Care, Nurseries, Emergency Rooms, as well as Police, EMS and other emergency response units. It is a broad spectrum kill of both gram positive and gram negative microorganisms.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with the label. Use READY TO USE WEX-CIDE for all disinfection of hard, non-porous inanimate surfaces in care facilities, for working surfaces, surgical instruments, stainless steel cloth or brush, or by immersing equipment to be disinfected. All hard surfaces remain wet for at least 10 minutes before drying with a cloth. Gross filth must be removed before applying cleaning solution. A fresh solution should be prepared at each use and becomes diluted or soiled. DO NOT MIX WITH OTHER CLEANERS.

When disinfecting surgical instruments, they must be pre-cleaned to remove gross soil and rough dried prior to immersion in, or treatment with READY TO USE WEX-CIDE. Rinse hollow instruments before filling with solution.

*READY TO USE WEX-CIDE KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in all settings in which blood or body fluids are soiled. Soiling of inanimate surfaces/objects with blood or body fluids, and inanimate surfaces soiled with blood or body fluids can be associated with the potential for transmission of HIV virus type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus): Wear disposable rubber or latex gloves, gowns, face masks, or eye coverings during cleaning of these surfaces. CONTACT TIME: Leave surfaces wet for 10 minutes necessary to control other common types of viruses and bacteria. CLEANING: Use solution to be used to remove and clean blood and other body fluids from surfaces. For lightly soiled areas wipe off visual soil/clean with solution, and disinfect with fresh solution. For large spills the contaminated area should be flooded with solution and decontaminated with fresh solution. DISPOSAL OF INFECTIOUS MATERIAL: Material should be autoclaved and disposed of according to federal, state, and local disposal.

BEST AVAILABLE COPY