

NOTICE OF PESTICIDE:

(Under the Federal Insecticide, Fungicide,
and Rodenticide Act, as amended)

REGISTRATION
 REREGISTRATION

NAME OF PESTICIDE PRODUCT

ULTRA

NAME AND ADDRESS OF REGISTRANT (Include ZIP code)

Wexford Labs, Inc.
325 Leffingwell
Kirkwood, MO 63122

NOTE: Changes in labeling formula differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above U.S. EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby Registered/Reregistered under the Federal Insecticide, Fungicide, and Rodenticide Act.

A copy of the labeling accepted in connection with this Registration/Reregistration is returned herewith.

Registration is in no way to be construed as an indorsement or approval of this product by this Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA section 3(c)(7)(A) provided that you:

1. Submit/submit all data required for registration/reregistration of your product under FIFRA section 3(c)(7) and section 4 when the Agency requires all registrants of similar products to submit such data.
2. Add the phrase "EPA Registration No. 4810-28" to your label before you release the product for shipment.
3. Submit five (5) copies of your final printed labeling before you release the product for shipment. Refer to the A-79 enclosure for a further description of final printed labeling.

If these conditions are not complied with the registration will be subject to cancellation in accordance with FIFRA section 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records.

Walter J. Pringle
Acting Director, Registration (32)
and Compliance Division
Registration Division, EPA
Washington, D.C. 20460

ATTACHMENT IS APPLICABLE

SIGNATURE OF APPROVING OFFICIAL

DATE

Pringle 10-29



ULTRA

**CONCENTRATED GERMICIDAL DETERGENT
CLEANS --- DISINFECTS --- DEODORIZES
GERMICIDAL, PSEUDOMONACIDAL, TUBERCULOCIDAL
FUNGICIDAL AND VIRUCIDAL***

in soft or hard water (A.O.A.C. 400 ppm Hardness as Calcium Carbonate) and
organic soil (5% serum).

| | |
|--------------------------------------|----------------|
| ACTIVE INGREDIENTS: | |
| ortho-Phenylphenol | 10.10% |
| ortho-Benzyl-para-Chlorophenol | 9.10% |
| INERT INGREDIENTS: | 80.80% |
| | 100.00% |

EPA REGISTRATION NO. 34810-

**KEEP OUT OF REACH OF CHILDREN
DANGER**

See side panel for additional precautionary statements.

Wexford Labs, Inc.

325 Leffingwell Avenue

Kirkwood, Missouri 63122

FOR INDUSTRIAL AND INSTITUTIONAL USE
Contents: One U.S. Gallon (128 ounces)

GERMICIDAL QUALIFICATIONS:

(AOAC Use Dilution Confirmation) in the presence of 400 ppm hardness
and 5% serum at

(AOAC Ps. aeruginosa, S. aureus, S. choleraesuis; including S. aureus
and [PR] and S. Epidemedia [MRSE]).

TUBERCULOCIDAL ACTIVITY TEST:

(AOAC M. Tuberculosis) in the presence of 400 ppm hardness as Ca
5% serum at 20°C for 10 minutes)

FUNGICIDAL TEST:

(AOAC T. mentagrophytes) in the presence of 400 ppm hardness as Ca
5% serum at

VIRUCIDAL ACTIVITY TEST:

(*Influenza A3 [Hong Kong], Adenovirus Type 2, Herpes simplex Types
in the presence of 400 ppm hardness and 5% serum, and HIV-1 [AIDS v
400 ppm hardness as Calcium Carbonate and 10% serum at

PRECAUTIONARY STATEMENTS

HAZARD TO HUMANS AND DOMESTIC ANIMALS. DANGER. Keep
Corrosive. Causes eye damage or skin irritation. Avoid breathing spray
on skin, or on clothing. Wear goggles or face shield and rubber gloves
concentrate. Harmful if swallowed.

STATEMENT OF PRACTICAL TREATMENT: In case of contact immerse
with plenty of water for at least 15 minutes. For eyes call a physician
contaminated clothing before reuse. If swallowed drink promptly a large
alcohol and call a physician.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate
Measures against circulatory shock, respiratory depression and convulsions.

STORAGE AND DISPOSAL

Store at room temperature. Do not leave containers of this product open
thoroughly. Do not reuse empty containers, wrap and put in trash.

Made in U.S.A.

**ACCEPTED
with COMMENTS
in EPA Letter Dated:**

OCT 6 0 1988

Under the Federal Insecticide,
Fungicide, and Rodenticide Act
as amended, for the pesticide
registration under EPA Reg. No.

34810-28

ULTRA

CONCENTRATED GERMICIDAL DETERGENT

ULTRA is used for one-step cleaning, disinfecting and deodorizing floors, walls, and other similar hard surfaces. Especially suited for use in Operating Rooms, Intensive Care, Nurseries, Emergency Areas, Patient Rooms and Rest Rooms. ULTRA is a synthetic germicidal detergent with broad spectrum kill of both gram positive and gram negative micro-organisms. It is an excellent cleaner, contains biodegradable detergents and is phosphate free. Do not mix with other cleaners.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. For all disinfection of hard surfaces in hospitals, health-care facilities, schools, and nursing homes, use ULTRA at a dilution of 1/3 ounce per gallon of water for walls, floors, bathroom fixtures, surgical instruments, stainless steel, chrome, glass and vinyl. ULTRA solution may be applied with a cloth, mop or brush, or by immersing equipment to be disinfected. Wet all surfaces thoroughly. The surface must remain wet for at least 10 minutes before drying with a cloth, or a mop in the case of floors. Gross filth and heavy soil must be removed before applying cleaning solution. Discard solution when it becomes dirty and replace with fresh solution. ULTRA is suitable for cleaning floors and walls in all departments of food processing plants at 1:384 use dilution. **DISINFECTING SURGICAL INSTRUMENTS:** These instruments must be thoroughly pre-cleaned to remove blood and serous exudates, rinsed and roughed dried prior to immersion in, or treatment with ULTRA.

***ULTRA KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1(AIDS virus). Wear rubber or latex gloves when performing manual cleaning of these surfaces. **CONTACT TIME:** Leave surfaces wet for 10 minute.

CLEANING PROCEDURE: ULTRA solution may be used to remove blood and other body fluids from surfaces and objects. For lightly soiled areas wipe off visual soil with ULTRA solution, and disinfect surface with a second application of the ULTRA solution. For large spills the contaminated area should be flooded with the ULTRA solution before cleaning, then decontaminated with fresh ULTRA solution.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other body fluids should be autoclaved and disposed of according to Federal, State and local regulations for infectious waste disposal.

Wexford Labs, Inc.

EPA ESTABLISHMENT NUMBER 34810-MO.1