

34810-37

9/14/2012

1/5



U.S. ENVIRONMENTAL PROTECTION AGENCY

Office of Pesticide Programs  
Antimicrobials Division (7510P)  
1200 Pennsylvania Avenue NW  
Washington, D.C. 20460

EPA Reg.

Date of Issuance:

Number:

9/14/2012

34810-37

Term of Issuance:

Unconditional

Name of Pesticide Product:

Wexford Disinfectant Wipes

NOTICE OF PESTICIDE:

Registration  
 Reregistration

(under FIFRA, as amended)

Name and Address of Registrant (include ZIP Code):

WEXFORD LABS, INC  
325 Leffingwell Avenue  
Kirkwood, MO 63122

Note: Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Antimicrobials Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby registered/reregistered under the Federal Insecticide, Fungicide and Rodenticide Act.

Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product (OPP Decision No. D-459534) is unconditionally registered in accordance with FIFRA sec 3(c)(5) provided that you:

1. Submit and/or cite all data required for registration of your product under FIFRA sec. 3(c)(5) when the Agency requires all registrants of similar products to submit such data; and submit acceptable responses required for re-registration of your product under FIFRA section 4.

2. Make the labeling changes listed below before you release the product for shipment:

a. Revise the "EPA Registration Symbol to read, "EPA Reg. No. 34810-37".

Signature of Approving Official:

Jacqueline Campbell, PM Team 34  
Regulatory Management Branch II  
Antimicrobials Division (7510P)

Date:

SEP 14 2012

Secretary


b.

3. Submit three (3) copies of your final printed labeling before distributing or selling the product bearing the revised labeling.

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records. Should you have any questions regarding this letter, please contact me by telephone (703) 308-6416 or by email at [Campbell-mcfarlane.jacqueline@epa.gov](mailto:Campbell-mcfarlane.jacqueline@epa.gov) or Stacey Grigsby by telephone at (703) 305-6440 or by email at [grigsby.stacey@epa.gov](mailto:grigsby.stacey@epa.gov). When you are submitting information or data in response to this letter, send a copy of this letter to accompany the submission in order to facilitate processing.

Sincerely,



Jacqueline Campbell  
Product Manager Team (34)  
Regulatory Management Branch II  
Antimicrobials Division (7510P)

Enclosure: Stamped Label

•  
•  
•  
•  
•  
•  
•  
•


L

L

315

**[Wexford Disinfectant Wipes]  
Cleans, Disinfects, Deodorizes  
Germicidal, Fungicidal, and Virucidal\***

Active Ingredient: Citric Acid .....0.6%  
Inert Ingredients: .....99.4%  
Includes biodegradable detergents  
Total: .....100.0%



EPA Reg. No. 34810-xx  
EPA Est. No. (insert Est. No. here)  
NET CONTENTS (insert container size here)  
[Country of origin (insert country)]  
[Manufactured in (insert country)]  
[Barcode]

**FLORAL / CITRUS / FRESH / CLEAN / LEMON / PINE / MINT [scent] [fragrance] [odor] [scented]**

**KEEP OUT OF REACH OF CHILDREN  
CAUTION**

**Effective in 5% serum against:**

**GERMICIDAL:** *Pseudomonas aeruginosa*<sup>1</sup> [ATCC 15442], *Staphylococcus aureus*<sup>1</sup> [ATCC 6538], *Salmonella enterica*<sup>1</sup> (formerly *Salmonella Choleraesuis*) [ATCC 10708], *Listeria monocytogenes*<sup>1</sup> [ATCC 19117], *Escherichia Coli O157:H7*<sup>1</sup> [ATCC 43888], *Staphylococcus aureus - MRSA*<sup>1</sup> [ATCC 33592], *Staphylococcus epidermis - MRSE*<sup>1</sup> [ATCC 51625], *Enterococcus faecalis Vancomycin Resistant*<sup>1</sup> (VRE) [ATCC 51575], *Acinetobacter baumannii*<sup>2</sup> [ATCC 19606] **TUBERCULOCIDAL:** *Mycobacterium bovis* (BCG)<sup>2</sup> at 20°C **VIrucIDAL:** Poliovirus<sup>2</sup> Type 1, Canine Parvovirus<sup>2</sup> (CPV), Respiratory Syncytial Virus<sup>1</sup>, Rotavirus<sup>1</sup>, Human Immunodeficiency Virus<sup>1</sup> Type 1 (Strain HTLV-III<sub>B</sub>), Herpes Simplex Virus<sup>1</sup> Type 1, Herpes Simplex Virus<sup>1</sup> Type 2 (Strain G), Vaccinia Virus<sup>1</sup>, Influenza A Virus<sup>1</sup> (Strain Hong Kong), kills Pandemic 2009 H1N1 influenza A virus<sup>1</sup>, Adenovirus<sup>1</sup> Type 2, Duck Hepatitis B Virus (DHBV)<sup>2</sup>, Hepatitis C (HCV)<sup>2</sup> (utilizing Bovine Viral Diarrhea Virus as a surrogate for Human Hepatitis C), Feline Calicivirus<sup>1</sup> (surrogate for Norovirus) **FUNGICIDAL:** *Trichophyton mentagrophytes*<sup>2</sup> [ATCC 9533]

- 1 - Five minutes Contact Time
- 2 - Ten minutes Contact Time

**PRECAUTIONARY STATEMENTS  
HAZARD TO HUMANS AND DOMESTIC ANIMALS**

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

**FIRST AID**

**IF IN EYES:** Hold eyelids open and flush with a steady, gentle stream of water for 15 minutes. Remove contact lenses, if present, after the first five minutes, then continue rinsing. Call a Poison Control Center or doctor for further treatment advice.

**Wexford Labs, Inc. 325 Leffingwell Avenue Kirkwood, Missouri 63122**

NET CONTENTS: X lbs, X oz.

X Pre-moistened Towelettes (X)" x (X)"  
(e.g. 100 Pre- moistened Towelettes 7"x 8")

Revised 01/05/12

[Note to Reviewer: Language in [ ] is considered optional or interchangeable.]

**DIRECTIONS FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Product is ready to use.

**DISPENSER DIRECTIONS:** To start feed, remove large cover. Pull up the towelette corner from the center of the towelette roll, twist it into a point and thread it through the hole in the container cover. Pull through about one inch. Replace cover. Pull out first towelette and snap off at a 90° angle. Remaining towelettes feed automatically ready for their next use. When through using, keep small center cap closed to prevent moisture loss.

**TO CLEAN/DEODORIZE:** Unfold premoistened cloth and wipe area to be cleaned for 30 seconds or until clean.

**TO DISINFECT:** Unfold premoistened cloth and wipe area to be cleaned and disinfected. Use enough wipes for treated surface to remain wet for 10 minutes. Wipe or let air dry. For gross filth or heavy soil, use a wipe to clean and then a second wipe to disinfect the surface. **DO NOT MIX WITH OTHER CLEANERS.**

**[AREAS OF USE:]**

[Homes [households]

[Kitchens]

[Hospitals]

[Health Care Facilities]

[Medical Offices]

[Dental offices]

[Day care centers]

[Nurseries]

[Schools]

[Nursing homes]

[Laboratories or other small animal facilities]

[Bathrooms]

[Institutions]

[Institutional facilities]

[Veterinary Clinics]

[Kennels]

[Grooming Facilities]

[Zoos]

[Pet Shops]

**[TYPES OF SURFACES:]**

[On washable hard, nonporous surfaces of:]

[Bathroom fixtures]

[Floors]

[Walls]

[Vinyl]

[Stainless Steel]

[Chrome]

[CuVerro™ Antimicrobial Copper Surfaces]

[Glass]

**DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION:**

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

**\*CLEAN-CIDE WIPES KILLS HBV, HCV and HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Hepatitis B Virus, Hepatitis C Virus and Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HBV, HCV AND HIV-1 (AIDS virus):**

**Personal Protection:** When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask, and eye coverings.

**Cleaning Procedure:** All blood and body fluids must be cleaned from surfaces before disinfection by the germicidal cloth. Use one germicidal cloth to remove gross soil.

**Contact Time:** Use a second germicidal cloth to thoroughly wet surface. Allow to remain wet 5 minutes for HIV-1, and 10 minutes for other organisms listed on the label, let air dry.

**Disposal of Infectious Material:** Dispose of used towelette in accordance with local regulations for infectious waste disposal.

**Veterinary Practice / Animal Care / Animal Laboratory / Zoos / Pet Shop / Kennels Disinfection Directions:**


For cleaning and disinfecting hard nonporous surfaces: equipment used for animal food or water, utensils, instruments, cages, kennels, stables, catteries, etc. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings and manure from floors, walls and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Unfold premoistened cloth and wipe area to be cleaned and disinfected. Use enough wipes for treated surface to remain visibly wet for 10 minutes. Wipe or allow to air dry. Thoroughly scrub all treated surfaces that come into contact with food, including equipment used for feeding or watering, with soap or detergent and rinse with potable water before reuse. Do not house animals or employ equipment until treatment has been absorbed, set or dried.

**STORAGE AND DISPOSAL:** Do not leave containers open during storage. No special handling is required. **Towelette:** Do not reuse towelette. Dispose of used towelette in trash. **Dispenser:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available.

**Made in the USA**

MSDS #: 3130-XX

Reorder #: 3130-XX



...  
...  
... of USA

... to ...  
... on the ...  
... se of used